



LAB	CODE	TEST NAME	CPT4 CODE
ACL	BDGLUC	1,3-BETA-D-GLUCAN Container :One gold gel 3.5 mL (also acceptable: plain red 4.0 mL) Stability :Refrigerated: 2 Weeks / Frozen: 2 Weeks	87449
ACL	CRJK	14-3-3 PROTEIN CSF Container :Spinal fluid (CSF) Stability :Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 2 Months	0035U 84182 86317
ACL	AFBPCR	16S rRNA - AFB Container :Tissue: Fresh (optimal) or paraffin embedded (see NOTE below) Fluid: (non-blood) from normally sterile sites collected in Vacutainer tubes without preservative or Lavender top (EDTA) tube Respiratory specimens can be tested for Mycobacterium tuberculosis complex and/or Mycobacterium avium complex only.	87551 87556
ACL	BCTPCR	16S rRNA - BACTERIAL Container :Tissue: Fresh (optimal) or paraffin embedded (see NOTE below) Fluid: (non-blood) from normally sterile sites collected in Vacutainer tubes without preservative or Lavender top (EDTA) tube Respiratory specimens can be tested for Mycobacterium tuberculosis complex and/or Mycobacterium avium complex only.	87801
ACL	FUNPCR	16S rRNA - FUNGAL Container :Tissue: Fresh (optimal) or paraffin embedded (see NOTE below) Fluid: (non-blood) from normally sterile sites collected in Vacutainer tubes without preservative or Lavender top (EDTA) tube Respiratory specimens can be tested for Mycobacterium tuberculosis complex and/or Mycobacterium avium complex only.	87801
ACL	17HYDR	17 OH PROGESTERONE Container :One gold gel 3.5 mL (also acceptable: plain red 4.0 mL or pink (EDTA) 6.0 mL or green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 6 Months	83498
ACL	17OHPG	17-HYDROXYPREGNENOLONE Container :One plain red 4.0 mL (also acceptable: gold gel 5.0 mL and Two lavender (EDTA) 3.0 mL and green (lithium heparin no gel) 6.0 mL and green (sodium heparin no gel) 4.0 mL and pink (K2EDTA) 6.0 mL) Stability :Frozen: 6 months	84143
ACL	21OHAB	21-HYDROXYLASE AB Container :One plain red 4.0 mL (also acceptable: One gold gel 3.5 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Month Refrigerated: 1 Week After separation from cells Frozen: 1 Month After separation from cells	83516
ACL	5NUCL	5' Nucleotidase Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 4 Hours / Refrigerated: 1 Week / Frozen: 2 Weeks Avoid rpeated freeze/thaw cycles	83915
ACL	5ADTSR	5A-DIHYDROTESTOSTERO Container :One plain red 4.0 mL (ALSO ACCEPTABLE: ONE GOLD GEL 5.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 5 days / Frozen: 6 Months	82642
PHC	ABGC	ABG (Arterial Blood Gases), Capillary, Routine Stability :Room temperature: 30 minutes Refrigerated: transport on ice if specimen will not be in lab within 30 minutes of draw Do Not Freeze	82803
PHC	ABGPLS	ABG (Arterial Blood Gases), Glucose, Hematocrit, Hemoglobin, Potassium, Sodium, and Oximetry Stability :Room temperature: 30 minutes Refrigerated: transport on ice if specimen will not be in lab within 30 minutes of draw Do Not Freeze	82805 83050 82947 84295 84132 82375 85014
PHC	ABG	ABG (Arterial Blood Gases), Routine Stability :Room temperature: 30 minutes Refrigerated: transport on ice if specimen will not be in lab within 30 minutes of draw Do Not Freeze	82803
PHC	TYPSRN	ABO & Rh Type and Antibody Screen, Blood Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	86850 86900 86901
PHC	TYPE	ABO & Rh Type, Blood Stability :Specimen stability:	86900 86901



LAB	CODE	TEST NAME	CPT4 CODE
		Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	
PHC	TYPE2	ABO/RH FOR BLOOD BANK USE ONLY, TYPE2 Stability :Ambient: 8 hours, Refrigerated: 4 Days, Frozen: Unacceptable	
PHC	ACETA	Acetaminophen, Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 14 days Frozen: 45 days	80329
PHC	ACET	Acetone, Plasma or Serum, Qualitative	
ACL	ACRB	ACETYLCHOL BINDING AB Container :One gold gel 3.5 mL Stability :Ambient: 2 weeks / Refrigerated: 2 Weeks / Frozen: 30 days	83519
ACL	ARBA	Acetylcholine Receptor Blocking Antibody Container :One gold gel 3.5 mL blood Stability :Ambient: 48 Hours after separation from cells Refrigerated: 2 Weeks after separation from cells Frozen: 1 Year after separation from cells (Avoid repeated freeze/thaw cycles)	83516
ACL	ACETREC	Acetylcholine Receptor Modulating Antibody Container :One gold gel 3.5 mL blood Stability :Ambient: 48 Hours after separation from cells Refrigerated: 2 Weeks after separation from cells Frozen: 1 Year after separation from cells (Avoid repeated freeze/thaw cycles)	83516
ACL	ACTHB	ACT Container :One lavender (EDTA) 3.0 mL Collect into a chilled EDTA lavender tube Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	82024
ACL	APCRAR	ACT PROT C RESISTANC Container :Two light blue (3.2% sodium citrate) 2.7 mL and: lavender (EDTA) 3.0 mL Stability :Ambient: Unacceptable for citrated plasma, 72 Hours for whole blood / Refrigerated: Unacceptable for citrated plasma, 14 Days for whole blood / Frozen: 30 Days for citrated plasma, Unacceptable for whole blood	85307
ACL	ADALIM	Adalimumab Activity and Neutralizing Antibody Container :One gold gel 5.0 mL blood Stability :AMBIENT: 48 hours/ REFRIGERATED:4 weeks / FROZEN: 1 year	80145 82397
PHC	ADAMTS	ADAM TS13 Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 2 weeks	85397
ACL	ADENO	ADENOVIRUS ANTIBODY Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 7 Days / Refrigerated: 14 Days / Frozen: 30 Days	86603
PHC	ADRSTI	Adrenal Stimulation Test Stability :Room temperature: 8 hours Refrigerated: 5 days Frozen: 30 days	82533 x3
PHC	AEPN	Adult Environmental Panel Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003 x25 82785
PHC	AFOODP	Adult Food Panel Stability :Ambient: 48 Hours / Refrigerated: 1 Weeks / Frozen: 1 Year	86003 x21 82785
ACL	AEPLC	Aeromonas and Plesiomonas Stool Culture Container :STOOL IN LIQUID C&S (CARY - BLAIR) TRANSPORT MEDIA RAW STOOL IN CLEAN LEAKPROOF CONTAINER. Stability :AMBIENT: 1 HOUR, NO PRESERVATIVE. YIELDS ARE OPTIMIZED IF STOOL IS IMMEDIATELY PLACED IN C&S (CARY-BLAIR) TRANSPORT MEDIA. AMBIENT: 4 DAYS WITH PRESERVATIVE / REFRIGERATED: UNACCEPTABLE / FROZEN: UNACCEPTABLE	87046
ACL	TXAAFP	AFP AMNIOTIC FLUID	82106
ACL	AFPMA	AFP OPEN NTD SCREEN Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	82105
ACL	AFPTM	AFP TUMOR MARKER Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 14 Days	82105



LAB	CODE	TEST NAME	CPT4 CODE
PHC	FLALB	Albumin, Body Fluid	82040
PHC	ALB	Albumin, Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 2 days Frozen: 6 months	82040
PHC	ETOH	Alcohol, Plasma or Serum Stability :Room temperature: 2 days Refrigerated: 3 days	80320
ACL	ALDO	ALDOLASE Container :One gold gel 3.5 mL Stability :Ambient 8 hours/Refrigerated: 5 Days / Frozen: 6 Months	82085
ACL	ALDOS	ALDOSTERONE Container :One gold gel 3.5 MI (also acceptable: plain red 4.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 5 Days / Frozen: 1 Month	82088
ACL	UALDOS	ALDOSTERONE,UR 24HR Container :One Urine Stability :Ambient: 4 days / 7 days / Frozen: 90 days	82088
ACL	ALRENR	Aldosterone/Renin Activity Ratio Container :One gold gel 5.0 mL blood and Two lavender (EDTA) 3.0 mL blood Stability :Ambient: Serum 8 hours, plasma 6 hours. / Refrigerated: Serum: 5 days, Plasma: Unacceptable / Frozen: 1 Month	82088 84244
ACL	CEDAR	ALGRN:CEDAR/RED TREE Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	ALKI	ALK PHOS ISOENZYMES Container :Two gold gel 3.5 mL Stability :Ambient: 1 Days / Refrigerated:7 Days / Frozen: 14 Days	84075 84080
PHC	ALP	Alkaline Phosphatase, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 2 months	84075
ACL	BCIGE	ALLERGEN BOTRYTIS CI Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	POPYR	ALLERGEN POPPY SEED Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	ANNATT	Allergen, Annatto Seed IgE Container :Gold gel 3.5 mL or Plain red 4.0 mL Stability :Ambient: 1 Month/ Refrigerated: 1 Month / Frozen: 1 Year	86003
ACL	CNDIDA	Allergen, Candida albicans IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 30 Days / Frozen: 30 Days	86003
ACL	COWMILK	Allergen, Cow Milk Components IgE Container :One gold gel 5.0 mL blood Also acceptable: One light green (lithium heparin) gel 5.0 mL or One lavender EDTA 3.0 mL Stability :Ambient: 1 Day/Refrigerated: 30 Days/Frozen: 30 Days	86008 x3
ACL	DPTER	Allergen, Dust Mites, D. Pteronyssinus IgE ImmunoCAP Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	OVALB	Allergen, Food, Ovalbumin IgE Container :One gold gel 3.5 mL Multiple specimen tubes should be avoided. Stability :Ambient: 48 Hours - after separation from cells Refrigerated: 2 Weeks - after separation from cells Frozen: 1 Year - after separation from cells	86008
ACL	EPCOC	Allergen, Fungi & Molds, Epicoccum purpurascens IgE ImmunoCAP Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	CHERB	Allergen, Fungi & Molds, Hormodendrum IgE ImmunoCAP Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MUCOR	Allergen, Fungi & Molds, Mucor racemosus IgE ImmunoCAP Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PNOTA	Allergen, Fungi & Molds, Penicillium chrysogenum IgE ImmunoCAP Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	AFUMG	Allergen, Fungi and Molds, A. fumigatus IgE ImmunoCAP Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	DFARN	Allergen, Mites, Dermatophagoides farinae IgE ImmunoCAP () Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ALMOND	ALLERGEN: ALMOND IgE Container :One gold gel 3.5 mL. Also acceptable: One: Green gel 5.0 mL or lavender 3.0 mL Stability :Ambient: 1 Day / Refrigerated: 30 Days / Frozen: 30 Days	86003
ACL	ASPEN	ALLERGEN: ASPEN IgE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 1 Month / Refrigerated: 1 Month / Frozen: 1 Year	86003
ACL	BARLEY	ALLERGEN: BARLEY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BEEF	ALLERGEN: BEEF IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CACAO	ALLERGEN: CACAO IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CARRO	ALLERGEN: CARROT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CASHW	ALLERGEN: CASHEW IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CELERY	ALLERGEN: CELERY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CLAMS	ALLERGEN: CLAM IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 30 Days / Frozen: 30 Days	86003
ACL	COFFEE	ALLERGEN: COFFEE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CORNS	ALLERGEN: CORN IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CRABA	ALLERGEN: CRAB IgE	86003



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	
ACL	GARLIC	ALLERGEN: GARLIC IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	GLUTN	ALLERGEN: GLUTEN IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	GRAPES	ALLERGEN: GRAPE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	KIWI	ALLERGEN: KIWI IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LEMON	ALLERGEN: LEMON IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LENTIL	ALLERGEN: LENTIL IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LIME	ALLERGEN: LIME IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MLTIGE	ALLERGEN: MALT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MANGO	ALLERGEN: MANGO IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MELON	ALLERGEN: MELON IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	NETTLE	ALLERGEN: NETTLE IgE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	OATIGE	ALLERGEN: OAT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ONIONS	ALLERGEN: ONION IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ORANG	ALLERGEN: ORANGE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PEASB	ALLERGEN: PEA IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	PEACHS	ALLERGEN: PEACH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PENUT	ALLERGEN: PEANUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PEARS	ALLERGEN: PEAR IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PERCHS	ALLERGEN: PERCH IgE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 1 Month / Refrigerated: 1 Month / Frozen: 1 Year	86003
ACL	PORK	ALLERGEN: PORK IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	RICE	ALLERGEN: RICE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	RYEA	ALLERGEN: RYE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SHRMPA	ALLERGEN: SHRIMP IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SQUID	ALLERGEN: SQUID IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	TROUT	ALLERGEN: TROUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	TUNAS	ALLERGEN: TUNA IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	TURKEY	ALLERGEN: TURKEY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WHEATS	ALLERGEN: WHEAT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WHEY	ALLERGEN: WHEY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BYEAST	ALLERGEN: YEAST IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	AVOCAD	ALLERGEN:AVOCADO IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	GOLDEN	ALLERGEN:GOLDENROD IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	HALIBT	ALLERGEN:HALIBUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MUSTD	ALLERGEN:MUSTARD IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ANTRY	ALPHA 1 ANTITRYPSIN Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Refrigerated: 3 Days / Frozen: 14 Days	82103
ACL	ALPSUB	ALPHA SUBUNIT Container :One plain red 4.0 mL Stability :Refrigerated: 7 Days / Frozen: 3 Months	83520
ACL	A1AT STL	Alpha-1-Antitrypsin Stool Container :One Stool, random Stability :Ambient: 24 Hours Refrigerated: 7 Days Frozen: 3 Months	82103
ACL	FANTR	ALRGN FIRE ANT Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	TRMULR	ALRGN RED MULBERRY Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	WYDKR	ALRGN YELLOW DOCK Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	TILAPI	ALRGN: TILAPIA IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ASPNI	ALRGN: A. NIGER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ATENS	ALRGN: A. TENUIS IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	APPLES	ALRGN: APPLE IgE Container :One gold gel 5.0 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	BNANA	ALRGN: BANANA IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BASIL	ALRGN: BASIL IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BLAPEP	ALRGN: BLACK PEPPER Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BLUMUS	ALRGN: BLUE MUSSEL Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL)	86003



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: 24 Hours / Refrigerated: 30 Days / Frozen: 30 Days	
ACL	BLUBRY	ALRGN: BLUEBERRY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 30 Days / Frozen: 30 Days	86003
ACL	BOXELD	ALRGN: BOX ELDER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BROCCI	ALRGN: BROCCOLI IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BUKWHT	ALRGN: BUCKWHEAT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CABBAG	ALRGN: CABBAGE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CTFISH	ALRGN: CATFISH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MLDCHS	ALRGN: CHEESE MOLD Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CHERRY	ALRGN: CHERRY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CHESNT	ALRGN: CHESTNUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CHCKF	ALRGN: CHICK FEATHER Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CHICK	ALRGN: CHICKEN IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CPEA	ALRGN: CHICKPEA IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CINAMN	ALRGN: CINNAMON IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	COKBUR	ALRGN: COCKLEBUR IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ROACH	ALRGN: COCKROACH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	COCONT	ALRGN: COCONUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL)	86003



LAB	CODE	TEST NAME	CPT4 CODE
		mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	
ACL	CODIGE	ALRGN: CODFISH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MILKC	ALRGN: COW MILK IgE Container :One blood gold gel 3.5 mL (also acceptable: One blood plain red 4.0 mL and One blood light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	CRANBY	ALRGN: CRANBERRY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CUCUM	ALRGN: CUCUMBER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	DUCKFR	ALRGN: DUCK FEATHERS Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	EGGWA	ALRGN: EGG WHITE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	EGGYOL	ALRGN: EGG YOLK IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ELMTR	ALRGN: ELM TREE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	ENGPLT	ALRGN: ENG PLANTAIN Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	FIRBSH	ALRGN: FIREBUSH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	FLOUND	ALRGN: FLOUNDER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	GRAGWD	ALRGN: GIANT RAGWEED Container :One gold gel 5.0 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	DUSTG	ALRGN: H DUST GREER Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	DUSTS	ALRGN: H DUST STIER Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	HELMN	ALRGN: H HALODES IgE Container :One gold gel 5.0 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	HAZNUT	ALRGN: HAZELNUT IgE	86003



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	
ACL	HIKPEC	ALRGN: HICKORY/PECAN Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	HONBEA	ALRGN: HONEYBEE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	HORSEA	ALRGN: HORSE DANDER Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	INSCTP	ALRGN: INSECT PANEL Panel IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL or lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 days	86003 x5
ACL	JONSON	ALRGN: JOHNSON GRASS Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	KDBEAN	ALRGN: KID BEAN IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LAMBQU	ALRGN: LAMBS QUARTER Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LATEXB	ALRGN: LATEX IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LETUCE	ALRGN: LETTUCE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LINSED	ALRGN: LINSEED IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	LBSTER	ALRGN: LOBSTER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MACADA	ALRGN: MACADAMIA IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MACKEL	ALRGN: MACKEREL IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MFESCU	ALRGN: MEADOW FESCUE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MOUEPI	ALRGN: MOUSE Epithelium IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 30 Days / Frozen: 30 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	MUGWRT	ALRGN: MUGWORT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MULBRY	ALRGN: MULBERRY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MUSHRM	ALRGN: MUSHROOM IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	OAKA	ALRGN: OAK TREE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	OYSIGE	ALRGN: OYSTER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PARSLY	ALRGN: PARSLEY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PECAN	ALRGN: PECAN NUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PENG	ALRGN: PENICILLIN G Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	PENV	ALRGN: PENICILLIN V Container :One blood gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	PENCLS	ALRGN: PENICILLINS Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003 x2
ACL	PGWEED	ALRGN: PIGWEED IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PINENT	ALRGN: PINE NUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PINAPP	ALRGN: PINEAPPLE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PISTAC	ALRGN: PISTACHIO IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	POTATO	ALRGN: POTATO IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	RHIZOP	ALRGN: R NIGRICANS Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	RASPB	ALRGN: RASPBERRY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	REDDY	ALRGN: RED DYE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	REDTOP	ALRGN: REDTOP GRASS Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
PHC	RESP7	ALRGN: RESP REGION 7 Container :Two gold gel 6.0 mL (also acceptable: Two plain red and Three light green (lithium heparin gel)) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003 x23 82285
ACL	RYEGRS	ALRGN: RYE GRASS IgE Container :One gold gel 3.5 mL (also acceptable: One light green (lithium heparin gel) 3.0 mL and One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	SALMN	ALRGN: SALMON IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SCALOP	ALRGN: SCALLOP IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SESAME	ALRGN: SESAME SEED Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SORR	ALRGN: SHEEP SORREL Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SRAGWD	ALRGN: SHORT RAGWEED Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SOLE	ALRGN: SOLE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SOYIGE	ALRGN: SOYBEAN IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	TOMATO	ALRGN: TOMATO IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WALNT	ALRGN: WALNUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WNUTTR	ALRGN: WALNUT TREE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WASHT	ALRGN: WHITE ASH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	WHIFSH	ALRGN: WHITEFISH IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WILLO	ALRGN: WILLOW TR IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 1 Day / Refrigerated: 30 Days / Frozen: 30 Days	86003
ACL	BASSR	ALRGN:BASS Container :One blood gold gel 3.5 mL (also acceptable: One blood plain red 4.0 mL) Stability :Ambient: 1 Month / Refrigerated: 1 Month / Frozen: 1 Year	86003
ACL	BECH	ALRGN:BEECH TREE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BRMUDA	ALRGN:BERMUDA GR IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BIRCHT	ALRGN:BIRCH TREE IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	BZ Nut	ALRGN:BRAZIL NUT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CATDND	ALRGN:CAT DANDER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CACREM	ALRGN:CEPHALOSPORIUM Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	CHEDCH	ALRGN:CHEDDAR CHEESE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	COTTWD	ALRGN:COTTONWOOD TR Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	COWHDR	ALRGN:COW HAIR&DANDE Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	K9DND	ALRGN:DOG DANDER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	FUSAR	ALRGN:F PROLIFERATUM Container :One gold gel 5.0 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	FLRAGR	ALRGN:FALSE RAGWEED Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	ALLFM	ALRGN:FEATHER MIX Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	GOOSE	ALRGN:GOOSE FEATHERS Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL)	86003



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	
ACL	GRAFRU	ALRGN:GRAPEFRUIT IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	GRBEAN	ALRGN:GREEN BEAN IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	GPIG	ALRGN:GUINEA PIG IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	JUNEGR	ALRGN:JUNE GRASS IgE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One light green (lithium heparin gel) 3.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	RMELDR	ALRGN:MARSHELDER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	MTJUNI	ALRGN:MT JUNIPER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	NVYBNR	ALRGN:NAVY BEAN Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	ORCH	ALRGN:ORCHARD GR IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PAWASP	ALRGN:PAPER WASP IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	PHOBET	ALRGN:PHOMA BETAE Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	RABEPI	ALRGN:RABBIT EPI IgE Container :One gold gel 3.5 mL (also acceptable: One light green (lithium heparin gel) 3.0 mL and One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86003
ACL	THISTL	ALRGN:RUSSIAN THISTL Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WORMR	ALRGN:SAGEBRSH/WORM Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	STRBRY	ALRGN:STRAWBERRY IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SUNFLR	ALRGN:SUNFLOWER IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	SWVERN	ALRGN:SW VERNAL GRS Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003



LAB	CODE	TEST NAME	CPT4 CODE
ACL	SYCTRE	ALRGN:SYCAMORE TREE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	TIMTHY	ALRGN:TIMOTHY GR IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WATMEL	ALRGN:WATERMELON IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	WSRAGR	ALRGN:WESTERN RAGWEE Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86003
ACL	WHORNT	ALRGN:WH FACE HORNET Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	YLHORN	ALRGN:YEL HORNET IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
ACL	YELJAK	ALRGN:YEL JACKET IgE Container :One gold gel 3.5 mL (also acceptable: light green (lithium heparin gel) 3.0 mL and lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86003
PHC	ALT	ALT/SGPT (Alanine Aminotransferase), Plasma or Serum Stability :Room temperature: 3 days Refrigerated: 7 days Frozen: 60 days at -40o C	84460
ACL	ALUM	ALUMINUM, SERUM Container :One royal blue (no additive) 6.0 mL Stability :Ambient: Indefinitely / Refrigerated: Indefinitely / Frozen: Indefinitely	82108
ACL	AMIKA	AMIKACIN, PEAK Container :One plain red 4.0 mL Stability :Ambient: 4 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80150
ACL	AMIKAT	AMIKACIN, TROUGH Container :One plain red 4.0 mL Stability :Ambient: 4 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80150
ACL	AMINOP	AMINO ACIDS QUANTITATIVE PLASMA Container :One green (sodium or lithium heparin no gel) 6.0 mL blood Stability :Refrigerated: 24 hours/Frozen: 1 Month	82139
ACL	UAMA	AMINO ACIDSQUANTITATIVE RANDOM URINE Container :One urine - random Stability :Refrigerated: 24 hours; Frozen: 1 month	82139
ACL	UALA	Aminolevulinic Acid (ALA), Urine Stability :Refrigerated: 4 Days / Frozen: 1 Month	82135
ACL	AMIOD	AMIODARONE Container :One plain red 4.0 mL (also acceptable: lavender (EDTA) 3.0 mL and pink (K2EDTA) 6.0 mL) Stability :Ambient: 1 month / Refrigerated: 6 weeks / Frozen: 6 Week	80299
ACL	AMTRP	AMITRIPTYLINE & NORT Container :One plain red 4.0 mL Separate serum from cells within 2 hours of collection Stability :Ambient: 5 Days / Refrigerated: 2 Weeks / Frozen: 6 Months	80335
PHC	AMON	Ammonia, Blood Stability :Room temperature: not accptable Refrigerated: 2 hours if spun and decanted Frozen: 3 weeks	82140
ACL	AMPCNF	Amphetamine Confirmation Stability :Ambient: 1 Week / Refrigerated: 1 month / Frozen: 3 Years	80325 80359
PHC	AMPHE	Amphetamines, Urine	80307



LAB	CODE	TEST NAME	CPT4 CODE
PHC	FLAMY	Amylase, Body Fluid Stability :Room temperature: 2 days Refrigerated: 7 days Frozen: 3 weeks	82150
PHC	AMY	Amylase, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 1 year	82150
ACL	UAMYT	AMYLASE, TIMED URINE Container :One Urine Stability :Ambient: 7 Days / Refrigerated: 7 Days / Frozen: 14 Days	82150
ACL	UAMYR	AMYLASE, URINE Container :One Urine, random Stability :Ambient: 7 Days / Refrigerated: 7 Days / Frozen: 14 Days	82150
ACL	ANACP	ANA COMPREHENSIVE PA Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86225 86235 x10
ACL	ANABL	ANA SCREEN WITH REFLEX Container :One gold gel 5.0 mL (also acceptable: plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86038
ACL	ANA	ANA WITH REFLEX IFA Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86038
ACL	ANDRO	ANDROSTENEDIONE Container :One gold gel 5.0 mL Also Acceptable: One green (sodium or lithium heparin) 6.0 mL Stability :Ambient: 24 hours / Refrigerated: 7 Days / Frozen: 6 Months	82157
ACL	CACE	ANG TENS CON ENZ,CSF Container :One Spinal fluid (CSF) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 6 Months	82164
ACL	ANGCON	ANGIOTENS. CONV ENZ Container :One gold gel 3.5 mL Stability :Ambient: 1 week / Refrigerated: 1 Week / Frozen: 6 Months	82164
ACL	Viracor: 2105	ANTI IgE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Weeks / Frozen: 2 Years	83516
ACL	MITOAB	ANTI MITOCHONDRIA AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Frozen: 14 Days	83516
ACL	AMULH	ANTI MULLERIAN HORMONE Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL and One green (lithium heparin no gel) 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 1 Week / Frozen: 3 Weeks	83520
ACL	THYAU	ANTI THYROID AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	86376 86800
ACL	CENAB	ANTI-CENTROMERE AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	ADNA	ANTI-dsDNA Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86225
ACL	ANCA Titer	ANTI-NEUTROPHIL CYTOPLASMIC AB IgG Container :One gold gel 3.5 mL Stability :Ambient: 48 hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86255
PHC	HEPA	Anti-Xa Level, Heparin, Plasma Stability :Specimen stability is 2 hours at room temperature pored off, 2 weeks at -20 degrees C, 12 months at -70 degrees C.	85520
PHC	ABIDB	Antibody Identification, Blood Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	86870
PHC	ABSNB	Antibody Screen, Blood	86850



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	
ACL	TTR	ANTIBODY TITER Container :One lavender (EDTA) 6.0 mL and One plain red 6.0 mL Stability :Ambient: 1 Day / Refrigerated: 3 Days	86886
ACL	APLPNL	Antiphospholipid Antibody Panel Container :Two light blue (3.2% sodium citrate) 2.7 mL and Two gold gel 5.0 mL (also acceptable for the serum sample: plain red 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 14 Days	85390 85613 85670 85732 86147 x3 86146 x3
ACL	ASO	ANTISTREPTOLYSIN O ANTIBODY Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 8 Days / Frozen: 14 Days	86060
PHC	ATIII	Antithrombin III Activity, Plasma Stability :Room temperature: 4 hours (unspun) Refrigerated: 4 hours (plasma) Frozen: 1 month at -30° C	85300
ACL	APOA1	APOLIPOPROTEIN A-1 Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 8 Days / Frozen: 3 Months	82172
ACL	APOLB	APOLIPOPROTEIN B Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 8 days / Frozen: 90 Days	82172
ACL	AQUA4 IGG	Aquaporin-4 Receptor Ab IgG w/ Reflex Container :Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.15 mL) Stability :Ambient: 48 Hours Refrigerated: 2 Weeks Frozen: 1 Year	86052 if reflex add 86256
ACL	ADH	ARGININE VASOPRESSIN Container :Three lavender (EDTA) 3.0 mL (also acceptable: pink (EDTA) 6.0 mL) Stability :Frozen: 1 Month	84588
ACL	ARIPSP	Aripiprazole Stability :Ambient: 2 Weeks Refrigerated: 2 Weeks Frozen: 2 Weeks	80342
ACL	ARSN	ARSENIC, BLOOD Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82175
ACL	UAS	ARSENIC, URINE Stability :Refrigerated: 14 Days / Frozen: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82175
PHC	AARTHP	Arthritis Panel; Rheumatoid Panel	84550 85652 86038 86060 86431
ACL	ARTE	Arthropod Exam Stability :Ambient: 1 month preserved / Refrigerated: unacceptable	87168
ACL	ASPIGG	Aspergillus fumigatus Antibody, IgG Container :One plain red 4.0 mL (also acceptable: One gold gel 3.5 mL) Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid repeated freeze/thaw cycles	86317
ACL	ASAG	ASPERGILLUS GALACTOMANNAN; ASPERGILLUS AG Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 7 Days	87305
PHC	AST	AST/SGOT (Aspartate Aminotransferase), Plasma or Serum Stability :Room temperature: 3 days Refrigerated: 7 days Frozen: 60 days	84450
ACL	BCGR	B CELL CLONALITY PCR Stability :Ambient: Blood or bone marrow, 24 hours / FFPE tissue, indefinitely or bone marrow, 14 days; Refrigerated: Blood	81261
ACL	BART	B. HENSELAE IgG IgM Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year avoid repeated freeze/thaw cycles	86611 x2



LAB	CODE	TEST NAME	CPT4 CODE
ACL	B2GLY	B2 GLYCOPROTEIN PANE Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 48 Hours / Frozen: 14 Days	86146 x3
ACL	BABES AB	Babesia microti Antibodies, IgG and IgM by IFA Container :Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.2 mL) Stability :Ambient: 48 Hours Refrigerated: 2 Weeks Frozen: 1 Year	86753 x2
ACL	REFID	Bacterial Isolate ID, Sensitivity Container :SPECIAL ACTIVELY GROWING BACTERIAL ISOLATE, IN PURE CULTURE, ON AGAR SLANTS, PLATES OR IN TRANSPORT MEDIA.	
ACL	BARCNF	Barbiturate Confirmation, Urine Container :3.5 ml (Min: 1.5 ml) urine refrigerated Transfer urine with no additives or preservatives to an ARUP Standard Transport Tube Stability :Ambient: 1 Week / Refrigerated: 1 Month / Frozen: 3 Years	80345
PHC	BARB	Barbiturate Screen, Urine Stability :Room temperature: 1 day Refrigerated: 5 days Frozen: indefinite at -200 C	80307
ACL	BARBQT	BARBITURATES,QUANT Container :3.5 mL (min: 1.5 mL) plasma or serum refrigerated Stability :Ambient: 1 Week After separation from cells. / Refrigerated: 2 Weeks After separation from cells. / Frozen: 3 Years After separation from cells.	80345
PHC	BMP	Basic Metabolic Panel (BMP)	80048
ACL	BCRQNT	BCR-ABL QUANTITATIVE Container :One lavender (EDTA) 6.0 mL (also acceptable: Two light blue (3.2% sodium citrate) 2.7 mL or bone marrow {1.0 mL in EDTA or Sodium or Lithium Heparin}) Stability :Refrigerated: Specimens should be stored and transported to the lab under refrigerated conditions. Testing must be performed within 40 hours of collection time. / Frozen: Unacceptable	81206
PHC	BENZO	Benzodiazepine Screen, Urine Stability :Room temperature: 1 day Refrigerated: 1 days Frozen: indefinite at -20 C	80101
ACL	BENCNF	Benzodiazepines Confirmation, Urine Container :One urine - random Stability :Ambient: 1 Week / Refrigerated: 1 Month / Frozen: 3 Years	80346
ACL	MICROB	BETA 2 MICROGLOBULIN Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82232
PHC	BHB	BETA HYDROXYBUTYRATE Stability :Ambient: 2 Hours	82010
ACL	B2T	Beta-2 Transferrin Container :Aural or nasal fluid in a sterile container. Specimen kept refrigerated may be collected if there is slow or intermittent flow or fluid. Stability :Ambient: 2 Weeks/Refrigerated: 2 Weeks/Frozen: 2 Weeks	86335
PHC	BICARB	Bicarbonate, Blood Stability :Room temperature: 2 hours Refrigerated: 4 hours	82374
ACL	BACID	BILE ACIDS, TOTAL Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 2 Weeks / Frozen: 3 Months	82239
PHC	ABIL	Bilirubin Fractionated, Blood Stability :Refrigerated: 5 days	82247 82248
PHC	DBIL	Bilirubin, Direct, Plasma or Serum Stability :Refrigerated: 5 days Frozen: 3 months at -20o C	82248
PHC	BILI CORD	Bilirubin, Total, Cord Blood Stability :Refrigerated: 2 days	82247
PHC	TBIL	Bilirubin, Total, Plasma or Serum	82247



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Refrigerated: 5 days	
ACL	BKPCR	BKV QUANT,URINE Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Month	87799
ACL	BKVQN	BKV QUANTITATIVE Container :One lavender (EDTA) 6.0 mL (also acceptable: One blue (sodium citrate)) Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Month	87799
ACL	BLAS	BLASTOMYCES ANTIGEN Container :One Urine, random (also acceptable: Gold Gel 3.5 mL, Plain Red 4.0 mL, Random Urine, Cerebral Spinal Fluid, Bronchoalveolar Lavage/Bronchial Alveolar Lavage (BRLA), Lavender (EDTA) 3.0 mL, Green (Sodium/Lithium Heparin) 4.0 mL, Light Blue (3.2% Sodium Citrate) 2.7 mL) Stability :Ambient: 2 Weeks / Refrigerated: 2 Weeks / Frozen: indefinitely	87449
PHC	BGCORD	Blood Gas, Cord Blood, Arterial Stability :Room temperature: 30 minutes Refrigerated: transport on ice if specimen will not be in lab within 30 minutes of draw Do Not Freeze	82803
PHC	VCBG	Blood Gas, Cord Blood, Venous Stability :Room temperature: 30 minutes Refrigerated: transport on ice if specimen will not be in lab within 30 minutes of draw Do Not Freeze	82803
ACL	BLFC	BLOOD,FUNGAL CULTURE Container :ONE YELLOW (SPS) 7.0 ML (ALSO ACCEPTABLE: BONE MARROW) ONE BACTEC MYCO/F LYTIC BOTTLE IS ALSO ACCEPTABLE, BUT NOT AVAILABLE FROM ACL LABORATORIES. Stability :AMBIENT: 72 HOURS	87103
ACL	APBONE	BONE SPECIFIC ALK PH Container :One gold gel 5.0 mL (Also acceptable: One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 2 Hours / Refrigerated: 2 Days / Frozen: 2 Months	84080
ACL	PERDNA	BORDETELLA PCR Container :ESwab - refrigerated Respiratory specimen in sterile tube - refrigerated FFPE - ambient Pleural fluid - refrigerated Biopsy Lung Tissue (in Eswab or UTM media) - Frozen. Stability :Ambient: 3 Days / Refrigerated: 14 Days / Frozen: Unacceptable (except Biopsy Lung Tissue)	87801
PHC	BRAF	BRAF Mutaion	81210 88381
ACL	BRONCI	Bronchial, Culture and Gram Stain Container :LOWER RESPIRATORY FLUIDS IN A STERILE LEAK PROOF CONTAINER REFRIGERATED. Stability :REFRIGERATED: 3 DAYS	87070 87205
ACL	BRUCAB	BRUCELLA AB TOTAL Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 6 Months	86622
PHC	BUNPO	BUN, POST-DIALYSIS Stability :Room temperature: 3 days Refrigerated: 7 days	84520
ACL	BUPM	Buprenorphine and Metabolites, Urine, Quantitative Container :One urine - random Stability :Ambient: 1 Week / Refrigerated: 2 Month / Frozen: 3 Years	80348
PHC	BURFLA	Bursal Fluid Analysis	89051
ACL	CDFDNA	C DIFFICILE TOXIN, PCR Container :One Stool Unpreserved Stability :Ambient: 48 Hours / Refrigerated: 5 Days / Frozen: Unacceptable	87493
ACL	CPEP	C PEPTIDE Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 24 Hours / Frozen: 14 Days	84681



LAB	CODE	TEST NAME	CPT4 CODE
PHC	CRP	C-Reactive Protein (CRP), Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 3 days for dimension. 7 days for Vista Frozen: 6 months at -20o C	86140
ACL	BCTX	C-TELOPEPTIDE BETA CROSS LINKED,SERUM Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL, light green (lithium heparin gel) 3.0 mL or pink (K2EDTA) 6.0 mL) Stability :Ambient: 4 Hours after separation from cells / Refrigerated: 8 Hours after separation from cells / Frozen: 3 Months after separation from cells	82523
ACL	CTRACH	C. TRACHOMATIS CULT Stability :Ambient: 1 Hour / Refrigerated: 2 Days / Frozen: 1 Month at -70 Celsius	87110 87140
ACL	C1ES	C1 ESTERASE INHIBIT Container :One gold gel 5.0 mL Stability :Ambient: 2 Hours / Refrigerated: Unacceptable / Frozen: 2 Weeks	86160
ACL	C1EFUN	C1 Esterase Inhibitor Functional Container :One gold gel 3.5 mL Stability :Ambient: 2 Hours / Frozen: 2 Weeks	86161
ACL	CC1Q	C1q Complement Protein Container :One lavender (EDTA) 3.0 mL Stability :Ambient: Unacceptable / Refrigerated: 48 Hours / Frozen: 1 Month	86160
ACL	C3	C3 Complement Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 6 Hours / Refrigerated: 24 Hours / Frozen: 14 Days	86160
ACL	C4	C4 Complement Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 6 Hours / Refrigerated: 24 Hours / Frozen: 14 Days	86160
ACL	CA2729	CA 27.29 Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 7 Days	86300
ACL	CA153	CA15-3 CANCER AG Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 24 Hours / Frozen: 14 Days	86300
ACL	CA199B	CA19-9 CANCER AG Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 14 Days	86301
ACL	CADM	CADMIUM, BLOOD Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82300
ACL	CALCIT	CALCITONIN Container :One gold gel 3.5 mL (also acceptable: One green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 3 Months	82308
PHC	IONCA	Calcium, Ionized, Blood	82330
PHC	CA	Calcium, Total, Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: 6 months	82310
ACL	UCAR	CALCIUM,URINE Container :Urine, random Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 14 Days	82340
ACL	UCAT	CALCIUM,URINE 24 HR Container :Urine Stability :Ambient: 3 Days / Refrigerated: 5 Days / Frozen: 14 Days	82340
ACL	CSA	CALCULI ANALYSIS Container :One Calculi Specimen Stability :Ambient: 1 Year / Refrigerated: 1 Year / Frozen: 1 Year	82365
ACL	CALPR	CALPROTECTIN, FECAL Container :One Stool, random Stability :Refrigerated: 10 Days	83993
ACL	CALR	CALR EXON 9 MUTAION Container :Lavender (EDTA) or green (sodium heparin) whole blood or bone marrow.	81219



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: 24 hours; Refrigerated: 5 days; Frozen: Unacceptable	
PHC	CA125B	Cancer Antigen 125 (CA 125), Serum Stability :Room temperature: 8 hours Refrigerated: 1 days	86304
ACL	CAAB	CANDIDA ANTIBODIES Container :One blood gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 6 Months	86628 x3
PHC	CANN	Cannabinoids (THC) Screen, Urine Stability :Refrigerated: 1 days	80101
ACL	CANMEC	Cannabinoids, Meconium Container :0.5 g meconium ambient Stability :Ambient: 1 Week/Refrigerated: 3 Months/Frozen: 1 Year	80349
PHC	TGRTOL	Carbamazepine, Serum Stability :Room temperature: 8 hours Refrigerated: 1 days	80156
ACL	CRESC	Carbapenem Resistant Culture Container :One rectal swab specimen submitted in ESwab (also acceptable: aerobic transport media) Stability :AMBIENT: 48 HOURS / REFRIGERATED: 48 HOURS / FROZEN: UNACCEPTABLE	87081
PHC	TCO2	Carbon Dioxide, Whole Blood	82374
PHC	CARBMO	Carbon Monoxide, Blood Stability :Refrigerated: 48 hours if unopened. Assay as soon as possible	82375
PHC	CEAMR	Carcinoembryonic Antigen (CEA), Serum Stability :Room temperature: 8 hours Refrigerated: 2 days	82378
ACL	CARDIO	CARDIOLIPIN Container :One gold gel 5.0 mL (also acceptable: One 10mL Red Top) Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 14 Days	86147 x3
ACL	CARN	CARNITINE FR/TOT, PL Container :One lavender (EDTA) 6.0 mL (also acceptable: One plain red 4.0 mL and One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 2 Hours / Refrigerated: 12 Hours / Frozen: 1 Month	82379
ACL	CAROT	CAROTENE Container :Two gold gel 3.5 mL (also acceptable: Plain red 6.0 mL) Stability :Ambient: unacceptable/ Refrigerated: 4 hours / Frozen: 1 Months	82380
ACL	CATAPL	CATECHOLAMINE FRA PL Container :Two light green (lithium heparin gel) 4.5 mL Stability :Frozen: 1 Month (-20 degree c) or 1 Year (-70 degree C)	82384
ACL	UCATA	Catecholamines, Urine Free Container :One 24 hour urine. Sample must be refrigerated during collection. Stability :Ambient: Unacceptable / Refrigerated: 1 Week Preserved / Frozen: 6 Months Preserved	82384
ACL	CCPEP	CCP PEPT AB Container :One gold gel 3.5 mL Stability :Refrigerated: 48 Hours / Frozen: 14 Days	86200
ACL	CEAFL	CEA FLUID Container :One body fluid in a clean sterile container. Acceptable body fluids include Abdominal (ABDFL), Ascites (ASFL), Cyst (CYST), Gastric (GASD), Pancreatic (PAFL), Peritoneal Dialysate (PDFL), Pericardial (PECFL), Peritoneal (Paracentesis) Fluid (PERFL), Ple Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82378
PHC	CELRX	Celiac Disease Screen Additional Antibodies	83516 x2
ACL	CELSCR	CELIAC SCR OVER 2YR Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	82784 83516
ACL	CELSR2	CELIAC SCR UNDER 2YR Stability :Ambient: Unacceptable / Refrigerated: 3 Days / Frozen: 14 Days	82784 83516 x3
ACL	CERUP	CERULOPLASMIN Container :One gold gel 3.5 mL (also acceptable: One red top) Stability :Ambient: 8 Hours / Refrigerated: 72 Hours / Frozen: 14 Days	82390
ACL	CHIKUN AB	CHIKUNGUNYA ANTIBODIES, IGG AND IGM Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid repeated freeze/thaw cycles	86790 x2



LAB	CODE	TEST NAME	CPT4 CODE
ACL	CHIKUNGUNYA PCR	Chikungunya virus by PCR Container :One gold gel 5.0 mL blood Stability :Ambient: 24 hours / Refrigerated: 5 days / Frozen: 2 weeks	87798
ACL	CGPT	CHLAM/GC AMP. PROBE Container :Aptima Collection Container Stability :Ambient: 30 Days Urines, 60 Days Swab / Refrigerated: 30 Days Urines, 60 Days Swab	87491 87591
PHC	CHLAMA	CHLAMYDIA AB IGG/IGM Container :One plain red 6.0 mL (also acceptable: One gold gel 3.5 mL) Stability :Ambient: 2 Days / Refrigerated: 14 Days / Frozen: 1 Year	86631 x3 86632 x3
ACL	CTPT	Chlamydia by Nucleic Acid Amplification Container :Aptima Collection Container Stability :Ambient: 30 Days Urines, 60 Days Swab / Refrigerated: 30 Days Urines, 60 Days Swab	87491
ACL	CHLAMG	CHLAMYDIA IgG PANEL Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86631 x3
ACL	CHLAMM	CHLAMYDIA IgM PANEL Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86632 x3
PHC	CL	Chloride, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 1 year at -20 o C	82435
PHC	UCLR	CHLORIDE, URINE RANDOM Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82436
ACL	UCLT	CHLORIDE,URINE 24 HR Container :Urine Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82436
PHC	HDL	Cholesterol, HDL (Alpha-1 Fraction), Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 7 days Frozen: 3 months at -20o C	83718
PHC	DLDL	Cholesterol, LDL, Direct, Plasma or Serum Stability :Room temperature: 1 days Refrigerated: 3 days Frozen: 3 months at -20o C	83721
PHC	FLCHOL	Cholesterol, Total, Body Fluid	82465
PHC	CHOL	Cholesterol, Total, Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: WMH 5 days: OMH/DNG/PAC 2 days Frozen: 3 months at -20o C	82465
ACL	TXCHMR	CHROM ANALY W REFLEX Container :One green (sodium heparin no gel) 4.0 mL and One lavender (EDTA) 6.0 mL . Both specimens necessary for testing. Stability :Ambient: 72 Hours / Refrigerated: 7 Days / Frozen: Unacceptable	81229 88230 88262 88291
ACL	TXCHCV	CHROM ANALYSIS CVS Stability :Ambient: 24 Hours / Frozen: Unacceptable Note: Specimen viability decreases with time. Send all samples to the cytogenetics lab for viability determination as CVS is an irretrievable sample.	88235 88267 88280 88285 88291
ACL	TXCHBM	CHROM BONE MARROW Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample. / Frozen: Unacceptable	88237 88264 88291
ACL	TXCHLK	CHROM LEUKEMIC BLD Container :One whole leukemic blood green (sodium heparin no gel) 4.0 mL Stability :Ambient: 48 Hours / Frozen: Unacceptable	88237 88264 88291
ACL	CROMAB	CHROMATIN AB, IgG Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
PHC	CHCOP	CHROMIUM & COBALT, SERUM	82495 83018
ACL	CHROMB	Chromium Blood Stability :Ambient: 1 Month / Refrigerated: 1 Month / Frozen: Unacceptable	82495
ACL	CHROM	CHROMIUM, SERUM	82495



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One royal blue (no additive) 6.0 mL Stability :Ambient: Indefinitely / Refrigerated: Indefinitely / Frozen: Indefinitely	
ACL	TXCHPS	CHROMO PARTIAL BLOOD Container :One green (sodium heparin no gel) 4.0 mL (also acceptable: One royal blue (sodium heparin)) Stability :Ambient: 72 Hours / Frozen: Unacceptable	88230 88261 88291
ACL	CHF10	CHROMOGENIC F10 ACT Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85260
ACL	CHROMA	CHROMOGRANIN A Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 2 Week / Frozen: 6 Weeks	86316
ACL	TXCHAM	CHROMOSOME AMNIOTIC Stability :Ambient: 24 Hours / Frozen: Unacceptable Note: Specimen viability decreases with time. Send all samples to the cytogenetics lab for viability determination as amniotic fluid is an irretrievable sample.	88235 88269 88280 88285 88291
ACL	TXCHBL	Chromosome Analysis Blood Container :One green (sodium heparin no gel) 4.0 mL (also acceptable: One royal blue (sodium heparin)) Stability :Ambient: 72 Hours / Frozen: Unacceptable	88230 88262
ACL	MAPOCR	Chromosome Analysis POC with Microarray Reflex AND Maternal Contamination Hold Stability :AMBIENT: 72 HOURS / REFRIGERATED: ACCEPTABLE, NOT PREFERRED. / FROZEN: UNACCEPTABLE	81229 88233 88262 88291
ACL	TXCHRT	CHROMOSOME TISSUE/POC Stability :Ambient: 72 Hours / Refrigerated: Acceptable / Frozen: Unacceptable Note: Specimen viability decreases with time. However this is an irretrievable specimen. Send specimen to cytogenetics for viability determination.	88233 88262 88291
ACL	UCIT	CITRATE, URINE Container :One urine - 24 hour Stability :Ambient: 8 Hours / Refrigerated: 1 Week / Frozen: 1 Month	82507
ACL	Clobazam	Clobazam Container :One plain red 6.0 mL blood Also acceptable: One lavender (EDTA) Stability :Ambient: 3 Days Refrigerated: 2 Weeks Frozen: 2 Months	80339
ACL	CLOMPR	Clomipramine and Metabolite Stability :Ambient: 5 Days after separation from cells / Refrigerated: 2 Weeks after separation from cells / Frozen: 6 Months after separation from cells	80335
PHC	CDFDNA	Clostridioides difficile Toxin by Real-Time PCR, Feces - STAT Order Stability :Room temperature: 24 hours, Refrigerated: 5 Days Refrigerated: 5 days Frozen: Do Not Freeze	87493
PHC	CLZ	Clozapine, Serum Stability :Ambient: 5 Weeks / Refrigerated: 2 Months / Frozen: 2 Months	80159
PHC	CMVAB	CMV IgG AND IgM Antibodies	86644 86645
ACL	CMVAVI	CMV IgG AVIDITY Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 1 Week / Refrigerated: 2 Weeks / Frozen: 1 Month	86644 x2
ACL	CMVQL	CMV QUALITATIVE Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Month	87496
ACL	CMVPLA	CMV QUANTITATIVE Container :One lavender (EDTA) 6.0 mL (also acceptable: One blue (sodium citrate)) Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Month	87497
PHC	COAGPO	Coagulation Heart Panel (Post-Op) Stability :Room temperature: Lavender tube - 8 hours/Blue tube - 4 hours (1 hour if patient is on heparin) Refrigerated: Lavender tube - 24 hours/Blue tube - Do Not Refrigerate Frozen: Lavender tube - Do Not freeze/Blue tube (plasma) - 6 months	85049 85384 85610 85730
ACL	COBALT	COBALT, SERUM Container :One royal blue (no additive) 6.0 mL (also acceptable: One royal blue (EDTA) 6.0 mL) Stability :Ambient: Indefinitely / Refrigerated: Indefinitely / Frozen: Indefinitely	83018
ACL	COBALB	COBALT,BLOOD Stability :Ambient: If the specimen is drawn and stored in the appropriate container, the trace element values do not change with time. / Refrigerated: Not preferred, but acceptable / Frozen: Unacceptable	83018



LAB	CODE	TEST NAME	CPT4 CODE
ACL	COCCNF	Cocaine Metabolite, Urine, Quantitative Container :One urine - random Stability :Ambient: 1 Week / Refrigerated: 1 Month / Frozen: 3 Years	80353
PHC	COCA	Cocaine Metabolites, Urine Stability :Refrigerated: 1 days Frozen: >5 days at -20o C	80307
ACL	COCCPAN	Coccidioides Antibody Reflexive Panel Container :One gold gel 5.0 mL blood. Stability :Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)	86635 x2
ACL	COQ10	COENZYME Q10 TOTAL Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL or One light green (lithium heparin gel) 4.5 mL) Stability :Ambient: Unacceptable / Refrigerated: 3 Weeks / Frozen: 1 Month	82491
ACL	CAGG	COLD AGGLUTININS Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 14 Days / Frozen: 1 Year	86157
ACL	COMPD	Complement Activity, Total Container :One plain red 6.0 mL Stability :Ambient: 2 Hours / Refrigerated: Unacceptable / Frozen: 2 Weeks	86162
ACL	C2	COMPLEMENT C2 Container :One gold gel 3.5 mL Stability :Ambient: 2 Hours / Frozen: 2 Weeks	86160
ACL	C5	COMPLEMENT C5 Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Frozen: 2 Weeks, after separation from cells. Avoid repeated freeze/thaw cycles.	86160
ACL	COMP3A	Complement Component 3A Container :One lavender (EDTA) 3.0 mL blood Stability :Frozen: 1 Year at -70C	86160
ACL	COMP4A	Complement Component 4A Container :One lavender (EDTA) 3.0 mL blood Stability :Frozen: 1 Year at -70C	86160
PHC	CBCC	Complete Blood Count (CBC) with Differential, Blood Stability :Room temperature: <24 hours Refrigerated: 36 hours: in some cases 48 hours - specimen dependent (to be determined by testing personnel) Frozen: Not acceptable	85025
PHC	CMP	Comprehensive Metabolic Panel (CMP)	80053
PHC	DATB	Coombs, Direct Antihuman Globulin (DAT), Blood Stability :Room temperature: 8 hours Refrigerated: 24 hours Frozen: Do Not Freeze	86880
ACL	COP	COPPER Container :One royal/dark blue (EDTA) 6.0 mL or royal blue (no additive) 6.0 MI Stability :Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82525
ACL	UCOP	COPPER,UR 24HR Container :One Urine Stability :Ambient: 30 Days / Refrigerated: 30 Days/ Frozen: 30 Days	82525
PHC	COVID	Coronavirus COVID-19 by PCR Container :VTM/UTM Stability :Ambient: Hologic Direct Load Tube (Huckleberry) or Multi-test Tube: 6 days Refrigerated: Hologic Direct Load Tube (Huckleberry) or Multi-test Tube: 6 days Frozen: Not Acceptable - All collection devices Refrigerated: UTM and Eswab: 4 days	87635 U0002 U0003 U0004
ACL	CORT120	CORTISOL 120 MIN Container :One gold gel 3.5 mL (also acceptable: One gold gel microtainer 600 mcL (clear or amber) or plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82533
ACL	CORT30	CORTISOL 30 MIN Container :One gold gel 3.5 mL (also acceptable: One gold gel microtainer 600 mcL (clear or amber) or plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82533



LAB	CODE	TEST NAME	CPT4 CODE
ACL	CORTAM	CORTISOL AM Container :One gold gel 3.5 mL (also acceptable: plain red 4.0 mL or gold gel microtainer 600 mcL (clear or amber)) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 28 Days	82533
ACL	CORTB	CORTISOL BASELINE Container :One gold gel 3.5 mL (also acceptable: One gold gel microtainer 600 mcL (clear or amber) or plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82533
ACL	CORTPM	CORTISOL PM Container :One gold gel 3.5 mL (also acceptable: plain red 4.0 mL or gold gel microtainer 600 mcL (clear or amber)) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 28 Days	82533
ACL	CORTSA	Cortisol Saliva Container :Salivette collection kit Stability :Ambient: 1 Week / Refrigerated: 3 Weeks / Frozen: 6 Months	82533
ACL	CORTPDEX	CORTISOL, POST DEXAMETHASONE SUPPRESSION Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours Refrigerated: 7 Days Frozen: 14 Days	82533
ACL	COVSPIKEIGG	COVID VACCINE IMMUNE RESPONSE IGG Container :One gold gel 5.0 mL (also acceptable: gold gel 3.5 mL) Stability :Ambient: 7 Days / Refrigerated: 14 Days / Frozen: 14 Days	86769
ACL	COXSA	COXSACKIE A9 AB Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86658
ACL	COXSB	COXSACKIE B VIRAL AB Container :One gold gel 3.5 mL (also acceptable: One Spinal fluid (CSF)) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86658 x6
PHC	CKNRMB	Creatine Kinase (CK), Total, No Reflex MB, Plasma or Serum Stability :Refrigerated: 7 days Frozen: 30 days	82550
PHC	CREACL	Creatinine Clearance, Serum and Urine Stability :Refrigerated: Serum - 2 days, Urine 4 days	81050 82575
PHC	GFRE	Creatinine with Estimated GFR Stability :Refrigerated: 2 days	82565
PHC	UCREAT	Creatinine, 24 Hour, Urine Stability :Refrigerated: 4 days	81050 82575
PHC	FCRET	Creatinine, Body Fluid Stability :Refrigerated: 2 days	82570
PHC	UCREAR	Creatinine, Random, Urine Stability :Refrigerated: 4 days	82570
ACL	CRODIS	Crohn Disease Prognostic Panel Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours after separation from cells / Refrigerated: 2 Weeks after separation from cells / Frozen: 1 Year after separation from cells, avoid repeated freeze/thaw cycles	83516 x3 86671
ACL	NTXB	Cross-Linked N-telopeptide, Serum Container :One gold gel 3.5 mL Stability :Ambient: 1 Day Refrigerated: 7 Days Frozen: 14 Days (Avoid repeated freeze/thaw cycles)	82523
PHC	XMTCH	Crossmatch, Blood Products	86920
ACL	HSCRIP	CRP HIGH SENSITIVITY Container :One gold gel 3.5 mL (also acceptable: One plain red) Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 14 Days	86141
ACL	CRYOFI	Cryofibrinogen Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: 1 Week after separation from cells / Refrigerated: 1 Week after separation from cells	82585
ACL	CRYO	CRYOGLOBULIN Container :6.0 mL red tube warmed Stability :Ambient: 7 Days / Refrigerated: Unacceptable / Frozen: Unacceptable	82595
PHC	CRYOPS	Cryoglobulin and Cryofibrinogen Panel, Serum and Plasma	82585 82595



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :After separation from cells: Ambient: 1 week: Refrigerated: 1 week for CRYOFIBRINOGEN only: Frozen: Unacceptable	
ACL	CPCL	Cryptococcus Antigen, Serum Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86403
PHC	SCRYST	Crystal Identification, Synovial Fluid Stability :Room temperature: 24 hours Refrigerated: 24 hours Frozen: Do Not Freeze	89060
PHC	CSFAH	CSF Cell Count & Differential	89051
ACL	CSFIG	CSF IgG Container :One Cerebrospinal fluid (CSF) Stability :Refrigerated: 48 Hours / Frozen: 14 Days	82784
ACL	CSFMCR	CSF VDRL Container :Spinal fluid (CSF) Stability :Refrigerated: 7 Days	86592
ACL	EVDCLS	CSF, EVD/Shunt Culture and Gram Stain Container :STERILE LEAKPROOF CONTAINER Stability :AMBIENT: 24 HOURS	87015 87070 87205
ACL	CSFCS	CULTURE AND SMEAR, CSF Container :STERILE LEAKPROOF CONTAINER Stability :AMBIENT: 24 HOURS	87015 87070 87205
ACL	ROCS	Culture, Aerobic with Smear 1 Container :ESWAB (ALSO ACCEPTABLE: AEROBIC CULTURE SWAB) Stability :AMBIENT: 48 HOURS	87070 87205
ACL	CULAAN	Culture, Anaerobic/Aerobic w Smear Container :ANAEROBIC TRANSPORT MEDIA, AMBIENT FLUID/TISSUES IN STERILE CONTAINER, AMBIENT Stability :AMBIENT: 48 HOURS IN ANAEROBIC TRANSPORT MEDIA.	87070 87075 87205
ACL	BLC	CULTURE, BLOOD Container :BACT/ALERT BLOOD CULTURE BOTTLES. (OTHER BLOOD CULTURE BOTTLE TYPES ARE ACCEPTABLE. HOWEVER, THESE BOTTLES WILL BE PROCESSED MANUALLY AND MAY LEAD TO INCREASED TIMED TO INITIAL DETECTION OF POSITIVE CULTURES.) Stability :AMBIENT: 24 HOURS	87040
ACL	BLMY	Culture, Blood for Mycobacteria Container :ONE YELLOW (SPS) 7.0 ML (ALSO ACCEPTABLE: BONE MARROW) ONE BACTEC MYCO/F LYTIC BOTTLE IS ALSO ACCEPTABLE, BUT NOT AVAILABLE FROM ACL LABORATORIES. Stability :AMBIENT: 72 HOURS	87116
ACL	CATHCL	Culture, Cath tip Aerobic Container :CATHETER TIP IN A STERILE CONTAINER, AMBIENT (DO NOT SUBMIT ASSOCIATED TUBING) Stability :AMBIENT: 24 HOURS	87070
ACL	ENVC	CULTURE, ENVIROMENTAL Container :ESWAB (ALSO ACCEPTABLE: AEROBIC CULTURE SWAB, STERILE CONTAINER, INOCULATED MEDIA INCLUDING HETEROTROPHIC MEDIA PLATE SAMPLER (HPS), RODAC PLATES, ETC.) ESWAB SPECIMEN COLLECTION DEVICE - WHITE CA Stability :AMBIENT: 24 HOURS / REFRIGERATED: 24 HOURS	87070
ACL	FNC	Culture, Fungal Hair, Skin and Nails Container :CLEAN DRY CONTAINER AT AMBIENT TEMPERATURE FOR HAIR, SKIN SCRAPINGS OR NAILS. Stability :AMBIENT: 1 WEEKS	87101
ACL	GCC	Culture, GC Container :FLUID OBTAINED FROM ANY BODY SITE IN STERILE LEAK-PROOF CONTAINER OR ESWAB (PREFERRED) OR AEROBIC CULTURE SWAB OBTAINED FROM ANY BODY SITE. Stability :AMBIENT: 12 HOURS	87081
ACL	SINGCL	Culture, Single Bacterial Identification Container :Eswab Stability :AMBIENT: 24 HOURS	87070
ACL	ENPC	CULTURE,ENTERIC PATHOGEN Container :STOOL IN LIQUID C&S (CARY - BLAIR) TRANSPORT MEDIA RAW STOOL IN CLEAN LEAKPROOF CONTAINER. Stability :AMBIENT: 1 HOUR, NO PRESERVATIVE. YIELDS ARE OPTIMIZED IF STOOL IS IMMEDIATELY PLACED IN C&S (CARY-BLAIR) TRANSPORT MEDIA. AMBIENT: 4 DAYS WITH PRESERVATIVE /	87045 87046 87427 87449



LAB	CODE	TEST NAME	CPT4 CODE
		REFRIGERATED: UNACCEPTABLE / FROZEN: UNACCEPTABLE	
ACL	CMRSA	CULTURE,MRSA Container :ESWAB COLLECTION KIT, AMBIENT AEROBIC CULTURE SWAB IN TRANSPORT MEDIA, AMBIENT. FLUID IN A STERILE LEAK-PROOF CONTAINER, REFRIGERATED. URINE IN A STERILE LEAK-PROOF CONTAINER, REFRIGERATED OR URINE IN A GRAY TOP BORIC ACID PRESERVATIVE TUBE, AMBIENT. Stability :AMBIENT: 48 HOURS FOR ESWAB, AEROBIC CULTURE SWAB, BORIC ACID URINE TUBE / REFRIGERATED: 24 HOURS FOR FLUIDS AND UNPRESERVED URINE	87081
ACL	FNCS	CULTURE/SMEAR, FUNGAL Stability :Ambient: 72 Hours / Refrigerated: 72 Hours	87102 88312
ACL	CAFBS	Culture/Smear, Mycobacteria Container :STERILE LEAK-PROOF CONTAINER Stability :AMBIENT: 24 HOURS FOR SPUTUM, BODY FLUID, TISSUE, BRONCHIAL SPECIMEN, OR GASTRIC FLUID / REFRIGERATED: 7 DAYS FOR SPUTUM, BODY FLUID, TISSUE, OR BRONCHIAL SPECIMEN. REFRIGERATED: 24 HOURS FOR GASTRIC FLUID	87015 87116 87206
ACL	CYCS	Cyclospora/Cystoisospora Screen Container :STOOL OR DUODENAL ASPIRATES PRESERVED IN SAF OR 10% FORMALIN Stability :AMBIENT: 1 MONTH (PRESERVED)	87177
PHC	CYCLOB	Cyclosporine, Blood (Draw Only)	80158
ACL	CSAPAR	CYCLOSPORINE-WHL BLD Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80158
ACL	CYSTC	CYSTATIN C Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL or green (lithium or sodium heparin) 3.0 mL) Stability :Ambient: unacceptable / Refrigerated: 1 Week / Frozen: 3 Months	82610
ACL	CYSTIC	CYSTIC FIBROSIS Container :One lavender (EDTA) 3.0 mL (also acceptable: One light blue (3.2% sodium citrate) 2.7 mL or 1 yellow (ACD) 10 mL) Stability :Ambient: 72 Hours / Refrigerated: 14 Days	81220
ACL	CNFNA	CYTOLOGY FINE NEEDLE ASPIRATE	88173
ACL	CMVGAB	CYTOMEGALOVIRUS IGG Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 14 Days	86644
ACL	CMVMAB	CYTOMEGALOVIRUS IGM Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: 3 Days / Frozen: 14 Days	86645
PHC	DIMER	D-Dimer, Plasma Stability :Room temperature: 8 hours Refrigerated: Not Acceptable Frozen: 1 month at -20° C	85379
ACL	DENGUE	Dengue virus (Type 1-4) PCR Container :One gold gel 5.0 mL blood Stability :Ambient: 24 hours / Refrigerated: 5 days / Frozen: 2 weeks	87798 x4
ACL	DENGUE AB	DENGUE VIRUS IgG/IgM Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid repeated freeze/thaw cycles	86790 x2
ACL	DEXA	DEXAMETHASONE Container :One gold gel 3.5 mL (also acceptable: One lavender (EDTA) 3.0 mL and One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 1 Week / Frozen: 6 Months	80299
ACL	DHEAS	DHEA SULFATE Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	82627
PHC	DIGX	Digoxin, Serum Stability :Room temperature: 8 hours Refrigerated: 7 days Frozen: 6 months	80162
ACL	DIPTET	DIPHTher/TETANUS ABS Container :One gold gel 3.5 mL Stability :Ambient: 7 Days / Refrigerated: 14 Days / Frozen: 2 Months	86648 86774



LAB	CODE	TEST NAME	CPT4 CODE
PHC	DICPNL	Disseminated Intravascular Coagulation (DIC) Panel Stability :Room temperature: 4 hours (Blue tubes) Refrigerated: Do Not Store Refrigerated Frozen: 6 months (decanted plasma from Blue tubes) / Do Not Freeze Lavender tube	85300 85384 85610 85730 85049 85379
ACL	DNASE	DNASE-B ANTIBODY Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86215
ACL	DRGSCM	Drug Detection Panel, Meconium, Qualitative Container :One 4.0 g (min: 2.0 g) meconium refrigerated. Meconium cannot be sent in a diaper, or with any foreign material in collection cup. Stability :Ambient: 1 Week / Refrigerated: 3 Weeks / Frozen: 1 Year	80326 80347 80364 80355
PHC	DRGSCU	Drug Screen (Qualitative), Urine Stability :Refrigerated: 24 hours Frozen: No maximum indicated	80307
PHC	DRGSCB	Drug Screen (Quantitative), Serum Stability :Room temperature: 8 hours Refrigerated: 14 days Frozen: 45 days	80320 80329
PHC	DRGURX	Drug Screen, Urine, with Reflex Confirmation Stability :Refrigerated: 1 days / Confirmation: 14 days	80307
ACL	EHISTO	E HISTOLYTICA AB Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid multiple freeze/thaw cycles	86753
ACL	EHRAB	E. CHAFFEENSIS ABS Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid repeated freeze/thaw cycles.	86666 x2
ACL	EARCLS	Ear, Culture, Aerobic Culture with Gram Stain Container :ESWAB (ALSO ACCEPTABLE: AEROBIC CULTURE SWAB) Stability :AMBIENT: 48 HOURS	87070 87205
ACL	EBVSRN	EBV AB SCREEN Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86665 x2
ACL	EB	EBV CHRONIC PNL Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Years	86663 86664 86665 x2
ACL	EBVQL	EBV QUAL Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Month	87798
ACL	EOSPRO	ECP Container :One gold gel 5.0 mL Stability :Ambient: 48 Hours / Refrigerated: 1 Week / Frozen: 2 Weeks	83520
ACL	EHRL	Ehrlichia/Anaplasma PCR, Blood Container :LAVENDER (EDTA) Stability :REFRIGERATED: 7 DAYS / FROZEN: 7 DAYS	87468 87484 87798 x2
ACL	STLYTE	ELECTROLYTES, FECAL Stability :Ambient: 30 minutes / Refrigerated: 7 Days / Frozen: 14 Days	82438 84999 x2
PHC	LYTES	Electrolytes, Plasma or Serum Stability :Room temperature: 1 day Refrigerated: 7 days Frozen: 14 days	80051
ACL	ENDO	ENDOMYSIAL IgA AB Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 1 Year	86256
ACL	EOSM	EOSINOPHIL SMEAR Stability :Ambient: Swab or Sputum 4 hours. Slides 7 days / Refrigerated: Swab or Sputum 24 hours. Slides 7 days	89190
ACL	EBVQN	EPSTEIN-BARR VIRUS QUANTITATIVE Container :ONE LAVENDER (EDTA) 6.0 ML (ALSO ACCEPTABLE: TWO LAVENDER (EDTA) 3.0 ML AND OR LIGHT BLUE (3.2% SODIUM CITRATE) 2.7 ML) Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Month	87799
ACL	ERYTHR	ERYTHROPOIETIN Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82668



LAB	CODE	TEST NAME	CPT4 CODE
PHC	ESR	ESR (Erythrocyte Sedimentation Rate), Blood	85652
PHC	ESTRD	Estradiol, Serum Stability :Room temperature: 1 day Refrigerated: 2 days Frozen: 6 months	82670
ACL	ESTRIOL	ESTRIOL, SERUM Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	82677
ACL	ESTRO	ESTROGEN SERUM FRAC Container :One gold gel 3.5 mL (also acceptable: One lavender (EDTA) 3.0 mL and One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 2 Days / Refrigerated: 1 Week / Frozen: 1 Month	82671
ACL	ESTRONE	ESTRONE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 1 Week / Frozen: 1 Month	82679
ACL	ETHOS	ETHOSUXIMIDE Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Ambient: 5 Days / Refrigerated: 2 Weeks / Frozen: 2 Months	80168
PHC	ETGLY	Ethylene Glycol, Serum Stability :Room temperature: Do Not Store at Room Temperature Refrigerated: OK Frozen: OK	82693
PHC	EXPOA	Exposure Panel Stability :Room temperature: 24 hours Refrigerated: 6 days Frozen: 14 days	84460 86803 86705 86706 86806 87340
PHC	EXRAS	EXRAS Mutation Analysis	81275 81403 81404 88381
PHC	ENA	Extractable Nuclear Antigen AB Stability :REFRIGERATED: 7 DAYS / FROZEN: 14 DAYS	86235 x6
ACL	EYECLS	Eye, Culture, Aerobic Culture with Gram Stain Container :ESWAB (ALSO ACCEPTABLE: AEROBIC CULTURE SWAB) Stability :AMBIENT: 48 HOURS	87070 87205
ACL	LEIDEN	FAC V LEIDEN GNTYP Container :One lavender (EDTA) 6.0 mL (also acceptable: 6.0 mL Pink (K3EDTA)) Stability :Ambient: 72 Hours / Refrigerated: 14 Days	81241
ACL	FACX	FACTOR 10 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85260
ACL	FACXI	FACTOR 11 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85270
ACL	FACXII	FACTOR 12 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: Acceptable up to 30 days.	85280
ACL	FCXIII	FACTOR 13 ACTIVITY Container :Two blood light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85290
ACL	FACII	FACTOR 2 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85210
ACL	FACV	FACTOR 5 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85220
ACL	FACVII	FACTOR 7 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85230
ACL	FACTIX	FACTOR 9 ASSAY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85250



LAB	CODE	TEST NAME	CPT4 CODE
PHC	FCVIII	Factor VIII Assay	85240
PHC	WBCSTL	Fecal Leukocytes Stability :Room temperature: 2 weeks Refrigerated: 2 weeks Frozen: 1 month	83630
ACL	FEL	FELBAMATE Container :One light green (lithium heparin gel) 4.5 mL (also acceptable: One green (lithium heparin no gel) 6.0 mL and One green (sodium heparin no gel) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Month	80299
ACL	PROTO	FEP/ZPP Container :One lavender (EDTA) 3.0 mL (also acceptable: royal/dark blue EDTA 6.0 mL) Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	84202
PHC	FERITN	Ferritin, Serum Stability :Refrigerated: 7 days Frozen: 6 months	82728
PHC	FTSCB	Fetal Cell Screen, Blood Stability :Room temperature: 8 hours Refrigerated: 24 hours Frozen: Do Not Freeze	85461
PHC	FFN	Fetal Fibronectin (fFN) Stability :Ambient: 8 Hours / Refrigerated: 72 Hours	82731
ACL	APTS	Fetal Hemoglobin Detection, Stool Stability :Ambient: 48 Hours	83033
PHC	FIB	Fibrinogen, Plasma Stability :Room temperature: 8 hours Refrigerated: Do Not Refrigerate Frozen:	85384
ACL	1STSCR	First Trimester Prenatal Screen Container :3.0 mL (Min: 1.0 mL) serum Refrigerated Separate cells within 2 hours of collection. Stability :Ambient: 72 Hours - after separation from cells Refrigerated: 2 Weeks - after separation from cells Frozen: 1 Year - after separation from cells. Avoid repeated freeze/thaw cycles.	81508
ACL	FISHAG	FISH AGGRESIVE LYMPHOMA PANEL Container :3.0 mL (min: 1.0 mL) non-diluted bone marrow aspirate in green sodium heparin tube or 5.0 mL (min: 2.0 mL) whole leukemic blood in green sodium heparin tube. Also acceptable - Lavender EDTA tubes acceptable but not preferred specimen type. Tissue Paraffin block is also acceptable but must include a marked H & E stained slide. Stability :Ambient: 48 Hours NOTE: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x2 88275
ACL	FISHBM	FISH Aquired Chromosome Abnormalities Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x2 88275 88291
ACL	FISHBA	FISH BCR/ABL Stability :Ambient: 48 Hours --- Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x2 88275 88291
ACL	FISHCL	FISH CLL PANEL Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x5 88275 x2 88291
ACL	FISHHE	FISH HYPEROSINOPHILI Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x9 88275 x4 88291
ACL	FISHMDS	FISH MDS PANEL Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x6 88275 x4
ACL	FISHMM	FISH Multiple Myeloma with reflex testing for IGH rearrangements Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x9 88275 x3 88291
ACL	FISHAL	FISH Panel for ALL Container :3.0 mL (min: 1.0 mL) non-diluted bone marrow aspirate in green sodium heparin tube or 5.0 mL (min: 2.0	88237 88271 x8 88275 x4



LAB	CODE	TEST NAME	CPT4 CODE
		mL) whole leukemic blood in green sodium heparin tube.	
		Also acceptable - Lavender EDTA tubes acceptable but not preferred specimen type. Tissue Paraffin block is also acceptable but must include a marked H & E stained slide. Stability :Ambient: 48 Hours NOTE: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	
ACL	FISHMF	FISH PARENTAL MA STU Container :One green (sodium heparin no gel) 6.0 mL Stability :Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	
ACL	FISHPR	FISH PML/RARA Stability :Ambient: 48 Hours Note: Send all samples to the cytogenetics lab for viability determination as bone marrow is an irretrievable sample.	88237 88271 x2 88275 88291
ACL	FISHPO	FISH POC Stability :Ambient: Note: Paraffin embedded tissue is stable indefinitely.	88271 x8 88275 x4 88291
ACL	FISHAN	FISH PROBE ANGLEMAN Container :One green (sodium heparin no gel) 6.0 mL Stability :Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	88230 88271 88273 88291
ACL	FISHDG	FISH PROBE DI GEORGE Container :One green (sodium heparin no gel) 6.0 mL Stability :Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	88230 88271 88273 88291
ACL	FISHPW	FISH PROBE PRADER WI Container :One green (sodium heparin no gel) 6.0 mL Stability :Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	88230 88271 88273 88291
ACL	FISHPN	FISH PROBE PRENATAL Stability :Ambient: 24 Hours Note: Send all samples to the cytogenetic laboratory for viability determination since specimen samples are irretrievable	88271 x5 88275 x2 88291
ACL	FISHBL	FISH/CONGENITAL MICD Container :One green (sodium heparin no gel) 6.0 mL Stability :Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	88230 88271 88273 88291
ACL	FLECA	FLECAINIDE Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL and One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 6 Weeks / Refrigerated: 6 Weeks / Frozen: 6 Weeks	80299
ACL	FLOW AP	FLOW CYTOMETRY AP SPECIMENS Container :Tissue (including FNA), Body Fluids, CSF: Specimens in RPMI (includes FNA and core biopsies) Stability :Ambient: 3 Days Refrigerated: 3 Days Frozen: Unacceptable	88184
ACL	FLOW BAL/BRONCH	FLOW Cytometry BAL/Broch CD4/CD8 Container :Collect in RPMI Stability :Ambient: 3 Days Refrigerated: Unacceptable Frozen: Unacceptable	88184
ACL	FLOW	FLOW Cytometry BM Container :4.0 mL (min: 2.0 mL) bone marrow, ambient.. Add bone marrow to RPMI. Stability :Ambient: 3 Days Refrigerated: Unacceptable Frozen: Unacceptable	88184
ACL	FLOW PB	FLOW Cytometry PB Container :3.0 mL (min: 3.0 mL) EDTA whole blood, ambient. Stability :Ambient: 3 Days Refrigerated: Unacceptable Frozen: Unacceptable	88184
ACL	FLTA	FLT3 ITD/TKD MUTATION ANALYSIS Container :One lavender (EDTA) 3.0 mL blood OR One ACD yellow tube 8.5 mL blood OR One lavender (EDTA) 3.0 mL bone marrow OR One ACD yellow tube 8.5 mL bone marrow Stability :Ambient: 7 Days / Refrigerated: 7 Days	85060
PHC	BALFLA	Fluid Analysis, BAL Fluid	89051
PHC	EFFLA	Fluid Analysis, Effusion-Body Fluid	89051 58007
PHC	CELPER	Fluid Analysis, Effusion-Peritoneal Fluid	89051
PHC	CELLPL	Fluid Analysis, Effusion-Pleural Fluid	89051
ACL	FLQRSC	Fluoroquinolone Resistant Enterobacteriaceae, Bacterial Culture Container :One Rectal swab specimen submitted in ESwab (also acceptable: aerobic transport media)	87081



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	
ACL	FLUOX	Fluoxetine/Norfluoxetine Container :Two plain red 6.0 mL (also acceptable: Three lavender (EDTA) 3.0 mL and Two pink (K2EDTA) 6.0 mL) Stability :Ambient: 2 Weeks / Refrigerated: 2 Weeks / Frozen: 9 Months	80332
ACL	FOL	FOLATE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days Folate is sensitive to light.	82746
ACL	ARXTRA	FONDAPARINUX ASSAY Container :Two blood light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85130
ACL	FRAGX	Fragile X Diagnosis Container :One lavender (EDTA) 6.0 mL (also acceptable: One light blue (3.2% sodium citrate) 2.7 mL and One pink (EDTA) 6.0 mL) Stability :Ambient: 72 Hours / Refrigerated: 14 Days	81243
ACL	CORTFR	Free Cortisol, Serum Container :Transfer 1 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 0.6 mL) Stability :Ambient: unacceptable Refrigerated: 7 Days after separation from cells Frozen: 6 Months after separation from cells	82530
ACL	UCORT	FREE CORTISOL, URINE Container :Urine 24 hour urine or random urine Stability :Refrigerated: 2 Weeks / Frozen: 6 Months	82530
ACL	FT3B	FREE T3 Container :One gold gel 3.5 mL (also acceptable: One gold gel microtainer 600 mcL (clear or amber) or plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	84481
ACL	FRUCT	FRUCTOSAMINE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One lavender (EDTA) 3.0 mL and One green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: 72 Hours / Refrigerated: 2 Weeks / Frozen: 2 Months	82985
ACL	FSHB	FSH Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	83001
ACL	T4HPLC	FT4 BY EQ DIALYSIS Container :One plain red 4.0 mL (also acceptable: gold gel 3.5 mL) Stability :Ambient: 4 Days / Refrigerated: 2 Weeks / Frozen: 1 Month	84439
ACL	FUNID	FUNGAL ANTIBODIES-ID Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 weeks / Frozen: 1 year Avoid repeated freeze/thaw cycles	86606 86612 86635 86698
ACL	KOHHSN	Fungal KOH Hair, Skin, Nails Container :HAIR, SKIN SCRAPINGS, OR NAILS IN A DRY, STERILE CONTAINER. SKIN SCRAPINGS FROM PEDIATRICS PATIENTS COLLECTED WITH ESWABS ARE ACCEPTABLE. Stability :AMBIENT 7 DAYS	87220
ACL	KOHOV	Fungal KOH Prep Oral/Vag Container :ESWAB AMBIENT STERILE CONTAINER FOR ESOPHAGEAL BRUSH AMBIENT Stability :AMBIENT 3 DAYS	87210
ACL	G6PDB	G-6-PD-QUANTITATIVE Container :One lavender (EDTA) 3.0 mL (also acceptable: One green (sodium or lithium heparin) 6.0 mL and One yellow (ACD) 10 mL) Stability :Ambient: 8 Hours / Refrigerated: 1 Week	82955
ACL	GAB	GABAPENTIN Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Ambient: 5 Weeks / Refrigerated: 6 Weeks / Frozen: 2 Months	80171
ACL	GABAU	Gabapentin Urine Stability :Ambient: 1 Month Refrigerated: 1 Month Frozen: 1 Month	80355
PHC	GFRR	Gadolinium Creatinine w/ eGFR Stability :Refrigerated: 2 days	82565
ACL	GAL3	GALECTIN-3 Container :One plain red 4.0 mL	82777



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: 48 Hours after separation from cells Refrigerated: 2 Weeks after separation from cells Frozen: 4 Months after separation from cells (Avoid repeated freeze/thaw cycles)	
PHC	GT	Gamma-Glutamyltransferase (GGT), Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 6 months	82977
ACL	UGHB	GAMMA-HYDROXYBUTYRIC ACID, URINE Container :One urine - random Stability :Ambient: 1 week / Refrigerated: 1 Weeks / Frozen: 3 weeks	80307
ACL	GHB	GAMMA-HYDROXYBUTYRIC ACID Container :One plain red 6.0 mL Stability :Ambient: 7 days / Refrigerated: 7 Days / Frozen: 2 months	80307
ACL	GASTR	GASTRIN Container :One gold gel 3.5 mL Stability :Refrigerated: 1 Week / Frozen: 1 Month	82941
ACL	GPARX	Gastrointestinal Parasites Panel by PCR Container :Stool in liquid Cary-Blair (enteric pathogen). Raw stool in clean leak-proof container. Stability :Ambient: 1 Hour, no preservative. 2 days in Cary-Blair./ Refrigerated: 4 Days in Cary-Blair / Frozen: Unacceptable	87505
ACL	GPPNL	Gastrointestinal Pathogen Panel Container :STOOL IN STERILE CONTAINER OR IN C&S (CARY- BLAIR) MEDIA Stability :AMBIENT: 1 HOUR, NO PRESERVATIVE. 2 DAYS IN CARY-BLAIR / REFRIGERATED: 4 DAYS IN CARY-BLAIR/ FROZEN: UNACCEPTABLE	87507
ACL	GLOMBM	GBM ANTIBODY Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	83516
ACL	NGPT	GC by Nucleic Acid Amplification Container :Aptima Collection Container Stability :Ambient: 30 Days Urines, 60 Days Swab / Refrigerated: 30 Days Urines, 60 Days Swab	87591
ACL	GENCLS	Genital Bacterial Culture with Gram Stain Container :ESWAB (ALSO ACCEPTABLE : AEROBIC CULTURE SWAB) Stability :AMBIENT: 24 HOURS	87070 87205
ACL	TXCHMA	GENOMIC MICROARRAY Container :One green (sodium heparin no gel) 4.0 mL and One lavender (EDTA) 6.0 mL , Both specimens necessary for testing. Stability :Ambient: 72 Hours / Refrigerated: 7 Days / Frozen: Unacceptable	81229 88291
ACL	MAPCO	Genomic Microarray, POC W/ Maternal Contamination hold Stability :AMBIENT: 72 HOURS / REFRIGERATED: ACCEPTABLE, NOT PREFERRED. / FROZEN: UNACCEPTABLE	81229 88291
ACL	TXPNMA	Genomic Microarray, Prenatal w/ Maternal Contamination Hold Stability :AMBIENT: 24 HOURS NOTE: SEND ALL SAMPLES TO THE CYTOGENETIC LABORATORY FOR VIABILITY DETERMINATION SINCE SPECIMEN SAMPLES ARE IRRETRIEVABLE.	81229 88291
PHC	GENTP	Gentamicin, Peak, Serum Stability :Refrigerated: 1 days	80170
PHC	GENTR	Gentamicin, Random, Serum Stability :Refrigerated: 1 days	80170
PHC	GENTT	Gentamicin, Trough, Serum Stability :Refrigerated: 1 days	80170
ACL	GLIAB	GLIADIN A&G DEAMIDAT Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 6 Months	83516 x2
ACL	GLUCA	GLUCAGON Container :Collect using one Protease Inhibitor tube (PPACK; Phe- Pro-Arg-chloromethylketone). Protease Inhibitor Stability :Refrigerated: 48 Hours after separation from cells / Frozen: 3 Months after separation from cells	82943
PHC	OSULLY	Glucose Tolerance Screen (O'Sullivan), Plasma or Serum	82950
PHC	GLUTOL	Glucose Tolerance, 2 Hour, Post Glucola®, Plasma or Serum	82947 82950
PHC	GESTT	Glucose Tolerance, 3 Hour, Gestational, Plasma or Serum	82951 x3 82952
PHC	FLG	Glucose, Body Fluid	82945



LAB	CODE	TEST NAME	CPT4 CODE
PHC	GLU	Glucose, Plasma or Serum Stability :Separated from cells: Room Temp 8 hours: refrigerated 3 days	82947
PHC	CSFG	Glucose, Spinal Fluid	82945
ACL	UGLR	GLUCOSE, URINE Container :Urine Stability :Refrigerated: 7 Days	82945
ACL	UGLT	GLUCOSE, URINE 24 HR Container :Urine Stability :Refrigerated: 7 Days	82945
ACL	GAD65	GLUTAMIC AC DECAR AB Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 11 Days / Frozen: 1 Month	86341
ACL	GANG	GM1 ANTIBODY PANEL Container :One gold gel 3.5 mL Stability :Refrigerated: 2 Weeks / Frozen: 1 Year	83516 x2
ACL	GRAS	GRAM SMEAR Container :SWAB, GLASS SLIDE IN FOLDER Stability :AMBIENT: 7 DAYS / REFRIGERATED: 7 DAYS	87205
ACL	GRAMBV	GRAM SMEAR, BACTERIAL VAGINOSIS Container :ESWAB Stability :AMBIENT: 3 DAYS	87205
ACL	GBSPCS	Group B Strep Screen by PCR with Sensitivity if indicated Container :ESWAB SPECIMEN COLLECTION DEVICE Stability :AMBIENT: 3 DAYS	87653
PHC	GBSDNA	GROUP B STREP SCREEN BY PCR-INTRAPARTUM ONLY Stability :Room temperature: 24 hours Refrigerated: 6 days Frozen: Do Not Freeze	87653
ACL	HPA	H PYLORI AB, IGG Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86677
ACL	HPSA	H.PYLORI STOOL AG Container :Clean, screw capped container Stability :Refrigerated: 3 Days / Frozen: 7 Days	87338
ACL	ORGCUL	Haemophilus ducreyi Culture Container :One ESwab Stability :Refrigerated: 72 Hours	87070
ACL	HALO	HALOPERIDOL Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Ambient: 4 Hours / Refrigerated: 1 Week / Frozen: 1 Month - Avoid repeated freeze/thaw cycles	80173
ACL	HPT	HAPTOGLOBIN Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 14 Days	83010
ACL	HPBDNA	HBV DNA Quantitation Container :One pink (EDTA) 6.0 mL blood (also acceptable: lavender (EDTA) 6.0 mL) Stability :Frozen: 1 Month	87517
ACL	HCGOUT	HCG (TUMOR MARKER) Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 2 Days / Frozen: 14 Days	84702
ACL	HEPCG	HCV GENOTYPE Container :One lavender (EDTA) 6.0 mL (also acceptable: One gold gel 5.0 mL) Stability :Frozen: 2 Months	87902
ACL	HEPCPC	HCV RNA QUANTITATIVE Container :One lavender (EDTA) 6.0 mL (also acceptable: One gold gel 5.0 mL) Stability :Ambient: 6 Hours-Whole blood / Refrigerated: 3 days-Plasma or serum / Frozen: 2 Months-Plasma or serum	87522
PHC	HCTA	Hematocrit, Blood Stability :Room temperature: 24 hours Refrigerated: 24 hours Frozen: Do Not Freeze	85014



LAB	CODE	TEST NAME	CPT4 CODE
PHC	HCTPCD	Hematocrit, Fluid Stability :Room temperature:24 hours Refrigerated: 24 hours Frozen: Do Not Freeze	85014
ACL	HGBA1	HEMOGLOBIN A1C Container :One lavender (EDTA) 3.0 mL (min 0.5) Stability :Ambient: 24 Hours / Refrigerated: 7 Days	83036
PHC	HGBA	Hemoglobin, Blood Stability :Room temperature: 24 hours Refrigerated: 24 hours Frozen: Do Not Freeze	85018
ACL	PLHGB	HEMOGLOBIN, PLASMA Container :One light green (lithium heparin gel) 3.0 mL (also acceptable: green (lithium heparin no gel) 2.0 mL) Separate plasma from cells ASAP. Stability :Refrigerated: 7 Days	83051
PHC	HGM	Hemogram, Blood Stability :Room temperature: 24 hours Refrigerated: 24 hours Frozen: Do Not Freeze	85027
ACL	HEMSID	HEMOSIDERIN, URINE Container :Urine Marble top Stability :Ambient: 2 Hours unpreserved, 48 hours in red/yellow top preservative tube / Frozen: Unacceptable	83070
ACL	HBSRAB	HEP B SURFACE AB Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86706
ACL	BSURAG	HEP B SURFACE AG Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	87340
ACL	SRAUH	Heparin Dependent Platelet Antibody-SRA Unfractionated Heparin Container :One plain red 6.0 mL Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	86022
PHC	HFP	Hepatic Function Panel	80076
ACL	HAVIGM	Hepatitis A Antibody, IgM Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86709
ACL	HEPAT	Hepatitis A IgG and IgM Screen Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86708
ACL	BCORAB	Hepatitis B Core Antibody IgM Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86705
ACL	HEPBT	Hepatitis B Core Total Antibody Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86704
ACL	HSABQR	Hepatitis B Surface Antibody, Quantitative Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86706
ACL	HBEAG	HEPATITIS BE ANTIGEN Container :One gold gel 3.5 mL (also acceptable: One lavender (EDTA) 3.0 mL or One plain red 4.0 mL) Stability :Refrigerated: 6 days / Frozen: indefinitely	87350
ACL	HCV	Hepatitis C Antibody with Reflex Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86803
ACL	HEPPNL	Hepatitis Panel Acute With Reflex HCV PCR Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80074
ACL	HSVM	HERPES SIMPLEX IGM Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 1 Year	86694
ACL	HSVPCR	HERPES SIMPLEX VIRUS 1 AND 2 BY PCR Container :Varies	87529 x2



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: Aptima device - 7 days, Other collection tubes - 8 hours / Refrigerated: Aptima device - 7 days, Other collection tubes - 72 hours / Frozen: Aptima device and Other collection tubes - 7 days	
ACL	HHV6	HERPES VIRUS 6, QUAN Container :One lavender (EDTA) 3.0 mL or One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 5 Days / Frozen: 3 Months	87533
ACL	HFE	HFE GENOTYPE Container :One lavender (EDTA) 6.0 mL (also acceptable: One light blue (3.2% sodium citrate) 2.7 mL or 6.0 ml Pink (K3EDTA)) Stability :Ambient: 72 Hours / Refrigerated: 14 Days	81256
ACL	HBELE	HGB ELECTROPHORESIS Container :One lavender (EDTA) 3.0 mL (also acceptable: Two lavender (EDTA) microtainer 500 mcL) Stability :Refrigerated: 7 Days	83020
PHC	TROP-HS	High Sensitivity Troponin I Stability :Room Temp: 8 hrs/ Refrigerated: 24 hrs/ Frozen: 14 days	84484
ACL	BHISTA	HISTAMINE, WHOLE BLOOD Container :One green (lithium heparin no gel) 6.0 mL (also acceptable: One green (sodium heparin no gel) 4.0 mL. Stability :Ambient: Unacceptable/ Refrigerated: Unacceptable / Frozen: 6 Months	83088
ACL	ANTHIS	HISTONE IgG AB Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 14 Days / Frozen: 1 Year	83516
ACL	HISTOU	HISTOPLASMA AG URINE Container :One urine - random Stability :Ambient: unacceptable/ Refrigerated: 14 Days / Frozen: 14 days	87385
PHC	HIVSCR	HIV AG/AB COMBO SCREEN Stability :Room temperature 2 days:Refrigerated 7 days	87806
ACL	HIVAGB	HIV Antigen/Antibody Screen Container :One gold gel 5.0 mL (also acceptable: Two lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	87389
ACL	HIVPCR	HIV RNA QUANTITATIVE Container :One blood lavender (EDTA) 6.0 mL Stability :Frozen: 1 Months	87536
ACL	HIVCF	HIV-1 RNA CONF Container :One blood lavender (EDTA) 6.0 mL (also acceptable: One blood gold gel 5.0 mL) Stability :Frozen: 1 Months	87535
ACL	HIV12	HIV1/HIV2 AB TEST Container :One gold gel 5.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) ***This test will be ordered by Core Labs ONLY and is a reflex from test code HIVSCR. For routine rapid HIV testing, order test code HIVSCR. Stability :Ambient: 2 Days / Refrigerated: 7 Days / Frozen: 14 Days	86703
ACL	HLAS	HLA B27-MOLECULAR Container :One ACD yellow tube Stability :Ambient: 72 Hours	81374
ACL	HLACEL	HLA CELIAC MOLECULAR Container :ACD yellow tube 8.5 mL Stability :Ambient: 5 Days	81376
PHC	PHLADR	HLA-ABC and DR, Phenotype Stability :Room temperature: 12 hours Refrigerated: Do Not Refrigerate Frozen: Do Not Freeze	81372 001 81375
PHC	ABCB	HLA-ABC, Phenotype Stability :Room temperature: 12 hours Refrigerated: Do Not Refrigerate Frozen: Do Not Freeze	81372
PHC	PHLASA	HLA-ABC, Single Antigen Stability :Room temperature: 12 hours Refrigerated: Do Not Refrigerate Frozen: Do Not Freeze	81372
ACL	CD57	HNK1 (CD57) Profile Stability :Ambient: 48 Hours	86356 86357
ACL	HOMOCY	HOMOCYSTEINE, TOTAL Container :One lavender (EDTA) 3.0 mL (also acceptable: One gold gel 3.5 mL) (Pediatric: lavender (EDTA) microtainer 500 mcL or gold gel microtainer 600 mcL (clear or amber))	83090



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Ambient: Unacceptable / Refrigerated: 48 Hours / Frozen: 14 Days	
ACL	HPVGT	HPV, Genotype 16, 18/45 Stability :Ambient: ThinPrep Pap Test specimen transport medium: 3 weeks.	87625
ACL	HPVHR	HPV, High Risk Stability :Ambient: ThinPrep Pap Test specimen transport medium: 3 weeks.	87624
ACL	HSV1GR	HSV1 GLYCOPROTEIN AB Container :One gold gel 3.5 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days (avoid repeated freeze/thaw cycles.)	86695
ACL	HSV2GR	HSV2 GLYCOPROTEIN AB Container :One gold gel 3.5 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days (avoid repeated freeze/thaw cycles)	86696
ACL	HTLV1C	HTLV I/II CONFIRMATI	
ACL	HGRH	HUMAN GROWTH HORMONE Container :One gold gel 3.5 mL. Separate serum from cells ASAP Stability :Ambient: 8 Hours/Refrigerated: 24 Hours / Frozen: 14 Days	83003
ACL	HUM2PR	HUMORAL IMMUNITY II Container :One gold gel 5.0 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86317 x16
ACL	HUMIPR	Humoral Immunity Panel 1 Container :Two Gold Gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 14 Days after separation from cells / Frozen: 6 Months after separation from cells	82784 x3 82787 x4 86317 x16
ACL	HYPNE1	HYPERSENS PNEUMO I Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86331 x4 86606 x2
ACL	HYPNE2	HYPERSENS PNEUMO II Container :Two gold gel 3.5 mL Stability :Ambient: 1 Day / Refrigerated: 1 Week / Frozen: 1 Month	86001 x8 86606 86609 x5
PHC	SPESCN	IFE with Immunoglobulins	86334 82734 x3
ACL	SIFEMX	IFE,SERUM w/IGG,A,M Container :One gold gel 5.0 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days	82784 x3 83883 x2 84155 86334
ACL	IG1234	IgG SUBCLASS 1,2,3,4 Container :One gold gel 3.5 mL Stability :Ambient: 2 Hours / Refrigerated: 8 Days / Frozen: 6 Months	82787 x4
PHC	IGSUBA	IgG SUBCLASSES WITH TOTAL IgG Container :One 2.0 mL (min: 1.0 mL) serum refrigerated Stability :Ambient: 7 Days / Refrigerated: 6 Weeks / Frozen: 12 Weeks	82787 x4 82784
PHC	IGHV	IGHV Mutation Analysis by Sequencing Stability :Ambient: 1 hour: Refrigerated: 7 days: Frozen: Unacceptable	81263
ACL	IMPDES	IMIPRAMINE & METAB Container :One plain red 4.0 mL (also acceptable: lavender (EDTA) 3.0 mL or pink (EDTA) 6.0 mL) Stability :Ambient: 5 Days / Refrigerated: 2 Weeks / Frozen: 6 Months	80335
PHC	IPF	Immature Platelet Fraction Stability :Room temperature: Preferred specimen should be less than 8 hours old. If stored at 4-8°C, allow it to come to room temperature before testing.	85055
ACL	UIFXMX	Immunofixation Electrophoresis, Urine Container :One Urine, random Stability :Refrigerated: 7 Days / Frozen: 14 Days	86335
ACL	IGA	Immunoglobulin A, IgA Container :One gold gel 3.5 mL Stability :Refrigerated: 3 Days / Frozen: 14 Days	82784
ACL	IGEB	Immunoglobulin E, IgE Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 4 Hours serum / Refrigerated: 2 Days serum / Frozen: 14 Days serum	82785
ACL	IGGS4	Immunoglobulin G Subclass 4 Container :One gold gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 14 Days / Frozen: 6 Months	82787
ACL	IGG	Immunoglobulin G, IgG	82784



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One gold gel 3.5 mL Stability :Refrigerated: 3 Days / Frozen: 14 Days	
ACL	IGM	Immunoglobulin M, IgM Container :One gold gel 3.5 mL Stability :Refrigerated: 3 Days / Frozen: 14 Days	82784
ACL	IMMGLB	Immunoglobulin, Quantitative Container :One gold gel 3.5 mL (also acceptable: One gold gel microtainer 600 mcL (clear or amber)) Stability :Refrigerated: 3 Days / Frozen: 14 Days	82784 x3
PHC	ITDNR	Infant Type and Coombs Non-Cord Stability :Specimen stability: Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	86880 86900 86901
PHC	ITD	Infant Type and Coombs- Cord Blood Stability :Specimen stability: Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	86880 86900 86901
ACL	IBDPNL	Inflammatory Bowel Disease Panel Container :One gold gel 3.5 mL Separate from cells ASAP or within 2 hours of collection Stability :Ambient: 48 Hours after separation from cells / Refrigerated: 2 Weeks after separation from cells / Frozen: 1 Year after separation from cells, avoid repeated freeze/thaw cycles	86255 86671 x2
ACL	INFLIX	Infliximab or Biosimilar Activity and Neutralizing Antibody Container :One gold gel 5.0 mL blood Stability :Ambient: 48 hours / Refrigerated: 4 weeks / Frozen: 1 Year	80230 82397
PHC	INFLAB	Influenza A&B, Rapid Antigen Test, Nasopharyngeal Stability :Room temperature: No Refrigerated: 72 hours Frozen: Do Not Freeze	87804 x2
ACL	INHA	INHIBIN A Container :One blood gold gel 5.0 mL (also acceptable: One blood plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	86336
ACL	INHIBB	INHIBIN B Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 72 Hours / Frozen: 1 Month	83520
ACL	INHSCR	INHIBITOR SCREEN Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: 4 Hours / Refrigerated: Unacceptable / Frozen: 30 Days	85610 85611 x2 85730 85732 x2
ACL	INSLN	INSULIN Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 14 Days	83525
ACL	INSABT	INSULIN ANTIBODY Container :One gold gel 3.5 mL (also acceptable: plain red 4.0 mL) Stability :Ambient: 1 Day/ Refrigerated: 7 days / Frozen: 2 months	86337
PHC	IRSP	Insulin Resistance Screen, Plasma and Serum	82947 83525
ACL	ILGFI	Insulin-Like Growth Factor 1 Container :One gold gel 3.5 mL Stability :Ambient: 24 Hours / Refrigerated: 2 Days / Frozen: 14 Days	84305
ACL	GFBP3	Insulin-Like Growth Factor Binding Protein 3 Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 14 Days	82397
ACL	IA2AB	Insulinoma Associated Antibody 2 Container :One gold gel 3.5 mL Stability :Ambient: 1 Day / Refrigerated: 7 Days / Frozen: 2 Month	86341
ACL	INTGAM	INTERFERON GAMMA Container :One gold gel 5.0 mL blood (also acceptable: plain red 4.0 mL blood and green (lithium heparin no gel) 6.0 mL blood) Stability :Ambient: 30 minutes after separation from cells / Refrigerated: Unacceptable / Frozen: 1 Year after separation from cells	83520
ACL	INTLEU	INTERLEUKIN 12 Container :One gold gel 5.0 mL blood (also acceptable: plain red 4.0 mL blood and green (lithium heparin no gel) 6.0	83520



LAB	CODE	TEST NAME	CPT4 CODE
		mL blood) Stability :Ambient: 30 minutes after separation from cells / Refrigerated: Unacceptable / Frozen: 1 Year after separation from cells	
ACL	IL6	INTERLEUKIN 6 Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 6 Hours / Refrigerated: 2 Days / Frozen: 14 Days	83520
ACL	IL2R	Interleukin-2 Receptor (CD25), Soluble (SIL2R) Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 30 Minutes after separation from cells Frozen: 1 Year after separation from cells	83520
ACL	INTRIN	INTRINSIC FACTOR AB Container :One gold gel 3.5 mL (ALSO ACCEPTABLE: PLAIN RED 4.0 ML) Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 1 Month	86340
ACL	BIODIN	IODIDE Container :One royal blue (no additive) 6.0 mL Stability :Ambient: 30 Days / Refrigerated: 30 Days / Frozen: 30 Days	83018
ACL	ICAOUT	IONIZED CALCIUM Container :Special Collection Required Stability :Ambient: 24 Hours / Refrigerated: 3 Days	82330
PHC	IIBC	Iron and Iron-Binding Capacity, Serum Stability :Refrigerated: 7 days Frozen: 6 months	83540 83550
PHC	IRON	Iron, Serum Stability :Refrigerated: 7 days when separated from cells	83540
ACL	ISLET	ISLET CELL ANTIBODY Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86341
PHC	ISALCO	Isopropyl Alcohol Stability :Room temperature: 8 hours Refrigerated: 7 days Frozen: 1 month	82491
ACL	ITRAC	ITRACONAZOLE Container :One plain red 4.0 mL (ALSO ACCEPTABLE: GREEN (LITHIUM HEPARIN NO GEL) 4.0 ML) Stability :Ambient: 15 days / Refrigerated: 15 days / Frozen: 60 days	80299
ACL	JAK2EXON	JAK2 (V617F) Mutation by ddPCR, Qualitative with Reflex to JAK2 Exon 12 Mutation Analysis by PCR JAK2 Exon 12 Container :Whole Blood: 5 mL whole blood. (Min: 1 mL) Bone Marrow: 3 mL bone marrow. (Min: 1 mL) Stability :Ambient: Unacceptable; Refrigerated: 4 days; Frozen: Unacceptable	81270
ACL	JAK2QN	JAK2 QUANTITATIVE Container :One lavender (EDTA) 6.0 mL (Also acceptable: pink (K2EDTA), yellow (ACD), light blue (sodium citrate), green (sodium heparin)) Stability :Ambient: 72 Hours / Refrigerated: 14 Days / Frozen: unacceptable	81270
ACL	JAK2	JAK2 Tyrosine Kinase V617 Mutation by PCR Container :3 ml (Min: 1 ml) whole blood refrigerated or 3 ml (Min: 1 ml) bone marrow refrigerate Stability :Ambient: 72 Hours / Refrigerated: 14 Days / Frozen unacceptable	81270
ACL	JAK2CALR	JAK2 V617 Mutation by ddPCR, Qualitative With Reflex to CALR/MPL Container :Whole Blood: 5 mL whole blood. (Min: 1 mL) Bone Marrow: 3 mL bone marrow. (Min: 1 mL) Stability :Ambient: Unacceptable; Refrigerated: 5 days; Frozen: Unacceptable	81270
ACL	JO1AB	JO 1 ANTIBODY Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	UKLFR	Kappa/Lambda Light Chains, Free With Ratio, Urine Stability :Ambient: 2 hours / Refrigerated: 7 Days / Frozen: Unacceptable	83883 x2
ACL	KAPLAM	KAPPA/LAMBDA RATIO Container :One blood gold gel 3.5 mL Stability :Refrigerated: 1 Week / Frozen: 2 Weeks	83521 x2
PHC	UKET	Ketones, Urine Stability :Refrigerated: 24 hours	82009
PHC	KLBT	Kleihauer-Betke (FetalDex®), Blood Stability :Room temperature: 8 hours Refrigerated: 24 hours	85460



LAB	CODE	TEST NAME	CPT4 CODE
		Frozen: Do Not Freeze	
PHC	KRAS	Kras Mutation Analysis	81275 81403 88381
ACL	LACOS	LACOSAMIDE Container :One plain red 4.0 mL Stability :Ambient: 72 Hours / Refrigerated: 2 Weeks / Frozen: 6 Months	80235
PHC	FLLDH	Lactate Dehydrogenase (LDH), Total, Body Fluid Stability :Room temperature: 3 days Refrigerated: unacceptable Frozen: unacceptable	83615
PHC	LDH	Lactate Dehydrogenase (LDH), Total, Serum Stability :AMBIENT: SERUM: 3 DAYS / REFRIGERATED: unacceptable / FROZEN: unacceptable	83615
PHC	CSFLDH	Lactate Dehydrogenase (LDH), Total, Spinal Fluid Stability :Room temperature: 3 days Refrigerated: unacceptable Frozen: unacceptable	83615
ACL	LACC	LACTATE, VENOUS Container :One gray (sodium fluoride/potassium oxalate) 2.0 mL Collection tube should be filled at least three-fourths, mixed gently, placed on ice and centrifuged immediately. Separate plasma from cells within 30 minutes of collection. Stability :Ambient: 4 Hours / Refrigerated: 7 Days / Frozen: 14 Days	83605
PHC	LACTA	Lactic Acid, Arterial; In-Hospital	83605
PHC	LACTV	Lactic Acid, Venous; In-Hospital	83605
ACL	FLM	LAMELLAR BODY COUNT Container :Amniotic fluid Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: Unacceptable	83664
ACL	LAM	LAMOTRIGINE Container :One plain red 6.0 mL (also acceptable: green (lithium heparin no gel) 6.0 mL and green (sodium heparin no gel) 6.0 mL and lavender (EDTA) 3.0 mL and pink (K2EDTA) 6.0 mL) Stability :Ambient: 7 days / Refrigerated: 1 Week / Frozen: 1 Month	80175
ACL	PBCAP	LEAD, BLOOD/CAP Container :One lavender (EDTA) microtainer 500 mcL (also acceptable: One green (lithium heparin) 1.0 mL microtainer) Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	83655
ACL	LEADB	LEAD, BLOOD/VEN Container :One lavender (EDTA) 3.0 mL (also acceptable: royal/dark blue sodium heparin or EDTA 6.0 mL or green (sodium or lithium heparin) 6.0 mL). Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	83655
ACL	LEGAB	LEGIONELLA AB Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	86713
ACL	LEGEIA	Legionella pneumophila Antigen, EIA Urine Container :STERILE LEAK PROOF CONTAINER Stability :Ambient: 24 Hours / Refrigerated: 14 Days / FROZEN: INDEFINITELY	87899
ACL	LEPTIN	LEPTIN Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 1 Day / Refrigerated: 2 Days / Frozen: 2 Months	83520
ACL	LEPTO	LEPTOSPIRA ANTIBODY Container :One gold gel 3.5 mL (also acceptable: green (sodium or lithium heparin) Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86720
ACL	LEV	LEVETIRACETAM Container :One plain red 6.0 mL Stability :Ambient: 2 Days / Refrigerated: 7 Days / Frozen: 14 Days	80177
ACL	LIDO	LIDOCAINE Container :One plain red 4.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80176
PHC	FLIP	Lipase, Body Fluid Stability :Room temperature: 24 Hours Refrigerated: 7 days	83690
PHC	LIPASE	Lipase, Plasma or Serum Stability :Room temperature: 24 Hours Refrigerated: 7 days	83690



LAB	CODE	TEST NAME	CPT4 CODE
PHC	LIPR	Lipid Panel Stability :Cholesterol: 5 days refrigerated triglycerides: 7 days refrigerated HDL: 5 days refrigerated Direct ldl: 3 day refrigerated if reflexed.	80061
ACL	LAPA2	LIPO ASSOC PHOS A2 Container :One gold gel 3.5 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Ambient: 7 Days /Refrigerated: 28 Days / Frozen: 28 Days	83698
ACL	LIPOA	LIPOPROTEIN (A) Container :One gold gel 3.5 mL (also acceptable: One green (lithium heparin no gel) 4.0 mL or One lavender (EDTA) 3.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 2 days / Frozen: 3 Months	83695
ACL	LIPELE	LIPOPROTEIN ELECTRO Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 1 Day / Refrigerated: 10 Days	80061 83700
PHC	LITH	Lithium, Serum Stability :Room temperature: 24 hours Refrigerated: 7 days Frozen: 6 months	80178
ACL	LKM	LKM ANTIBODY Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Frozen: 14 Days	86376
ACL	LUPPRO	LUPUS ANTICOAGULANT Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85390 85613 85670 85732
PHC	LHB	Luteinizing Hormone (LH), Serum Stability :Room temperature: 8 hours Refrigerated: 2 days	83002
ACL	LYMEWB	Lyme Antibodies IgG & IgM by Western Blot Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 2 weeks / Frozen: 1 year Avoid repeated freeze/thaw cycles; For long term storage, store at -70 degrees C	86617 x2
ACL	CSLYME	Lyme Antibodies, IgG and IgM by Immunoblot (CSF) Container :One Spinal fluid (CSF) Stability :Ambient: 8 Hours / Refrigerated: 14 Days / Frozen: 1 Year	86618
ACL	LYPCR	LYME DISEASE by PCR Container :One blood lavender (EDTA) 3.0 mL (also acceptable: One blood plain red 4.0 mL and One body fluid Synovial Fluid and One Cerebrospinal fluid (CSF)) Stability :Ambient: 8 Hours / Refrigerated: 3 Days / Frozen: 1 Year	87476
ACL	LYME	LYME IGG/M AB SCREEN Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: unacceptable/Refrigerated: 7 Days / Frozen: 14 Days	86618
ACL	LYSO	LYSOZYME Container :One gold gel 3.5 mL Stability :Ambient: unacceptable / Refrigerated: 5 days / Frozen: 1 Months	85549
ACL	MYCSER	M PNEUMONIAE AB IgA Container :One plain red 4.0 MI (ALSO ACCEPTABLE: GOLD GEL 3.5 ML) Stability :Frozen: 120 Days	86738
ACL	MACPRO	MACROPROLACTIN Container :One gold gel 5.0 mL blood. Also acceptable: One green lithium heparin 6.0 mL or One lavender (EDTA) 6.0 mL Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 1 Month	84146 x2
ACL	MAGRBC	MAGNESIUM RBC Container :One royal blue (EDTA) 6.0 mL Stability :Ambient: 7 Days / Refrigerated: 1 month	83735
PHC	MG	Magnesium, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 1 year	83735
ACL	UMGR	MAGNESIUM, URINE Container :One Urine, random Stability :Refrigerated: 7 Days / Frozen: 14 Days	83735
ACL	UMGT	MAGNESIUM, URINE 24HR Container :Urine	83735



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :Refrigerated: 7 Days / Frozen: 14 Days	
ACL	RMALP	MALARIA RAPID SCREEN Container :One lavender (EDTA) 3.0 mL blood Stability :Ambient: 24 Hours / Refrigerated: 24 Hours / Frozen: Unacceptable	
PHC	SS02120	Measles IgM and IgG Stability :Refrigerated: 5 days	86765 x2
PHC	VR01713	Measles Virus PCR Container :UTM Stability :Refrigerated: 48 hours	87798
ACL	MERCB	MERCURY, BLOOD Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	83825
ACL	URHMS	METALS PNL, UR/OCCU Stability :Refrigerated: 5 Days / Frozen: 14 Days Stability based upon creatinine requirements.	82175 82300 82570 83655 83825
ACL	HMSB	METALS PNL1, BLOOD Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82175 82300 83655 83825
ACL	UMET	METANEPHRINE,UR 24HR Container :One Urine Stability :Ambient: unacceptable / Refrigerated: 2 Week / Frozen: 1 Month	83835
ACL	PMET	Metanephrines, Plasma Container :One lavender (EDTA) 6.0 mL (also acceptable: One pink K2EDTA and Green Sodium or Lithium heparin, NO GEL) Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 1 Month	83835
ACL	METHQNT	Methadone and Metabolite Quantitation- Serum/Plasma Container :One plain red 4.0 mL blood Also acceptable: Two gray (sodium fluoride/potassium oxalate) 2.0 mL, One green (sodium heparin) 4.0 mL, One lavender (EDTA) 3.0 mL, or One pink (K2EDTA) 6.0mL Stability :Ambient: 1 Week after separation from cells Refrigerated: 2 Weeks after separation from cells Frozen: 3 Years after separation from cells	80358
ACL	MTHCNF	Methadone and Metabolite, Urine, Quantitative Container :1.0 mL (Min: 0.5 mL) urine refrigerated Stability :Ambient: 1 Week / Refrigerated: 1 Month/ Frozen: 3 Years	80358
PHC	METHA	Methadone, Urine Stability :Refrigerated: 24 Hours	80307
PHC	METHB	Methemoglobin, Blood Stability :Room temperature: 30 minutes	83050
ACL	MRSPCR	Methicillin-Resistant Staphylococcus aureus (MRSA) by Rapid PCR, Screening Container :Eswab Stability :Ambient: 48 Hours / Refrigerated: 5 Days / Frozen: Unacceptable	87641
ACL	METHO	METHOTREXATE Container :One plain red 6.0 mL Stability :Ambient: 4 Hours / Refrigerated: 7 days / Frozen: 2 Weeks	80320
ACL	MMA	METHYLMALONIC ACID Container :One green (sodium or lithium heparin) 6.0 mL (also acceptable: One plain red 4.0 mL and One gold gel 3.5 mL) Stability :Ambient: unacceptable / Refrigerated: 1 week / Frozen: 1 month	83921
ACL	MEX	MEXILETINE Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Ambient: 2 Days / Refrigerated: 5 Days / Frozen: 2 Months	80299
ACL	MALB	MICROALBUMIN, 24HR Container :24 hour Urine Stability :Refrigerated: 3 Days	82043
ACL	MACRE	MICROALBUMIN,RANDOM Container :Urine. Stability :Refrigerated: 3 Days	82043 82570
PHC	MVO2	MIXED VENOUS O2 Stability :Rm Temp: 30 minutes	82810
ACL	MGFLWP	Monoclonal Gammopathy Follow Up Reflex IFESE	83883 x3 84155 84165



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One blood gold gel 3.5 mL Stability :Refrigerated: 1 Week / Frozen: 2 Weeks	
ACL	MGSCRN	Monoclonal Gammopathy Screen Reflex IFESE, QIGS Container :One blood gold gel 3.5 mL Stability :Refrigerated: 1 Week / Frozen: 2 Weeks	83883 x2 84155 84165
PHC	MONO	Mononucleosis Screen Stability :Refrigerated: Whole blood, 24 hours	86308
ACL	MPLCODON515	MPL CODON 515 Container :Lavender (EDTA) or green (sodium heparin) whole blood or bone marrow. Stability :Ambient: 24 hours; Refrigerated: 5 days; Frozen: Unacceptable	81338
ACL	MYLPER	MPO/PR-3 AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	83516 x2
ACL	MRMSSA CX	MRSA/MSSA CULTURE Container :ESWAB Stability :AMBIENT: 24 HOURS / REFRIGERATED: 48 HOURS / FROZEN: UNACCEPTABLE	87081
ACL	MTHFR	MTHFR GENOTYPE Container :One lavender (EDTA) 6.0 mL (also acceptable: One light blue (3.2% sodium citrate) 2.7 mL or 6.0 ml Pink (K3EDTA)) Stability :Ambient: 72 Hours / Refrigerated: 14 Days	81291
ACL	MSPNLB	Multiple Sclerosis Panel with CSF Container :Spinal fluid (CSF) and One gold gel 5.0 mL Stability :Refrigerated: 48 Hours / Frozen: 14 Days	82040 82042 82784 x2 83916 84157
ACL	MUMPIS	MUMPS IGG AB Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 1 Year (avoid repeated freeze/thaw cycles)	86735
ACL	MUMPSM	MUMPS IgM ANTIBODY Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 14 Days / Frozen: 1 Year	86735
PHC	MUMPA	Mumps Virus Antibodies Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 1 Year (avoid repeated freeze/thaw cycles)	86735 x2
ACL	RFMUSK	MuSK ANTIBODY TEST Container :One plain red 4.0 mL (also acceptable: gold gel 3.5 mL) Stability :Ambient: 72 Hours / Refrigerated: 42 Days / Frozen:1 year	83519
ACL	MGPNL2	MYASTHENIA GRAVIS PANEL 2 W/REFLEX Container :TWO PLAIN RED 4.0 ml (ALSO ACCEPTABLE: TWO GOLD GEL 3.5 ml) Stability :Ambient: 72 Hours / Refrigerated: 14 Days / Frozen: 30 Days	83519 x3
ACL	AFBU	MYCOBACTERIAL CULTURE,URINE Container :Refrigerated, in sterile leakproof container. Stability :Ambient: 24 Hours Refrigerated: 7 Days	87015 87116
PHC	MTBADT	Mycobacterium tuberculosis Complex, Molecular Detection, PCR Container :Preferred: Sterile container Acceptable: Biopsy specimen of tissue fixed with formalin in a paraffin block Stability :Unprocessed: Ambient: 3 days: Refrigerated: 1 week: Frozen: 1 month: Processed: Ambient: Unacceptable: Refrigerated: 1 week: Frozen: 1 month	87556 87798
ACL	MAG	MYELIN ASSOC GLYCOPR Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	83516
ACL	MYLPET	MYELOPEROXIDASE AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	83516
ACL	SMYO	MYOGLOBIN Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 8 Hours serum or plasma / Refrigerated: 7 Days serum(plasma,refrigerated: 3 days) / Frozen: 14 Days serum or plasma	83874
ACL	UMYO	MYOGLOBIN, URINE Container :ARUP Standard Transport Tube prefilled with Sodium Carbonate Stability :Ambient: 4 Hours / Refrigerated: 72 Hours / Frozen: 1 Month Stability is based on samples that have been transferred to ARUP Standard Transport Tube prefilled with Sodium Carbonate or pH to 8-9 by adding 10 percent Na2CO3.	83874



LAB	CODE	TEST NAME	CPT4 CODE
PHC	NTPBNP	N-Terminal Pro Brain Natriuretic Peptide Stability :Ambient: 3 Days / Refrigerated: 3 Days	83880
ACL	JNT04	NATIVE JOINT CULTURE Stability :Ambient: 48 Hours	87070 87075 87205
ACL	NEOHSV	Neonatal HSV Screen Stability :Ambient: 8 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	87529 x2
ACL	NSE	Neuron-Specific Enolase, Serum Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: Unacceptable Refrigerated: 1 Week Frozen: 1 Year (avoid repeated freeze/thaw cycles)	86316
PHC	NBSCR	Newborn Screening Panel, Blood Spots	82017 82261 82775 83020 83498 83516 83788 84030 84443 84510
ACL	NICKEL	NICKEL, SERUM Container :One royal blue (no additive) 6.0 mL Stability :Ambient: Indefinitely / Refrigerated: Indefinitely / Frozen: Indefinitely	83885
ACL	NICOU	NICOTINE & METABOLITES, URINE Container :Urine Stability :Ambient: 10 days / Refrigerated: 10 days / Frozen: 8 months	80323
ACL	NICOT	Nicotine and Metabolites, Serum or Plasma, Quantitative Container :Two gold gels 3.5 mL (also acceptable: Two plain red 4.0 mL, two lavender 3.0 mL, one green sodium heparin 6.0 mL) Stability :Ambient: 1 week / Refrigerated: 2 weeks / Frozen: 3 years	80323
ACL	AQUA4 CSF	NMO AQ4 IgG CSF Container :One Spinal fluid (CSF) Stability :Refrigerated: 2 Weeks/Frozen: 1 Year	86255
ACL	NMR	NMR LIPOPROFILE WITH LIPID PANEL Container :One plain red 6.0 mL Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: Unacceptable	80061 83704
ACL	CYTOLOGY	NON-GYNECOLOGICAL CYTOLOGY	VARIES
ACL	NORVRX	NOROVIRUS Stability :Ambient: 1 Hour, no preservative. 2 days in Cary-Blair / Refrigerated: 4 Days in Cary-Blair / Frozen: Unacceptable	87798 x2
ACL	NORTRI	NORTRIPTYLINE Container :(ALSO ACCEPTABLE: LAVENDER (EDTA) 3.0 ML, PINK (EDTA) 6.0 ML) Stability :Ambient: 5 Days / Refrigerated: 2 Weeks / Frozen: 6 Months	80335
PHC	SVCO2	O2 SATURATION CENTRAL VENOUS LINE Stability :Rm Temp: 30 minutes	82810
PHC	OBRH	OB RH Container :Collect: Full, large lavender-top (EDTA) tube. 3.5mL min. Serum gel tube is not acceptable. Fill tube completely. Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	86901
PHC	OBP	Obstetrics Panel Stability :Refer to individual test listings.	80055
PHC	OBPAN	Obstetrics Panel with HIV Stability :Refer to individual test listings.	80081
ACL	IFOB	OCCULT BLD,TUBE TEST Container :Stool Stability :Ambient: 7 Days sample collection device / Refrigerated: 14 Days sample collection device	82274
ACL	OLAN	OLANZAPINE Container :One plain red 4.0 mL (ALSO ACCEPTABLE: LAVENDER (EDTA) 3.0 ML, PINK (K2EDTA) 6.0 ML) SEPARATE SERUM OR PLASMA FROM CELLS WITHIN 2 HOURS OF COLLECTION. Stability :Frozen: 2 Weeks	80342
ACL	OPAT	Opiates Urine Quantitative Container :One urine - random Stability :Ambient: 1 Week / Refrigerated: 1 Month / Frozen: 3 Years	80361 80365
PHC	OPIAT	Opiates, Urine Stability :Refrigerated: 24 Hours Freeze if > 24 Hours	80307



LAB	CODE	TEST NAME	CPT4 CODE
ACL	UORAC	ORGANIC ACIDS SCREEN URINE Container :ONE URINE - RANDOM Stability :REFRIGERATED: 14 DAYS / FROZEN: 416 DAYS	83919
PHC	OSMOL	Osmolality, Serum Stability :Room temperature: 24 hours Refrigerated: 7 days Frozen: 14 days	83930
PHC	UOSMOL	Osmolality, Urine Stability :Room temperature: 24 hours Refrigerated: 7 days Frozen: 14 days	83935
ACL	FRAGSA	OSMOTIC FRAGILITY, RBC Container :Two green (sodium or lithium heparin) 4.0 mL (also acceptable: Two lavender (EDTA) 3.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 72 Hours / Frozen: Unacceptable	85555
ACL	OSTEO	OSTEOCALCIN Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) (LAVENDER (EDTA) 3.0 ML, PINK (K2EDTA) 6.0 ML, GREEN (SODIUM OR LITHIUM HEPARIN) 4.0 ML) Stability :Ambient: 8 hours/ Refrigerated: 72 Hours / Frozen: 3 Months	83937
ACL	OPUR	Ova and Parasite Exam, Body Fluid or Urine Stability :Body Fluid: Ambient: Unacceptable; Refrigerated: 72 hours; Frozen: Unacceptable Urine: Ambient: 72 hours; Refrigerated: 72 hours; Frozen: Unacceptable CSF: Ambient: 72 hours; Refrigerated: Unacceptable; Frozen: Unacceptable	87177 87209
ACL	OVAPF	Ova and Parasite Exam, Fecal (Immunocompromised or Travel History) Container :5g (Minimum: 5g) of stool preserved in SAF. Transport at ambient temperature. Stability :Ambient: 9 Months Refrigerated: 9 Months Frozen: Unacceptable	87177 87209
ACL	UOXAL	OXALATE,UR 24HR Container :One Urine 24 hour urine. Sample must be refrigerated during collection. Stability :Refrigerated: 1 Week / Frozen: 1 Week	83945
ACL	OMHC	OXCARBAZEPINE Container :One plain red 4.0 mL (also acceptable: lavender (EDTA) 3.0 mL and pink (K2EDTA) 6.0 mL). Stability :Ambient: 6 weeks / Refrigerated: 6 Week / Frozen: 3 Month	80183
PHC	OXY	Oxycodone, Urine Stability :Refrigerated: 24 hours Freeze if >24 hours	80305
ACL	PANCEF	PANC. ELASTASE,FECAL Container :One stool Stability :Ambient: 1 week / Refrigerated: 1 Week / Frozen: 1 Year	83520
ACL	HERPAN	Pancreatitis Panel Container :One lavender (EDTA) 3.0 mL blood (also acceptable: ACD yellow tube blood) Stability :Ambient: 72 Hours / Refrigerated: 1 Week / Frozen: Unacceptable	81223 81404 x2 81479
ACL	PAP	Pap Test, Thin Prep Stability :Ambient: 3 Weeks	88175
ACL	PNA	PARANEOPLASTIC AUTOB Container :Two gold gel 5.0 mL (also acceptable: plain red 6.0 mL blood) Stability :Ambient: 48 Hours / Refrigerated: 48 Hours / Frozen: 21 Days	83519 x5 86255 x10
ACL	PARID	PARASITE EXAM Container :70% ALCOHOL (RUBBING ALCOHOL) OR OTHER FIXATIVE IN CONTAINER WITH A TIGHT FITTING LID. Stability :AMBIENT: 1 MONTH (PRESERVED)	87169
PHC	PTHBL	Parathyroid Hormone (PTH), Baseline Stability :Room temperature: 2 days Refrigerated: 3 days Frozen: 6 months	83970
PHC	PTHP	Parathyroid Hormone (PTH), Post Excision Stability :Room temperature: 2 days Refrigerated: 3 days Frozen: 6 months	83970
PHC	PTHW	Parathyroid Hormone (PTH), Whole, Plasma Stability :Room temperature: 2 days Refrigerated: 3 days Frozen: 6 months	83970



LAB	CODE	TEST NAME	CPT4 CODE
PHC	PTHB	Parathyroid Hormone (PTH), with Minerals, Plasma	82310 82565 83970 84100
ACL	PARAB	PARIETAL CELL AB Container :One gold gel 3.5 mL Stability :Refrigerated: 48 Hours / Frozen: 14 Days	86255
ACL	PRVB19	PARVO B19 BY PCR Container :One pink (K2EDTA) 6.0 mL (also acceptable: lavender (EDTA) 3.0 mL and gold gel 3.5 mL and plain red 6.0 mL and amniotic fluid and cerebrospinal fluid and synovial fluid and tissue and FFPE (Formalin fixed, paraffin embedded tissue)) SPECIMEN SOURCE IS MANDATORY Stability :Ambient: 24 Hours (Bone Marrow: 1 Week, Tissue: Unacceptable, FFPE: Indefinitely) / Refrigerated: 5 Days (Bone Marrow: 1 Week, Tissue: Unacceptable, FFPE: Indefinitely) / Frozen: 6 Months (Bone Marrow: 1 Week, Tissue: 6 Months, FFPE: Indefinitely)	87798
ACL	PARVO	PARVOVIRUS B19 AB Container :ONE PINK (K2EDTA) 6.0 ML (ALSO ACCEPTABLE: LAVENDER (EDTA) 3.0 ML AND GOLD GEL 3.5 ML AND PLAIN RED 6.0 ML AND AMNIOTIC FLUID AND CEREBROSPINAL FLUID AND SYNOVIAL FLUID AND TISSUE AND FFPE (FORMALIN FIXED, PARAFFIN EMBEDDED TISSUE)) Stability :AMBIENT: 24 HOURS (BONE MARROW: 1 WEEK, TISSUE: UNACCEPTABLE, FFPE: INDEFINITELY) / REFRIGERATED: 5 DAYS (BONE MARROW: 1 WEEK, TISSUE: UNACCEPTABLE, FFPE: INDEFINITELY) / FROZEN: 6 MONTHS (BONE MARROW: 1 WEEK, TISSUE: 6 MONTHS, FFPE: INDEFINITELY)	86747 x2
ACL	PNCOMP	PEANUT COMPONENTS Container :One gold gel 5.0 mL Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003 x6
PHC	PEPN	Pediatric Environmental Panel Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	86003 x14 82785
PHC	PFOODP	Pediatric Food Panel 10 Stability :Ambient: 48 Hours / Refrigerated: 1 Weeks / Frozen: 1 Year	86003 x10 82785
ACL	PEMPNL	PEMPHIGOID AB PANEL Container :One plain red 6.0 mL (also acceptable: One gold gel 5.0 mL) Stability :Ambient: 1 Week / Refrigerated: 2 Weeks / Frozen: Indefinitely	83516 x2 88346 88350 x3
ACL	PGSPNL	PEMPHIGUS SCREEN Container :One plain red 6.0 mL (also acceptable: One gold gel 5.0 mL) Stability :Ambient: 1 week / Refrigerated: 2 weeks / Frozen: Indefinitely	83516 x2 88346 88350
ACL	MALSM	PERIPHERAL BLOOD PARASITES Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 24 Hours. If longer than 24 hours before shipment, please make 2 thick and 4 thin slides and send unstained.	87207
PHC	RSMR	Peripheral Blood Smear, Path Review Stability :Room temperature: >24 hours Refrigerated: 36 hours: in some cases 48 hours - specimen dependent (to be determined by testing personnel) Frozen: Not acceptable	85060
ACL	PF4IGR	PF4 Hep IGG Reflex Container :Two blood plain red 6.0 mL Clot activator blood tube acceptable, Gel blood tubes should not be used. Stability :Ambient: 4 Hours / Refrigerated: 48 Hours / Frozen: 30 Days	86022
PHC	PHWLB	pH, Blood	82800
PHC	FPH	pH, Body Fluid	81003
PHC	UPH2	pH, Urine	81003
ACL	UPHEN	Phencyclidine (PCP), Urine, Quantitative Container :urine - random Stability :Ambient: 1 week / Refrigerated: 14 Days / Frozen: 1 year	83992
PHC	PHENO	Phenobarbital, Serum Stability :Refrigerated: 24 Hours	80184
ACL	PHEALA	PHENYLALANINE + TYRO Container :One light green (lithium heparin gel) 4.5 mL (also acceptable: One green (lithium heparin no gel) 6.0 mL and One green (sodium heparin no gel) 6.0 mL) Stability :Refrigerated: 24 Hours / Frozen: 2 Weeks	84030 84510
PHC	PHNY	Phenytoin, Serum Stability :Refrigerated: 1 day	80185
ACL	UPHOST	PHOS, URINE 24 HR	84105



LAB	CODE	TEST NAME	CPT4 CODE
		Container :Urine Stability :Ambient: 3 Days / Refrigerated: 7 Days / Frozen: 2 Months	
ACL	PETH	Phosphatidylethanol (PEth), Whole Blood, Quantitative Container :lavender (EDTA) 3.0 mL or pink (K2EDTA) or green (lithium heparin no gel) 4.0 mL or gray (sodium fluoride/potassium oxalate) 2.0 mL Stability :Refrigerated: 2 Weeks Frozen: 1 Month Ambient: Unacceptable	80321
ACL	PLA2R	Phospholipase A2 Receptor (PLA2R) Antibody, IgG with Reflex to Titer Container :1.0 ml (Min: 0.2 ml) serum refrigerated Stability :Refrigerated: 2 Weeks Frozen: 1 Year	86255
PHC	PHOS	Phosphorus, Plasma or Serum Stability :Room temperature: 2 day Refrigerated: 4 days Frozen: 3 months	84100
ACL	UPHOSR	PHOSPHORUS,URINE Container :Urine Stability :Ambient: unacceptable / Refrigerated: 7 Days / Frozen: 2 Months	84105
ACL	PINS	PINWORM EXAM Container :SPECIAL PINWORM SPECIMEN COLLECTION DEVICE. Stability :AMBIENT: 7 DAYS	87172
ACL	PCMX	PLASMA CELL PROFILE Container :Two gold gel 3.5 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days	82784 x3 83521 x2 84155 84165 x2 86334
ACL	PLASM	PLASMINOGEN ACTIVITY Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: 4 Hours/ Frozen: 30 Days	85420
ACL	PLTABS	PLATELET AB SCREEN Container :Two plain red 6.0 mL Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: Samples older than 7 days are acceptable if they have been frozen	86022
PHC	PLTAGG	PLATELET AGGREGATION	
PHC	PLTDI	Platelet Antibody, Drug Induced Stability :Room temperature: N/A Refrigerated: 7 days Frozen: > 7 days	86022
ACL	PLTAUT	Platelet Autoantibodies Container :Five yellow (ACD A ONLY) 8.5 mL Stability :Refrigerated: 5 Days	86022 86023
PHC	PLTA	Platelet Count, Blood	85595
PHC	PLFS	Platelet Function Screen, Blood Stability :Room temperature: 4 hours Refrigerated: Do Not Refrigerate Frozen: Do Not Freeze	85576
ACL	P2Y12	PLATELET P2Y12 Stability :Ambient: 4 Hours	85576
ACL	PML	PML RARA Translocation Container :One lavender (EDTA) 6.0 mL Stability :Ambient: 1 Hour / Refrigerated: 48 Hours / Frozen: Unacceptable	81315
ACL	PDFA	Pneumocystis Direct Fluorescent Antibody Container :INDUCED SPUTUM, BRONCHIAL WASH, BRONCHOALVEOLAR LAVAGE, TRACHEAL ASPIRATE. BRONCHIAL BRUSHINGS OR LUNG TISSUE MAY ALSO BE SUBMITTED IN WISCONSIN. Stability :AMBIENT: 2 HOURS / REFRIGERATED: 7 DAYS	87821
ACL	PNHMK	PNH ANALYSIS Container :Three lavender (EDTA) 3.0 mL. Collect one lavender (EDTA) 3.0 mL for the CBC order. Stability :Ambient: 48 Hours / Refrigerated: 48 Hours / Frozen: Unacceptable	88184 88185 x8 88188
ACL	UPOR	PORPHYRINS UR FRACT Container :One urine - 24 hour Stability :Refrigerated: 4 Days / Frozen: 1 Month	84120



LAB	CODE	TEST NAME	CPT4 CODE
ACL	PORTOT	PORPHYRINS, SERUM Container :One plain red 6.0 mL Stability :Ambient: Unacceptable / Refrigerated: 4 Days / Frozen: 1 Month	84311
ACL	POSACN	POSACONAZOLE Container :One plain red 6.0 mL (also acceptable: One green (lithium heparin no gel) 4.0 mL) Stability :Ambient: 15 Days/Refrigerated: 15 Days/Frozen: 60 Days	80187
PHC	K	Potassium, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 1 year	84132
PHC	UKR	Potassium, Random Urine Stability :Refrigerated: 24 hours	84133
ACL	UKT	POTASSIUM, URINE 24HR Container :Urine Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 14 Days	84133
ACL	PAB	PREALBUMIN Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 8 hours/ Refrigerated: 3 Days / Frozen: 14 Days	84134
PHC	PREGS	Pregnancy Test (Qualitative), Chorionic Gonadotropin (hCG), Serum Stability :Refrigerated: 3 days Frozen: 3 months	84703
PHC	PREGU	Pregnancy Test (Qualitative), Chorionic Gonadotropin (hCG), Urine Stability :Refrigerated: 3 days	81025
PHC	HCGQN	Pregnancy Test (Quantitative), Chorionic Gonadotropin (hCG), Beta-Subunit, Plasma or Serum Stability :Room temperature: 24 hours Refrigerated: 7 days Frozen: 1 year	84702
ACL	PREGN	PREGNENOLONE Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL and One lavender (EDTA) 3.0 mL and One green (sodium or lithium heparin) 4.0 MI Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 6 Months after separation of cells	84140
ACL	FISHRP	Prenatal FISH with Reflex Maternal Contamination Hold Stability :AMBIENT: 24 HOURS NOTE: SEND ALL SAMPLES TO THE CYTOGENETIC LABORATORY FOR VIABILITY DETERMINATION SINCE SPECIMEN SAMPLES ARE IRRETRIEVABLE.	88291 88271 x5 88275
ACL	PRIM	PRIMIDONE/PHENOBARB Container :One plain red 4.0 mL Stability :Ambient: 2 Days / Refrigerated: 7 Days / Frozen: 14 Days	80184 80188
ACL	PANAPA	PROCAINAMIDE/NAPA Container :One plain red 4.0 mL (also acceptable: One green (sodium heparin no gel) 4.0 mL) Stability :Ambient: 1 Week / Refrigerated: 2 Weeks / Frozen: 6 Months	80192
PHC	PCT	Procalcitonin Stability :Room temperature: 24 Hours Refrigerated: 48 hours, if separated from the clot Frozen: 3 months (1 freeze/thaw cycle is acceptable)	84145
PHC	PRGST	Progesterone, Serum Stability :Room temperature: 8 hours Refrigerated: 1 days Frozen: 14 days	84144
ACL	PROINS	PROINSULIN, INTACT Container :One gold gel 3.5 mL (also acceptable: One lavender (EDTA) 3.0 mL or plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 48 Hours / Frozen: 2 Months (avoid repeated freeze/thaw cycles)	84206
ACL	PROLIN	PROLACTIN Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	84146
ACL	PROLM	PROLACTIN MACROADENOMA Container :One gold gel 3.5 mL (also acceptable: One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 8 Hours/ Refrigerated: 48 Hours / Frozen: 3 Months	84146
ACL	PPXCNF	Propoxyphene and Metabolite, Urine Stability :Ambient: 72 Hours / Refrigerated: 1 Week / Frozen: 3 Months	80367
ACL	PCA3TM	Prostate Cancer Biomarker Stability :Ambient: 5 Days / Refrigerated: 5 Days / Frozen: 3 Months	81313



LAB	CODE	TEST NAME	CPT4 CODE
PHC	PSAD	Prostate-Specific Antigen (PSA), Diagnostic, Plasma or Serum Stability :Refrigerated: 3 day	84153
PHC	PSASCR	Prostate-Specific Antigen (PSA), Screening, Plasma or Serum Stability :Refrigerated: 3 day	84153
PHC	PSAS	Prostate-Specific Antigen (PSA), Total and Free, Plasma or Serum Stability :Refrigerated: 8 hours	84153 84154
PHC	PSARF	Prostate-Specific Antigen (PSA), Total with Reflex PSA Free, Plasma or Serum Stability :Refrigerated: 8 hours	84153
ACL	JNT14	PROSTHESIS JOINT CULTURE Stability :Ambient: 48 Hours	87070 87075 87205
ACL	UPETMX	PROT ELECT,URIN 24HR Container :One urine - 24 hour Stability :Refrigerated: 5 Days/ Frozen: 14 days	84156 84166
ACL	UPERMX	PROT ELECT,URIN RAND Container :One urine - random Stability :Refrigerated: 5 Days/ Frozen: 14 days	84156 84166
ACL	PROTC	Protein C Activity, Reflex Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85303
ACL	SPRMX	Protein Electrophoresis Reflex Immunofixation Electrophoresis/IgG Quantitative Container :One gold gel 3.5 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days	84155 84165
ACL	SPEMX	Protein Electrophoresis, Serum Container :One gold gel 5.0 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days	84155 84165
ACL	PROTS	Protein S Activity, Reflex Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85306
PHC	UPROR	Protein, Random, Urine Stability :Refrigerated: 3 days	84156
PHC	FLPRO	Protein, Total, Body Fluid Stability :Refrigerated: 3 days	84157
PHC	TP	Protein, Total, Plasma or Serum Stability :Refrigerated: 3 days	84155
PHC	CSFP	Protein, Total, Spinal Fluid Stability :Refrigerated: 3 days	84157
PHC	UPROT	PROTEIN, URINE 24 HR Stability :Ambient: 4 Hours / Refrigerated: 3 Days	84156
PHC	UPROCR	Protein/Creatinine Ratio, Random Urine Stability :Refrigerated: 3 days	84156 82570
ACL	PROGEN	PROTHROMBIN GENOTYPE Container :One lavender (EDTA) 6.0 mL (also acceptable: One light blue (3.2% sodium citrate) 2.7 mL or 6.0 mL Pink (K3EDTA)) Stability :Ambient: 72 Hours / Refrigerated: 14 Days	81240
ACL	CHOLIN	PSEUDOCHOLINESTERASE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 7 Days	82480
PHC	PT	PT (Prothrombin Time), Plasma Stability :Room temperature: 24 hours Refrigerated: Do Not Refrigerate Frozen: 6 months	85610
ACL	PTHRP	PTH RELATED PEPTIDE Container :One Protease Inhibitor tube Stability :Frozen: 90 Days	82542
PHC	PTT	PTT (Activated Partial Thromboplastin Time), Plasma Stability :Room temperature: 2 hour if patient is on Heparin, 4 hours if not on Heparin Refrigerated: Do Not Refrigerate Frozen: 6 months	85730
ACL	AFP4MA	QUAD MATERNAL SCREEN	82105 82677 84702



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	86336
ACL	QUANTP	QUANTIFERON TB PLUS Container :TB Collection Tube Set (4 - 1mL Test Specific Tubes) Stability :Ambient: 14 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	86480
PHC	CDDA	Quantitative Duodenal Aspirate Culture Container :Sterile Syringe with air removed (remove needle and cap syringe), Red Top Tube Stability :Ambient: 24 Hours	87071 87073 87102
PHC	CQNT	Quantitative Tissue Culture Stability :Ambient: 24 Hours	87071
PHC	SARSRAP AG	Quidel Sofia SARS Rapid Antigen FIA Container :Preferred: Nasal Swab Sample: Use the Puritan nylon flocked swab supplied by the Laboratory. Also acceptable: Nasopharyngeal Swab Sample: Use a nylon flocked nasopharyngeal swab provided by the Laboratory. Stability :Ambient: 48 hours / Refrigerated: 48 hours	87426
ACL	QUINID	QUINIDINE Container :One plain red 4.0 mL (also acceptable: One green (sodium heparin no gel) 4.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 1 Week / Frozen: 2 Weeks	80194
ACL	RABIES AB	Rabies Antibody Screen by RFFIT, Serum Container :Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 1 mL) Stability :Ambient: 1 Week Refrigerated: 2 Weeks Frozen: 1 Month	86382
ACL	RMEPNL	RAPID MENINGITIS/ENCEPHALITIS CSF Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 3 Months	87483
PHC	RFP	Renal Function Panel	80069
ACL	RENIN	RENIN ACTIVITY Container :Two lavender (EDTA) 3.0 mL Stability :Ambient: 6 Hours after separated from cells / Refrigerated: Unacceptable / Frozen: 1 Month	84244
ACL	RPPNL	Respiratory Pathogen Panel Container :The ACL Molecular Pathology Department has validated the assay for use with nasopharyngeal swabs, nasopharyngeal aspirates, bronchial lavage, and bronchial washings in a sterile tube (minimum volume: 2 mL) Stability :Refrigerated: 3 Days / Frozen: 1 Months @ -70 deg C	87633
PHC	RSVRAP	Respiratory Syncytial Virus (RSV) Rapid Antigen, Nasopharyngeal Stability :Room temperature: No Refrigerated: 72 hours Frozen: Do Not Freeze	87807
PHC	RESP4	RESPIRATORY VIRUS 4-PLEX TEST BY PCR Container :VTM/UTM Stability :Refrigerated: 3 days/ Frozen (-70): 1 Month	0241U
PHC	ARETIC	Reticulocyte Count, Blood Stability :Room temperature: 24 hours Refrigerated: 48 hours Frozen: Do Not Freeze	85046
PHC	RHE	Reticulocyte Hemoglobin Stability :Room temperature: 12 Hours Refrigerated: 24 Hours Frozen: Reject	85046
ACL	RHF	RHEUMATOID FACTOR Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 7 Days	86431
ACL	RNPAB	RIBONUCLEOPROTEIN AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	RIBSOM	RIBOSOMAL P PROTEIN Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	ROCKY	Rickettsia rickettsii IgG and IgM Antibodies Container :One gold gel 5.0 mL blood Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year avoid repeated freeze/thaw cycles	86757 x2
ACL	TYPHUS	Rickettsia Typhi IgG & IgM Antibodies	86757 x2



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One gold gel 3.5 mL Remove serum from cells within 2 hours of collection. Stability :Ambient: 48 Hours / Refrigerated: 2 Weeks / Frozen: 1 Year	
ACL	RNAIII	RNA POLYMERASE 3 IgG Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	83516
ACL	ROTA VI	ROTAVIRUS Ag Container :Transport stool in a leak-proof container. Stability :Ambient: 4 Hours / Refrigerated: 24 Hours / Frozen: 1 Week	87425
ACL	RPRQ	RPR QUANT Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL or lavender (EDTA) 3.0 mL) Stability :Refrigerated: 5 Days or plasma 2 days / Frozen: 14 Days	86592
ACL	RUBIS	RUBELLA ANTIBODY IGG Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86762
ACL	MEASL	RUBEOLA IGG (IMMUNE) Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86765
ACL	MEAMAB	RUBEOLA IGM AB Container :One gold gel 3.5 mL Stability :Ambient: Unacceptable / Refrigerated: 48 hours / Frozen: 1 Year	86765
PHC	SAL	Salicylate, Serum Stability :Refrigerated: 7 days	80329
ACL	COVIGG	SARS-CoV-2 IgG Container :One gold gel 5.0 mL (also acceptable: gold gel 3.5 mL) Stability :Ambient: 2 Days / Refrigerated: 7 Days / Frozen: 14 Days	86769
ACL	COV FLU PCR	SARS-COV-2/Influenza by PCR Stability :Ambient: Hologic Direct Load Tube (Huckleberry): 6 days Refrigerated: Hologic Direct Load Tube (Huckleberry) : 6 days Frozen: Unacceptable Refrigerated: UTM and Eswab: 7 days	87636
ACL	SCL70	SCLERODERMA AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	SEL	SELENIUM, PLASMA Container :One royal blue (EDTA) 6.0 mL (also acceptable: One royal blue (no additive) 6.0 mL) Stability :Ambient: Indefinitely / Refrigerated: Indefinitely / Frozen: Indefinitely	84255
PHC	SEMEN	Semen Analysis, Complete (Fertility Work-Up)	89320
PHC	SEMPV	Semen Analysis, Post Vasectomy	89310
ACL	SPR3	SERINE PROTEASE 3 Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	83516
ACL	SER	SEROTONIN, SERUM Container :One gold gel 3.5 MI Stability :Refrigerated: 24 Hours / Frozen: 30 Days	84260
ACL	VISCTY	SERUM VISCOSITY Container :One plain red 6.0 mL Stability :Ambient: 8 Hours / Refrigerated: 4 Days / Frozen: 4 Days	85810
ACL	SHB	Sex Hormone Binding Globulin Container :One blood gold gel 3.5 mL Stability :Ambient: 4 Hours / Refrigerated: 6 Days / Frozen: 14 Days	84270
ACL	HGBS	Sickle Cell Anemia Screen without Reflex Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 7 Days / Refrigerated: 14 Days	85660
ACL	SINCLS	Sinus and Nasopharynx Culture with Gram Stain Container :ESWAB Stability :AMBIENT: 48 HOURS	87070 87205
ACL	SIROLR	SIROLIMUS Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80195



LAB	CODE	TEST NAME	CPT4 CODE
ACL	SSARO	SJOGRENS SYNDROME A Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	SSB	SJOGRENS SYNDROME B Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	SMRNP	SM/RNP ANTIBODY IgG Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	SMAAB	SMITH AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86235
ACL	SMAB	Smooth Muscle Antibody (F Actin) (SMA) Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Frozen: 14 Days	83516
PHC	NA	Sodium, Plasma or Serum Stability :Room temperature: 1 week Refrigerated: 1 week Frozen: 1 year	84295
PHC	UNAR	Sodium, Random, Urine Stability :Refrigerated: 1 day	84300
ACL	UNAT	SODIUM, URINE 24HR Container :Urine Stability :Refrigerated: 7 Days / Frozen: 14 Days	84300
PHC	SOFIA FLUSAR	Sofia flu & SARS Antigen FIA Container :Preferred: Nasal Swab Sample: Use the Puritan nylon flocked swab supplied by the Laboratory. Also acceptable: Nasopharyngeal Swab Sample: Use a nylon flocked nasopharyngeal swab provided by the Laboratory. Stability :Ambient: 48 hours / Refrigerated: 48 hours	87428
ACL	SOTA	SOTALOL Container :One plain red 4.0 mL Stability :Ambient: 8 Hours / Refrigerated: 1 Month / Frozen: 1 Month	80299
PHC	SPEFUP	SPE with reflex IFE	81455 84165
PHC	USPG	Specific Gravity, Urine	81003
PHC	CSFC	Spinal Fluid Analysis	82945 84157 89051
ACL	SMACNA	Spinal Muscular Atrophy (SMA) Copy Number Analysis Stability :Ambient: 72 Hours / Refrigerated: 2 Week / Frozen: Unacceptable	81401
ACL	RTCS	Sputum Culture and Smear Container :SPUTUMS IN STERILE CONTAINER Stability :REFRIGERATED: 3 DAYS	87070 87205
ACL	MRMSSA PCR	Staphylococcus Aureus Methicillin Sensitive (MSSA) and Methicillin Resistant (MRSA) PCR Container :ESwab from both nasal and groin sources Stability :Ambient: 48 Hours Refrigerated: 5 Days Frozen: Unacceptable	87641 87640
ACL	SSAT	Stonerisk Diagnostic Profile Container :One urine - 24 hour Stone risk container Stability :Ambient: acceptable for specimen vials in transport box / Frozen: unacceptable	81050 82140 82340 82507 82570 83735 83945 83986 84105 84133 84300 84392 84560
ACL	FATQLT	STOOL FOR FAT STAIN Container :Stool Stability :Ambient: 6 Hours / Refrigerated: 24 Hours / Frozen: 1 Month	89125
ACL	JCVIDX	Stratify JCV Antibody and Index with Reflex to Inhibition Assay Container :One gold gel 3.5 mL (also acceptable: One lavender (EDTA) 3.0 mL or plain red 3.0 mL) Stability :Ambient: 7 Days / Refrigerated: 2 Weeks / Frozen: 90 Days	86711
ACL	STREPU	Strep Pneumoniae Antigen Detection Container :10.0 ML (MIN: 2.0 ML) URINE AMBIENT IF NOT PERFORMED ON SITE, SHOULD BE TRANSPORTED REFRIGERATED Stability :AMBIENT: 24 HOURS / REFRIGERATED: 14 DAYS / FROZEN: 14 DAYS	87899



LAB	CODE	TEST NAME	CPT4 CODE
PHC	RAPSTC	Strep Screen, Rapid with Reflex Culture, Throat Stability :Room temperature: 72 hours Refrigerated: 72 hours Frozen: Do Not Freeze	87880
PHC	RAPSTR	Strep Screen, Rapid, Throat Stability :Room temperature: 72 hours Refrigerated: 72 hours Frozen: Do Not Freeze	87880
PHC	STRABA	Streptococcal Antibodies Profile Container :Preferred: Gold Top Acceptable: Red Top Stability :Refrigerated: 8 Days / Frozen: 14 Days for ASO. Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year for dnase	86060 86215
ACL	STTH	Streptococcus Group A (Streptococcus pyogenes), Bacterial Culture Container :THROAT SWAB SUBMITTED IN ESWAB [ALSO ACCEPTABLE: THROAT SUBMITTED ON AEROBIC CULTURE SWAB (RAYON OR DACRON SWAB WITH A PLASTIC SHAFT)] ESWAB SPECIMEN COLLECTION DEVICE COLLECT SPECIMEN FROM THE TONSILS, TONSILLAR CRYPT, AND/OR BACK OF THROAT. AVOID TOUCHING TEETH, GUMS, TONGUE OR CHEEK SURFACES. PERIANAL, RECTAL, VAGINAL OR VAGINAL/RECTAL ESWAB MAY ALSO BE COLLECTED IF STREPTOCOCCUS PERIANAL OR VAGINAL DERMATITIS IS SUSPECTED. Stability :AMBIENT: 48 HOURS	87081
ACL	POLABR	Streptococcus Pneumoniae AB, IgG (14 Serotypes) Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid repeated freeze/thaw cycles	86317 x14
ACL	STRP23	Streptococcus pneumoniae Antibodies, IgG Container :One Gold Gel 3.5 mL Stability :Ambient: 48 Hours (After separation from cells) / Refrigerated: 2 Weeks (After separation from cells) / Frozen: 1 Year (After separation from cells)	86317 x23
ACL	STRSER	STRONGYLOIDES AB IgG Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86682
ACL	SULFON	Sulfonylurea Hypoglycemic Panel Container :One plain red 4.0 mL (also acceptable: Two gray sodium fluoride/potassium oxalate 2.0 mL) Stability :Ambient: 2 Days / Refrigerated: 28 Days / Frozen: 24 Months	80377
ACL	SWOCN	SwabOne Candida/Trichomonas Panel by NAA Container :APTIMA Multi-Test Swab Specimen Collection Kit Stability :Ambient: 30 Days Refrigerated: 30 Days	87481 x2 87661
ACL	SWOMG	SwabOne Mycoplasma Genitalium Container :ONE APTIMA UNISEX OR MULTI-TEST SPECIMEN COLLECTION KIT Stability :Ambient: 7 Days / Refrigerated: 7 Days / Frozen: 7 Days	87563
ACL	SWOPNL	SwabOne Vaginitis Panel Container :ONE APTIMA UNISEX OR MULTI-TEST SPECIMEN COLLECTION KIT Stability :AMBIENT: 48 HOURS / REFRIGERATED: 7 DAYS	87481 x2 87563 87661 87801
PHC	SYNFLA	Synovial Fluid Analysis	89051 89060
PHC	SYNAH	Synovial Fluid Cell Count and Diff Stability :Ambient 2 hours; Refrigerator 24 hours in EDTA	89051
ACL	MHATP	T PALLIDUM IgG TPPA Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Refrigerated: 5 Days / Frozen: 14 Days	86780
ACL	TRVAG	T VAGINALIS NAA Stability :Ambient: 30 days urines (in APTIMA urine transport tube), 60 days swab (in APTIMA collection kit) / Refrigerated: 30 days urines (in APTIMA urine transport tube), 60 days swab (in APTIMA collection kit)	87661
ACL	TOXPCR	T. GONDII BY PCR Container :One lavender (EDTA) 3.0 mL (also acceptable: pink (K2EDTA) and gold gel 3.5 mL and red 6.0 mL and amniotic fluid and cerebrospinal fluid and ocular fluid and tissue) Stability :Ambient: 8 Hours (Tissue: Unacceptable) / Refrigerated: 5 Days (Tissue: Unacceptable) / Frozen: 3 Months (Tissue: 3 Months)	87798
ACL	SYPT	T. Pallidum Total IgG/IgM, Reverse Syphilis Screen Algorithm Container :One gold gel 5.0 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days	86780



LAB	CODE	TEST NAME	CPT4 CODE
ACL	T3U	T3 UPTAKE Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	84479
ACL	T3REV	T3, REVERSE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One lavender (EDTA) 3.0 mL) Stability :Ambient: 1 Day / Refrigerated: 1 Week / Frozen: 3 Months	84482
ACL	T4	T4 (Thyroxine) Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	84436
PHC	FT4	T4 (Thyroxine), Free, Plasma or Serum Stability :Refrigerated: 7 days	84439
ACL	FK506S	TACROLIMUS Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80197
PHC	FK506	Tacrolimus (FK506/Prograf®), Blood, Draw Only Stability :Room temperature: 8 days Refrigerated: 8 days Frozen: 30 days	80197
ACL	TCELLB	TCELL T4/T8 RATIO Container :One green (sodium heparin no gel) 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Ambient: 36 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	86355 86359 86360 86357
ACL	TSFFR	Testosterone Free, Females or Children Container :One gold gel 5.0 mL (also acceptable: green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 6 Months	84402
ACL	TSSHFR	Testosterone, Bioavailable & Sex Binding Globulin(w/Total Testost),Females or Children Container :One gold gel 5.0 mL (also acceptable: green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 6 Months	84270 84403
ACL	TESTFR	Testosterone, Female or Children Container :One gold gel 5.0 mL (also acceptable: green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 6 Months	84403
ACL	TSTFFR	Testosterone, Free and Total(Includes Sex Hormone Binding Globulin), Female or Children Container :One gold gel 5.0 mL (also acceptable: green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 6 Months	84270 84403
ACL	TSTFPB	Testosterone, Free, Total, Bioavailable & SHBG, Adult Male Container :One gold gel 3.5 mL blood Stability :Ambient: 4 Hours After separation from cells / Refrigerated: 6 days After separation from cells / Frozen: 14 days After separation from cells	84270 84403
ACL	TESTG	Testosterone, Free/Total, Adult Males by ED/LC-MS/MS Container :One gold gel 3.5 mL (also acceptable: One green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 2 Months	84402 84403
ACL	TESTO	TESTOSTERONE, MALE Container :One gold gel 3.5 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 4 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	84403
ACL	TETTOX	TETANUS ANTIBODY IgG Container :One gold gel 3.5 mL Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year	86317
ACL	UTETCA	THC Metabolite, Urine, Quantitative Container :One urine - random Stability :Ambient: 1 Week / Refrigerated:1 Month / Frozen: 1 Month	80349
PHC	THEOPH	Theophylline, Serum Stability :Refrigerated: 1 day	80198
ACL	THIOMT	Thiopurine Metabolites Container :Two lavender (EDTA) 3.0 mL (also acceptable: pink (EDTA) 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 5 Days / Frozen: unacceptable	80299
ACL	TPMTRB	Thiopurine Methyltransferase RBC Container :Two lavender (EDTA) 3.0 mL (also acceptable: pink (EDTA) 6.0 mL and Two green (sodium heparin no gel) 4.0 mL and Two green (lithium heparin no gel) 4.0 mL) Stability :Ambient: 3 Hours / Refrigerated: 6 Days / Frozen: unacceptable	84433
ACL	THTCL	Throat, Bacterial Culture Container :THROAT SAMPLES SHOULD BE COLLECTED ON THE WHITE CAP ESWAB. NASOPHARYNGEAL SAMPLE SHOULD BE COLLECTED ON THE GREEN CAP (MINI-TIP) ESWAB.	87070



LAB	CODE	TEST NAME	CPT4 CODE
		Stability :AMBIENT: 48 HOURS	
PHC	TT	Thrombin Time, Plasma Stability :Room temperature: 4 hours Refrigerated: Do Not Refrigerate Frozen: 6 months	85670
PHC	THBHRA	Thrombosis Panel, High Risk	86147 x3 86146 x3 81240 81241 85613 85670 85732 83090 85303 85305 85300
PHC	THBLRA	Thrombosis Panel, Low Risk	86147 x3 86146 x3 81240 81241 85613 85670 85732
ACL	ATHYR	Thyroglobulin Antibody Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 8 Hours/ Refrigerated: 48 Hours / Frozen: 14 Days	86800
ACL	THYRO	THYROGLOBULIN REFLEX Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	84432 86800
ACL	TPO	THYROID PEROXIDASE A Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 14 Days	86376
ACL	TSG	THYROID STIMULATE IG Container :One gold gel 3.5 mL Stability :Ambient: 24 Hours / Refrigerated: 6 Days / Frozen: 3 Months	84445
PHC	TSHR	Thyroid-Stimulating Hormone (TSH), Ultra-Sensitive with Reflex Free T4, Plasma or Serum Stability :Refrigerated: 7 days	84443
PHC	HSTSH	Thyroid-Stimulating Hormone (TSH), Ultra-Sensitive, Plasma or Serum Stability :Refrigerated: 7 days	84443
ACL	TBG	THYROXINE BIND GLOB Container :One gold gel 5.0 mL Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Month	84442
PHC	TOBRP	Tobramycin, Serum-Peak Stability :Refrigerated: 1 days	80200
PHC	TOBRR	Tobramycin, Serum-Random Stability :Refrigerated: 1 days	80200
PHC	TOBRT	Tobramycin, Serum-Trough Stability :Refrigerated: 1 days	80200
ACL	TOPIRB	TOPIRAMATE Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL and One green (sodium or lithium heparin) 6.0 mL) Stability :Ambient: unacceptable / Refrigerated: 1 Week / Frozen: 4 Weeks	80201
ACL	TT3	TOTAL T3 Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 5 Days / Frozen: 14 Days	84480
ACL	TOXIG	TOXOPLASMA IgG AB Container :One gold gel 5.0 mL Stability :Refrigerated: 7 Days / Frozen: 14 Days	86777
PHC	TOXOPL	Toxoplasma IgG and IgM Antibodies	86777 86778
ACL	TOXMAB	TOXOPLASMA IGM Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 14 Days	86778
ACL	TPMTGN	TPMT Genotyping 4 Variants Container :One lavender (EDTA) 3.0 mL (also acceptable:lavender (EDTA) 6.0 mL or ACD yellow tube) Stability :Ambient: 72 Hours / Refrigerated: 1 Weeks / Frozen: 1 Month	81306 81335
ACL	TRAMUSR	Tramadol, Urine Screen with Reflex to Quantitation Stability :Ambient: 1 Week Refrigerated: 1 Month Frozen: 1 Month	80307
ACL	TRANS	TRANSFERRIN	84466



LAB	CODE	TEST NAME	CPT4 CODE
		Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	
ACL	TGFBPL	Transforming Growth Factor Beta, Plasma Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 30 Minutes - after separation from cells Refrigerated: Unacceptable Frozen: 60 Days - after separation from cells	83520
ACL	TGFBSE	Transforming Growth Factor Beta, Serum Container :One gold gel 3.5 mL (also acceptable: plain red 6.0 mL) Stability :Ambient: 30 Minutes - after separation from cells Refrigerated: Unacceptable Frozen: 60 Days - after separation from cells	83520
ACL	TRAZ	TRAZODONE Container :One plain red 4.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Separate serum or plasma from cells within 2 hours of collection. Stability :Ambient: 24 Hours - after separation from cells Refrigerated: 2 Weeks - after separation from cells Frozen: 4 Months - after separation from cells	80338
PHC	TRICHO	Trichomonas Rapid Test Stability :Room temperature:24 Hours Refrigerated: 36 Hours Frozen: unacceptable	87808
PHC	TRIG	Triglycerides, Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 7 days at WMH Pewaukee, Oconomowoc and Mukwonago Lab refrigerated stability is 2 Days	84478
PHC	FTG	Triglycerides, Pleural or Peritoneal Fluid Stability :Refrigerated: 2 days	84478
ACL	TRYP	TRYPTASE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 72 hours / Frozen: 1 Month	83520
ACL	TRAB	TSH RECEPTOR AB Container :One gold gel 3.5 mL Stability :Ambient: 24 Hours / Refrigerated: 6 Days / Frozen: 12 Months	83520 84445
ACL	TTGAB	TTG IGA Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 6 months	83516
ACL	TTGIGG	TTG IGG Container :One gold gel 3.5 mL Stability :Ambient: 8 Hours / Refrigerated: 2 Days / Frozen: 6 months	83516
ACL	U5HIAA	U 5 HIAA Container :24 hour urine or random urine. Sample must be refrigerated during collection. Refrigeration is the most important aspect of sample preservation. Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 2 Weeks	83497
ACL	ETOHUR	U-ETHANOL Container :Urine Stability :Ambient: 24 Hours / Refrigerated: 14 Days / Frozen: 14 Days	80301
ACL	UMETMX	UPE MONOCLONAL, 24HR Container :Urine Stability :Refrigerated: 7 Days / Frozen: 14 Days	84156 84166 86335
ACL	UMERMX	UPE MONOCLONAL, RAND Container :One urine - random Stability :Refrigerated: 7 Days / Frozen: 14 Days	84156 84166 86335
ACL	UPETRX	UPE w/ IFE, Urine 24 hour Container :One urine - 24 hour Stability :Refrigerated: 5 Days	84156 84166
ACL	UPERRX	UPE w/ IFE, Urine Random Container :One urine - random Stability :REFRIGERATED: 7 DAYS / FROZEN: 14 DAYS	84156 84166
ACL	UUNT	UREA NITRO,URIN 24HR Container :One Urine Stability :Refrigerated: 3 Days	84540
PHC	BUN	Urea Nitrogen (BUN), Plasma or Serum Stability :Room temperature: 3 days	84520



LAB	CODE	TEST NAME	CPT4 CODE
		Refrigerated: 7 days	
PHC	FLUN	Urea Nitrogen, Body Fluid Stability :Room temperature: 3 days Refrigerated: 7 days	84520
ACL	UUNR	UREA NITROGEN,URINE Container :Urine Stability :Refrigerated: 3 Days	84540
PHC	URICR	Uric Acid Rasburicase Stability :Ambient: Unacceptable Refrigerated: 4 Hours (if plasma removed from cells), 2 hours in original Li Hep tube Frozen: Unacceptable	84550
PHC	URIC	Uric Acid, Plasma or Serum Stability :Room temperature: 1 days Refrigerated: 3 days	84550
ACL	URURCR	URIC ACID, URINE RANDOM Container :Urine Stability :Refrigerated: 3 Days	84560
ACL	UURICT	URIC ACID,URINE 24HR Container :Urine Stability :Refrigerated: 3 Days	84560
PHC	UACNRX	Urinalysis w/ Microscopic Analysis, No Reflex Culture	81001
PHC	UACOM	Urinalysis w/ Microscopic Analysis, Reflex Culture	81001
PHC	UDPMIC	Urinalysis with Reflex Microscopic Analysis	81003
PHC	UDPST	Urinalysis, Dipstick	81003
PHC	MICRO	Urinalysis, Microscopic Stability :Room temperature: 2 hours Refrigerated: 24 hours Frozen: Do Not Freeze	81015
ACL	URC	Urine Culture Stability :AMBIENT: 2 HOURS UNPRESERVED; 48 HOURS IN GRAY TOP PRESERVATIVE TUBE. / REFRIGERATED: 24 HOURS UNPRESERVED; 48 HOURS IN A GRAY TOP PRESERVATIVE TUBE.	87086
ACL	UREOS	URINE EOSINOPHILS Container :One 10.0 mL (min: 5.0 mL) urine ambient in red/yellow top urine preservative tube. Stability :Ambient: 2 Hours unpreserved, 48 hours in red/yellow top preservative tube / Refrigerated: 2 Hours unpreserved, 48 Hours in a red/yellow top preservative tube.	87205
ACL	UroVysion	Urovysion Fish Container :THINPREP® UROCYTE URINE COLLECTION KIT, Stability :REFRIGERATED: 3 DAYS(MUST BE FIXED WITH PRESERVCYT®)	88210
ACL	CUIND	Urticaria-Inducing Activity Container :One plain red 4.0 mL Stability :Ambient: 30 minutes / Refrigerated: Unacceptable / Frozen: 1 Year	86352
ACL	USTEKQNT	Ustekinumab Quantitation with Antibodies, Serum Container :0.5 mL (Min: 0.35 mL) serum refrigerated in ARUP tube Stability :Refrigerated: 21 Days Frozen: 21 Days	80299 83520
PHC	VALP	Valproic Acid, Serum Stability :Refrigerated: 1 day	80164
ACL	CULVSC	Vancomycin Resistant Enterococcus Culture Container :ESWAB COLLECTION KIT, AMBIENT AEROBIC CULTURE SWAB IN TRANSPORT MEDIA, AMBIENT FLUIDS IN A STERILE LEAK-PROOF CONTAINER, REFRIGERATED URINE IN BORIC ACID PRESERVATIVE TUBE, AMBIENT URINE IN A STERILE LEAK-PROOF CONTAINER, REFRIGERATED STOOL IN CARY-BLAIR, AMBIENT Stability :AMBIENT: 48 HOURS FOR URINE IN BORIC ACID PRESERVATIVE TUBES, ESWAB AND AEROBIC CULTURE SWAB, 4 DAYS FOR CARY-BLAIR / REFRIGERATED: 24 HOURS FOR FLUIDS AND UNPRESERVED URINE	87081
PHC	VANCP	Vancomycin, Serum-Peak Stability :Refrigerated: 1 day	80202



LAB	CODE	TEST NAME	CPT4 CODE
PHC	VANCR	Vancomycin, Serum-Random Stability :Refrigerated:1 day	80202
PHC	VANCT	Vancomycin, Serum-Trough Stability :Refrigerated:1 day	80202
ACL	UVMA	Vanillylmandelic Acid, Urine Stability :Ambient: Unacceptable / Refrigerated: 1 Week / Frozen: 2 Weeks	84585
ACL	VIST	VARIC. ZOSTER IGG AB Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	86787
PHC	VZABSA	VARICELLA ZOSTER IGG/IGM Stability :VIST: Refrigerated: 7 Days / Frozen: 14 Days VIGM: Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Year	86787 x2
ACL	VIGM	VARICELLA ZOSTER IGM Container :One gold gel 3.5 mL Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Year	86787
ACL	VZPCR	VARICELLA ZOSTER PCR Container :Varies Stability :Ambient: 8 Hours / Refrigerated: 72 Hours / Frozen: 1 Month	87798
ACL	VEGF	Vascular Endothelial Growth Factor Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 4 Hours / Refrigerated: 6 Hours / Frozen: 6 Months	83520
ACL	VASSRN	VASCULITIS AB SCREEN Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Refrigerated: 7 Days / Frozen: 14 Days	83516 x3
ACL	VIPB	VASOACTIVE INT PEP Container :One Protease Inhibitor tube (PPACK; Phe-Pro-Arg-chloromethylketone). Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 90 Days	84586
ACL	VEDOQNT	Vedolizumab Quantitation with Antibodies, Serum Container :Separate from cells ASAP or within 2 hours of collection. Transfer 1.5 mL serum to an ARUP Standard Transport Tube. (Min: 0.75 mL) Stability :Ambient: Unacceptable Refrigerated: 28 Days Frozen: 28 Days	80280 82397
ACL	VIBC	Vibrio Culture Container :STOOL transferred into CARY-BLAIR (ENTERIC PATHOGEN) Stability :AMBIENT: 1 HOUR UNPRESERVED - 4 DAYS PRESERVED IN CARY-BLAIR (C&S) TRANSPORT MEDIA	87046
ACL	VITDO	VIT D, 25-HYDROXY Container :One gold gel 5.0 mL (also acceptable: plain red) Stability :Ambient: 24 hours after separation from cells. / Refrigerated: 7 Days after separation from cells. / Frozen: 14 Days after separation from cells.	82306
ACL	VITAT	VITAMIN A Container :One pink (EDTA) 6.0 mL (also acceptable: gold gel 3.5 mL or lavender (EDTA) 3.0 mL or light green (lithium heparin gel) 3.0 mL or green (lithium heparin no gel) 4.0 mL or green (sodium heparin no gel) 4.0 mL) Stability :Refrigerated: 1 month / Frozen: 1 Year	84590
ACL	VB1	VITAMIN B1,WHOLE BLD Container :One green (sodium or lithium heparin no gel) 6.0 mL (also acceptable: One lavender (EDTA) 3.0 mL) Stability :Refrigerated: 4 Hours / Frozen: 6 Months	84425
ACL	B12	VITAMIN B12 Container :One gold gel 5.0 mL (also acceptable: One plain red 6.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82607
ACL	B12FOL	VITAMIN B12/FOLATE Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	82607 82746
ACL	RIBO	VITAMIN B2 Container :One green (lithium heparin no gel) 4.0 mL (also acceptable: one green (lithium heparin no gel) 2.0 mL and one light green (lithium heparin gel) 3.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 1 Month	84252
ACL	VITVB6	VITAMIN B6, PLASMA Container :One green (sodium or lithium heparin) 6.0 mL blood Also acceptable: One lavender (EDTA) 3.0 mL, One pink (K2EDTA) 6.0 mL, One gold gel 5.0 mL or One plain red 6.0 mL	84207



LAB	CODE	TEST NAME	CPT4 CODE
		Also acceptable: One lavender (EDTA) 3.0 mL, One pink (K2EDTA) 6.0 mL, One gold gel 5.0 mL or One plain red 6.0 mL Stability :Ambient: unacceptable / Refrigerated: 1 Week after separation from cells / Frozen: 2 Months after separation from cells	
ACL	VITCB	VITAMIN C Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 1 Month and protected from light	82180
ACL	VITAD	VITAMIN D,1,25-DIHYD Container :One gold gel 3.5 mL Stability :Ambient: 7 Days / Refrigerated: 2 Week / Frozen: 2 weeks	82652
ACL	VITE	VITAMIN E Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Year	84446
ACL	VITK	VITAMIN K Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL and One lavender (EDTA) 3.0 mL) Stability :Refrigerated: 1 Month / Frozen: 6 Months	84597
PHC	VS	Volatile Screen, Serum Stability :Room temperature: Do Not Store at Room Temperature Refrigerated: Forward STAT Frozen: Do Not Freeze	80320
ACL	VWFR	VON WILLEBRAND ACT Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: 4 Hours / Refrigerated: Unacceptable / Frozen: 30 Days	85245
ACL	VWAG	VON WILLEBRAND AG Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85246
ACL	VOR	VORICONAZOLE Container :One plain red 4.0 mL Stability :Ambient: 28 Days / Refrigerated: 28 Days / Frozen: 28 Days	80285
ACL	VWM	VWF MULTIMER Container :One light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 6 Months	85247
PHC	WBCA	WBC Count, Blood	85048
ACL	WNVAB	WEST NILE AB PANEL Container :One gold gel 3.5 mL Stability :Ambient: 48 Hours / Refrigerated: 2 weeks / Frozen: 1 year	86788 86789
ACL	WNVCSF	WEST NILE VIRUS CSF Container :One Spinal fluid (CSF) Stability :Ambient: 8 Hours / Refrigerated: 2 weeks / Frozen: 1 Year Avoid repeated freeze thaw cycles	86788 86789
PHC	WETPRE	Wet Preparation, Vaginal Stability :Room temperature: 24 hours Refrigerated: 36 hours Frozen: Do Not Freeze	87210 87808
ACL	NTX	X-LINK N-TELOPEPTIDE Container :One urine - random Stability :Ambient: 24 Hours / Refrigerated: 7 days / Frozen: 6 Months	82523
ACL	YSTC	Yeast Culture Container :ESWAB COLLECTION KIT - AMBIENT OR AEROBIC CULTURE SWAB IN TRANSPORT MEDIA - AMBIENT OR 4.0 ML (MIN: 4.0 ML) URINE (MUST BE FILLED TO THE MARK ON THE SIDE OF THE GRAY TOP TUBE. IF INSUFFICIENT VOLUME OF URINE FOR GRAY TOP TUBE, URINE MUST BE REFRIGERATED AT 4 DEGREES C WITHIN 2 HOURS OF COLLECTION AND SETUP WITHIN 24 HOURS) - REFRIGERATED Stability :AMBIENT: 3 DAYS ESWAB OR AEROBIC CULTURE SWAB; 48 HOURS IN GRAY TOP PRESERVATIVE TUBE / REFRIGERATED: 48 HOURS IN A GRAY TOP PRESERVATIVE TUBE / Refrigerated: 24 Hours low volume urine in sterile cup/FROZEN: UNACCEPTABLE	87102
ACL	YERC	Yersinia Culture Container :RECTAL ESWAB OR ESWAB OF ANY BODY SITE (ALSO ACCEPTABLE: AEROBIC CULTURE SWAB). STOOL CULTURE CARY-BLAIR (ENTERIC PATHOGEN Stability :AMBIENT: 48 HOURS FOR URINE IN BORIC ACID PRESERVATIVE TUBES, ESWAB AND AEROBIC CULTURE SWAB, 4 DAYS FOR CARY-BLAIR / REFRIGERATED: 24 HOURS FOR FLUIDS AND UNPRESERVED URINE	87046



LAB	CODE	TEST NAME	CPT4 CODE
ACL	ZIKAPCR	Zika Virus by PCR, Blood Container :One gold gel 5.0 mL blood Stability :Ambient: Unacceptable / Refrigerated: 5 days / Frozen:6 weeks	87662
ACL	ZIKA AB	ZIKA VIRUS IGM ANTIBODY CAPTURE (MAC) BY ELISA Container :One gold gel 3.5 mL (also acceptable: One plain red 4.0 mL) Stability :Ambient: 2 Days / Refrigerated: 2 Weeks / Frozen: 1 Year Avoid repeated freeze/thaw cycles	86794
ACL	ZINC 8	Zinc Transporter 8 Antibody Container :One gold gel 3.5 mL blood Stability :Ambient: 5 Days after separation from cells Refrigerated: 1 Week after separation from cells Frozen: 3 Months after separation from cells	86341
ACL	ZNRBC	Zinc, RBC Container :One royal blue (EDTA) 6.0 mL Stability :Ambient: 2 Weeks after separation from plasma Refrigerated: 2 Weeks after separation from plasma Frozen: Unacceptable	84630
ACL	ZIN	ZINC,BLOOD Stability :Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	84630
ACL	ZON	ZONISAMIDE Container :One plain red 6.0 mL blood Also acceptable: One lavender (EDTA) 3.0 mL Stability :Ambient: 6 Hours / Refrigerated: 1 Week / Frozen: 3 Month	80203