



DEGREE OF HYDROXOCOBALAMIN INTERFERENCE WITH COMMON LABORATORY TESTS

Due to the intense red color of hydroxocobalamin (OHCob) and its absorbance at 352 and 525 nm, interference studies were performed on common Chemistry and Hematology tests at Legacy Laboratories. The following tables summarize these findings. If the test is not listed in the tables below, then the degree of interference is unknown. Delta is an estimate of the amount of change of that analyte at 1000 mg/L OHCob. Healthcare providers should interpret patient results with caution since these results are only an estimate of OHCob interference for each of the analytes. The actual *in vivo* concentration of OHCob in the patient's blood or urine is based on many factors including, but not limited to, dose, dose rate, and time since administration, and genetics. The tables are organized as described below:

| TABLE | DEPARTMENT | MATRIX |
|-------|---------------------|-------------------------|
| 1 | Chemistry | Serum or plasma |
| 2 | Chemistry | Urine |
| 3 | Immunochemistry | Serum or plasma |
| 4 | Respiratory Therapy | Heparinized whole blood |
| 5 | Hematology | EDTA whole blood |
| 6 | Coagulation | Citrated plasma |
| 7 | Urinalysis | Urine |

Table 1: OHCob Interference on Blood Chemistry Tests¹ – Various concentrations of OHCob was added to pooled Serum and Heparinized Plasma specimens.

| INTERFERENCE | ASSAY | METHOD | DELTA |
|--------------|---------------------------------------|------------|-------|
| None | Albumin | Beckman AU | <6% |
| | Alkaline Phosphatase (ALP) | Beckman AU | <6% |
| | CO ₂ , Total (Bicarbonate) | Beckman AU | <6% |
| | Calcium (Ca) | Beckman AU | <6% |
| | Chloride (Cl) | Beckman AU | <6% |
| | GGT | Beckman AU | <6% |
| | Glucose | Beckman AU | <6% |
| | Lactate | Beckman AU | <6% |
| | Lactate Dehydrogenase (LDH) | Beckman AU | <6% |
| | Lipase | Beckman AU | <6% |
| | Potassium (K) | Beckman AU | <6% |
| | Protein, Total (TP) | Beckman AU | <6% |
| | Sodium (Na) | Beckman AU | <6% |
| | Urea Nitrogen (BUN) | Beckman AU | <6% |
| | Triglyceride | Beckman AU | <6% |



DEGREE OF HYDROXOCOBALAMIN INTERFERENCE WITH COMMON LABORATORY TESTS

| INTERFERENCE | ASSAY | METHOD | DELTA |
|--------------|----------------------|----------------------------------|--------------------|
| Increase | Cholesterol, Total | Beckman AU | 51.0% |
| | Ethanol | Beckman AU/ Siemens Syva Emit | 11.0% |
| | Magnesium | Beckman AU | 34.0% |
| | Phosphorus | Beckman AU | 8.0% |
| | Prealbumin | Beckman AU | 11.0% |
| Decrease | ALT | Beckman AU | -5 to -73% |
| | Amylase | Beckman AU | -32.0% |
| | AST | Beckman AU | -1 to -51% |
| | Bilirubin, Direct | Beckman AU | -0.1 to -0.2 mg/dL |
| | Bilirubin, Total | Beckman AU | -0.1 to -0.2 mg/dL |
| | Creatinine | Beckman AU | -0.1 to -0.2 mg/dL |
| | Creatine Kinase (CK) | Beckman AU | -15.0% |
| | Iron | Beckman AU | -8.0% |
| | Uric Acid | Beckman AU | -28.0% |

Table 2: OHCob Interference on Urine Chemistry Tests¹ – Various concentrations of OHCob was added to pooled urine specimens.

| INTERFERENCE | ASSAY | METHOD | DELTA |
|--------------|----------------------|------------|-----------|
| None | Calcium (U Ca) | Beckman AU | <6% |
| | Chloride (U Cl) | Beckman AU | <6% |
| | Creatinine | Beckman AU | <6% |
| | Glucose | Beckman AU | <6% |
| | Phosphorus | Beckman AU | <6% |
| | Potassium (U K) | Beckman AU | <6% |
| | Sodium (U Na) | Beckman AU | <6% |
| | Urea Nitrogen (UUN) | Beckman AU | <6% |
| | Uric Acid | Beckman AU | <6% |
| Increase | Magnesium (U Mg) | Beckman AU | 0.6 mg/dL |
| | Protein, Total (UTP) | Beckman AU | 8.0% |

Table 3: OHCob Interference on Blood Immunochemistry Tests¹ – Various concentrations of OHCob was added to pooled Serum and Heparinized Plasma specimens.

| INTERFERENCE | ASSAY | METHOD | DELTA |
|----------------|------------------|------------------|--------|
| Decreased | CK-MB, Blood | Beckman Access 2 | -15.0% |
| Unpredictable* | Myoglobin, Blood | Beckman Access 2 | na |

*Unpredictable = the concentration fluxuated in the both the positive and negative direction



DEGREE OF HYDROXOCOBALAMIN INTERFERENCE WITH COMMON LABORATORY TESTS

Table 4: OHCob Interference on Respiratory Therapy Tests, Blood Gases, and Co-Oximetry Tests^{1,2} - Various concentrations of OHCob was added to pooled venous whole blood.

| INTERFERENCE | ASSAY | METHOD | DELTA |
|----------------|-------------------------------|--------------------|-------------------------|
| None | Base Excess | Siemens RapidPoint | <6% ¹ |
| | Calcium, Ionized | Siemens RapidPoint | <6% ¹ |
| | Chloride | Siemens RapidPoint | <6% ¹ |
| | HCO ₃ ⁻ | Siemens RapidPoint | <6% ¹ |
| | Glucose | Siemens RapidPoint | <6% ¹ |
| | Sodium | Siemens RapidPoint | <6% ¹ |
| | pH | Siemens RapidPoint | <6% ¹ |
| | Potassium | Siemens RapidPoint | <6% ¹ |
| Decrease | pO ₂ | Siemens RapidPoint | -7.0% ¹ |
| | Carboxyhemoglobin (COHb) | Siemens RapidPoint | Up to -10% ² |
| | Methemoglobin (MetHb) | Siemens RapidPoint | Up to -7% ² |
| | Total Hemoglobin | Siemens RapidPoint | Up to -4% ² |
| | Deoxyhemoglobin | Siemens RapidPoint | Up to -1% ² |
| Increase | Oxyhemoglobin | Siemens RapidPoint | Up to 8% ² |
| Unpredictable* | pCO ₂ | Siemens RapidPoint | na ¹ |

*Unpredictable = the concentration fluxuated in the both the positive and negative direction

Table 5: OHCob Interference on Hematology Tests¹ - Various concentrations of OHCob was added to EDTA whole blood.

| INTERFERENCE | ASSAY | METHOD | DELTA |
|--------------|---------------|--------|-------|
| None | WBC | Sysmex | <2% |
| | RBC | Sysmex | <1% |
| | Hematocrit | Sysmex | <1% |
| | MCV | Sysmex | <1% |
| | RDW-CV | Sysmex | <1% |
| | Platelets | Sysmex | <2% |
| | MPV | Sysmex | <1% |
| | Neutrophils % | Sysmex | <1% |
| | Neutrophils # | Sysmex | <2% |
| | Lymphocytes % | Sysmex | <1% |
| | Lymphocytes # | Sysmex | <3% |
| | Monocytes % | Sysmex | <2% |
| | Monocytes # | Sysmex | <1% |
| | Eosinophils % | Sysmex | <4% |
| | Eosinophils # | Sysmex | <1% |
| | Basophils % | Sysmex | <5% |
| | Basophils # | Sysmex | <1% |



DEGREE OF HYDROXOCOBALAMIN INTERFERENCE WITH COMMON LABORATORY TESTS

| INTERFERENCE | ASSAY | METHOD | DELTA |
|--------------|-----------------|--------|-----------|
| Increase | Hemoglobin | Sysmex | 6.0% |
| | Reticulocytes % | Sysmex | 5% to 15% |

Table 6: OHCob Interference on Coagulation Tests¹ – Interference studies were not determined using Legacy’s current instrumentation. Based on previous studies using a different method, OHCob causes significant changes in PT-INR, PTT, and Fibrinogen results.¹ Coagulation testing is not recommended on patients who have been administered this drug.

| INTERFERENCE | ASSAY |
|--|------------|
| Results change significantly. Coagulation testing not recommended while the patient is taking this drug. | PT-INR |
| | PTT |
| | Fibrinogen |

Table 7: OHCob Interference on Urinalysis³ - Various concentrations of OHCob was added to Urine.

| INTERFERENCE | ASSAY | DELTA |
|--------------|--------------|-------|
| Increase | pH | >10% |
| | Glucose | >10% |
| | Protein | >10% |
| | Erythrocytes | >10% |
| | Leukocytes | >10% |
| | Ketones | >10% |
| | Bilirubin | >10% |
| | Urobilinogen | >10% |
| | Nitrate | >10% |

REFERENCES:

1. Beaudoin D, Peck J, Schulmerich M, Bettger M, Kong L, Chan B, Smith C, Morishita C, Westfall E, Legacy Laboratory Services in-house studies: *Hydroxocobalamin Interferes with Common Laboratory Tests [Abstract]*. Clin Chem, 2008.
2. *Urgent Field Safety Notice - Hydroxocobalamin Interference*, Siemens Healthcare Diagnostics, Inc., POC 18-010.A.US, July 2018.
3. *Cyanokit® Package Insert* (drug’s manufacturer package insert), Dey, L.P. Napa, CA 94558. Last revision: 12/2006.