



| LAB | CODE      | TEST NAME   | CPT4 CODE |
|-----|-----------|---|-----------|
| ACL | 17OHPG    | 17-Hydroxypregnenolone<br>Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL or two lavender (EDTA) 3.0 mL or one green (lithium heparin, no gel) 6.0 mL or one green (sodium heparin, no gel) 6.0 mL)<br>Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 6 Months | 84143     |
| ACL | 17HYDR    | 17-Hydroxyprogesterone<br>Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL or one lavender (EDTA) 6.0 mL or one green (sodium or lithium heparin, no gel) 6.0 mL)<br>Stability :Ambient: 3 Days / Refrigerated: 7 Days / Frozen: 6 Months  | 83498     |
| ACL | 17OHP Neo | 17-Hydroxyprogesterone, Neonatal/Infant<br>Container :One plain red (also acceptable: one lavender (EDTA) 3.0 mL or one green (sodium heparin or lithium heparin) 4.0 mL)<br>Stability :Ambient: 48 Hours / Refrigerated: 7 Days / Frozen: 2 Years  | 83498     |
| ACL | 21OHAB    | 21 Hydroxylase Antibody<br>Container :One plain red 6.0 mL (also acceptable: one gold gel 5.0 mL)<br>Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 1 Month  | 83516     |
| ACL | U5HIAA    | 5 HIAA, Urine 24 Hour<br>Container :One Urine<br>Stability :Ambient: Unacceptable / Refrigerated: 1 Week / Frozen: 2 Months   | 83497     |
| ACL | 5NUCL     | 5' Nucleotidase<br>Container :One gold gel 5.0 mL<br>Stability :Ambient: 4 Hours / Refrigerated: 1 Week / Frozen: 2 Weeks (avoid repeated freeze/thaw cycles)   | 83915     |
| ACL | 5ADTSR    | 5-A-Dihydrotestosterone<br>Container :One plain red 6.0 mL (also acceptable: one gold gel 5.0 mL)<br>Stability :Ambient: 48 Hours / Refrigerated: 5 days / Frozen: 6 Months   | 82642     |