



LAB	CODE	TEST NAME	CPT4 CODE
ACL	CPEP	C Peptide Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 24 Hours / Frozen: 14 Days	84681
PHC	CRP	C-Reactive Protein (CRP), Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 3 days for dimension. 7 days for Vista Frozen: 6 months at -20o C	86140
ACL	BCTX	C-Telopeptide, Beta Cross Linked Container :One gold gel 5.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 48 Hours / Frozen: 3 Months	82523
ACL	C1ES	C1 Esterase Inhibitor Antigen Container :One gold gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 14 Days / Frozen: 1 Month	86160
ACL	C1EFUN	C1 Esterase Inhibitor Functional Container :One gold gel 5.0 mL Stability :Ambient: 2 Hours / Refrigerated: Unacceptable / Frozen: 2 Weeks	86161
ACL	C3	C3 Complement Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	86160
ACL	C4	C4 Complement Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	86160
ACL	CAD	Cadmium, Blood Container :One royal/dark blue (EDTA) 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82300
ACL	CALCIT	Calcitonin, Blood Container :One gold gel 5.0 mL (also acceptable: one green (sodium or lithium heparin, with or without gel) 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 3 Months	82308
ACL	ICAOUT	Calcium, Ionized (CLINIC ONLY) Container :One gold gel 3.5 mL NOTE: Tube must be at least half full Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: Unacceptable	82330
PHC	IONCA	Calcium, Ionized, Blood	82330
PHC	CA	Calcium, Total, Plasma or Serum Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: 6 months	82310
ACL	UCAT	Calcium, Urine 24 Hour Container :One 24 hour urine Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 14 Days	82340
ACL	UCAR	Calcium, Urine Random Container :One random urine Stability :Ambient: Unacceptable / Refrigerated: 5 Days / Frozen: 14 Days	82340
ACL	CSA	Calculi Analysis Stability :Ambient: 1 Year / Refrigerated: 1 Year / Frozen: 1 Year	82365