



LAB	CODE	TEST NAME	CPT4 CODE
ACL	FACX	Factor 10 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85260
ACL	FACXI	Factor 11 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85270
ACL	FACXII	Factor 12 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85280
ACL	FACII	Factor 2 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 30 Days	85210
ACL	FACV	Factor 5 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85220
ACL	LEIDEN	Factor 5 Leiden Container :One lavender (EDTA) 6.0 mL (also acceptable FOR EXTENUATING CIRCUMSTANCES ONLY: one light blue (3.2% sodium citrate) 2.7 mL) Stability :Ambient: 72 Hours / Refrigerated: 14 Days / Frozen: Unacceptable	81241
ACL	FACVII	Factor 7 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85230
PHC	FCVIII	Factor 8 Assay Container :One FULL light blue-top (citrate) tube Stability :Frozen: 30 Days	85240
ACL	FACTIX	Factor 9 Assay Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Frozen: 30 Days	85250
ACL	FXIIIAG	Factor XIII Antigen Level Container :Two light blue (3.2% sodium citrate) 2.7 mL Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 14 Days	85290
ACL	FACID	Fatty Acids Profile, Essential Container :One green (sodium or lithium heparin) 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL or one plain red 6.0 mL or one gold gel 5.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 1 Week / Frozen: 3 Months	82542
ACL	FATQLT	Fecal Fat, Qualitative Container :Stool in sterile container Stability :Ambient: 6 Hours / Refrigerated: 24 Hours / Frozen: 1 Month	89125
PHC	WBCSTL	Fecal Leukocytes Stability :Room temperature: 2 weeks Refrigerated: 2 weeks Frozen: 1 month	83630
ACL	IFOB	Fecal Occult Blood, Tube Test Container :One stool specimen Stability :Ambient: 15 Days / Refrigerated: 30 Days / Frozen: Unacceptable	82274
ACL	FEL	Felbamate Level Container :One plain red 4.0 mL (also acceptable: one lavender (EDTA) 3.0 mL or one green (sodium heparin no gel) 4.0 mL) Stability :Ambient: 48 Hours / Refrigerated: 1 Month / Frozen: 6 Months	80167