



LAB	CODE	TEST NAME	CPT4 CODE
PHC	SPESCN	IFE with Immunoglobulins	86334 82734 x3
ACL	IG1234	IgG SUBCLASS 1,2,3,4 Container :One gold gel 3.5 mL Stability :Ambient: 2 Hours / Refrigerated: 8 Days / Frozen: 6 Months	82787 x4
ACL	IGGSUB	IgG Subclasses 1,2,3,4 Container :One gold gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 2 Weeks / Frozen: 6 Months	82787 x4
PHC	IGSUBA	IgG with Subclasses Container :IGG: One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) IGGSUB: One gold gel 5.0 mL	82787 x4 82784
PHC	IGHV	IGHV Mutation Analysis by Sequencing Stability :Ambient: 1 hour: Refrigerated: 7 days: Frozen: Unacceptable	81263
PHC	IPF	Immature Platelet Fraction Stability :Room temperature: Preferred specimen should be less than 8 hours old. If stored at 4-8°C, allow it to come to room temperature before testing.	85055
ACL	SIFEMX	Immunofixation Electrophoresis with IgG, IgA, IgM, and Kappa/Lambda Free Light Chains, Serum Container :One gold gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	82784 x3 83521 x2 86334
ACL	UIFXMX	Immunofixation Electrophoresis, Urine Container :One Urine Stability :Refrigerated: 7 Days / Frozen: 14 Days	86335
ACL	IGA	Immunoglobulin A, IgA Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	82784
ACL	IGDQNT	Immunoglobulin D, Quantitative Container :One gold gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 2 Weeks / Frozen: 6 Months	82784
ACL	IGEB	Immunoglobulin E, IgE Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL) Stability :Ambient: 4 Hours / Refrigerated: 2 Days / Frozen: 14 Days	82785
ACL	IGG	Immunoglobulin G (IgG) Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	82784
ACL	IGGS4	Immunoglobulin G Subclass 4 Container :One gold gel 5.0 mL Stability :Ambient: Unacceptable / Refrigerated: 14 Days / Frozen: 6 Months	82787
ACL	IGM	Immunoglobulin M, IgM Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	82784
ACL	IMMGLB	Immunoglobulin Profile (IgG, IgM, IgA) Container :One gold gel 3.5 mL (also acceptable: one gold gel microtainer 600 mcL (clear or amber) or one plain red 4.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	82784 x3