



LAB	CODE	TEST NAME	CPT4 CODE
ACL	MACPRO	Macroprolactin Container :One gold gel 5.0 mL (also acceptable: one green lithium heparin 6.0 mL or one lavender (EDTA) 3.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 1 Month	84146 x2
ACL	MAGRBC	Magnesium RBC Container :One royal blue (EDTA) 6.0 mL Stability :Ambient: 2 Weeks / Refrigerated: 2 Weeks / Frozen: Unacceptable	83735
PHC	MG	Magnesium, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 1 year	83735
ACL	UMGT	Magnesium, Urine 24 Hour Container :One 24 hour urine Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	83735
ACL	UMGR	Magnesium, Urine Random Container :One random urine Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	83735
ACL	RMALP	MALARIA RAPID SCREEN Container :One lavender (EDTA) 3.0 mL blood Stability :Ambient: 24 Hours / Refrigerated: 24 Hours / Frozen: Unacceptable	
PHC	VR01713	Measles Virus PCR Container :UTM Stability :Refrigerated: 48 hours	87798
ACL	DRGSCM	Meconium Drug Screen Container :Collect all meconium (blackish material) excreted until milk/formula-based stool (yellow-green) appears. Transfer meconium into a sterile container. Stability :Ambient: 1 Week / Refrigerated: 3 Weeks / Frozen: 1 Year	80346 80348 80353 80361 80356 80365 80373 80354 80362 80355 80359 80325 80360 80345 80372 80358 83992 80323 80368
ACL	MSHA	Melanocyte Stimulating Hormone Container :One lavender (EDTA) 3.0 mL Stability :Ambient: Unacceptable / Refrigerated: 24 Hours / Frozen: 1 Month	83519
ACL	MELAT	Melatonin Container :One plain red 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: 2 Weeks / Refrigerated: 1 Month / Frozen: 1 Year	80375
ACL	MERCB	Mercury, Blood Container :One royal blue (EDTA) 6.0 mL (also acceptable: one lavender (EDTA) 6.0 mL) Stability :Ambient: 14 Days / Refrigerated: 14 Days / Frozen: Unacceptable	83825
ACL	PMET	Metanephrines, Frac/Free, Blood Container :One green (sodium or lithium heparin, no gel) 6.0 mL (also acceptable: one lavender (EDTA) 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 10 Days / Frozen: 1 Month	83835
ACL	UMET	Metanephrines, Urine Container :One Urine Stability :Ambient: Unacceptable / Refrigerated: 2 Weeks / Frozen: 1 Month	83835
ACL	METHQNT	Methadone and Metabolite, Quantitative Serum or Plasma Container :One plain red 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL or two gray (sodium fluoride/potassium oxalate) 2.0 mL or one green (sodium heparin no gel) 4.0 mL) Stability :Ambient: 1 Week / Refrigerated: 2 Weeks / Frozen: 3 Years (avoid repeated freeze/thaw cycles)	80358
ACL	MTHCNF	Methadone Confirmation, Urine Container :One random urine Stability :Ambient: 1 Week / Refrigerated: 1 Month / Frozen: 3 Years	80358