



LAB	CODE	TEST NAME	CPT4 CODE
ACL	NTXU	N-Telopeptide, Urine Container :One urine Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 4 Weeks (avoid repeated freeze/thaw cycles)	82523
PHC	NTPBNP	N-Terminal Pro Brain Natriuretic Peptide Stability :Ambient: 3 Days / Refrigerated: 3 Days	83880
ACL	GCC	N. Gonorrhoeae (GC) Culture Container :Fluid obtained from any appropriate site in an ESwab (preferred) or sterile leak-proof container. Stability :Ambient: 24 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	87081
ACL	JNT04	Native Joint Culture and Gram Stain (4 Day Hold) Stability :Ambient: 48 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	87070 87075 87205
ACL	NEOHSV	Neonatal HSV Screen Stability :Ambient: 8 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable	87529 x2
ACL	NSE	Neuron Specific Enolase, Serum Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 1 Week / Frozen: 1 Year (avoid repeated freeze/thaw cycles)	86316
PHC	NEUTAB	Neutrophil Antibody Screen Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: > 7 Days	86021
PHC	NBSCR	Newborn Screening Panel, Blood Spots	82017 82261 82775 83020 83498 83516 83788 84030 84443 84510
ACL	NEXTGENMY	Next Generation Sequencing (NGS) for Acute Myeloid Leukemia Container :Blood: Two lavender (EDTA) 3.0 mL (also acceptable FOR EXTENUATING CIRCUMSTANCES ONLY: two green (sodium heparin, no gel) 3.0 mL or two green (lithium heparin, no gel) 3.0 mL) Bone Marrow: One lavender (EDTA) 3.0 mL (also acceptable FOR EXTENUATING CIRCUMSTANCES ONLY: one green (sodium heparin, no gel) 3.0 mL or one green (lithium heparin, no gel) 3.0 mL) Stability :Ambient: 72 Hours / Refrigerated: 2 Weeks / Frozen: Unacceptable	81455
ACL	NEXTGENCF	Next Generation Sequencing (NGS) for Cystic Fibrosis Carrier Screening Container :One lavender (EDTA) 6.0 mL Stability :Ambient: 72 Hours / Refrigerated: 14 Days / Frozen: Unacceptable	81220
ACL	NICKEL	Nickel, Serum Container :One royal blue (no additive) 6.0 mL Stability :Ambient: Indefinitely / Refrigerated: Indefinitely / Frozen: Indefinitely	83885
ACL	NICOT	Nicotine and Metabolites Level, Quantitative Container :Two plain red 6.0 mL (also acceptable: three lavender (EDTA) 3.0 mL or two green (sodium heparin) 6.0 mL) Stability :Ambient: 1 Week / Refrigerated: 2 Weeks / Frozen: 3 Years (avoid repeated freeze/thaw cycles)	80323
ACL	NICOU	Nicotine and Metabolites, Urine Container :One random urine Stability :Ambient: 10 Days / Refrigerated: 10 Days / Frozen: 8 Months	80323
ACL	AQUA4 CSF	NMO AQ4 IgG CSF Container :One Spinal fluid (CSF) Stability :Refrigerated: 2 Weeks/Frozen: 1 Year	86255
ACL	CYTOLOGY	Non-Gynecological Cytology Container :Varies	Varies