



LAB	CODE	TEST NAME	CPT4 CODE
PHC	SVCO2	O2 SATURATION CENTRAL VENOUS LINE Stability :Rm Temp: 30 minutes	82810
ACL	TTR	OB RBC Antibody Titer Monitoring Container :One lavender (EDTA) 6.0 mL AND one plain red 6.0 mL Stability :Ambient: 1 Day / Refrigerated: 3 Days	86886
PHC	OBRH	OB RH Container :Collect: Full, large lavender-top (EDTA) tube. 3.5mL min. Serum gel tube is not acceptable. Fill tube completely. Stability :Room temperature: 8 hours Refrigerated: 4 days Frozen: Do Not Freeze	86901
PHC	OBPAN	Obstetrics Panel with HIV Stability :Refer to individual test listings.	80081
ACL	OLAN	Olanzapine, Blood Container :One plain red 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 1 Month	80342
ACL	OPAT	Opiates Confirmation, Urine Container :One random urine Stability :Ambient: 1 Week / Refrigerated: 1 Month / Frozen: 3 Years	80361 80365
PHC	OPIAT	Opiates, Urine Stability :Refrigerated: 24 Hours Freeze if > 24 Hours	80307
ACL	UORAC	Organic Acids Screen, Urine Container :One random urine Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 1 Month	83918
PHC	OSMOL	Osmolality, Serum Stability :Room temperature: 24 hours Refrigerated: 7 days Frozen: 14 days	83930
PHC	UOSMOL	Osmolality, Urine Stability :Room temperature: 24 hours Refrigerated: 7 days Frozen: 14 days	83935
ACL	FRAGSA	Osmotic Fragility, RBC Container :Two green (sodium or lithium heparin) 4.0 mL (also acceptable: two lavender (EDTA) 3.0 mL) Stability :Ambient: Unacceptable / Refrigerated: 72 Hours / Frozen: Unacceptable	85555
ACL	OSTEO	Osteocalcin, Blood Container :One gold gel 3.5 mL (also acceptable: one lavender (EDTA) 3.0 mL or one green (lithium heparin) 4.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 72 Hours / Frozen: 3 Months	83937
ACL	OPUR	Ova and Parasite Exam, Body Fluid or Urine Container :One body fluid, CSF, or urine Stability :Body Fluid: Ambient: Unacceptable / Refrigerated: 72 Hours / Frozen: Unacceptable Urine: Ambient: 72 Hours / Refrigerated: 72 Hours / Frozen: Unacceptable CSF: Ambient: 72 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	87177 87209
ACL	OVAPF	Ova and Parasite Exam, Fecal (Immunocompromised or Travel History) Container :One stool specimen IMPORTANT: Stool must be transferred to either EcoFix or BOTH Pink & Gray PVA/10% Formalin vials within one hour of collection. Stability :Stool in EcoFix: Ambient: 9 Months / Refrigerated: 9 Months / Frozen: Unacceptable Stool in Pink & Gray PVA/10% Formalin: Ambient: 9 Months / Refrigerated: 9 Months / Frozen: Unacceptable	87177 87209
ACL	UOXAL	Oxalate, Urine 24 Hour Container :One 24 hour urine Stability :Refrigerated: 1 Week / Frozen: 1 Week	83945