



LAB	CODE	TEST NAME	CPT4 CODE
ACL	AFP4MA	QUAD MATERNAL SCREEN Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	82105 82677 84702 86336
ACL	QUANTP	Quantiferon TB Plus Stability :Ambient: 14 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	86480
PHC	CDDA	Quantitative Duodenal Aspirate Culture Container :Sterile Syringe with air removed (remove needle and cap syringe), Red Top Tube Stability :Ambient: 24 Hours	87071 87073 87102
PHC	CQNT	Quantitative Tissue Culture Stability :Ambient: 24 Hours	87071
PHC	SARSRAP AG	Quidel Sofia SARS Rapid Antigen FIA Container :Preferred: Nasal Swab Sample: Use the Puritan nylon flocked swab supplied by the Laboratory. Also acceptable: Nasopharyngeal Swab Sample: Use a nylon flocked nasopharyngeal swab provided by the Laboratory. Stability :Ambient: 48 hours / Refrigerated: 48 hours	87426