



| LAB | CODE | TEST NAME | CPT4 CODE |
|-----|------------|---|----------------------------|
| ACL | AFP4MA | QUAD MATERNAL SCREEN Container :One gold gel 5.0 mL Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days | 82105 82677 84702 86336 |
| ACL | QUANTP | Quantiferon TB Plus Stability :Ambient: 14 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable | 86480 |
| PHC | CDDA | Quantitative Duodenal Aspirate Culture Container :Sterile Syringe with air removed (remove needle and cap syringe), Red Top Tube Stability :Ambient: 24 Hours | 87071 87073 87102 |
| PHC | CQNT | Quantitative Tissue Culture Stability :Ambient: 24 Hours | 87071 |
| PHC | SARSRAP AG | Quidel Sofia SARS Rapid Antigen FIA Container :Preferred: Nasal Swab Sample: Use the Puritan nylon flocked swab supplied by the Laboratory. Also acceptable: Nasopharyngeal Swab Sample: Use a nylon flocked nasopharyngeal swab provided by the Laboratory. Stability :Ambient: 48 hours / Refrigerated: 48 hours | 87426 |