



LAB	CODE	TEST NAME	CPT4 CODE
ACL	YSTC	Yeast Culture Container :Collect one of the following: One eSwab or aerobic culture swab (acceptable sources: mouth, throat, vagina, or penis) Random urine in sterile container Stability :Eswab or Aerobic Culture Swab: Ambient: 72 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable Urine in Gray Boric Acid Tube: Ambient: Unacceptable / Refrigerated: 48 Hours / Frozen: Unacceptable Urine in Sterile Container: Ambient: Unacceptable / Refrigerated: 24 Hours / Frozen: Unacceptable	87102
ACL	YERC	Yersinia Culture, Stool Container :One stool specimen IMPORTANT: Stool must be transferred to Cary-Blair media within one hour of collection. Stability :Stool in Sterile Container: Ambient: 1 Hour Stool in Cary-Blair: Ambient: 4 Days / Refrigerated: Unacceptable / Frozen: Unacceptable	87046