


BD Vacutainer® Plus Plastic Citrate Tube

BD Vacutainer® Plus Plastic Citrate Tube Draw Volume Guide

*Sufficient volume achieved if blood drawn falls above minimum fill indicator.
For blood transfer, **do not** fill above illustrated dashed maximum line.*

Note: The quantity of blood drawn into evacuated tubes varies with altitude, ambient temperature, barometric pressure, tube age, venous pressure and filling technique.



*According to CLSI guideline, Dec. 2003, Doc. H1-A5, Vol. 23, No. 33.