


Legacy Hospitals / Labcorp Microbiology and Molecular Test Quick Reference Guide – Inpatient

Please refer to Lab Test Table on MyLegacy Intranet for most up to date specimen collection information

Test Name	Container	Infor / Lawson Ordering Information	
Genital			
Chlamydia & Gonorrhea by PCR, urine (CGAMD)	Urine Aptima	428206	
Chlamydia & Gonorrhea by PCR (CGAMD)	Aptima Swab	428207	
	Aptima Swab	428207	
	Aptima Swab	428207	
Vaginitis Panel, DNA probe	BD Affirm	283791	
Vaginitis Panel with STDs (Vaginitis Plus (VG+), NuSwab)	Aptima swab	428207	
Vaginitis Plus With Candida, NuSwab	Aptima swab or Aptima Multitest swab	428207	
		428208 (multitest swab)	
Vaginitis Plus With HSV, NuSwab	Aptima swab or Aptima Multitest swab	428207	
		428208 (multitest swab)	
BK Virus Quantitation by PCR, urine	Cobas PCR Media Urine Sample Kit	396120	
Herpes Simplex Virus 1/2 DNA by PCR	Transsystem double swab in viral transport media (M4rt)	274146 (swab)	
		142081 (M4rt)	

Legacy Hospitals / Labcorp Microbiology and Molecular Test Quick Reference Guide – Inpatient

Please refer to Lab Test Table on MyLegacy Intranet for most up to date specimen collection information

Test Name	Container	Infor / Lawson Ordering Information	
Xpert Xpress Multiplex Vaginal Panel (MVP)	Xpert MVP Swab Collection Kit	372282	
Trichomonas vaginalis by PCR (T VAG PCR)	T Vag GeneXpert Urine Tube		
Wet Mount	Culture sterile swab AND sterile saline	140758 (swab) Call Lab for sterile saline	
Respiratory			
B pertussis DNA by PCR	NP swab in viral transport media (M4rt) or NP wash/aspirate in sterile container	248997 (soft aluminum wire swab) 142081 (M4rt)	
Beta Strep Rapid Screen	Double swab	274146	
Coronavirus SARS-CoV-2 with Influenza AB RSV Xpress test (COV19 INFL RSV XPS)	NP or flocked swab in viral transport media (M4rt)	428998 (flocked swab) 142081 (M4rt)	
MRSA PCR Nasal Screen MRSA PCR SCN)	Transystem double swab	274146	

Legacy Hospitals / Labcorp Microbiology and Molecular Test Quick Reference Guide – Inpatient

Please refer to Lab Test Table on MyLegacy Intranet for most up to date specimen collection information

Test Name	Container	Infor / Lawson Ordering Information	
Nasal Complete MRSA/SA PCR	Transsystem double swab	274146	
Respiratory Panel with Coronavirus SARS-CoV-2 by PCR (RESP PAN COV19 PCR)	NP or flocked swab in viral transport media (M4rt)	428998 (flocked swab) 142081 (M4rt)	
Other			
Routine Bacterial Culture (aerobic only)	Single or double swab	140758 (single swab) 274146 (double swab)	
Routine Bacterial Culture (anaerobic and aerobic)	Port-A-Cul tube	105886	
Culture Fungus (not for skin, hair, or nail specimens)	Single or double swab	140758 (single swab) 274146 (double swab)	
MRSA SA Skin PCR Screen (MRSA SA PCR SKIN SCN)	Transsystem double swab	274146	
Herpes Simplex Virus 1/2 DNA by PCR	Transsystem double swab in viral transport media (M4rt)	274146 (swab) 142081 (M4rt)	

Legacy Hospitals / Labcorp Microbiology and Molecular Test Quick Reference Guide – Inpatient

Please refer to Lab Test Table on MyLegacy Intranet for most up to date specimen collection information

Test Name	Container	Infor / Lawson Ordering Information
-----------	-----------	-------------------------------------

