

Frequently Ordered Tests with Tube Types

For an updated listing of tests, please refer to our online directory of services at <https://cdos.halfpenny.com/labcorp/sacredheartlabs>

For questions, please contact the Laboratory at **(850) 416-7796**
or email **gcascensionoutreach@labcorp.onmicrosoft.com**

Green PST "Gel" Tube (Plasma) or Red SST "Gel" Tube (Serum Separator)

(Spin within 2 hours of collection)

Profiles

Comprehensive Metabolic Panel
Basic Metabolic Panel
Hepatic Function Profile
Electrolytes Panel
Renal Function Panel
Iron profile
Lipid profile



Single Tests

AFP
Albumin
Alcohol
Alk Phos
ALT
Amikacin
Amylase
AST
Beta HCG Quant
BUN
CA 125
Calcium
Carbamazepine Level (Tegretol)
Cardiac Enzymes
CEA
Chloride
Cholesterol HDL
CKMB
C02
Cortisol

Creatinine
Creatine Kinase
CRP
Digoxin
Dilantin
Direct Bilirubin
Estradiol
Ferritin
Folate
Free T3
Free T4
FSH
Gentamicin
GGT
Glucose
High Sens CRP
Homocysteine
IGE
Iron
LDL Cholesterol



Single Tests Continued

LH
Lidocaine
Lipase
Magnesium
Phenobarbital
Phosphorous
Potassium
Prealbumin
Progesterone
Prolactin
Sodium
Testosterone
Theophylline
Tobramycin
Total Bilirubin

Total IBC
Total Protein
Total T3
Total T4
Transferrin
Triglyceride
Troponin Quant
TSH
Uric Acid
Valproic Acid
Vancomycin
Vitamin B12
Vitamin D-25
Hydroxy

Red SST "Gel" Tube (Serum Separator)

(Spin within 2 hours of collection)

Serum Drug Screen
Mononucleosis Screen
Free PSA
Hepatitis Panel
ANA
ASO
PSA
Toxoplasma IGG
Toxoplasma IGM
Vitamin D 125 (send out)



Plain Red Tube (NO GEL)

Lithium



Lavender Top Tube

CBC
Hemogram
Platelets
Retics
Hemoglobin A1C (GHGB)

ESR (2 pediatric bullets for children)
BTNP (24 hr stability)
Insulin
Vitamin B1 (whole blood)



Blue Top Tube (Must be filled to black line on tube and GENTLY mixed well)

Protime (24 hr stability)
Fibrinogen

PTT (4 hr stability)
D-Dimer (8 hr stability)



Blood Bank Specimens

Pink Top Tube (Please ensure to fill to the line)

Required Information: Full Name, MRN, Date, Time of Collection, Employee # of phlebotomist(s)

Infants <4 Months 2 Micro Lavender Pedi

Microtainers must contain the full patient label.

NO FOOTIES. Samples with footies will be rejected!



Less Frequently Ordered Individual Tests with Tube Types

Short Plain Green Top (Lithium Heparin) NO GEL

- Ammonia (**on ice**)



Short Gray (Sodium Fluoride/Potassium Oxalate)

- Lactic Acid (**on ice**)



Tall Plain Green Top (Lithium Heparin)

NO GEL

- Quantiferon TB Gold

