





TEST	CONTAINER TYPE
Stool Culture (Routine), includes E.coli 0157:H7  Vibrio Stool Culture  Yersinia Stool Culture	<b>C&amp;S Fecal Swab</b>  
Stool Culture (Routine), includes E.coli 0157:H7  Vibrio Stool Culture  Yersinia Stool Culture	<b>C&amp;S PARA PAK</b>  
Parasite Identification  Ova, Cyst, and Parasite Exam w/o WBC  Giardia Antigen Stool by DFA  Giardia/Cryptosporidium Antigen Stool by DFA	<b>PVA and Formalin Vials</b>  
Parasite Identification WBCs Stool Occult Blood - Stool Stool Fat, Semi-Quantitative Stool pH STRP SRS Clostridium difficile Antigen / Toxin EIA Helicobacter Pylori Stool Antigen Rotavirus	<b>Stool Cup</b>  

<p>Chlamydia Trachomatis by NAAT N. Gonorrhoeae by NAAT Trichomonas Vaginalis NAAT</p>	<p><b>Xpert CT/NG Vaginal Endocervical Specimen Collection Kit</b></p> 
<p>Fluid Culture w/ Anaerobes and Smear (Body Fluid Culture)</p> <p>Wound Culture w/Anaerobes and Smear (Any wound Aspirate)</p>	<p><b>Sterile Syringe</b></p> 
<p>CSF Culture w/Smear Meningitis/Encephalitis Pathogen Panel by NAAT CSF Herpes Simplex Virus 1 / 2 DNA by NAAT (CSF) Varicella Zoster Culture CMV Culture</p>	<p><b>CSF Vials</b></p> 
<p>2019 Coronavirus SARS-CoV-2 Influenza A and B by PCR Influenza A and B by EIA Herpes Simplex Virus 1 / 2 DNA by PCR (NOT CSF) Herpes Simplex Virus "HSV" Culture Varicella Zoster Culture RSV "Respiratory Syncytial Virus" Respiratory Pathogens Panel by NAAT Bordetella Pertussis by NAAT Mycoplasma Pneumonia by NAAT</p>	<p><b>Viral-Universal Transport Media</b></p> 

<p>Chlamydia Trachomatis by NAAT Urine  N. Gonorrhoeae by NAAT  Trichomonas Vaginalis NAAT (Urine)  Influenza A and B by PCR  Influenza A and B by EIA  Urine Culture  Respiratory Culture Sputum  CR Respiratory Culture (Cystic Fibrosis Patients Only)  Bronchial Culture  Tissue Culture  Catheter Tip Culture  RSV "Respiratory Syncytial Virus"  CMV Culture (Urine)</p>	<p style="text-align: center;"><b>Sterile Container</b></p> 
<p>Wound Culture  Genital Culture  GC Culture Screen  Ear, Eye, Or Skin Culture  Throat/NP Culture  CR Respiratory Culture (Cystic Fibrosis Patients Only)  Strep Screen Confirmation (Group A Strep Culture)  Group A Streptococcus by NAAT  Group A Streptococcus by NAAT  MRSA Screen by NATT  VRE by NATT (Vancomycin Resistant Enterococcus)  Nasal MSSA and MRSA by NAAT</p>	<p style="text-align: center;"><b>E-Swab</b></p> 
<p>Cytomegalovirus (CMV) Qualitative NAAT (Saliva)</p>	