

## OnLine Test Menu



CODE	TEST NAME	CPT4 CODE
LAB5706	A. FUMIGATUS, ALLERGEN Performing Laboratory:LABCORP BURLINGTON	
LAB3185	ABORH TYPE Container:One 6 mL EDTA (pink) or One 4 mL EDTA (lavender) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	86900 86901
LAB2113	ABSOLUTE NEUTROPHIL COUNT WITH WBC Container:4 mL EDTA (lavender) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	85048
LAB43	ACETAMINOPHEN LEVEL Container:One 5 mL serum separator tube (gold) or one 3 mL plasma separator tube (light green) (Min: 1.0 mL whole blood) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	80329
LAB505	ACETONE, QUANTITATIVE, SERUM Container:5 ml serum separator tube or 3 ml plasma separator tube Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	82010
LAB836	ACETYLCHOLINE RECEPTOR, BINDING Container:One 5 ml serum separator tube (gold) Temperature:Refrigerated Performing Laboratory:ARUP Laboratories	83519
LAB838	ACETYLCHOLINE RECEPTOR, BLOCKING Container:One 5 ml serum separator tube (gold) Temperature:Refrigerated Performing Laboratory:ARUP Laboratories	83516
LAB837	ACETYLCHOLINE RECEPTOR, MODULATING Container:One 5 ml serum separator tube (gold) Temperature:Refrigerated Performing Laboratory:ARUP Laboratories	83516
LAB511	ACTH Performing Laboratory:LABCORP DENVER	
LAB325	ACTIVATED PARTIAL THROMBOPLASTIN TIME Container:One 2.7 mL 3.2% sodium citrate (lt. blue) (Min: 2.7 mL). Temperature:Frozen Performing Laboratory:CHI Health Laboratory	85730
LAB2336	ADAMTS13 (VWF CLEAV PROTEASE) Performing Laboratory:LABCORP BURLINGTON	
LAB2317	ADENOVIRUS DETECTION BY PCR Performing Laboratory:VIRACOR EUROFINS	
LAB877	AFB CULTURE Container:Sterile, non-leaking container. No swabs Temperature:All specimens other than CSF: Refrigerate Performing Laboratory:CHI Health Laboratory	87116
LAB246	AFB CULTURE, BLOOD Container:SPS Vacutainer, Heparin tube (dark green, no gel) Temperature:Room temperature Performing Laboratory:CHI Health Laboratory	87116
LAB559	AFP TUMOR MARKER Container:One 5 mL serum separator tube (gold) (Min: 1.0 mL whole blood) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	82105