

CODE	TEST NAME	CPT4 CODE
LAB5516	DAT POLY Container:One 6 mL EDTA tube or One 4 mL EDTA tube (lavender) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	86880
LAB313	D-DIMER, QUANTITATIVE Container:One 2.7 mL sodium citrate (lt. blue) (Min: 2.7 mL). Tube must be full to frosted line. Temperature:Frozen Performing Laboratory:CHI Health Laboratory	85379
LAB7350	DEAMINATED GLIADIN PEPTIDE-IGA (DGP-IGA) Container:Serum separator tube (gold) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	83516
LAB7351	DEAMINATED GLIADIN PEPTIDE-IGG (DGP-IGG) Container:Serum separator tube (gold) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	83516
LAB522	DHEA Performing Laboratory:LABCORP BURLINGTON	
LAB524	DHEA-SULFATE Performing Laboratory:LABCORP LINCOLN	
LAB23	DIGOXIN LEVEL Container:One 5 mL serum separator tube (gold) or one 3 mL plasma separator tube (light green) (Min: 1.0 mL whole blood) Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	80162
LAB792	DIPHTHERIA / TETANUS ANTIBODY PANEL Performing Laboratory:LABCORP BURLINGTON	
LAB2201	DIPHTHERIA AB IGG Performing Laboratory:LABCORP BURLINGTON	
LAB2129	DRUG SCREEN 5 URINE Container:Sterile container Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	80307
LAB2295	DRUG SCREEN 5 WITH CONFIRMATION Container:sterile container Temperature:Refrigerated Performing Laboratory:CHI Health Laboratory	80307
LAB2128	DRUG SCREEN 8 URINE Container:Sterile container Performing Laboratory:CHI Health Laboratory	80307
LAB2204	DRUG SCREEN 8 WITH CONFIRMATION Container:Sterile container Performing Laboratory:CHI Health Laboratory	80307