

LAB	CODE	TEST NAME
Evergreen	17OH-LC	17-OH Progesterone
Evergreen	UQHIAA-LC	5 Hydroxyindoleacetic Acid Ur Qnt
Evergreen	5NT-LC	5 Nucleotidase
LABCORP	791000	1,3-Dimethylamylamine, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for eight days at room temperature, refrigerated, or frozen.
LABCORP	761490	10+Oxycodone-Scr, Oral Fluid Container:Quantisal(TM) Oral Fluid collection device, contact lab for collection kit. Temperature:Maintain specimen at room temperature for transport to laboratory; refrigerate if transport is delayed.
LABCORP	500171	11-Deoxycortisol: HPLC/MS-MS Container:Gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube Temperature:Freeze
LABCORP	500550	11-Desoxycortisol HPLC/MS-MS Container:Red-top tube or gel-barrier tube Temperature:FREEZE
LABCORP	500560	11-Desoxycortisol, 2 Specs Container:Red-top tube, gel-barrier tube, lavender-top (EDTA) tube or green-top (heparin) tube Temperature:FREEZE
LABCORP	504677	11-Hydroxyandrostenedione Container:Serum from Red Top Tube OR EDTA Plasma tube OR Heparin plasma tube Temperature:Freeze Stability: Room temperature: 14 days Refrigerated: 14 days Frozen: 3 months Freeze/thaw: Stable x6
LABCORP	504680	11-Hydroxytestosterone Container:Serum from Red Top Tube OR EDTA Plasma tube OR Heprin Plasma Tube Temperature:Frozen Stability: Room Temperature: 14 days Refrigerated: 14 days Frozen: 3 months Freeze/Thaw: Stable x6
LABCORP	504674	11-Ketotestosterone Container:Serum from Red Top Tube OR EDTA Plasma tube OR Heprin Plasma Tube Temperature:Freeze Stability: Room temperature: 14 days Refrigerated: 14 days Frozen: 3 months Freeze/thaw: Stable x6
LABCORP	504683	11-oxo-Androgens Panel Container:Serum from Red Top Tube, EDTA Plasma tube, or Heprin Plasma tube Temperature:Freeze
LABCORP	504550	14.3.3 ETA, Rheumatoid Arthritis Container:Gel-barrier tube OR red-top tube Temperature:FREEZE Stability: Room temperature up to 3 days Refrigerated up to 3 days Frozen up to 7 days Freeze/thaw cycles: Stable x6
LABCORP	831163	17 Panel Nail Test Container:Foil in secure envelope Temperature:Maintain specimen at room temperature.
LABCORP	140715	17-Hydroxypregnenolone, Mass Spectrometry Container:Red-top tube, gel-barrier tube, or lavender top (EDTA) tube Temperature:FREEZE