

LAB	CODE	TEST NAME
Evergreen	C1Q-LC	C1Q Binding Assay
Evergreen	C3-LC	C3 Complement
Evergreen	C4-LC	C4 Complement
Evergreen	CAFF	CAFFEINE LEVEL (UW) Container:Gold SST Temperature:Refrigerated
Evergreen	CALC-LC	Calcitonin Level
Evergreen	LAB53	CALCIUM Container:Gold/Tiger (serum); Alternate: Mint green (lithium heparin); Red (serum)
Evergreen	ICA	Calcium, Ionized Container:mint green (lithium heparin) tube. Gold/Tiger (serum) tubes acceptable if allowed to clot completely first, and must NOT be uncapped (Opened).
Evergreen	ICAD	Calcium, Ionized (Dialysate) Container:Mint green (lithium heparin) tube. Gold/Tiger (serum) tubes acceptable if allowed to clot completely first, and must NOT be uncapped (Opened).
Evergreen	CAUQ	Calcium, Urine (Quant, Timed) Container:24-hour dark plastic urine container. Temperature:Refrigerated.
Evergreen	CALUR	Calcium, Urine (Random) Container:Sterile leak-proof plastic urine container.
Evergreen	CALPFC-LC	Calprotectin, Fecal
Evergreen	CA125-LC	Cancer Antigen 125
Evergreen	CA153-LC	Cancer Antigen 15 3
Evergreen	CA19 9-LC	Cancer Antigen 19 9
Evergreen	CA2729-LC	Cancer Antigen 27.29