

LAB	CODE	TEST NAME
Evergreen	OBP	OB Panel (BWNW) Container:1- 7mL Lavender, 1-m Lavender, 2 - 10 mL SST
Evergreen	OCBGU	Occult Blood, Gastric Container: Sterile Plastic Container with Lid
Evergreen	OLIGO(MS) -LC	Oligoclonal Band Profile (MS Panel)
Evergreen	OLIGO-LC	Oligoclonal Bands (IgG), CSF/Serum
Evergreen	OPI-LC	Opiate Screen (Reflex)
Evergreen	LAB107	OSMOLALITY Container: SST Temperature: Serum specimens should be tested as soon as possible after collection in order to obtain accurate osmolality results in clinical laboratories. When immediate processing of serum specimens is not possible, samples can be refrigerated or stored at room temperature within twenty-four (24) hours of collection. Prior to analysis, serum specimens must be warmed to room temperature.
Evergreen	OSMUR	Osmolality, Urine (Random) Container: Leakproof plastic urine container.
Evergreen	OSTEO-LC	Osteocalcin
Evergreen	OP	Ova and Parasites Container: Clean leakproof sterile container OR Totalfix Temperature: Store and transport refrigerated.
Evergreen	OXCAR-LC	Oxcarbazepine Metabolite Level
Evergreen	OXYMET-LC	Oxycodone and Metabolite Free
LABCORP	791128	o,m,p-Xylene (as Methylhippuric Acid), Occupational Exposure, Urine Container: Plastic urine container without preservative Temperature: Refrigerate; specimens may also be frozen.
LABCORP	790532	o,m,p-Xylene, Occupational Exposure, Whole Blood Container: Green-top (Heparin) tubes Temperature: Refrigerate; specimen may also be frozen.
LABCORP	602673	O001-IgE Cotton Lint Container: Red-top tube or gel-barrier tube Temperature: Room temperature
LABCORP	602674	O070-IgE Seminal Fluid Container: Red-top tube or gel-barrier tube Temperature: Room temperature