

LAB	CODE	TEST NAME
Evergreen	UFH	Unfractionated Heparin (anti-Xa) Container:Light Blue (Buffered 3.2% Sodium Citrate)
Evergreen	MYCO	Ureaplasma/Mycoplasma hominis Container:Viral, Chlamydia, or Mycoplasma culture transport provided by LabCorp, or other appropriate transport medium
Evergreen	URAC	Uric Acid Container:Gold/Tiger (serum)    Alternate: Mint green (lithium heparin); Red (serum)
Evergreen	URICBF	Uric Acid, Body Fluid Container:Collect 1 mL body fluid, centrifuge and separate to remove cellular material. Send frozen in plastic container
Evergreen	URICUQ	Uric Acid, Urine (Quant) Container:24 hour dark plastic urine container Temperature:Ambient (room temperature)
Evergreen	URICUR	Uric Acid, Urine (Random) Container:Sterile leakproof plastic urine container. Temperature:Ambient (room temperature)
Evergreen	LAB10123	URINALYSIS W/MICROSCOPIC, REFLEX TO CULTURE IF INDICATED Container:For Urinalysis: 1 Urine Yellow Top (No Additive) OR 1 Urine/Red/Yellow Tiger (w/Preservative) Top for Outpatient For Culture: 1 Urine Grey/Gray Top Acceptable sample types: Sterile cup, conical tube, or yellow BD round bottom tube Gray top tubes are not acceptable for Urinalysis' Chemistries, only for urine cultures Temperature:Refrigerated
Evergreen	UAM	Urinalysis with Microscopic Container:Sterile leakproof plastic container. Temperature:Store and transport refrigerated.
Evergreen	LAB347	URINALYSIS WITH REFLEX MICROSCOPIC, REFLEX CULTURE IF INDICATED Container:For Urinalysis: 1 Urine Yellow Top (No Additive) OR 1 Urine/Red/Yellow Tiger (w/Preservative) Top for Outpatient For Culture: 1 Urine Grey/Gray Top Acceptable sample types: Sterile cup, conical tube, or yellow BD round bottom tube Gray top tubes are not acceptable for Urinalysis' Chemistries, only for urine cultures Temperature:Refrigerated
Evergreen	LAB393	Urine Fentanyl Screen Container:Urine Cup Temperature:Room Temp
Evergreen	OBUR24	Urine OB U24 Container:See Individual Tests
Evergreen	UMO	Urine, Microscopic Only Container:Sterile leakproof plastic container. Temperature:Store and transport refrigerated.
LABCORP	511200	UGT1A1 Irinotecan Toxicity Container:Lavender-top (EDTA) tube, yellow-top (ACD) tube OR LabCorp buccal swab kit Temperature:Maintain whole blood specimen at room temperature for 7 days or refrigerated for 1 month, frozen for 2 years. Maintain buccal swabs at room temperature for 6 months. **Do NOT Freeze glass tubes**
LABCORP	470115	Uniparental Disomy (UPD) Parent, DNA Analysis Container:Lavender-Top (EDTA) Tube, Yellow-top (ACD) tube Temperature:Room temperature. Maintain specimen at room temperature. Do not freeze.
LABCORP	470074	Uniparental Disomy (UPD) Proband, DNA Analysis Container:Lavender-Top (EDTA) Tube, yellow-top (ACD) tube; sterile plastic conical tube or two confluent T-25 flasks for fetal testing Temperature:Room temperature. Maintain specimen at room temperature. Do not freeze.