



Online Test Menu

Laboratory Catalog

LAB	CODE	TEST NAME	CPT4 CODE
Mission Health	ZINCLC	Zinc, Plasma or Serum LC Place Of Service:LabCorp	
Mission Health	ZONILC	Zonisamide(Zonegran), Serum LC Place Of Service:LabCorp	
LABCORP	790660	Zaleplon, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80299
LABCORP	790710	Zaleplon, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80368
LABCORP	489000	ZAP-70 in B-CLL Container:Green-top (heparin) tube OR lavender-top (EDTA) tube Temperature:Store and ship whole blood at ambient temperature. Ship bone marrow at ambient temperature.	NOCPT
LABCORP	163084	Zika Virus Antibody, IgM Container:Gel-barrier tube OR Red-Top tube OR Serum Transfer tube Temperature:Refrigerate	86794
LABCORP	010170	Zinc Protoporphyrin (ZPP) Container:Royal blue-top (EDTA) tube OR tan-top lead-free tube; submit original tube. Temperature:Room temperature. Stable for five days at room temperature or for seven days when refrigerated.	84202
LABCORP	010180	Zinc Protoporphyrin (ZPP), Workplace Container:Royal blue-top (EDTA) tube OR tan-top lead-free tube; submit original tube. Temperature:Room temperature. Stable for five days at room temperature or for seven days when refrigerated.	84202
LABCORP	070111	Zinc Random Urine	84630 82570
LABCORP	070029	Zinc, RBC Container:Royal blue-top (EDTA) tube Temperature:Room Temperature	84630
LABCORP	003434	Zinc, Urine Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84630 82570
LABCORP	070032	Zinc, Whole Blood Container:Royal blue-top (EDTA) tube Temperature:Room Temperature	84630
LABCORP	790384	Ziprasidone, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:Maintain specimen at room temperature; stable for three days. Stable refrigerated 14 days	80299
LABCORP	503995	ZNT8 Antibodies Container:Red-top tube, gel-barrier tube or green-top(heparin) tube Temperature:Frozen Stability: Room temperature up to 7 days Refrigerated up to 7 days Frozen up to 7 days Freeze/Thaw Cycle: Stable x6	86341
LABCORP	808135	Zolpidem, Serum or Plasma Container:Red-top tube or green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for three days at room temperature, for three days refrigerated, and for two weeks frozen.	80299