

TEST CODE	TEST NAME	CPT4 CODE
LAB2297	JAK2 V617F Mutation Lab Mnemonic:JAK2M Container:Lavender Top (EDTA) Place Of Service:NorthShore University HealthSystem Department of Pathology	81270
LAB1147	JC VIRUS DETECT BY PCR, SCF Alert:** NEW LABCORP TEST ** Lab Mnemonic:JCVDX Place Of Service:LabCorp	
115101	JAK2 Exon 12, 13, 14 and 15 Mutation Analysis Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube OR yellow-top (ACD) tube Temperature:Refrigerate	81279
481020	JAK2 V617F Mutation Analysis, Quantitative Container:Lavender-top (EDTA) tube OR green-top (sodium heparin) tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
489422	JAK2 V617F Mutation Analysis, Quantitative with Reflex toCALR Mutation Analysis, JAK2 Exon 12-15 Mutation Analysisand MPL Mutation Analysis Container:Lavender-top (EDTA) tube OR green-top (sodium heparin) tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
489418	JAK2 V617F Mutation Analysis, Quantitative, With Reflex toCALR Mutation Analysis and MPL Mutation Analysis Container:Lavender-top (EDTA) tube OR green-top (sodium heparin) tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
489417	JAK2 V617F Mutation Analysis, Quantitative, With Reflex toJAK2 Exon 12-15 Mutation Analysis Container:Lavender-top (EDTA) tube OR green-top (sodium heparin) tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81270
826359	Jalapeno Pepper, IgE Container:Gel-barrier tube Temperature:Maintain specimen at room temperature	86003
139380	JC Virus DNA, PCR, Urine Container:Sterile urine container, no preservative Temperature:FROZEN (preferred) OR refrigerate	87798
139370	JC Virus DNA, PCR, Whole Blood Container:Lavender-top (EDTA) tube	87798
820115	JCV qPCR (Plasma) Container:Lavender-top (EDTA) tube Temperature:Room temperature or frozen	87799
820121	JCV qPCR (Serum) Container:Red-top tube Temperature:Room temperature or frozen	87799
511295	Jewish Ancestry Reproductive Profile Container:Lavender-top (EDTA) tube AND yellow-top (ACD) tube Temperature:Refrigerate yellow-top tube. Maintain lavender-top tube at room temperature.	83080 81200 81220
006030	Jka/Jkb Red Blood Cell (RBC) Antigen Typing Container:Lavender-Top (EDTA) tube Temperature:Specimens should be tested as soon as possible or within 72 hours of collection. Samples that cannot be tested immediately should be stored at 2 degrees Celsius to 8 degrees Celsius. Do NOT freeze red cells.	86905 x2
511490	Joubert Syndrome Type II Container:Lavender-top (EDTA) tube, yellow-top (ACD) tube, sterile plastic conical tube or two confluent T-25 flasks for fetal testing, OR LabCorp buccal swab kit Temperature:Room temperature 4(degrees)C	81479