

TEST CODE	TEST NAME	CPT4 CODE
LAB2490	K079-IgE PHTHALIC ANHYDRIDE Alert:** NEW LABCORP TEST ** Lab Mnemonic:MALEZ Place Of Service:LabCorp	
032425	K+Na Container:Transfer tube(1) Temperature:Refrigerate	84295 84132
602660	K072-IgE Ispaghula Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
064675	K072-IgE Ispaghula Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
650267	K080-IgE Formaldehyde Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
602963	K081-IgE Ficus Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
605018	K082-IgG Latex Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86001
602747	K083-IgE Cottonseed Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
605102	K083-IgG Cottonseed Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86001
605112	K084-IgG Sunflower Seed Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86001
602671	K087-IgE Alpha-Amylase Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
630542	Kallmann Syndrome Genetic Panel Container:Whole blood: Lavender-Top (EDTA) Tube; Oral swab: OCD-100 DNA Genotek device ONLY; or Extracted DNA: Sterile screw capped vial	81405 81479 81407 81406
500128	Kaolin Clotting Time (KCT) Container:Blue-top (sodium citrate) tube Temperature:Freeze Stability: Room temperature up to 4 hours	85347
844809	Kennedy Disease Mutation Anal. Container:Lavender-top (EDTA) tube Temperature:Maintain specimen at room temperature	81401
764258	Ketamine Container:Use plastic urine drug bottle and tamper-evident seal for forensic specimen. Collection kits are available by request from the laboratory. Temperature:Maintain specimen at room temperature Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307