

LAB	CODE	TEST NAME	CPT4 CODE
PRL	LAB16	ELECTROLYTE PANEL Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	80051
PRL	LAB565	ELECTROLYTES, URINE RANDOM Place of Service:LabCorp	
PRL	LAB565 (PRL)	ELECTROLYTES, URINE RANDOM Container:Random Urine in a sterile, leak-proof container (urine cup) without preservative Place of Service:OMD, OPH, OSV	
PRL	LAB24025	Electrolytes, Whole Blood Container:3.0 mL Light Green Top Tube (Lithium Heparin); Also acceptable: Heparinized Blood gas syringe Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	80051
PRL	LAB13085	Electrophoresis Serum, Reflex IFE Container:5.0 mL Gold Top Tube (SST) Also acceptable: 6.0 mL Red Top Tube (Plain) Place of Service:LabCorp Halsey	84165 86334
PRL	LAB13089	Electrophoresis, CSF Place of Service:LabCorp	84166 84157
PRL	LAB13088	Electrophoresis, Urine, Reflex IFE Synonyms:Electrophoresis, Protein, Random Urine With Reflex to IFE, Urine; Urine Electrophoresis, Random Urine With Reflex to IFE, Urine; Urine Protein Electrophoresis, Random Urine With Reflex to IFE, Urine Container:Random Urine in a sterile, leak-proof container (urine cup) without preservative Also acceptable: 24-hour or timed urine in a sterile leak-proof urine container without preservative Place of Service:LabCorp Halsey	86334 86334
PRL	LAB25202	ENCEPHALOPATHY AUTOIMMUNE/PARANEOPLASTIC EVAL, CSF Synonyms:Behavioral Change; Confusion; Encephalitis; Psychosis; Limbic encephalitis; Encephalopathy-Autoimmune Eval, CSF; ENC2; Mayo Test Code: ENC2; Epic Name: Encephalopathy Autoimmune Eval, CSF Container:CSF Sterile Tube Place of Service:Mayo - Rochester	86255 x19 86341 if reflexed, 86255 x7; 86256 x8; 84182 x7
PRL	LAB25203	ENCEPHALOPATHY AUTOIMMUNE/PARANEOPLASTIC EVAL, SERUM Synonyms:Autoimmune Encephalitis Panel; Behavioral Change; Confusion; Encephalitis; Limbic encephalitis; Psychosis; Encephalopathy-Autoimmune Eval, S; ENS2; Mayo Test Code: ENS2; Epic Name: Encephalopathy Autoimmune Eval, Serum Container:Two (2) - 6.0 mL Red Top Tube (Plain) ; Also acceptable: ; Two (2) - 5.0 mL Gold Top Tube (SST) Place of Service:Mayo - Rochester	86255 x19 86341 if reflexed, 83519; 84182 x7; 86255 x7; 86256 x8
PRL	LAB774	ENDOMYSIAL AB, IGA Place of Service:LABCORP SPOKANE	
PRL	LAB2090	ENTAMOEBA HISTOLYTICA AB Synonyms:E histolytica, Amebic Antibodies Place of Service:LABCORP BURLINGTON	
PRL	LAB1333	Enterovirus RNA, NAAT Container:BD Universal Viral Transport System (UTM); Also acceptable: Viral Transport Media (VTM); Sterile, Leak-Proof Container; CSF Sterile Tube; 5.0 mL White Top Tube (PPT Pearl); 4.0 mL Lavender Top Tube (EDTA); 6.0 mL Pink Top Tube (EDTA) Place of Service:Performed at Labcorp Halsey, Portland Oregon	87498
PRL	LAB328	Eosinophil Smear Place of Service:LabCorp	
PRL	LAB1891	Epidermal Growth Factor Receptor Mutation Synonyms:EGFR MUT Container:Formalin Fixed, Paraffin Embed tissue (FFPE) in 10% neutral buffered formalin Place of Service:ORL	81235 88381 G0452
PRL	LAB654	Epstein-Barr Antibody to Early D Antigen (EA-D), IgG Synonyms:EBV Early Antigen Antibodies, IgG; EBV-EA Antibodies; Early antigen-diffuse; EA(D)	