

LAB	CODE	TEST NAME	CPT4 CODE
PRL	LAB494	RAPID PLASMA REAGIN, QUAL Synonyms:Nontreponemal Test; Rapid Plasma Reagin; RPR Container:5.0 mL Gold Top Tube (SST); Also acceptable: 6.0 mL Red Top Tube (Plain) Place of Service:LabCorp Halsey	86592 if reflexed add 86593
PRL	LAB2541	RAS/RAF Targeted Mutation Panel, NGS Synonyms:RAS/RAF Container:Formalin Fixed, Paraffin Embed tissue (FFPE) in 10% neutral buffered formalin Place of Service:ORL	81210 81275 81311 88381 G0452
PRL	MISC (3117-01)	RBC Antigen Type by Molecular Method - Single Blood Group Place of Service:Bloodworks Northwest	
PRL	MISC (RBC Antigen Type by Serological Method - Common Antigens)	RBC Antigen Type by Serological Method - Common Antigens Container:6.0 mL Pink Top Tube (EDTA); Also acceptable: 4.0 mL Lavender Top Tube (EDTA)	
PRL	LAB1842	RBC, COUNT ONLY Container:4.0 mL Lavender Top Tube (EDTA); Cord blood is an acceptable source for H&H testing at PSVMC Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF, OPH PCI FRANZ, OPH PCI NEWBERG, OPH PCI WESTSIDE, OWF PCI WILLAMETTE FALLS	85018 85014
PRL	MISC (3117-04)	Red Cell Genotyping for Multiple Blood Groups Place of Service:Bloodworks Northwest	
PRL	LAB19	RENAL FUNCTION PANEL Container:3.0 mL Light Green Top Tube (Lithium Heparin); 5.0 mL Gold Top Tube (SST); 6.0 mL Red Top Tube (Plain) Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	80069
PRL	LAB532	RENIN ACTIVITY Synonyms:Plasma Renin Activity (PRA); PRA Place of Service:LABCORP BURLINGTON	
PRL	LAB1307	RESPIRATORY PATHOGEN PANEL,NAAT Container:BD Universal Viral Transport System (UTM); Also acceptable: Viral Transport Media (VTM); Sterile Saline (0.9%); Sterile Leak-Proof Container Place of Service:OHR, OMD, OMW, ONB, OPH, OSS, OSV, OWF	
PRL	LAB296 (PRL)	RETIC COUNT Container:4.0 mL Lavender Top Tube (EDTA); Also acceptable: 6.0 mL Pink Top Tube (EDTA) Place of Service:OMD, OPH, OSV	85045
PRL	LAB296	Reticulocyte Count Container:4.0 mL Lavender Top Tube (EDTA); Also acceptable: 6.0 mL Pink Top Tube (EDTA) Place of Service:LabCorp Halsey	85045
PRL	MISC (3117-07)	RhD Evaluation Place of Service:Bloodworks Northwest	
PRL	MISC (3117-09)	RhD Variant - Weak D Types 1,2,3 Place of Service:Bloodworks Northwest	
PRL	MISC (3117-05)	RhD Zygosity (including CDE systems) Place of Service:Bloodworks Northwest	
PRL	LAB206	Rheumatoid Factor Place of Service:Labcorp	