

Online Test Menu Laboratory Catalog

ricolcu	Cerreer	200010101, 0010108	
LAB	CODE	TEST NAME	CPT4 CODE
PSMMC	LAB24180	1,3 BETA-D-GLUCAN Place Of Service:KSL Diagnostics Inc via LABCORP	
PSMMC	LAB975	11-DEOXYCORTISOL Place Of Service:ESOTERIX ENDOCRINOLOGY	
PSMMC	LAB1760	14-3-3 PROTEIN, CSF (PRION DISEASE) Place Of Service:NATL PRION DISEASE PATH SURV	
PSMMC	LAB351	17-HYDROXYCORTICOSTEROIDS Place Of Service:ARUP LABORATORIES INC	
PSMMC	LAB2360	17-HYDROXYPREGNENOLONE Place Of Service:ESOTERIX ENDOCRINOLOGY	
PSMMC	LAB720	17-HYDROXYPROGESTERONE Place Of Service:LABCORP SPOKANE	
PSMMC	LAB16042	21-HYDROXYLASE AB Place Of Service:ESOTERIX ENDOCRINOLOGY	
PSMMC	LAB353	5 HIAA, QUANTITATIVE, URINE, RANDOM Place Of Service:LABCORP BURLINGTON	
PSMMC	LAB352	5 HIAA, URINE, 24HR, QUANT Place Of Service:LABCORP BURLINGTON	
LABCORP	071910	071910 7 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	074666	074666 9 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	500600	1,25-Dihydroxy,Vitamin D by MS Container:Gel-barrier tube, lavender-top (EDTA) tube, or green-top (heparin) tube Temperature:FREEZE Stability: Room temperature up to 2 days Refrigerated up to 2 days Frozen up to 544 days Freeze/Thaw cycle = 6	82652
LABCORP	791000	1,3-Dimethylamylamine, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for eight days at room temperature, refrigerated, or frozen.	80307
LABCORP	761490	10+Oxycodone-Scr, Oral Fluid Container:Quantisal(TM) Oral Fluid collection device, contact lab for collection kit. Temperature:Maintain specimen at room temperature for transport to laboratory; refrigerate if transport is delayed.	80307
LABCORP	768561	11+Oxyco+Crt-Bund+1-Confirm Container:Urine container Temperature:Maintain specimen at room temperature Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80346 80307