

LAB	CODE	TEST NAME	CPT4 CODE
SAH	17AHY	17-OH PROGESTERONE Place Of Service:LabCorp	
SAH	UOSMO	24 HR URINE FOR OSMOLALITY Place Of Service:LabCorp	
SAH	ULY24	24 HR. URINE LYTES Container:URINE TIMED SPECIMEN Place Of Service:MLAB	
SAH	24NA	24HR URINE SODIUM Container:Plastic urine container without preservative Temperature:Room temperature Place Of Service:MLAB	
SAH	5NUCL	5' NUCLEOTIDASE Place Of Service:LabCorp	
SAH	5HIAU	5-HIAA QN,RANDOM URINE Place Of Service:LabCorp	
SAH	HIA24	5-HIAA,QUANT.,24 HR URINE Place Of Service:LabCorp	
LABCORP	071910	071910 7 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	074666	074666 9 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	500600	1,25-Dihydroxy,Vitamin D by MS Container:Gel-barrier tube, lavender-top (EDTA) tube, or green-top (heparin) tube Temperature:FREEZE Stability: Room temperature up to 2 days Refrigerated up to 2 days Frozen up to 544 days Freeze/Thaw cycle = 6	82652
LABCORP	791000	1,3-Dimethylamylamine, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for eight days at room temperature, refrigerated, or frozen.	80307
LABCORP	761490	10+Oxycodone-Scr, Oral Fluid Container:Quantisal(TM) Oral Fluid collection device, contact lab for collection kit. Temperature:Maintain specimen at room temperature for transport to laboratory; refrigerate if transport is delayed.	80307
LABCORP	768561	11+Oxyco+Crt-Bund+1-Confirm Container:Urine container Temperature:Maintain specimen at room temperature Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80346 80307
LABCORP	500171	11-Deoxycortisol: HPLC/MS-MS Container:Gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube Temperature:Freeze	82634
LABCORP	500550	11-Desoxycortisol HPLC/MS-MS Container:Red-top tube or gel-barrier tube Temperature:FREEZE	82634