

LAB	CODE	TEST NAME	CPT4 CODE
SAH	GAD65	GAD-65 Place Of Service:LabCorp	
SAH	GASTR	GASTRIN, SERUM Place Of Service:LabCorp	
SAH	GCCHL	GC CHLAMYDIA BY NA AMPLIFICATION(GENITAL) Place Of Service:LabCorp	
SAH	GCCHA	GC CHLAMYDIA BY NA AMPLIFICATION(URINE) Place Of Service:LabCorp	
SAH	ALLGP	GENERAL ALLERGEN PROFILE Place Of Service:LabCorp	
SAH	GENTR	GENTA RANDOM Container:Red top(RED) Temperature:Refrigerate Place Of Service:MLAB	80170
SAH	GENTP	GENTAMICIN PEAK Container:Red top(RED) Temperature:Room Temperature Place Of Service:MLAB	80170
SAH	GENTT	GENTAMICIN TROUGH Container:Red top(RED) Temperature:Maintain specimen at room temperature Place Of Service:MLAB	80170
SAH	GFR	GFR Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature.[1] Place Of Service:MLAB	82565
SAH	GGT	GGT Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature. Place Of Service:MLAB	82977
SAH	GIACR	GIARDIA & CRYPTOSPORIDIUM EIA Place Of Service:LabCorp	
SAH	GIARE	GIARDIA LAMBLIA AG Place Of Service:LabCorp	
SAH	GL03P	GLU 3HR (GEST) Container:Grey top(T4C) Temperature:Maintain specimen at room temperature.[3] Place Of Service:MLAB	82951
SAH	GL1PP	GLU 1HR PP Container:Grey top(GYT) Temperature:Maintain at room temperature or refrigerate up to 14 days Place Of Service:MLAB	82947
SAH	FLGLU	GLUCOSE (FLUID) Container:Sterile container(FL1) Temperature:Room temperature Place Of Service:MLAB	82945