

LAB	CODE	TEST NAME	CPT4 CODE
SAH	LACID	LACTIC ACID Container:Grey top(GYT) Temperature:Room temperature Place Of Service:MLAB	83605
SAH	LAMO	LAMOTRIGINE, SERUM Place Of Service:LabCorp	
SAH	LDH	LDH Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature.[1] Place Of Service:MLAB	83615
SAH	FLLDH	LDH (FLUID) Container:Sterile container(FL1) Temperature:Room temperature Place Of Service:MLAB	83615
SAH	LDPED	LEAD (PEDIATRIC) Place Of Service:LabCorp	
SAH	24LDU	LEAD 24 HRS URINE Place Of Service:LabCorp	
SAH	LEGUA	LEGIONELLA AG UR Place Of Service:LabCorp	
SAH	LEGAB	LEGIONNAIRES AB IGG/IGM Place Of Service:LabCorp	
SAH	LEVE	LEVETIRACETAM Place Of Service:LabCorp	
SAH	LIDO	LIDOCAINE Container:Red-top tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of a gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Temperature:Room Temperature Place Of Service:MLAB	80176
SAH	LIPAS	LIPASE Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature. Place Of Service:MLAB	83690
SAH	LIPID	Lipid Panel(Coronary Risk) Container:Mint green top (lithium heparin)(LGR) Temperature:Maintain specimen at room temperature. Place Of Service:MLAB	80061
SAH	LIPOA	LIPOPROTEIN (A) Place Of Service:LabCorp	
SAH	LITH	LITHIUM Container:Red top(RED) Temperature:Refrigerate Place Of Service:MLAB	80178
SAH	LUPUS	LUPUS ANTICOAGULANT REFLEX Place Of Service:LabCorp	