

LAB	CODE	TEST NAME	CPT4 CODE
SAH	OB2	OBSTETRIC PNL Container:1 Lav top, 1 Pink Top, 1 SST and 1 Red Top Temperature:Refrigerate serum tubes and one lavender-top (EDTA) tube. One lavender-top (EDTA) tube may be stored at room temperature for 24 hours.	80055
SAH	OCCBG	OCCULT BLOOD GASTRIC Place Of Service:LabCorp	
SAH	OLIGB	OLIGOCLONAL BANDS Place Of Service:LabCorp	
SAH	ORGUR	ORGANIC ACID ANALYSIS, URINE Place Of Service:LabCorp	
SAH	ORGID	ORGANISM ID (YEAST) Place Of Service:LabCorp	
SAH	OSMOS	OSMOLALITY SERUM Place Of Service:LabCorp	
SAH	CLONP	OVA & PARASITE Place Of Service:LabCorp	
SAH	UOXAL	OXALATE 24 HR URINE Place Of Service:LabCorp	
SAH	OXCAR	OXCARBAZEPINE Place Of Service:LabCorp	
LABCORP	791128	o,m,p-Xylene (as Methylhippuric Acid), OccupationalExposure, Urine Container:Plastic urine container without preservative Temperature:Refrigerate; specimens may also be frozen.	82570 83921
LABCORP	790532	o,m,p-Xylene, Occupational Exposure, Whole Blood Container:Green-top (Heparin) tubes Temperature:Refrigerate; specimen may also be frozen.	84600
LABCORP	180851	O+P Exam, Formalin Only Container:Para Pak-Pink Temperature:Maintain specimen at room temperature.	87177
LABCORP	602673	O001-IgE Cotton Lint Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602674	O070-IgE Seminal Fluid Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	502051	Obstetric Complications Profile (Consultative) Container:Blue-top (sodium citrate) tube citrated platelet-free plasma AND lavender-top (EDTA) tube OR yellow-top (ACD) tube AND red-top tube OR gel-barrier tube. Temperature:FREEZE plasma/serum; maintain whole blood at room temperature. Stability: Whole blood: Room temperature up to seven days.	85306 86146 x3 85290 86147 x3 85300 85303 85613 x2 81240 81241 85598 83090 85730 86148 x2 85415 85732 x2