

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SJMC	LAB1786	1 Hour Gestational Glucose Tolerance Test Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	82950
SJMC	LAB975	11-DEOXYCORTISOL Place Of Service:Labcorp	
SJMC	LAB720	17-HYDROXYPROGESTERONE Place Of Service:Labcorp	
SJMC	LAB24841	17-OH PROGESTERONE Place Of Service:Esoterix	
SJMC	LAB88	2 Hour Gestational Glucose Tolerance Test Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	82951
SJMC	LAB1775	2 Hour Non-Gestational Glucose Tolerance Test Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	82951
SJMC	LAB2598	3 Hour Gestational Glucose Tolerance Test Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	82952
SJMC	LAB353	5 HIAA, QUANTITATIVE, URINE, RANDOM Place Of Service:Labcorp	
SJMC	LAB352	5 HIAA, URINE, 24HR, QUANT Place Of Service:Labcorp	
LABCORP	071910	071910 7 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	074666	074666 9 Drug-Scr Container:Plastic urine container Temperature:Maintain specimen at room temperature. Stability: If arrival at lab will extend beyond 7 days then refrigerate.	80307
LABCORP	500600	1,25-Dihydroxy,Vitamin D by MS Container:Gel-barrier tube, lavender-top (EDTA) tube, or green-top (heparin) tube Temperature:FREEZE Stability: Room temperature up to 2 days Refrigerated up to 2 days Frozen up to 544 days Freeze/Thaw cycle = 6	82652
LABCORP	791000	1,3-Dimethylamylamine, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Specimen is stable for eight days at room temperature, refrigerated, or frozen.	80307
LABCORP	761490	10+Oxycodone-Scr, Oral Fluid Container:Quantisal(TM) Oral Fluid collection device, contact lab for collection kit. Temperature:Maintain specimen at room temperature for transport to laboratory; refrigerate if transport is delayed.	80307
LABCORP	500171	11-Deoxycortisol: HPLC/MS-MS Container:Gel-barrier tube, lavender-top (EDTA) tube OR green-top (heparin) tube Temperature:Freeze	82634