

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SJMC	LAB16	Electrolyte Panel Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	80051
SJMC	LAB565	ELECTROLYTES, URINE RANDOM Place Of Service:St. Patrick Hospital	
SJMC	LAB2215	EOSINOPHIL SMEAR, RESPIRATORY Place Of Service:St. Patrick Hospital	89190
SJMC	LAB2199	EOSINOPHIL SMEAR, URINE Place Of Service:St. Patrick Hospital	89190
LABCORP	006018	E/e Red Blood Cell (RBC) Antigen Typing Container:Lavender-Top (EDTA) tube Temperature:Specimens should be tested as soon as possible or within 72 hours of collection. Samples that cannot be tested immediately should be stored at 2 degrees Celsius to 8 degrees Celsius. Do NOT freeze red cells.	86905 x2
LABCORP	605024	E001-IgG Cat Dander Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86001
LABCORP	602684	E003-IgE Horse Dander Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602685	E004-IgE Cow Dander Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602456	E005-IgE Dog Dander Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602686	E006-IgE Guinea Pig Epithelium Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602688	E070-IgE Goose Feathers Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602994	E071-IgE Mouse Epithelium Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602689	E072-IgE Mouse Urine Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602690	E073-IgE Rat Epithelia Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602691	E074-IgE Rat Urine Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003