

Online Test Menu Laboratory Catalog



| LAB | CODE | TEST NAME | CPT4 CODE |
|---------|----------|---|-------------|
| SJMC | LAB96 | Lactate Dehydrogenase Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 83615 |
| SJMC | LAB2377 | LACTATE DEHYDROGENASE, BODY FLUID Place Of Service:St. Patrick Hospital | 83615 |
| SJMC | LAB95 | Lactic Acid Container:GRAY-K OXALATE/NA FLUORIDE Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 83605 |
| SJMC | LAB99 | Lipase Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 83690 |
| SJMC | LAB13233 | LIPASE, BODY FLUID Place Of Service:St. Patrick Hospital | 83690 |
| SJMC | LAB2478 | Lipid Panel Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 80061 |
| SJMC | LAB29 | LITHIUM LEVEL Place Of Service:St. Patrick Hospital | 80178 |
| LABCORP | 007012 | Lacosamide Container:Red-top tube, lavender-top (EDTA) tube, OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Temperature:Maintain specimen at room temperature. | 80235 |
| LABCORP | 620044 | Lactate (CSF) Container:Sterile screw capped vial Temperature:Freeze as soon as possible after collection | 83605 |
| LABCORP | 023804 | Lactate and Pyruvate Container:2 Gray-top (sodium fluoride) tube Temperature:Refrigerate | 84210 83605 |
| LABCORP | 001842 | Lactate Dehydrogenase (LD) Isoenzymes Container:Red-top tube or gel-barrier tube Temperature:Maintain specimen at room temperature. Do NOT freeze. Do NOT refrigerate for more than three days. | 83615 83625 |
| LABCORP | 100156 | Lactic Acid Dehydrogenase (LD), Body Fluid Container:Sterile body fluid container Temperature:Room temperature | 83615 |
| LABCORP | 812030 | Lactic Acid, Body Fluid Container:Sterile container Temperature:FREEZE | 83605 |
| LABCORP | 122408 | Lactic Acid, CSF Container:Clean sterile tube Temperature:Refrigerate for up to 15 days or freeze. Stable at room temperature up to 3 days. | 83605 |
| LABCORP | 120543 | Lactic/Pyruvic Ratio Container:Gray-Top (Sodium Fluoride) Tube | 84210 83605 |