

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SJMC	LAB695	Occult Blood Stool, For Colorectal Neoplasm Screening Container: Sterile Container/Urine Cup, Occult Blood Cards Temperature: Refrigerated Place Of Service: Providence St Joseph Medical Center Laboratory	82270
SJMC	LAB694	Occult Blood Stool, Not Colorectal Neoplasm Screen Container: Sterile Container/Urine Cup, Occult Blood Cards Temperature: Refrigerated Place Of Service: Providence St Joseph Medical Center Laboratory	82272
SJMC	LAB107	OSMOLALITY, SERUM Place Of Service: St. Patrick Hospital	83930
SJMC	LAB420	OSMOLALITY, URINE Place Of Service: St. Patrick Hospital	83935
LABCORP	791128	o,m,p-Xylene (as Methylhippuric Acid), Occupational Exposure, Urine Container: Plastic urine container without preservative Temperature: Refrigerate; specimens may also be frozen.	82570 83921
LABCORP	790532	o,m,p-Xylene, Occupational Exposure, Whole Blood Container: Green-top (Heparin) tubes Temperature: Refrigerate; specimen may also be frozen.	84600
LABCORP	180851	O+P Exam, Formalin Only Container: Para Pak-Pink Temperature: Maintain specimen at room temperature.	87177
LABCORP	602673	O001-IgE Cotton Lint Container: Red-top tube or gel-barrier tube Temperature: Room temperature	86003
LABCORP	602674	O070-IgE Seminal Fluid Container: Red-top tube or gel-barrier tube Temperature: Room temperature	86003
LABCORP	502051	Obstetric Complications Profile (Consultative) Container: Blue-top (sodium citrate) tube citrated platelet-free plasma AND lavender-top (EDTA) tube OR yellow-top (ACD) tube AND red-top tube OR gel-barrier tube. Temperature: FREEZE plasma/serum; maintain whole blood at room temperature. Stability: Whole blood: Room temperature up to seven days.	85306 86146 x3 85290 86147 x3 85300 85303 85613 x2 81240 81241 85598 83090 85730 86148 x2 85415 85732 x2
LABCORP	790690	o-Cresol, Occupational Exposure, Urine Container: Plastic urine container without preservative Temperature: Refrigerated; specimen may also be frozen	84600
LABCORP	832805	Octopus (Octopus vulgaris) IgE Container: Serum Separator (Gel) Tube Temperature: Room temperature	86003
LABCORP	819217	Octreotide (Serum/Plasma) Container: Red-top tube or lavender-top (EDTA) tube	80299
LABCORP	811513	Olanzapine, Serum or Plasma Container: Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature: SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80299
LABCORP	019216	Oligoclonal Banding, Serum and Cerebrospinal Fluid Container: Red-top tube or gel-barrier tube AND plastic (CSF) tube Temperature: Room temperature	83916