

Online Test Menu Laboratory Catalog



| LAB | CODE | TEST NAME | CPT4 CODE |
|---------|----------|--|-------------|
| SJMC | LAB34 | Salicylate Level Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 80307 |
| SJMC | LAB1788 | Sedimentation Rate Container:LAV-EDTA or PINK-K2EDTA Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 85652 |
| SJMC | LAB891 | Semen Analysis, Post-Vas Container:Clean Cup Place Of Service:Providence St Joseph Medical Center Laboratory | 89321 |
| SJMC | LAB2325 | Shigatoxin, 1&2 EIA Container:CARY-BLAIR MEDIA, alt: clean cup Temperature:Room Temp Place Of Service:Providence St Joseph Medical Center Laboratory | 87427 87015 |
| SJMC | LAB122 | Sodium Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 84295 |
| SJMC | LAB197 | SODIUM, BODY FLUID Place Of Service:St. Patrick Hospital | 84302 |
| SJMC | LAB2475 | SODIUM, URINE, 24HR Place Of Service:St. Patrick Hospital | 81050 84300 |
| SJMC | LAB444 | SODIUM, URINE, RANDOM Place Of Service:St. Patrick Hospital | 84300 |
| SJMC | LAB23996 | Source Patient Exposure Panel | |
| SJMC | LAB13048 | Specific Gravity, Urine Container:Clean Cup Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 81003 |
| SJMC | LAB885 | Strep A Screen, Rapid Container:Sterile Culture Swab with transport media, alt: Sterile Cotton/Poly Swab Temperature:Room Temp Place Of Service:Providence St Joseph Medical Center Laboratory | 87880 |
| SJMC | LAB1371 | Strep B DNA Probe, NAAT Container:Sterile Culture Swab with transport media Temperature:Room Temp Place Of Service:Providence St Joseph Medical Center Laboratory | 87653 |
| SJMC | LAB1856 | STREPTOCOCCUS GROUP A SCREEN, REFLEX CULTURE Container:Sterile Culture Swab with transport media, alt: Sterile Cotton/Poly Swab Temperature:Room Temp Place Of Service:Providence St Joseph Medical Center Laboratory | 87880 |
| LABCORP | 006040 | S/s Red Blood Cell (RBC) Antigen Typing Container:Lavender-Top (EDTA) tube Temperature:Specimens should be tested as soon as possible or within 72 hours of collection. Samples that cannot be tested immediately should be stored at 2 degrees Celsius to 8 degrees Celsius. Do NOT freeze red cells. | 86905 x2 |
| LABCORP | 140110 | S-100B Protein, Serum Container:Red-top tube or serum-separator tube Temperature:Refrigerate | 86316 |