

Online Test Menu Laboratory Catalog



| LAB | CODE | TEST NAME | CPT4 CODE |
|---------|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------|
| SJMC | LAB24 | Valproic Acid Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 80164 |
| SJMC | LAB40 | Vancomycin Level Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 80202 |
| SJMC | LAB41 | Vancomycin, Peak Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 80202 |
| SJMC | LAB39 | Vancomycin, Trough Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 80202 |
| SJMC | LAB2466 | VITAMIN B-12 Place Of Service:St. Patrick Hospital | 82607 |
| SJMC | LAB2301 | Vitamin D Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory | 82306 |
| LABCORP | 180042 | Vaginitis (VG) With Candida (Six Species), NuSwab(R) Container:APTIMA(R) vaginal or unisex swab Temperature:Room temperature | 87801 87661 87798 x3 |
| LABCORP | 180071 | Vaginitis (VG), Herpes Simplex Virus (HSV), NuSwab(R) Container:Aptima(R) Multitest Swab Specimen Collection Kit or Aptima(R) Unisex Swab Specimen Collection Kit Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C) | 87801 87661 87798 x3 87529 x2 |
| LABCORP | 180039 | Vaginitis (VG), NuSwab(SM) Container:APTIMA(R) vaginal or unisex swab Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C) | 87798 x3 87801 87661 |
| LABCORP | 180068 | Vaginitis Plus (VG+) With Candida (Six Species), NuSwab(R) Container:Aptima(R) Multitest Swab Specimen Collection Kit or Aptima(R) Unisex Swab Specimen Collection Kit Temperature:Room temperature | 87801 87591 87661 87798 x3 87491 |
| LABCORP | 180066 | Vaginitis Plus (VG+), Herpes Simplex Virus (HSV), NuSwab(R) Container:Aptima(R) Multitest Swab Specimen Collection Kit or Aptima(R) Unisex Swab Specimen Collection Kit Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C) | 87591 87661 87798 x3 87529 x2 87491 87801 |
| LABCORP | 180021 | Vaginitis Plus (VG+), NuSwab(SM) Container:APTIMA(R) vaginal or unisex swab Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C) | 87491 87801 87591 87661 87798 x3 |
| LABCORP | 180026 | Vaginitis/Vaginosis, DNA Probe Container:Affirm VPIII Ambient Temperature Transport System (ATTS) Temperature:Maintain specimen at room temperature for up to 72 hours. | 87660 87480 87510 |
| LABCORP | 706499 | Valproic Acid (Total+Free) Container:Red-stopper tube. Transfer separated serum to a plastic transport tube. Do not use a serum separator tube; the use of serum separator tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Temperature:Room Temperature | 80164 80165 |
| LABCORP | 070789 | Valproic Acid, Free, Serum or Plasma Container:Red-top tube OR lavender-top (EDTA) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Temperature:Room Temperature | 80165 |