

Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SJMC	LAB24	Valproic Acid Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	80164
SJMC	LAB40	Vancomycin Level Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	80202
SJMC	LAB41	Vancomycin, Peak Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	80202
SJMC	LAB39	Vancomycin, Trough Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	80202
SJMC	LAB2466	VITAMIN B-12 Place Of Service:St. Patrick Hospital	82607
SJMC	LAB2301	Vitamin D Container:GREEN-LI HEP, alt: GOLD-SST or RED TOP-PLAIN Temperature:Refrigerated Place Of Service:Providence St Joseph Medical Center Laboratory	82306
LABCORP	180042	Vaginitis (VG) With Candida (Six Species), NuSwab(R) Container:APTIMA(R) vaginal or unisex swab Temperature:Room temperature	87801 87661 87798 x3
LABCORP	180071	Vaginitis (VG), Herpes Simplex Virus (HSV), NuSwab(R) Container:Aptima(R) Multitest Swab Specimen Collection Kit or Aptima(R) Unisex Swab Specimen Collection Kit Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87801 87661 87798 x3 87529 x2
LABCORP	180039	Vaginitis (VG), NuSwab(SM) Container:APTIMA(R) vaginal or unisex swab Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87798 x3 87801 87661
LABCORP	180068	Vaginitis Plus (VG+) With Candida (Six Species), NuSwab(R) Container:Aptima(R) Multitest Swab Specimen Collection Kit or Aptima(R) Unisex Swab Specimen Collection Kit Temperature:Room temperature	87801 87591 87661 87798 x3 87491
LABCORP	180066	Vaginitis Plus (VG+), Herpes Simplex Virus (HSV), NuSwab(R) Container:Aptima(R) Multitest Swab Specimen Collection Kit or Aptima(R) Unisex Swab Specimen Collection Kit Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87591 87661 87798 x3 87529 x2 87491 87801
LABCORP	180021	Vaginitis Plus (VG+), NuSwab(SM) Container:APTIMA(R) vaginal or unisex swab Temperature:Maintain specimen at room temperature or refrigerate (2 degrees C to 30 degrees C)	87491 87801 87591 87661 87798 x3
LABCORP	180026	Vaginitis/Vaginosis, DNA Probe Container:Affirm VPIII Ambient Temperature Transport System (ATTS) Temperature:Maintain specimen at room temperature for up to 72 hours.	87660 87480 87510
LABCORP	706499	Valproic Acid (Total+Free) Container:Red-stopper tube. Transfer separated serum to a plastic transport tube. Do not use a serum separator tube; the use of serum separator tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Temperature:Room Temperature	80164 80165
LABCORP	070789	Valproic Acid, Free, Serum or Plasma Container:Red-top tube OR lavender-top (EDTA) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gel-barrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Temperature:Room Temperature	80165