

LAB	MNEMONIC	TEST NAME	CPT CODE
SPH	ZAP70	ZAP-70 in B-CLL	
SPH	ZINC	Zinc, Plasma or Serum	
SPH	ZIPRA	Ziprasidone	
SPH	ZONIS	Zonisamide(Zonegran), Serum	
LABCORP	790660	Zaleplon, Serum or Plasma Container:Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended. Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80299
LABCORP	790710	Zaleplon, Urine Container:Plastic urine container without preservative Temperature:SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.	80368
LABCORP	489435	ZAP-70 Panel (ASR) Container:Green-top (heparin) tube OR lavender-top (EDTA) tube	88187 88184 88185 x4
LABCORP	163084	Zika Virus Antibody, IgM Container:Gel-barrier tube OR Red-Top tube OR Serum Transfer tube Temperature:Refrigerate	86794
LABCORP	010170	Zinc Protoporphyrin (ZPP) Container:Royal blue-top (EDTA) tube OR tan-top lead-free tube; submit original tube. Temperature:Room temperature. Stable for five days at room temperature or for seven days when refrigerated.	84202
LABCORP	010180	Zinc Protoporphyrin (ZPP), Workplace Container:Royal blue-top (EDTA) tube OR tan-top lead-free tube; submit original tube. Temperature:Room temperature. Stable for five days at room temperature or for seven days when refrigerated.	84202
LABCORP	070111	Zinc Random Urine	84630 82570
LABCORP	070029	Zinc, RBC Container:Royal blue-top (EDTA) tube Temperature:Room Temperature	84630
LABCORP	003434	Zinc, Urine Container:Plastic urine container, no preservative Temperature:Maintain specimen at room temperature.	84630 82570
LABCORP	070032	Zinc, Whole Blood Container:Royal blue-top (EDTA) tube Temperature:Room Temperature	84630
LABCORP	503995	ZNT8 Antibodies Container:Red-top tube, gel-barrier tube or green-top(heparin) tube Temperature:Frozen Stability: Room temperature up to 7 days Refrigerated up to 7 days Frozen: 64 days Freeze/Thaw Cycle: Stable x6	86341