

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SPH	LAB2127	JAK2 MUTATION Place Of Service:GenPath	
SPH	LAB24732	JAK2 MUTATION QUAL WITH REFLEX TO CALR/MPL Place Of Service:GenPath	
SPH	LAB24787	JAK2 V617F MUTATION ANALYSIS, QUAL, W REFLEX TO JAK2 EXON 12-15 Place Of Service:GenPath	
SPH	LAB24206	JC VIRUS AB W/INDEX, QUAL, REFLEX INHIB ASSAY Place Of Service:LabCorp	
SPH	LAB24850	JC VIRUS DNA, PCR, QUALITATIVE, CSF Place Of Service:LabCorp	
SPH	LAB770	JO 1 AB, IGG Place Of Service:LabCorp	
LABCORP	489514	JAK2 E12-15 + CALR + MPL NGS Panel Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81219 81279 81338
LABCORP	489501	JAK2 Exon 12, 13, 14 and 15 Mutation Analysis by NGS Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81279
LABCORP	489615	JAK2 V617F + CALR + MPL NGS Panel Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81219 81270 81338
LABCORP	139370	JC Virus DNA, PCR, Whole Blood Container:Lavender-top (EDTA) tube	87798
LABCORP	820115	JCV qPCR (Plasma) Container:Lavender-top (EDTA) tube Temperature:Room temperature or frozen	87799
LABCORP	820121	JCV qPCR (Serum) Container:Red-top tube Temperature:Room temperature or frozen	87799