

# Online Test Menu Laboratory Catalog



LAB	CODE	TEST NAME	CPT4 CODE
SPH	LAB24708	UGT1A1 IRINOTECAN TOXICITY Place Of Service:LabCorp	
SPH	LAB12964	UREA NITROGEN, BODY FLUID Container:Green top, Red top, or Sterile Container from the Body Fluid collection kit Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	84520
SPH	LAB123748	UREA NITROGEN, URINE, 24HR Container:24 Hr Urine Container Temperature:24 hr Urine Container should be stored on ice or refrigerated. Specimen should be refrigerated for transport. IF you are not sending the entire collection, the Total Volume of the entire collection needs to be recorded on the specimen.	84540
SPH	LAB12963	UREA NITROGEN, URINE, RANDOM Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	84540
SPH	LAB23927	UREAPLASMA/MYCOPLASMA HOMINIS Place Of Service:LabCorp	
SPH	LAB141	URIC ACID Container:Gel-Barrier Tube is preferred (Send Entire Tube). Green tops preferable, with gold tops or Red Tops being acceptable. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	84550
SPH	LAB12960	URIC ACID, BODY FLUID Container:Green top, Red top, or Sterile Container from the Body Fluid collection kit Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	
SPH	LAB841	URIC ACID, URINE, 24 HOUR Place Of Service:LABCORP SEATTLE	
SPH	LAB347	URINALYSIS Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	
SPH	LAB348	URINALYSIS WITH MICROSCOPIC Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	81001
SPH	LAB2479	URINALYSIS WITH MICROSCOPIC IF INDICATED Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	
SPH	LAB2480	URINALYSIS WITH MICROSCOPIC WITH CULTURE IF INDICATED Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	
SPH	LAB1828	URINALYSIS, DIPSTICK ONLY Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	81003
SPH	LAB2209	URINALYSIS, MICROSCOPIC ONLY Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	81015
SPH	LAB12834	URINALYSIS, REFLEX MICROSCOPIC AND/OR CULTURE Container:Urine Specimen Cup or Sterile Container. Temperature:Ambient for immediate testing, or refrigerated for storage and transport to lab.	