

LAB	CODE	TEST NAME	CPT4 CODE
BSHSI-R	JAKL	JAK2 MUTATION ANALYSIS	
BSHSI-R	JAKRL	JAK2 V617F QUAL, RFX E12-15	
BSHSI-R	JAKRML	JAK2 V617F,RFLX CALR/MPL	
BSHSI-R	JCVL	JC VIRUS DNA BY PCR, QL	
BSHSI-R	JCVCL	JC VIRUS DNA, CSF, BY PCR	
BSHSI-R	JO1L	JO-1 AB	
LABCORP	489514	JAK2 E12-15 + CALR + MPL NGS Panel Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81219 81279 81338
LABCORP	489501	JAK2 Exon 12, 13, 14 and 15 Mutation Analysis by NGS Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81279
LABCORP	489615	JAK2 V617F + CALR + MPL NGS Panel Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C.	81219 81270 81338
LABCORP	826359	Jalapeno Pepper, IgE Container:Gel-barrier tube Temperature:Maintain specimen at room temperature	86003
LABCORP	139380	JC Virus DNA, PCR, Urine Container:Sterile urine container, no preservative Temperature:FROZEN (preferred) OR refrigerate	87798
LABCORP	820115	JCV qPCR (Plasma) Container:Lavender-top (EDTA) tube Temperature:Room temperature or frozen	87799
LABCORP	820121	JCV qPCR (Serum) Container:Red-top tube Temperature:Room temperature or frozen	87799
LABCORP	006030	Jka/Jkb Red Blood Cell (RBC) Antigen Typing Container:Lavender-Top (EDTA) tube Temperature:Specimens should be tested as soon as possible or within 72 hours of collection. Samples that cannot be tested immediately should be stored at 2 degrees Celsius to 8 degrees Celsius. Do NOT freeze red cells.	86905 x2
LABCORP	828292	Juniper Mix IgE Container:Red-top tube or gel-barrier tube Temperature:Maintain specimen at room temperature	86005