

LAB	CODE	TEST NAME	CPT4 CODE
BSHSI-R	KSTNSL	KIDNEY STONE PANEL	
BSHSI-R	FTD	Kleihauer-Betke Stain for Fetal Hemoglobin Container:One Pink EDTA tube 7-ml or one Lavender EDTA tube Temperature:Room temperature	85460
BSHSI-R	KOHP	KOH Prep Container:Swab in a sterile, additive free tube. Add saline to cover swab tip. Break off swab. Temperature:Room Temperature or Refrigerated (at 2-8°C)	87220
BSHSI-R	KOHS	KOH, Hair, Skin, or Nail Container:Sterile Collection Container Temperature:Room Temperature or Refrigerated (at 2-8°C)	87220
LABCORP	032425	K+Na Container:Transfer tube(1) Temperature:Refrigerate	84295 84132
LABCORP	602660	K072-IgE Ispaghula Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	064675	K072-IgE Ispaghula Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602662	K075-IgE Isocyanate TDI Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	650036	K075-IgE Isocyanate TDI Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	650044	K076-IgE Isocyanate MDI Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602663	K076-IgE Isocyanate MDI Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602665	K077-IgE Isocyanate HDI Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	650051	K077-IgE Isocyanate HDI Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	650069	K078-IgE Ethylene Oxide Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003
LABCORP	602666	K078-IgE Ethylene Oxide Container:Red-top tube or gel-barrier tube Temperature:Room temperature	86003