

LAB	CODE	TEST NAME	CPT4 CODE
BSSC	URCLRL	Urea clearance	
BSSC	UNUL	Urea Nitrogen, 24 hr Urine	
BSSC	URICA	Uric Acid Container:Plasma Separator Tube (PST) or Serum Separator Tube (SST) Temperature:Ambient	84550
BSSC	FUAL	Uric Acid, Fluid	
BSSC	UA	Urinalysis Container:No additive yellow top conical Urine Tube Temperature:Refrigerate	
BSSC	UC	Urine Culture Container:BD Gray Urine Tube Temperature:BD Gray Urine Tube transport room temp or refrigerate.	
BSSC	GUT1	Urine Glucose - 1 Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUT2	Urine Glucose - 2 Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUT3	Urine Glucose - 3 Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUT4	Urine Glucose - 4 Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUT5	Urine Glucose - 5 Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUT6	Urine Glucose - 6 Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUF	Urine Glucose - Fasting Container:No additive Urine tube Temperature:Refrigerate	
BSSC	GUH	Urine Glucose - Half Hour Container:No additive Urine tube Temperature:Refrigerate	
BSSC	UMIC	Urine Microscopic Container:No additive yellow top conical Urine Tube Temperature:Refrigerate	