

LABORATORY CATALOG

Where good health begins.

| LAB | CODE | TEST NAME | CPT4 CODE |
|---------|---------|---|-------------------|
| LRH | JAK2 | JAK2 V617F MUTATION | |
| LRH | JCPCR | JC POLYOMA VIRUS DNA | |
| LRH | CSFJC | JC POLYOMA VIRUS DNA, CSF | |
| LRH | JCVINDX | JCVirus AB with Index | |
| LRH | JO1 | JO-1 IGG AUTOANTIBODIES | |
| LABCORP | 489514 | JAK2 E12-15 + CALR + MPL NGS Panel Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81219 81279 81338 |
| LABCORP | 489501 | JAK2 Exon 12, 13, 14 and 15 Mutation Analysis by NGS Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81279 |
| LABCORP | 489615 | JAK2 V617F + CALR + MPL NGS Panel Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81219 81270 81338 |
| LABCORP | 481020 | JAK2 V617F Mutation Analysis, Quantitative Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube, Tan-top (K2-EDTA) tube, OR Pink-top (K2-EDTA) tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81270 |
| LABCORP | 489555 | JAK2V617F Mutation Analysis, Quantitative With Reflex toCALR Mutation Analysis, JAK2 Exon 12-15 Mutation Analysisand MPL Mutation Analysis by NGS Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81270 |
| LABCORP | 482960 | JAK2V617F Mutation Analysis, Quantitative, With Reflex toCALR Mutation Analysis and MPL Mutation Analysis by NGS Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81270 |
| LABCORP | 489540 | JAK2V617F Mutation Analysis, Quantitative, With Reflex toJAK2 Exon 12-15 Mutation Analysis by NGS Container:Lavender-top (EDTA) tube, green-top (sodium heparin) tube yellow-top (ACD) tube, Tan-top (K2-EDTA) Tube, or pink-top (K2-EDTA) Tube Temperature:Ship at room temperature. If specimen is to be stored prior to shipment, store at 2(degree)C to 8(degree)C. | 81270 |
| LABCORP | 826359 | Jalapeno Pepper, IgE Container:Gel-barrier tube Temperature:Maintain specimen at room temperature | 86003 |
| LABCORP | 139380 | JC Virus DNA, PCR, Urine Container:Sterile urine container, no preservative Temperature:FROZEN (preferred) OR refrigerate | 87798 |
| LABCORP | 820115 | JCV qPCR (Plasma) Container:Lavender-top (EDTA) tube Temperature:Room temperature or frozen | 87799 |