

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	DXY11	11-DESOXYCORTISOL Place Of Service:LABCORP	
UHEALTH Labs	ETAPRO	14.3.3 eta Protein Container:GoldTop Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	RA143B	14.3.3 ETA, RHEUM. ARTHRITIS Place Of Service:LABCORP	
UHEALTH Labs	SRA	14C-SEROTONIN RELEASE ASSAY, UFH Container:RED Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	16SR	16S rRNA Next Generation Sequencing Container:Collect tissue or fluid into a DNA free container labeled with at least two identifiers. Temperature:dry ice Place Of Service:UHEALTH LABS	87801
UHEALTH Labs	HYDX17	17-HYDROXYCORTICOSTEROIDS Place Of Service:LABCORP	
UHEALTH Labs	HPREG	17-HYDROXYPREGNENOLONE Container:No additive (red-top) (preferred),SST (red-top/glass) Place Of Service:QUEST	
UHEALTH Labs	UKETOS	17-KETOSTEROIDS W/CREATININE 24HR URN Container:PLASTIC URINE CONTAINER Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	17OHCR	17-OH P/CREAT RATIO, URINE Place Of Service:LABCORP	
UHEALTH Labs	17OHLC	17-OH PROGESTERONE LCMS Place Of Service:LABCORP	
UHEALTH Labs	17OPAS	17-OH Progesterone. LC/MS, ACTH Stimulation Place Of Service:LabCorp	
UHEALTH Labs	COVIDS	2019 Novel Coronavirus (COVID-19), NAA using Saliva Collection Place Of Service:LABCORP	
UHEALTH Labs	25HYD	25-HYDROXYVITAMIN D LCMS D2+D3 Place Of Service:LABCORP	
UHEALTH Labs	HMGCR	3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase Ab,IgG Container:Serum Separator Tube Temperature:Refrigerated Place Of Service:ARUP	83516
UHEALTH Labs	5NUCL	5' NUCLEOTIDASE Place Of Service:LABCORP	
UHEALTH Labs	FLUOCY	5-FLUOROCYTOSINE LEVEL,SERUM Container:Plain red Temperature:STORAGE/FROZEN Place Of Service:ARUP	
UHEALTH Labs	HIAA5Q	5-HIAA QN, RANDOM URINE Place Of Service:LABCORP	
UHEALTH Labs	5HIAAQ	5-HIAA,QUANT., 24 HR URINE Place Of Service:LABCORP	
UHEALTH Labs	ALPHA1	A-1-ANTITRYPSIN QN & PHENO Container:SERUM SEPARATOR TUBE (SST) Alt. Spec. Types: SERUM COLLECTED IN: SST (RED-TOP/PLASTIC) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	A1GLYO	A1C W/GLYCOMARK (R) REFLEX Place Of Service:LABCORP	
UHEALTH Labs	ABSCI	AB Screen Inst Container:Pink (EDTA) Tube Temperature:Transport at room temperature	86850

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	ABL1K	Place Of Service:UHealth Labs ABL1 KINASE DOMAIN MUTATION Container:EDTA (purple-top) tube Temperature:Refrigerated cold pack Place Of Service:UHEALTH LAB	81170
UHEALTH Labs	ABOTP	ABO GROUP AND RH TYPE Container:ACD-A or ACD-B (yellow-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ABORHO	ABO GROUPING AND RHO(D) TYPING Place Of Service:LABCORP	
UHEALTH Labs	ABRHI	ABO RH Inst Container:Pink (EDTA) Tube Temperature:Transport at room temperature Place Of Service:UHealth Labs	86900 86901
UHEALTH Labs	ACETA	ACETAMINOPHEN LEVEL Container:GOLD TOP/ RED TOP Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	
UHEALTH Labs	ACET	ACETONE Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	ACGAAB	ACETYLCHOLINE RECEPTOR GANGLIONIC Container:Red Top Temperature:Ship Frozen Place Of Service:Quest	83519
UHEALTH Labs	ACHRS	ACHR BINDING ABS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	ACBLAB	ACHR BLOCKING ABS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	AFBBC	Acid fast Bacilli Blood Culture Container:BACTEC® Myco® F/Lytic aerobic blood bottle Temperature:Room temperature Place Of Service:UHEALTH LABS	
UHEALTH Labs	AMALTD	ACID MALTASE DEF.,GSD TYPE II Container:Lavender/Specimen Place Of Service:DUKE	
UHEALTH Labs	ACTH2S	ACTH Stimulation, 2 Specimens Container:Serum Separator Tube (SST?) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ACTHP	ACTH, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	ACTH	ACTH, PLASMA Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	ACTHPL	ACTH,PLASMA (3 SPECIMENS) Place Of Service:LABCORP	
UHEALTH Labs	SMOAB	ACTIN (SMOOTH MUSCLE) AB (IGG) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MUSCLE	ACTIN (SMOOTH MUSCLE) ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	ACTMC	ACTINOMYCES CULTURE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	APCRES	ACTIVATED PROTEIN C RESISTANCE Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	AHEPP	ACUTE HEPATITIS PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80074
UHEALTH Labs	ACYLPL	ACYLCARNITINE PROFILE, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	ADAM	ADAMTS-13 ACTIVITY Container:BLUE Temperature:FROZEN Place Of Service:Uhealth Labs	
UHEALTH Labs	ADAS13	ADAMTS13 ACTIVITY REFLEX PANEL Place Of Service:LABCORP	
UHEALTH Labs	ADAMR	ADAMTS13 Activity with Reflex to Inhibitor Container:3.2% sodium citrate (light blue-top) tube    Transport Container: Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	ADACSF	ADENOSINE DEAMINASE, CSF Container:Sterile leak-proof container Temperature:FROZEN-PREFERRED    ROOM-UNACCEPTABLE    REFRIGERATED-UNACCEPTABLE Place Of Service:QUEST	
UHEALTH Labs	ADPERI	ADENOSINE DEAMINASE, PERITONEAL Place Of Service:LABCORP	
UHEALTH Labs	ADPLF	ADENOSINE DEAMINASE, PLEURAL Place Of Service:LABCORP	
UHEALTH Labs	ADEVQ	ADENOVIRUS DNA PCR (PLASMA) Place Of Service:VIRACOR	
UHEALTH Labs	ADEVQU	Adenovirus Quantitative Real-time PCR, Urine Container:Collect in a sterile urinalysis container then transfer to sterile, screw top tube for shipment. Temperature:Room Temperature or frozen temperature Monday through Friday. Specimens shipped at ambient temperature must be received within 96 hrs. of collection. Place Of Service:Eurofins Viracor	
UHEALTH Labs	ADEPCR	ADENOVIRUS, PCR QUANT Place Of Service:LABCORP	
UHEALTH Labs	ADHLC	ADH Place Of Service:LABCORP	
UHEALTH Labs	ADRHYD	Adrenal 21-Hydroxylase Autoantibodies Place Of Service:LABCORP	
UHEALTH Labs	RTCTG	AER/ANA TISSUE CULT WITH GS Container:STERILE CONTAINER Temperature:STORE AT RT, DO NOT REFRIGERATE. Place Of Service:UHealth Labs	
UHEALTH Labs	AFSO	AFB DIRECT SMEAR ONLY Temperature:REFRIGERATE IF SAMPLE IS TO BE SPLIT FOR OTHER TESTS. Place Of Service:UHealth Labs	
UHEALTH Labs	AFBID1	AFB IDENTIFICATION 1 Container:Double-walled container Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	AFBSO	AFB SMEAR ONLY Container:Sterile container with tight screw-cap seal or green-top (sodium heparin) tube or Isolator™ Temperature:Refrigerate Place Of Service:Uhealth	87206
UHEALTH Labs	AFBSMT	AFB SUCEPTIBILITY M. TURBERCULOSIS COMPLEX PRIMARY PANEL Container:Culture Temperature:Refrigerated Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	AFBPAN	AFB SUSCEPT, DRUG PANEL (1) Container:Double-walled container Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	MICRG	AFB,RAPID GROWER MIC Container:Double-walled safety container Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	AFPL3P	AFP WITH AFP-L3% Place Of Service:LABCORP	
UHEALTH Labs	DL24AL	ALA DELTA, 24-HOUR URINE Place Of Service:LABCORP	
UHEALTH Labs	DLALA	ALA DELTA, RANDOM URINE Place Of Service:LABCORP	
UHEALTH Labs	ALB	ALBUMIN LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82040
UHEALTH Labs	ALB24U	ALBUMIN, 24-HR URINE Place Of Service:LABCORP	
UHEALTH Labs	ALCEB	Alcohol, Ethyl, Blood Container:Fluoride oxalate (gray-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	AEQNU	Alcohol, Ethyl, Quantitative, Random Urine Container:Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALC	ALCOHOL, SERUM ETHANOL LEVEL Container:GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80320
UHEALTH Labs	ALDOL	ALDOLASE Place Of Service:LABCORP	
UHEALTH Labs	ALDST	ALDOSTERONE LCMS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	ALDUR	ALDOSTERONE, URINE Place Of Service:LABCORP	
UHEALTH Labs	ALDRAT	ALDOSTERONE/PLASMA RENIN ACTIVITY RATIO,LC/MS/MS Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALKPHO	ALK PHOS ISOENZYME Place Of Service:LABCORP	
UHEALTH Labs	ALKPBS	ALK PHOSPHATASE, BONE SPECIFIC Place Of Service:LABCORP	
UHEALTH Labs	ALKP	ALKALINE PHOSPHATASE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84075
UHEALTH Labs	ALLPR	ALL FISH PROFILE Container:sodium heparin (green-top) tube or RPMI for fluids and tissues Temperature:Use refrigerated cold pack for transport. Place Of Service:UHEALTH LABS	88368 88369 88377 x5
UHEALTH Labs	ALLFP	ALL FUSION PROFILE Container:sodium heparin (green-top) tube OR for Fixed cell pellets (3:1 methanol and acetic acid fixative) Temperature:Use refrigerated cold pack for transport Place Of Service:GENOPTIX	81450

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	ALMIG	ALLERGEN ALMOND IGG Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BEETS	ALLERGEN IGE BEET ROOT Temperature:TRANSPORT @ ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PEANG	ALLERGEN PEANUT IGG Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALLFOD	ALLERGEN PROFILE, BASIC FOOD Place Of Service:LABCORP	
UHEALTH Labs	DOGIGE	Allergen Profile, Dog, IgE With Component Reflex Place Of Service:LABCORP	
UHEALTH Labs	ALPRFF	Allergen Profile, Food-Fish Place Of Service:LabCorp	
UHEALTH Labs	SHELLF	Allergen Profile, Food-Shellfish IgE Place Of Service:LabCorp	
UHEALTH Labs	NUTIGE	Allergen Profile, Nut, IgE with Component Reflexes Place Of Service:LabCorp	
UHEALTH Labs	EELAL	ALLERGEN SPECIFIC IGE EEL * Container:GoldTop Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALLPR4	Allergen w/Comp Rflx Area 4 Place Of Service:LABCORP	
UHEALTH Labs	ALLF15	ALLERGENS (15) FOOD Place Of Service:LABCORP	
UHEALTH Labs	ALLR19	ALLERGENS W/TOTAL IGE AREA 19 Place Of Service:LABCORP	
UHEALTH Labs	ALLR2	ALLERGENS W/TOTAL IGE AREA 2 Place Of Service:LABCORP	
UHEALTH Labs	ALRGWT	ALLERGENS W/TOTAL IGE AREA 3 Place Of Service:LABCORP	
UHEALTH Labs	ALLER4	ALLERGENS W/TOTAL IGE AREA 4 Place Of Service:LABCORP	
UHEALTH Labs	AREA5	ALLERGENS W/TOTAL IGE AREA 5 Place Of Service:LABCORP	
UHEALTH Labs	ALLPER	ALLERGENS, PERENNIAL Place Of Service:LABCORP	
UHEALTH Labs	AALLP	ALLERGY PANEL 12, ANIMAL GROUP Container:Serum Separator Tube (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALL4E	ALLERGY PANEL, ANIMAL GROUP Place Of Service:LABCORP	
UHEALTH Labs	ALCERL	ALLERGY PANEL, CEREAL GROUP Place Of Service:LABCORP	
UHEALTH Labs	ALMOLD	ALLERGY PANEL, MOLD GROUP Place Of Service:LABCORP	
UHEALTH Labs	ALNUT	ALLERGY PANEL, NUT MIX GROUP Place Of Service:LABCORP	
UHEALTH Labs	ALSALD	ALLERGY PANEL, SALAD GROUP Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	ALSEAF	ALLERGY PANEL, SEAFOOD GROUP Place Of Service:LABCORP	
UHEALTH Labs	ALVEGE	ALLERGY PANEL, VEGETABLE GROUP Place Of Service:LABCORP	
UHEALTH Labs	AGHBA	Alpha Globin (HBA1&HBA2) Seq. & Del/Dup Container:Lavender (EDTA), pink (K2EDTA), or Yellow (ACD Solution A or B) Temperature:Refrigerated Place Of Service:ARUP	81259 81269
UHEALTH Labs	ALSBU	ALPHA SUBUNIT (FREE) Place Of Service:LABCORP	
UHEALTH Labs	A1ADEF	ALPHA-1-ANTITRYPSIN DEFICIENCY Place Of Service:LABCORP	
UHEALTH Labs	ALPPHE	ALPHA-1-ANTITRYPSIN PHENOTYP Place Of Service:LABCORP	
UHEALTH Labs	ALP1A	ALPHA-1-ANTITRYPSIN, SERUM Place Of Service:LABCORP	
UHEALTH Labs	AATFEC	ALPHA-1-ANTITRYPSIN,FECES Container:PLASTIC SCREW-CAP CONTAINER Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	ALP2M	ALPHA-2 MACROGLOBULINS, QN Place Of Service:LABCORP	
UHEALTH Labs	A2PLSM	ALPHA-2-ANTIPLASMIN Place Of Service:LABCORP	
UHEALTH Labs	AFP	ALPHA-FETOPROTEIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84460
UHEALTH Labs	AFPPF	Alpha-Fetoprotein, Pleural Fluid Container:sterile plastic screw-cap container Temperature:Frozen Place Of Service:QUEST	86316
UHEALTH Labs	ALPGAL	Alpha-Gal IgE Panel Place Of Service:LABCORP	
UHEALTH Labs	AGMAN	Alpha-Globin Common Mutation Analysis Container:EDTA (lavender-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALPRAZ	ALPRAZOLAM (XANAX), SERUM Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALT	ALT (SGPT) ALANINE AMINOTRANSFERASE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82107
UHEALTH Labs	AALTER	ALTERNARIA ALTERNATA (M6) IGE Container:Serum Separator Tube (SST®) Place Of Service:QUEST	
UHEALTH Labs	ALTPAC	ALTERNATIVE PATHWAY ACTIVITY Place Of Service:LABCORP	
UHEALTH Labs	AL24U	ALUMINUM, 24 HR URINE Container:ACID-WASHED PLASTIC URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALUMPB	ALUMINUM, PLASMA/SERUM Place Of Service:LABCORP	
UHEALTH Labs	AMEBAB	AMEBIASIS ANTIBODIES Place Of Service:LABCORP	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	AMIKQ	AMIKACIN Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	AMIPK	AMIKACIN PEAK, SERUM Place Of Service:LABCORP	
UHEALTH Labs	AMIRAN	AMIKACIN RANDOM, SERUM Place Of Service:LABCORP	
UHEALTH Labs	AMITR	AMIKACIN TROUGH, SERUM Place Of Service:LABCORP	
UHEALTH Labs	AMIKPQ	AMIKACIN, PEAK Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	AMIKTQ	AMIKACIN, TROUGH Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	AAPCSF	Amino Acid Profile, Qn, CSF Place Of Service:LabCorp	
UHEALTH Labs	AAPRQN	AMINO ACID PROFILE, QN, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	AAPRQU	AMINO ACID PROFILE, QN, URINE Place Of Service:LABCORP	
UHEALTH Labs	AMIOD	AMIODARONE (CORDARONE) Place Of Service:LABCORP	
UHEALTH Labs	AMITPS	AMITRIPTYLINE (ELAVIL), SERUM Place Of Service:LABCORP	
UHEALTH Labs	AMLFSH	AML FISH PANEL Place Of Service:LABCORP	
UHEALTH Labs	AMLNGS	AML MOLECULAR PROFILE,NGS Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	AMLPCR	AML MOLECULAR PROFILE,PCR Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	AMLPP	AML PROGNOSTIC PANEL, NORM KARYO Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	AMMON	AMMONIA Container:LAVENDER TOP TUBE (EDTA) Temperature:FROZEN Place Of Service:UHealth Labs	82140
UHEALTH Labs	AMYL	AMYLASE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82150
UHEALTH Labs	AMYISO	AMYLASE ISOENZYMES Place Of Service:LABCORP	
UHEALTH Labs	ANAIF	ANA IFA RFLX TITER SCLERODERMA AB Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ANAAQ	ANA IFA SCREEN W/REFL TITR/PATT Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ANAREF	ANA IFA, W/REFL TO TITER/PATTERN/CASCADE Container:Transport tube Temperature:Room temperature Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	ANARFX	ANA,MPLX IA,RFX 5-BIOMARKER PROFILE (ENA) Place Of Service:LABCORP	
UHEALTH Labs	AAALKM	ANA+AMA+ASMA+LKM AB Container:Red-top tube or gel-barrier tube Temperature:Room Temperature or Refrigerated Place Of Service:LABCORP	83516 x2 86038 001 86376 001
UHEALTH Labs	AAASJO	ANA+AMA+ASMA+LKM Ab+SjoSSA Container:Red-top tube or gel-barrier tube Temperature:Room Temperature Place Of Service:LABCORP	
UHEALTH Labs	ANAB1T	ANA+ANCA+MPO+PR 3 Container:Red-top tube or gel-barrier tube Temperature:Room Temperature or Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	AEDDS	ANA+ENA+DNA/DS+SCL70+SJOSSA/B Container:Red-top tube or gel-barrier tube Temperature:Room Temperature or Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	ANARAQ	ANA+RA QN+DNA/DS+A ENA+SCL Container:Red-top tube or gel-barrier tube Temperature:Room Temperature or Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	ANCAB	ANCA Container:GoldTop, Red Temperature:Refrigerated Place Of Service:JMH	
UHEALTH Labs	ANCAC	ANCA CONFIRMATORY (MPO/PR3) Container:GoldTop, Red1mL Temperature:Refrigerated Place Of Service:JMH	
UHEALTH Labs	ANCAP	ANCA SCREEN WITH MPO AND PR3 WITH REFLEX TO ANCA TITER Container:Plastic screw-cap vial Place Of Service:QUEST	
UHEALTH Labs	ANCARE	ANCA SCREEN WITH REFLEX TO TITER Container:Transport tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	ANCAQ	ANCA VASCULITIDES Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ANDRST	ANDROSTENEDIONE LCMS Place Of Service:LABCORP	
UHEALTH Labs	ANGDMA	ANGIOEDEMA PANEL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ACECSF	ANGIOTENSIN CONVERT ENZ.,CSF Place Of Service:LABCORP	
UHEALTH Labs	ACEPL	ANGIOTENSIN CONVERTING ENZYME(ACE)POLYMORPHISM Container:EDTA (lavender-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ACEB	ANGIOTENSIN-CONVERTING ENZYME Place Of Service:LABCORP	
UHEALTH Labs	AB2GPI	Anti B2 Glycoprotein I (IgG, IgM) Container:RED TOP TUBE (NO SERUM SEPARATOR) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	86146 x2
UHEALTH Labs	GMCARD	ANTI CARDIOLIPIN Ab (IgG,IgM) Container:RED TOP TUBE (NO SERUM SEPARATOR) Temperature:TRANSPORT AT ROOM TEMPERATURE	86147 x3



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CARDM	Place Of Service:UHealth Labs ANTI CARDIOLIPIN M Container:RED TOP TUBE (NO SERUM SEPARATOR) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86147 x3
UHEALTH Labs	AMERK	Anti Merkel Cell Panel Container:RED TOP or SST tube Temperature:Refrigerate Place Of Service:UHEALTH LABS	0058U
UHEALTH Labs	ANECSE	ANTI NEURONAL CELL AB, CSF Place Of Service:LabCorp	
UHEALTH Labs	ANNEX	Anti-Annexin V IgG, IgM Place Of Service:LabCorp	
UHEALTH Labs	ABIDTT	ANTIBODY ID TITER AND TYPING, RBC Container:N ACD-A OR ACD-B (YELLOW-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RBCREF	ANTIBODY SCREEN, RBC W/REFL ID, TITER AND AG Container:ACD-A or ACD-B (yellow-top) tube Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	ACARD	Anti-Cardiolipin Antibody Panel Container:RED TOP TUBE (NO SERUM SEPARATOR) Temperature:Transport Refrigerated Place Of Service:Uhealth Labs	86147 x3
UHEALTH Labs	CENTRO	ANTI-CENTROMERE B ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	ACHAB	ANTICHROMATIN ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	SADPNL	ANTIDEPRESSANT PANEL, QL, EXP. Place Of Service:LABCORP	
UHEALTH Labs	ADHP	Antidiuretic Hormone (ADH) Profile Place Of Service:LABCORP	
UHEALTH Labs	ADSDAB	ANTI-DSDNA ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	EPIDAB	ANTIEPIDERMAL ABS W/REFLEX TITER Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	AXNUC	ANTIEXTRACTABLE NUCLEAR AG Place Of Service:LABCORP	
UHEALTH Labs	YANST	ANTIGEN STIMULATION Temperature:SAMPLES SHOULD BE SENT AT ROOM TEMPERATURE, UNLESS INFORMED OTHERWISE. PACKAGE SECURELY TO AVOID BREAKAGE AND EXTREME WEATHER CONDITIONS. Place Of Service:CINCIN	
UHEALTH Labs	ANABIA	ANTIGLIADIN ABS, IGA Place Of Service:LABCORP	
UHEALTH Labs	GLIABS	ANTIGLIADIN ABS, IGG Place Of Service:LABCORP	
UHEALTH Labs	GLOMAB	ANTIGLOMERULAR BM AB Place Of Service:LABCORP	
UHEALTH Labs	AHISAB	ANTI HISTONE ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	ANTIGE	ANTI-IGE Place Of Service:LABCORP	
UHEALTH Labs	ANJO1B	ANTI-JO-1 Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	NTMACT	ANTIMICROBIAL SUSCEP (NTM/AERO ACTINOMYCETES) Place Of Service:NJH	
UHEALTH Labs	AMHB	ANTI-MULLERIAN HORMONE (AMH) Place Of Service:LABCORP	
UHEALTH Labs	MPOAB	Antimyeloperoxidase (MPO) Antibodies Place Of Service:LABCORP	
UHEALTH Labs	ANTCYT	ANTINEUTROPHIL CYTOPLASMIC AB Place Of Service:LABCORP	
UHEALTH Labs	ANCAPR	Antineutrophil Cytoplasmic Antibody (ANCA) Profile Place Of Service:LABCORP	
UHEALTH Labs	ANAAB	ANTINUCLEAR AB TITER AND PATTERN Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ANAD	ANTINUCLEAR ANTIBODIES DIRECT Place Of Service:LABCORP	
UHEALTH Labs	AAIFA	ANTINUCLEAR ANTIBODIES, IFA Place Of Service:LABCORP	
UHEALTH Labs	PANISL	ANTIPANCREATIC ISLET CELLS Place Of Service:LABCORP	
UHEALTH Labs	ANTCLA	ANTIPARIETAL CELL ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	PHOSIG	ANTIPHOSPHATIDYLSERINE IGG/M/A Place Of Service:LABCORP	
UHEALTH Labs	APSP	ANTI-PHOSPHOLIPID SYNDROME Place Of Service:LABCORP	
UHEALTH Labs	ANTP	Antiphospholipid Syndrome Comp Place Of Service:LabCorp	
UHEALTH Labs	ANTIPL	ANTIPLASMIN ACTIVITY Container:Citrated Light Blue Top Tube Temperature:FROZEN Place Of Service:Uhealth Labs	85410
UHEALTH Labs	PF4HEP	Anti-Platelet Factor 4 Temperature:FROZEN Place Of Service:Uhealth Labs	86022
UHEALTH Labs	PR3ABS	ANTIPROTEINASE 3 (PR-3) ABS Place Of Service:LABCORP	
UHEALTH Labs	PHIGIM	Anti-PS/PT Abs IgG, IgM Place Of Service:LABCORP	
UHEALTH Labs	ANTPAB	ANTIRIBOSOMAL P ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	ANRNA	Anti-RNA Polymerase III (RDL) Place Of Service:LabCorp	
UHEALTH Labs	AN70AB	ANTISCLERODERMA-70 ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	SMA	ANTI-SMOOTH MUSCLE AB Container:Transport tube Temperature:Frozen Place Of Service:UMHEP	
UHEALTH Labs	ASTREP	ANTISTREPTOLYSIN O AB Place Of Service:LABCORP	
UHEALTH Labs	AT	ANTITHROMBIN ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85300
UHEALTH Labs	ATAG	ANTITHROMBIN ANTIGEN Container:BLUE TOP TUBE (SODIUM CITRATE)	85301

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	
UHEALTH Labs	ANTITH	ANTITHROMBIN III ACTIVITY W/REFL Container:3.2% sodium citrate (light blue-top) tube Temperature:Frozen Place Of Service:QUEST	
UHEALTH Labs	ATPAN	ANTITHROMBIN III PANEL Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85301
UHEALTH Labs	ATSHAB	ANTI-TSH AUTOANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	AXALMW	ANTI-Xa LOW MOLEC.WEIGHT HEPARIN Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	AXASTD	ANTI-Xa STANDARD HEPARIN Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	APCSDD	APC Sequencing and Deletion/Duplication Container:Lav Top-EDTA or ACD tube (yellow top) Temperature:Room Temp Place Of Service:QUEST	
UHEALTH Labs	APIXA	Apixaban Place Of Service:LABCORP	
UHEALTH Labs	APLQNT	APL(PML-RARA) QUANT Place Of Service:LABCORP	
UHEALTH Labs	ALIPO	APOLIPOPROTEIN A1 Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	APOL	APOLIPOPROTEIN B Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	APOA1B	APOLIPOPROTEIN EVALUATION Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PTT	APTT Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:ROOM TEMPERATURE Place Of Service:UHealth Labs	85730
UHEALTH Labs	AQUAP	Aquaporin-4 Antibody (IgG), CBA Container:Red-top tube (no gel) Temperature:Room temperature Place Of Service:Quest	86255
UHEALTH Labs	ARBOGM	ARBO VIRUS IGG/IGM PANEL, SERUM Place Of Service:LABCORP	
UHEALTH Labs	ARBOM	Arbo Virus IgM Panel Serum Place Of Service:LabCorp	
UHEALTH Labs	UARS	ARSENIC, 24 HOUR URINE Container:ACID WASHED PLASTIC URINE CONTAINER Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	ARSBLD	ARSENIC, BLOOD Place Of Service:LABCORP	
UHEALTH Labs	ARSH	ARSENIC, HAIR Container:PLASTIC CONTAINER. Temperature:ROOM TEMPERATURE	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:QUEST	
UHEALTH Labs	ARSEN	ARSENIC, NAILS Container:PLASTIC CONTAINER (ACID WASHED OR TRACE METAL-FREE) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ARYLSA	Arylsulfatase A Deficiency, Leukocytes Place Of Service:LABCORP	
UHEALTH Labs	ASPRQN	ASPERGILLUS AB, QN, DID Place Of Service:LABCORP	
UHEALTH Labs	ASPBAL	ASPERGILLUS AG, BAL/SERUM Place Of Service:LABCORP	
UHEALTH Labs	YASPGA	ASPERGILLUS GALACTOMANNAN (SERUM) Container:serum separator tube (SST) Place Of Service:VIRACOR	
UHEALTH Labs	ASPMIR	Aspergillus Galactomannan EIA (MiraVista) Container:Serum: Collect serum specimens in serum separator or red top tube. Temperature:Refrigerated/ Frozen Place Of Service:UHEALTH LABS	87305
UHEALTH Labs	ASPGBW	Aspergillus Galactomannan EIA, BW Container:Sterile, screw top tube Temperature:Specimen should be stored at 2 to 8C or frozen in a non-self-defrosting freezer. Place Of Service:Viracor	87305
UHEALTH Labs	ASPGGA	ASPERGILLUS GALACTOMANNAN(BAL) Container:sterile, screw-cap tube Temperature:SPECIMEN SHOULD BE STORED AT 2 TO 8C OR FROZEN IN A NON-SELF-DEFROSTING FREEZER AND SHIPPED WITH FROZEN GEL PACKS OR DRY ICE FOR OVERNIGHT DELIVERY AT VIRACOR. Place Of Service:VIRACOR	
UHEALTH Labs	ASPGPN	ASPERGILLUS RT PCR PNL (BAL) Container:Sterile, leak-proof container Temperature:Frozen Place Of Service:VIRACOR	
UHEALTH Labs	AST	AST (SGOT) ASPARTATE AMINOTRANSFERASE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84450
UHEALTH Labs	APNUM	ATYPICAL PNEUMONIA PANEL(BAL) Container:Sterile, leak-proof container Temperature:Frozen Place Of Service:VIRACOR	
UHEALTH Labs	ASTPAL	AUSTRALIAN PINE (T73) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ACC	AUTOMATED CELL COUNT Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85025
UHEALTH Labs	BAFF	B Cell Activating Factor Level (BAFF) Container:Lavender Temperature:Room Temperature Place Of Service:Cincinnati	83520
UHEALTH Labs	YBCCELL	B CELL PANEL Container:EDTA (lavender-top) tube Place Of Service:CINCIN	
UHEALTH Labs	BSUBS	B Cell Subset Analysis Place Of Service:LabCorp	
UHEALTH Labs	BTNP	B TYPE NATRIURETIC PEPTIDE (BNP) Container:EDTA (lavender-top) tube Temperature:FROZEN Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	BABEAB	BABESIA MICROTI ANTIBODY PANEL Place Of Service:LABCORP	
UHEALTH Labs	BABMIC	Babesia microti, RT DNA PCR Place Of Service:LABCORP	
UHEALTH Labs	BMAG	BACTERIAL MENINGITIS AG PANEL Container:STERILE PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	BAHIA	Bahia Grass IgE (G017) Place Of Service:LabCorp	
UHEALTH Labs	BALLPR	B-ALL MRD PROFILE Container:sodium heparin (green-top) tube Temperature:refrigerated cold pack Place Of Service:UHEALTH LABS	88184 88185 x12 88188
UHEALTH Labs	BARBUR	BARBITURATE (GC/MS), URINE Place Of Service:LABCORP	
UHEALTH Labs	BARLAL	BARLEY (F6) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BARTAB	BARTONELLA ANTIBODY PANEL Place Of Service:LABCORP	
UHEALTH Labs	BMP	BASIC METABOLIC PROFILE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80048
UHEALTH Labs	BAFNA	BasoFunction HRT Naproxen Container:Green Top Sodium Heparin Tube Temperature:Ambient Place Of Service:Viracor	
UHEALTH Labs	BAYBER	Bayberry/Sweet Gale Place Of Service:LABCORP	
UHEALTH Labs	BCLPCR	BCL 1/JH T(11,14) RTPCR LEUMETA R Container:EDTA (LAVENDER-TOP) TUB Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	BCRMJ	BCR ABL MAJOR BY RT PCR Temperature:REFRIGERATED Place Of Service:ARUP	
UHEALTH Labs	GREA	BCR ABL1 GENE REARRANGEMENT,QNT PCR Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BCRKIN	BCR-ABL1 KINASE DOMAIN Place Of Service:LABCORP	
UHEALTH Labs	BCABLB	BCR-ABL1, CML/ALL, PCR, QUANT Place Of Service:LABCORP	
UHEALTH Labs	BCRRFX	BCR-ALB1 QL WITH RFX TO BCR-ALB1QN Temperature:REFRIDGERATED Place Of Service:ARUP	
UHEALTH Labs	BRCAR	BCRAVANTAGE,REARRANGEMENTS Container:EDTA (lavender-top) tube Place Of Service:QUEST	
UHEALTH Labs	BEEFS	BEEF Place Of Service:LABCORP	
UHEALTH Labs	BNZOUR	BENZODIAZEPINES CONFIRM, URINE Place Of Service:LABCORP	
UHEALTH Labs	BENSC	Benzodiazepines Screen with Confirmation, Urine Container:Plastic urine container Temperature:ROOM TEMPERATURE	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:QUEST	
UHEALTH Labs	BENQU	Benzodiazepines, Quantitative, Urine Container:Plastic urine container Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	BENZT	BENZTROPINE (COGENTIN) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BERMU	BERMUDA GRASS Place Of Service:LABCORP	
UHEALTH Labs	UR2BG	BETA 2 MICROGLOBULIN, RANDM UR Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	B2MIC	BETA 2 MICROGLOBULIN, SERUM Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82232
UHEALTH Labs	BAMYCF	Beta Amyloid 42/40 Ratio, CSF Place Of Service:LABCORP	
UHEALTH Labs	BTHBB	BETA THALASSEMIA: HBB Place Of Service:LABCORP	
UHEALTH Labs	BT2IGA	BETA-2 GLYCOPROTEIN I AB, IGA Place Of Service:LABCORP	
UHEALTH Labs	BT2IGG	BETA-2 GLYCOPROTEIN I AB, IGG Place Of Service:LABCORP	
UHEALTH Labs	BT2IGM	BETA-2 GLYCOPROTEIN I AB, IGM Place Of Service:LABCORP	
UHEALTH Labs	B2MS	BETA-2 MICROGLOBULIN, SERUM Place Of Service:LABCORP	
UHEALTH Labs	B2TRNS	BETA-2 TRANSFERRIN Place Of Service:LABCORP	
UHEALTH Labs	BTA2CD	BETA-2(IGG/A/M)+CARD(IGG/A/M) Place Of Service:LABCORP	
UHEALTH Labs	BGLC	Beta-Globin Complete Container:Preferred: EDTA (lavender-top) or ACD solution B (yellow-top) tube Alternate : EDTA (royal-blue) tube, sodium heparin (green-top) tube or ACD solution A (yellow-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BHCG	BETA-HCG, BLOOD Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84702
UHEALTH Labs	BTAHYB	BETA-HYDROXYBUTYRATE Place Of Service:LABCORP	
UHEALTH Labs	BAFRTP	BILE ACIDS, FRACT/TOTAL, PREG Container:SERUM SEPARATOR TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BILAS	BILE ACIDS, FRACTIONATED & TOTAL Container:No additive (red-top) No additive serum separator tube Temperature:REFRIGERATED PREFERRED ROOM TEMPERATURE ACCEPTABLE FROZEN ACCEPTABLE Place Of Service:QUEST	
UHEALTH Labs	BLCMS	BILE ACIDS, FRACTIONATED LCMS Place Of Service:LABCORP	
UHEALTH Labs	BILET	BILE ACIDS, TOTAL Container:PLASTIC SCREW-CAP VIAL	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	BILFR	BILIRUBIN, FRACTIONATED Container:Foil-wrapped spun SST® • Amber plastic screw-cap vial • Foil-wrapped plastic screw-cap vial Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82247
UHEALTH Labs	BKQNT	BK QUANT PCR (PLASMA/SERUM) Place Of Service:LABCORP	
UHEALTH Labs	BKQPCU	BK QUANT PCR (URINE) Place Of Service:LABCORP	
UHEALTH Labs	BKPCR	BK VIRUS BKV qPCR(PLASMA) Container:L-Lavender EDTA Temperature:Room Temperature Place Of Service:Viracor	
UHEALTH Labs	BKPCRU	BK VIRUS BKV qPCR(URINE) Container:sterile urinalysis container then transfer to sterile, screw top tube for shipment Temperature:ambient or frozen temperature Place Of Service:VIRACOR	
UHEALTH Labs	BKVQT	BK VIRUS DNA, QN RT PCR, PLASMA Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	BKDNAU	BK Virus DNA, Quantitative, Real-Time PCR, Urine Container:Plastic urine container Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	BKVIRG	BK VIRUS IGG ANTIBODY Container:EDTA (lavender top) tube Temperature:ROOM TEMPERATURE Place Of Service:VIRACOR	
UHEALTH Labs	BLBEAN	Black Bean IgE Container:Collect 1-2 mL whole blood in red top tube Temperature:Ship at ambient or frozen temperature Monday through Friday. Place Of Service:LABCORP	
UHEALTH Labs	BLADCA	BLADDER CANCER, FISH Container:Plastic urine container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	BLSTAB	BLASTOMYCES ABS, QN, DID Place Of Service:LABCORP	
UHEALTH Labs	BLASAB	Blastomyces Antibody by I.D. Container:Serum separated from the clot in a leak-proof container. Temperature:Frozen Place Of Service:MiraVista	86612
UHEALTH Labs	BQAEIA	Blastomyces Quantitative Antigen Container:Plastic urine container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	BCGPP	BLD CULT GRAM POS PANEL, NUC ACID PROBE Container:BD BACTEC™ Plus™ Aerobic/F (blue label and cap) and Lytic/10 Anaerobic/F bottles (purple label and cap), Temperature:Room temperature Place Of Service:UHEALTH LABS	
UHEALTH Labs	BDIAL	Bleeding Diathesis Profile, Limited Container:Light-blue top (citrate) Temperature:Frozen Place Of Service:UHEALTH LABS	85240 85246 85291 85366 85379 85384 85390 x26 85397 85610 85670 85730
UHEALTH Labs	BLMTRP	BLOMIA TROPICALIS (D201) IGE Container:GoldTop Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	BC	BLOOD CULTURE Container:AEROBIC BOTTLE (BLUE TOP) ANAEROBIC BOTTLE (PURPLE TOP) PEDIATRIC BOTTLE (PINK TOP) Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	BLDPAR	BLOOD PARASITES Container:EDTA TUBE (LAVENDER TOP TUBE) Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87207
UHEALTH Labs	BUN	BLOOD UREA NITROGEN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84520
UHEALTH Labs	BDYF	BODY FLUID ANALYSIS Container:RED TOP TUBE (NO SERUM SEPARATOR) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85025
UHEALTH Labs	FLDCUL	BODY FLUID CULTURE (EXCEPT BLOOD OR CSF) Temperature:STORE AT RT DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	BMCULT	BONE MARROW CULTURE Container:MYCO/F LYTIC (RED) AND PEDS BOTTLE (PINK) Temperature:STORE AT RT DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	LBR	BONE MARROW EVALUATION Place Of Service:UHEALTH LABS	
UHEALTH Labs	BPIGGM	Bordetella pertussis Antibodies, IgG, IgM Place Of Service:LABCORP	
UHEALTH Labs	BPIGM	Bordetella pertussis Antibodies, IgM Place Of Service:LabCorp	
UHEALTH Labs	BPRAG	BORDETELLA PERTUSSIS IGG/IGA Container:Red-top tube or gel-barrier tube Temperature:Refrigerate Place Of Service:LABCORP	86615 x2
UHEALTH Labs	BPPCRB	Bordetella pertussis Plus Real-time PCR Panel, BAL Container:Collect in a sterile, screw top tube. Temperature:Frozen (dry ice) Place Of Service:UHEALTH LABS	87798
UHEALTH Labs	BPPPCR	Bordetella pertussis/parapertussis, RT DNA PCR Place Of Service:LabCorp	
UHEALTH Labs	BRAFMT	BRAF Mutation Analysis, Papillary Thyroid Cancer Container:Fixative, slides or paraffin embedded tissue block Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BRAZAL	BRAZIL NUT (F18) IGE Container:GoldTop Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BRCASS	BRCASSURE COMPREHENSIVE TEST Place Of Service:LABCORP	
UHEALTH Labs	BRUCAB	BRUCELLA ANTIBODY IGG/IGM Place Of Service:LABCORP	
UHEALTH Labs	NATRI	B-TYPE NATRIURETIC PEPTIDE Place Of Service:LABCORP	
UHEALTH Labs	BUKWAL	BUCKWHEAT (F11) IGE Container:GoldTop Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BP230	Bullous Pemphigoid BP230 IgG Container:Transport tube	83520

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	BUNCR	BUN/CREATININE RATIO WITH eGFR Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82565
UHEALTH Labs	BUPRB	BUPROPION (WELLBUTRIN) Place Of Service:LABCORP	
UHEALTH Labs	BUSP	BUSPIRONE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	BUSUL	BUSULFAN LEVEL Container:Sodium Heparin green top tubes Temperature:SAMPLES MUST BE SHIPPED ON A MINIMUM OF 3 KG OF DRY ICE USING AN OVERNIGHT CARRIER. SAMPLES THAT ARRIVED THAWED WILL NOT BE ANALYZED. Place Of Service:SCCA	
UHEALTH Labs	BACVAG	BV+Ct/Ng/Tv NAA+Yeast Culture Place Of Service:LabCorp	
UHEALTH Labs	CDFNAA	C DIFFICILE TOXIN GENE NAA Place Of Service:LABCORP	
UHEALTH Labs	CDCYB	C DIFFICILE, CYTOTOXIN B Place Of Service:LABCORP	
UHEALTH Labs	CDCAN	C. DIFFICILE CYTOTOXIN ANTIBODY NEUTRALIZATION Container:Sterile, plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	C1ESTB	C1 ESTERASE INHIBITOR, FNC Place Of Service:LABCORP	
UHEALTH Labs	C1EST	C1 ESTERASE INHIBITOR, SERUM Place Of Service:LABCORP	
UHEALTH Labs	C1QF	C1q Function (C1QF) Container:Plain red top Temperature:frozen Place Of Service:NATIONAL JEWISH	86161
UHEALTH Labs	C3C4DN	C3, C4 and Anti-dsDNA (Double-stranded) Antibodies Container:Red-top tube or gel-barrier tube Temperature:Room Temperature or Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	C3C4CT	C3+C4+COMPT Place Of Service:LABCORP	
UHEALTH Labs	C4BND	C4 BINDING PROTEIN Place Of Service:LABCORP	
UHEALTH Labs	CFBPP	C4 BINDING PROTEIN, PLASMA Container:3.2% SODIUM CITRATE (LIGHBLUE-TO Temperature:FROZEN ON DRY ICE Place Of Service:QUEST	
UHEALTH Labs	C4DF	C4D FRAGMENT, EIA Container:EDTA (LAVENDER-TOP) TUB Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CA125	CA 125 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86304
UHEALTH Labs	CA153	CA 15-3 Container:Serum Separator Tube (SST®) Temperature:Room temperature Place Of Service:UHEALTH LABS	86300

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CA199	CA 19-9 Container:Red-top tube or gel-barrier tube Temperature:Refrigerate Place Of Service:UHealth	86301
UHEALTH Labs	YFCA19	CA 19-9, FLUID Container:Cup, Red plain tube Temperature:Refrigerate Place Of Service:ARUP	
UHEALTH Labs	CA2729	CA 27.29 Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	C2729B	CA 27.29 Place Of Service:LABCORP	
UHEALTH Labs	CAFFSR	CAFFEINE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	CALCIT	CALCITONIN, SERUM Place Of Service:LABCORP	
UHEALTH Labs	VD125	CALCITRIOL(1,25 DI-OH VIT D) Place Of Service:LABCORP	
UHEALTH Labs	CA	CALCIUM LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84295
UHEALTH Labs	UCACR	CALCIUM, 24 HOUR URINE (W/ CREATININE) Container:JUG- PLASTIC URINE CONTAINER Temperature:TRANSPORT REFRIGERATED (ON ICE) Place Of Service:UHealth Labs	82570 84300
UHEALTH Labs	CALFEC	CALPROTECTIN, FECAL Place Of Service:LABCORP	
UHEALTH Labs	CALRMU	CALR MUTATION ANALYSIS Place Of Service:LABCORP	
UHEALTH Labs	CAMPYC	CAMPYLOBACTER, CULTURE Container:Cary-Blair stool culture transport medium Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CNDAAB	CANDIDA ANTIBODIES IGG,IGA,IGM Place Of Service:LABCORP	
UHEALTH Labs	CANABS	CANNABINOIDS Container:A RED-TOP TUBE (NO GEL) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CARBX	CARBAMAZEPINE + CARB EPOXIDE Place Of Service:LABCORP	
UHEALTH Labs	CABAM	CARBAMAZEPINE AND METABOLITE Container:A RED-TOP TUBE (NO GEL) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TEG	CARBAMAZEPINE LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERRED), GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80156
UHEALTH Labs	CARBAM	Carbamazepine-10, 11 Epoxide Place Of Service:LABCORP	
UHEALTH Labs	CEA	CARCINOEMBRYONIC Ag Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82378

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CARDP	CARDIAC PROFILE Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CAPOE	Cardio IQ APOE Genotype Container:Whole blood collected in an EDTA (lavender-top) tube Temperature:Refrigerated (cold packs) Place Of Service:Quest	81401
UHEALTH Labs	CAPOB	CARDIO IQ APOLIPOPROTEIN B Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CIQLIP	CARDIO IQ LIPID PANEL Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CLIPA	Cardio IQ Lipoprotein (a) Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CARDIQ	Cardio IQ ST2, Soluble Container:Goldtop; Lavender; Red Temperature:Refrigerated (cold packs) Place Of Service:QUEST	
UHEALTH Labs	CARTF	CARNITINE, TOTAL AND FREE Place Of Service:LABCORP	
UHEALTH Labs	CRTENE	CAROTENE, BETA Place Of Service:LABCORP	
UHEALTH Labs	CTVMAU	CATECHOLAMINE+VMA, 24 HR Place Of Service:LABCORP	
UHEALTH Labs	CATEPL	CATECHOLAMINES, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	CATFRU	CATECHOLAMINES, UR., FREE, RANDOM Place Of Service:LABCORP	
UHEALTH Labs	CTE24U	CATECHOLAMINES,UR, FREE,24 HR Place Of Service:LABCORP	
UHEALTH Labs	CARCFV	CATECHOLAMN FRCT VMA 24H WO CREAT Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CATHC	CATHETER AND INDWELLING DEVICES CULTURE Container:STERILE CUP Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	CBCSMR	CBC (W PLT AND DIFF) AND SMEAR REVIEW Container:LAVENDER TOP TUBE(EDTA) AND YELLOW TOP TUBE(ACD A) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85007
UHEALTH Labs	CBCPR	CBC (WITH PLT AND DIFF),INCLUDES PATH REVIEW Container:LAVENDER TOP TUBE(EDTA) AND YELLOW TOP TUBE(ACD A) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85025
UHEALTH Labs	CYCPEP	CCP AB (IGG) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CCPIG	CCP ANTIBODIES IGG/IGA Place Of Service:LABCORP	
UHEALTH Labs	CD107A	CD107a Mobilization (NK Cell Degranulation) Container:Sodium Heparin (Dark Green) Temperature:AMBIENT	88184 88185 88187 86352

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CD20B	Place Of Service:CINCIN CD20, B-Cells Container:EDTA (lavender top) Temperature:Room temperature Place Of Service:UHEALTH LAB	86355 86356
UHEALTH Labs	CD3C	CD3 COUNT Container:Lavender Top Tube for peripheral blood / Yellow Top tube for bone Marrow Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CD34FC	CD34 Container:Lavender Top Tube for peripheral blood / Yellow Top tube for bone Marrow Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	YCD45	CD45R/CD45RO Container:Lavendr/MIC-LAV Place Of Service:CINCIN	
UHEALTH Labs	CEBPAM	CEBPA MUTATION DETECTION Container:Lavender(EDTA) Or Bone marrow(EDTA) Temperature:REFRIDGERATED Place Of Service:ARUP	
UHEALTH Labs	CELERY	CELERY IgE ALLERGEN Container:Serum Separator Tube (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CLACMP	CELIAC DISEASE COMPREHENSIVE Place Of Service:LABCORP	
UHEALTH Labs	CELC	CELIAC DISEASE PANEL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CELIAC	Celiac HLA DQ Association Place Of Service:LabCorp	
UHEALTH Labs	CERUL	CERULOPLASMIN Place Of Service:LABCORP	
UHEALTH Labs	CHCHEE	Cheese, Cheddar Type (F081) IgE Place Of Service:LABCORP	
UHEALTH Labs	MOCHEE	Cheese, Mold Type (F082) IgE Place Of Service:LABCORP	
UHEALTH Labs	CHSTAL	CHESTNUT (RF299) IGE Container:Serum Separator Tube (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CHICKP	Chick Pea IgE (F309) Container:Red-top tube or gel-barrier tube Temperature:Room temperature Place Of Service:LABCORP	86003
UHEALTH Labs	CHFTAL	CHICKEN FEATHERS (E85) IGE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CKVPCR	CHIKUNGUNYA VIRUS RNA, QUAL, RT PCR Container:Serum Separator Tube (SST) Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CHIKV	Chikungunya Virus, PCR, Serum Place Of Service:LABCORP	
UHEALTH Labs	ALLCLD	CHILDHOOD ALLERGY PROFILE+IGE Place Of Service:LABCORP	
UHEALTH Labs	CHLPNE	CHL.PNEUMONIAE (IGG/IGM/IGA) Place Of Service:LABCORP	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CHLIGG	CHLAMYDIA ANTIBODIES, IGG Place Of Service:LABCORP	
UHEALTH Labs	CHLMP	CHLAMYDIA PSITTACI Ab PANEL Container:No additive (red-top) Temperature:REFRIGERATED PREFERRED ROOM TEMPERATURE UNACCEPTABLE FROZEN ACCEPTABLE Place Of Service:QUEST	
UHEALTH Labs	CHPSIM	Chlamydia psittaci IgM Ab Place Of Service:LabCorp	
UHEALTH Labs	CHTRAC	CHLAMYDIA TRACHOMATIS CULTURE Place Of Service:LABCORP	
UHEALTH Labs	CTRACC	CHLAMYDIA TRACHOMATIS CULTURE Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CLMYPN	CHLAMYDIA TRACHOMATIS PANEL Place Of Service:LABCORP	
UHEALTH Labs	CLMNAA	CHLAMYDIA TRACHOMATIS, NAA Place Of Service:LABCORP	
UHEALTH Labs	CTNGTV	Chlamydia trachomatis,Neisseria Gonorrhoeae, and Trichom Place Of Service:LabCorp	
UHEALTH Labs	CHLGCA	CHLAMYDIA/GC AMPLIFICATION Place Of Service:LABCORP	
UHEALTH Labs	CTNG	CHLAMYDIA/NEISSERIA QUAL, DNA PCR Container:URINE CUP KIT/ VIGINAL SWAB Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	87491 87591
UHEALTH Labs	TRIETH	CHLORAL HYDRATE, BLOOD Container:GRAY TOP TUBE Temperature:CHILL AND SHIP WITH COLD PACK Place Of Service:QUEST	
UHEALTH Labs	DRUG16	CHLORDIAZEPOXIDE AND METABOLITES Container:RED TOP TUBE (NO GEL) Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CL	CHLORIDE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82435
UHEALTH Labs	CHLORQ	CHLORPROMAZINE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CHOL	CHOLESTEROL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82465
UHEALTH Labs	CHOLSE	CHOLINESTERASE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	CHRMPL	CHROMIUM, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	CHRMBC	CHROMIUM, RBC Place Of Service:LABCORP	
UHEALTH Labs	CROMS	CHROMIUM, SERUM Container:ACID WASHED OR METAL-FREE VIAL Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	CHGRA	CHROMOGRANIN A Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CHROAB	Chromosome Analysis, Blood Container:Sodium heparin (green-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CHRMP	Chromosome SNP Microarray-Pediatric (Reveal) Place Of Service:LabCorp	
UHEALTH Labs	CHROMO	CHROMOSOME, BLOOD, ROUTINE Place Of Service:LABCORP	
UHEALTH Labs	CHRFA	Chromosomes, DEB Assay for Fanconi Anemia Container:sodium heparin (green-top) tube Temperature:Room temperature Place Of Service:QUEST	88230 88249
UHEALTH Labs	URTIB	CHRONIC URTICARIA Place Of Service:LABCORP	
UHEALTH Labs	STTRYP	CHYMOTRYPSIN, STOOL Container:STERILE SCREW CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CTCC	CIRCULATING TUMOR CELLS, COLON Container:CELLSAVE(TM) PRESERVATIVE TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CTCP	CIRCULATING TUMOR CELLS-PROSTATE Container:CELLSAVE(TM) PRESERVATIVE TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CTRATU	CITRIC ACID (CITRATE), URINE Place Of Service:LABCORP	
UHEALTH Labs	CRECIT	Citric Acid (Citrate):Creatinine Ratio, 24 Hour Urine Container:Use random urine transport tube or urine Monovette® with pH stabilizer (LabCorp N° 56015). Monovette® instructions accompany units. Temperature:Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	CKTISO	CK, TOTAL+ISOENZYMES, SERUM Place Of Service:LABCORP	
UHEALTH Labs	CLDSPM	CLADOSPORIUM HERBARUM (M2) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CLLTP	CLL TARGETED PROFILE Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	CLOB	Clobazam(ONFI), Serum or Plasma Place Of Service:LabCorp	
UHEALTH Labs	CLOMIP	CLOMIPRAMINE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	CLNZPS	CLONAZEPAM (KLONOPIN(R)), SERUM Place Of Service:LABCORP	
UHEALTH Labs	CLOPI	Clopidogrel CYP2C19 Genotyping Place Of Service:LabCorp	
UHEALTH Labs	CLORAZ	CLORAZEPATE (TRANXENE(R)) Container:Plastic, polypropylene screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CLODPC	Clostridium difficile Toxin/GDH with Reflex to PCR Container:Clean plastic screw-cap, leak-proof container Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CDIFFP	CLOSTRIDIUM DIFFICILE,EIA W/RFLX TO PCR Container:CLEAN DRY CONTAINER	87449

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:REFRIGERATED, FOR TRANSPORTATION/DRIVING CONDITIONS ONLY Place Of Service:UHealth Labs	
UHEALTH Labs	CLOZA	CLOZAPINE (CLOZARIL), SERUM Place Of Service:LABCORP	
UHEALTH Labs	CMVAN	CMV ANTIVIRAL RESISTANCE (PLASMA) Temperature:ROOM TEMPERATURE Place Of Service:VIRACOR	
UHEALTH Labs	CMVDNA	CMV DNA, QUANT RT PCR, PLASMA Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	87496
UHEALTH Labs	CMVPRC	CMV PCR Place Of Service:LABCORP	
UHEALTH Labs	CMVANC	CMV QUANT REAL-TIME PCR (CSF) Container:sterile, screw top tube. Temperature:Frozen Place Of Service:VIRACOR	
UHEALTH Labs	CMVANL	CMV QUANT REAL-TIME PCR(BAL) Container:sterile, screw top tube. Temperature:ambient or frozen temperature Place Of Service:VIRACOR	
UHEALTH Labs	COBLT	COBALT, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	COCAM	Cocaine Metabolites Clinical Screen, with Confirmation, Urine Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	COCCF	Coccidioidal Complement Fixation (CF) Container:Serum separator tube or plain red Temperature:Refrigerate Place Of Service:UHEALTH LABS	86171
UHEALTH Labs	COCID	Coccidioidal Immunodiffusion (ID) Container:Serum separator tube or plain red Temperature:Refrigerate Place Of Service:UHEALTH LABS	86331
UHEALTH Labs	CIDCF	Coccidioidal Reflex ID to CF Container:Serum separator tube or plain red Temperature:Refrigerate Place Of Service:UHEALTH LABS	86171 86331
UHEALTH Labs	COCAB	COCCIDIOIDES ANTIBODY Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	COCCAG	Coccidioides QN Antigen, EIA, Serum Container:Collect serum specimens in serum separator or red top tube. Temperature:Room temperature Place Of Service:MiraVista	87449
UHEALTH Labs	COCKR	COCKROACH, AMERICAN Place Of Service:LABCORP	
UHEALTH Labs	CODFI	CODFISH (F3) IgE Container:No additive (red-top) Temperature:REFRIGERATED PREFERRED ROOM TEMPERATURE ACCEPTABLE FROZEN ACCEPTABLE Place Of Service:QUEST	
UHEALTH Labs	COQ10B	COENZYME Q10, TOTAL Place Of Service:LABCORP	
UHEALTH Labs	AGGLUT	COLD AGGLUTININ TITER, QUANT Place Of Service:LABCORP	
UHEALTH Labs	CHAG	COLD HEMAGGLUTININS Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	PFADP	COLLAGEN ADP Container:BLUE Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CCTX	COLLAGEN TYPE I C-TELOPEPTIDE Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CMC1Q	COMPLEMENT C1Q, QUANTITATIVE Place Of Service:LABCORP	
UHEALTH Labs	CMPC2B	COMPLEMENT C2 Place Of Service:LABCORP	
UHEALTH Labs	C3	COMPLEMENT C3 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86160
UHEALTH Labs	COMP34	COMPLEMENT C3 AND C4 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86160
UHEALTH Labs	C4	COMPLEMENT C4 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86160
UHEALTH Labs	C4A	Complement C4a Place Of Service:LabCorp	
UHEALTH Labs	CH50T	COMPLEMENT, TOTAL (CH50) Place Of Service:LABCORP	
UHEALTH Labs	CADP	Complete Atopic Dermatitis Panel Container:Gold or Red Temperature:Room Temperature Place Of Service:Viracor	83520 86003 82785
UHEALTH Labs	CBC	COMPLETE BLOOD COUNT (CBC) Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85025
UHEALTH Labs	UMACUM	COMPLETE URINALYSIS (W/ MICROSCOPY) Container:Urine Cup Urine Transport Tube Temperature:Urine cup: Room temperature up to 2 hours and refrigerated up to 48 hours. Urine transport tube: Refrigerated or room temperature for up to 72 hours. Place Of Service:UHealth Labs	81001
UHEALTH Labs	COMDRG	COMPLIANCE DRUG ANALYSIS, UR Place Of Service:LABCORP	
UHEALTH Labs	CMP	COMPREHENSIVE METABOLIC PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80053
UHEALTH Labs	COOMB	COOMBS', DIRECT Place Of Service:LABCORP	
UHEALTH Labs	CPRRBC	COPPER, PLASMA AND RBC Place Of Service:LABCORP	
UHEALTH Labs	COPPER	COPPER, SERUM Place Of Service:LABCORP	
UHEALTH Labs	COPUR	COPPER, URINE Place Of Service:LABCORP	
UHEALTH Labs	CORT	CORTISOL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82533

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CORTA	Cortisol - AM Place Of Service:LabCorp	
UHEALTH Labs	CTDXRF	CORTISOL DEXAMETHASONE REFLEX Place Of Service:LABCORP	
UHEALTH Labs	CRT2SP	CORTISOL, 2 SPECIMENS Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CORTFT	CORTISOL, FREE & TOTAL Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	COR24F	CORTISOL, FREE 24-HOUR URINE LC/MS/MS Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	COACTH	Cortisol, Free, Response to ACTH Stimulation Container:Plastic screw-cap vials Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CORTT	CORTISOL, TOTAL Container:Serum Separator Tube (SST) Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	CORTLC	Cortisol, Total, LC/MS/MS Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CNCLU	Cortisol:Cortisone Ratio, 24 Hour Urine Container:Plastic urine container Place Of Service:LABCORP	
UHEALTH Labs	CORTU	CORTISON, URINARY FREE Place Of Service:LABCORP	
UHEALTH Labs	CORTNE	CORTISONE, 24 HOUR URINE Container:Plastic urine container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CORTSN	CORTISONE, SERUM Place Of Service:QUEST	
UHEALTH Labs	CORTUR	CORTISONE, URINE Container:24 hour Urine Container Place Of Service:LABCORP	
UHEALTH Labs	COWDAL	COW DANDER (E4) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	COXBVA	Coxsackie B Virus Antibodies Place Of Service:LABCORP	
UHEALTH Labs	COXGM	COXSACKIE A IGG/IGM ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	CPEP	C-PEPTIDE Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CPEPS	C-PEPTIDE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	CPK	CPK Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82550

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CKMB	CPK MB ISOENZYME Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CRAB	CRAB IgE ALLERGEN Place Of Service:QUEST	
UHEALTH Labs	CRP	C-REACTIVE PROTEIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86140
UHEALTH Labs	CRPCAR	C-REACTIVE PROTEIN, CARDIAC Place Of Service:LABCORP	
UHEALTH Labs	URCREA	CREATINE, 24 HOUR URINE Container:PLASTIC URINE CONTAINER Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CREA	Creatine, Serum Place Of Service:LABCORP	
UHEALTH Labs	CREATS	CREATINE, SERUM Container:Transport Container: Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CRTUPH	CREATININE AND PHOSPHORUS, 24 HOUR Container:Plastic urine container, with or without 6N HCl Temperature:Room Temperature Place Of Service:LABCORP	
UHEALTH Labs	CRT24U	CREATININE, 24-HOUR URINE Place Of Service:LABCORP	
UHEALTH Labs	CREATB	CREATININE,BLOOD Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82565
UHEALTH Labs	CRE5HI	CREATU+5 HIAA Place Of Service:LABCORP	
UHEALTH Labs	CRECAU	CREATU+CATECU Temperature:CreatU Temp: Frozen CatecU Temp: Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	ORTAC	CROSSMATCH BLOOD BANK PRODUCT-OR ONLY Place Of Service:JMH	
UHEALTH Labs	CRYFGN	CRYOFIBRINOGEN Container:2 PLASTIC SCREW-CAP VIALS Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CRYRFX	CRYOGLOBULIN, QL, SERUM, RFLX Place Of Service:LABCORP	
UHEALTH Labs	CRYPQ	CRYPTOCOCCAL AG SCR W/REFL TITR Container:CSF: STERILE, SCREW-CAP TUBE; SERUM: PLASTIC SCREW-CAP VIAL Temperature:SERUM OR CSF: REFRIGERATED (COLD PACKS) OR FROZEN Place Of Service:QUEST	
UHEALTH Labs	CRYAGB	Cryptococcal Antigen in Blood only Container:Gold Top, Red Top Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CRYAGC	Cryptococcal Antigen in CSF only Container:Sterile container Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CRYDFA	CRYPTOSPORIDIUM AG, STOOL BY DFA Container:10% formalin or Total-Fix transport vial	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	GIACRP	CRYPTOSPORIDIUM/GIARDIA AG, STOOL Container:PARA PAK CONTAINER( ALL 3) Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87015
UHEALTH Labs	SYNCRY	CRYSTALS, SYNOVIAL FLUID Container:LAVENDER TOP TUBE (EDTA) or GREEN TOP TUBE (Sodium-Heparin) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	89060
UHEALTH Labs	CSFCNT	CSF CELL COUNT Container:Sterile container Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CSFC	CSF CULTURE Container:STERILE CUP Temperature:STORE AT RT, DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	CSFGLU	CSF GLUCOSE Container:CSF Liquid/ CSF transfer tube container Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CSFGP	CSF GLUCOSE, TOTAL PROTEIN Container:CSF Liquid/ CSF transfer tube container Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CSFTP	CSF TOTAL PROTEIN Container:CSF Liquid/ CSF transfer tube container Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CTGCRL	CT/GC NAA, RECTAL Place Of Service:LABCORP	
UHEALTH Labs	CTELOP	C-TELOPEPTIDE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	UCLS	CULTURE, URINE, SPECIAL Container:Urine collection kits (gray-top tube); Sterile container for specimens < 5 mL (for pediatric patients only) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	VICULT	Culture, Viral, Body Fluids, Tissues Container:Tissue or conjunctival swab: Submit in V-C-M medium (green-cap) tube or equivalent Respiratory, newborn urine, CSF, and body fluids: Submit in an equal volume of V-C-M medium (green-cap) tube or equivalent, or in a sterile leak-proof container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	IMCUTD	CUTANEOUS DIRECT IMMUNOFLUORESCEN Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CYAN	CYANIDE BLOOD Container:GREEN-TOP (SODIUM HEPARIN) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CCP	CYCLIC CITRULLINATED PEPTIDE IGG, SERUM Container:No additive (red-top) Temperature:REFRIGERATED PREFERRED, ROOM TEMPERATURE ACCEPTABLE, FROZEN ACCEPTABLE Place Of Service:JMH	
UHEALTH Labs	CYCLSP	Cyclospora Smear, Stool Place Of Service:LABCORP	
UHEALTH Labs	YCYCLR	CYCLOSPORIN, LC/MS/MS Place Of Service:IML	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CYCLOB	CYCLOSPORINE, BLOOD Place Of Service:LABCORP	
UHEALTH Labs	CYSCGF	Cystatin C w/Glomerular Filtration Rate,Est.(eGFR) Container:Transport tube Temperature:Room temperature Place Of Service:Quest	82610
UHEALTH Labs	CYSC	CYSTATIN-C Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	CF97	Cystic Fibrosis, 97 Variants Place Of Service:LABCORP	
UHEALTH Labs	RESPCF	CYSTIC FIBROSIS,RESPIRATORY CULT W/GS Container:STERILE CUP Temperature:STORE AT RT, DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	CYSTSS	CYSTICERCOSIS (TAENIA SOLIUM) Place Of Service:LABCORP	
UHEALTH Labs	CYSTGA	CYSTICERCUS ANTIBODY, ELISA Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CYSIGG	Cysticercus IgG ELISA, CSF Place Of Service:LABCORP	
UHEALTH Labs	CYQ24B	CYSTINE, QUANT. 24 HR URINE Place Of Service:LABCORP	
UHEALTH Labs	UCYSQN	CYSTINE, QUANTITATIVE, RANDOM URINE Container:Transport Container: Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	CYTP13	Cytokine Panel 13, Serum Container:RED Temperature:Frozen Place Of Service:QUEST	83520 x12
UHEALTH Labs	CYTO TESTING	CYTOLOGY TESTING Place Of Service:UHealth Labs	
UHEALTH Labs	CYTOIG	CYTOMEGALOVIRUS (CMV) AB, IGG Place Of Service:LABCORP	
UHEALTH Labs	CMVIGB	CYTOMEGALOVIRUS (CMV) AB, IGM Temperature:Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	YCTL	CYTOTOXIC T LYMPHOCYTES FUNCTION Temperature:SAMPLES SHOULD BE SENT AT ROOM TEMPERATURE, UNLESS INFORMED OTHERWISE. PACKAGE SECURELY TO AVOID BREAKAGE AND EXTREME WEATHER CONDITIONS. Place Of Service:CINCIN	
UHEALTH Labs	DXYLOS	D XYLOSE ABSORPTION TEST 2 HR Container:SODIUM FLUORIDE (GRAY-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DPTERE	D001-IG D PTERONYSSINUS Place Of Service:LABCORP	
UHEALTH Labs	DFARIE	D002-IGE D FARINAE Place Of Service:LABCORP	
UHEALTH Labs	BLOTRE	D201-IGE BLOMIA TROPICALIS Place Of Service:LABCORP	
UHEALTH Labs	DDIMER	D-DIMER Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN	85379

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	DWDC	Place Of Service:UHealth Labs DEEP WOUND CULTURE WITH GRAM STAIN Container:STERILE CUP OR CULTURETTE Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	80707
UHEALTH Labs	DENGUE	DENGUE FEVER ANTIBODIES (IgG,IgM) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DNGVGM	DENGUE VIRUS IGG AND IGM Place Of Service:LABCORP	
UHEALTH Labs	DNGVIR	DENGUE VIRUS IGG ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	DENGVR	DENGUE VIRUS IGM ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	DNGRNA	Dengue Virus RNA, Qualitative, Real-Time PCR Container:GOLD TOP/RED TOP Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	DOCSB	DEOXYCORTICOSTERONE (DOC),S Place Of Service:LABCORP	
UHEALTH Labs	DFAR	DERM. FARINAE (MITES) (D2) IgE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DRMPT	DERMATOPHAGOIDES PTERONYSSINUS IGE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DGMMAB	DES GAMMA CARBOXY PROTHROMBIN Place Of Service:LABCORP	
UHEALTH Labs	DESIP	DESIPRAMINE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	DESIGG	Desmoglein 1 and Desmoglein 3 (Pemphigus), IgG Abs Place Of Service:LABCORP	
UHEALTH Labs	DESMAB	DESMOGLEIN ANTIBODIES (1 AND 3) Container:Transport tube Temperature:Refrigerated (cold packs) Place Of Service:Quest	83520 x2 84999
UHEALTH Labs	DXASER	DEXAMETHASONE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	DHEAU	DHEA, LC/MS/MS Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DHEASR	DHEA, SERUM Place Of Service:LABCORP	
UHEALTH Labs	DHEAST	DHEA, SERUM Place Of Service:LABCORP	
UHEALTH Labs	DHTFRB	DHT,FREE,LCMS/DIALYSIS Place Of Service:LABCORP	
UHEALTH Labs	DRUG19	DIAZEPAM AND METABOLITE Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	DIG	DIGOXIN LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERRED), GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE	80162

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:UHealth Labs	
UHEALTH Labs	DPYD	Dihydropyrimidine Dehydrogenase (DPYD), 3 Variants Container:Lavender (EDTA), pink (K2EDTA), or yellow (ACD Solution A or B) Temperature:Refrigerated. Place Of Service:ARUP	81232 81346
UHEALTH Labs	DHYTST	DIHYDROTESTOSTERONE Place Of Service:LABCORP	
UHEALTH Labs	PTN	DILANTIN LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERRED), GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80185
UHEALTH Labs	DRVV	DILUTE RUSSEL VIPER VENON TEST Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85613
UHEALTH Labs	DIPAHB	DIPHTHERIA ANTITOXOID AB Place Of Service:LABCORP	
UHEALTH Labs	DATQ	DIRECT ANTIGLOBULIN TEST (DAT) Container:ACD-A or ACD-B (yellow-top) tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	ANTIGL	DIRECT ANTIGLOBULIN W/REFL ANTI C3,ANTI IGG Container:ACD-A or ACD-B (yellow-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DBIL	DIRECT BILIRUBIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:PROTECTED FROM LIGHT, ROOM TEMPERATURE. Place Of Service:UHealth Labs	
UHEALTH Labs	DAT	Direct Coombs Container:Pink (EDTA) Tube Temperature:Transport at room temperature Place Of Service:UHealth Labs	86880
UHEALTH Labs	DLDL	DIRECT LDL Container:Serum Separator Tube (SST) Temperature:Room Temperature Place Of Service:UHEALTH LABS	83721
UHEALTH Labs	DLDLQ	DIRECT LDL Container:SERUM SEPARATOR TUBE Alt. Spec. Types: PLASMA COLLECTED IN: SODIUM HEPARIN (GREEN-TOP) TUBE OR LITHIUM HEPARIN (GREEN-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DISAC	Disaccharidases, Tissue Place Of Service:LabCorp	
UHEALTH Labs	DISOP	DISOPYRAMIDE Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	DIUR	Diuretic Screen, Urine Place Of Service:LabCorp	
UHEALTH Labs	DLAC	D-Lactate, Plasma Container:Grey top (potassium oxalate/sodium fluoride) (T275). Submission Container/Tube: Plastic vial Temperature:Frozen Place Of Service:Mayo	83605
UHEALTH Labs	DNAQ	DNA (DS) ANTIBODY Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	DNAAB	DNA (ds) ANTIBODY, HIGH AVIDITY Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	DNAB	DNASE B ANTIBODY Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	WNVDON	DONOR West Nile Virus (NAT) Assay Place Of Service:LabCorp	
UHEALTH Labs	DHHHP	Donor,HIV-1/HBV/HCV NAT Procleix Rfx HIV1/HBV/HCV Container:EDTA (lavender-top) tube Temperature:Refrigerated (cold packs) Place Of Service:QUEST	87801
UHEALTH Labs	DOXE	DOXEPIN Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DPYDG	DPYDG Genotyping Place Of Service:LABCORP	
UHEALTH Labs	DAP10	DRUG ABUSE PANEL 10 Container:PLASTIC URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DRGP10	DRUG ABUSE PANEL 10 WITH CONFIRMATION,URINE Container:Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DAPU	Drug Abuse Panel 10-50 plus Ethanol, w/o Confirm., Ur Container:Plastic urine container Temperature:60 mL random urine Place Of Service:Quest	80307
UHEALTH Labs	DAP750	Drug Abuse Panel 7-50, with Confirmation, Urine Container:Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DAS9	DRUG ABUSE PANEL 9, SERUM Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ALCWCU	Drug Monitoring, Alcohol Metabolite, with Confirmation, Urine Container:Clinical drug test transport vial; Alternate :urine collection container Temperature:Room Temperature Place Of Service:Quest	
UHEALTH Labs	DAU	DRUG OF ABUSE SCREEN, URINE Container:URINE random collection CUP/ Yellow top tube no additive Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	MECON	DRUG SCREEN PANEL 9, MECONIUM Container:STERILE LEAKPROOF CONTAINER Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	DSCWU	Drug Screen, Clinical 1, without Confirmation, Urine Container:Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DSCR	DRUG SCREEN,COMP (SERUM/PLASMA) Container:SERUM: PLASTIC SCREW-CAP VIAL. OPTIONAL WHOLE BLOOD: UNOPENED GRAY-TOP TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	DRGURN	DRUG SCREEN,COMPREHENSIVE (URINE) Container:Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	DTM1OF	Drug Toxicology Monitoring 1, Quantitative, Oral Fluid Container:Oral-Eze collection device Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	DOASV	DRUGS OF ABUSE SCRIN SER AND VOLAT Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	CRITHD	DS DNA ANTIBODY CRITHIDIA IFA RFL Container:NO ADDITIVE (RED-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SHIG3	E coli O157 CULTURE, STOOL Container:Cary-Blair stool culture transport medium Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	CATDE	E001-IGE CAT DANDER Place Of Service:LABCORP	
UHEALTH Labs	HRSEDE	E003-IGE HORSE DANDER Place Of Service:LABCORP	
UHEALTH Labs	COWDE	E004-IGE COW DANDER Place Of Service:LABCORP	
UHEALTH Labs	DOGDE	E005-IGE DOG DANDER Place Of Service:LABCORP	
UHEALTH Labs	GOOSEE	E070-IGE GOOSE FEATHERS Place Of Service:LABCORP	
UHEALTH Labs	RAUIGE	E074-IgE Rat Urine Place Of Service:LabCorp	
UHEALTH Labs	CHIKFE	E085-IGE CHICKEN FEATHERS Place Of Service:LABCORP	
UHEALTH Labs	EARCUL	EAR CULTURE Container:STERILE CUP, SCREW-CAPPED, STERILE TUBE OR CONTAINER Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	EBVVCA	EBV AB VCA, IGG Place Of Service:LABCORP	
UHEALTH Labs	EBVIGM	EBV AB VCA, IGM Place Of Service:LABCORP	
UHEALTH Labs	EBVAP	EBV Antibody Profile Place Of Service:LabCorp	
UHEALTH Labs	YEBVPB	EBV DNA PCR (BAL) Place Of Service:VIRACOR	
UHEALTH Labs	EBVPC	EBV DNA PCR (WHBLD) Container:EDTA (lavender top) tube Place Of Service:VIRACOR	
UHEALTH Labs	EBVDP	EBV DNA QUAL REAL-TIME PCR Container:EDTA (LAVENDER-TOP) Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	EBVIGG	EBV EARLY ANTIGEN AB, IGG Place Of Service:LABCORP	
UHEALTH Labs	EBVG	EBV EARLY ANTIGEN D AB (IGG) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	EBVAGB	EBV NUCLEAR ANTIGEN AB, IGG Place Of Service:LABCORP	
UHEALTH Labs	EBVGM	EBVCA (IGG+IGM)+EBVNIG Place Of Service:LABCORP	
UHEALTH Labs	EBVEA	EBVEA+EBV Place Of Service:LabCorp	
UHEALTH Labs	ECHAB	ECHINOCOCCUS ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	ECHOP	ECHOVIRUS ANTIBODY PANEL, CSF Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ECHVIP	ECHOVIRUS SERUM AB PANEL Place Of Service:LABCORP	
UHEALTH Labs	EGGY	EGG YOLK (F75) IgE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	EGGPLA	Eggplant IgE (F262) Place Of Service:LABCORP	
UHEALTH Labs	EHRLIC	EHRlichia CHAFFEENSIS Ab(IgG,IgM) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LYTES	ELECTROLYTES, BLOOD Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80051
UHEALTH Labs	LYTFEC	ELECTROLYTES, FECAL Place Of Service:LABCORP	
UHEALTH Labs	ENCCSF	ENCEPHALITIS ANTIBODY PANEL CSF Container:Specimen Container: PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	ENCE	ENCEPHALITIS EVALUATION, COMP Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ENDIGA	ENDOMYSIAL ANTIBODY IGA Place Of Service:LABCORP	
UHEALTH Labs	ENTPAN	ENTERIC BACTERIAL PANEL,STOOL by PCR Container:Fresh Stool or Cary-Blair Temperature:Room Temperature or Refrigerated Place Of Service:UHEALTH LABS	87798 x8
UHEALTH Labs	ENTVCF	Enterovirus Qn Real-time RT-PCR, CSF Container:Sterile, screw top tube Temperature:Store frozen and ship on dry ice for overnight delivery. Place Of Service:Viracor	87799
UHEALTH Labs	YENTER	ENTEROVIRUS qPCR (PL) Container:EDTA (lavender-top) tube Place Of Service:VIRACOR	
UHEALTH Labs	NTROV	ENTEROVIRUS RT-PCR Place Of Service:LABCORP	
UHEALTH Labs	ENVCUL	ENVIRONMENTAL CULTURE Container:TRYPTICASE SOY AGAR WITH LECITHIN AND POLYSORBATE 80, ISOLATOR PACK (SURFACE SAMPLING) Temperature:STORE AT RT DO NOT REFRIGERATE Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	ENZBIO	ENZYME BIOTINIDASE DEF, SERUM Place Of Service:LABCORP	
UHEALTH Labs	EOSCNT	EOSINOPHIL COUNT Place Of Service:LABCORP	
UHEALTH Labs	EOSNA	EOSINOPHIL COUNT, NASAL Container:SLIDE IN SLIDE HOLDER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ESCTNA	EOSINOPHIL COUNT, NASAL Place Of Service:LABCORP	
UHEALTH Labs	EOSPUT	EOSINOPHIL COUNT, SPUTUM Container:Slide in slide holder; Alternate: Sputum submitted in a sterile cup Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	EPILEP	Epilepsy, Autoimmune Evaluation, Serum Container:Preferred: Red top; Acceptable: Serum gel Temperature:Refrigerated (preferred) Place Of Service:Mayo	83519 x5 86255 x14
UHEALTH Labs	EAESF	Epilepsy, Autoimmune Evaluation, Spinal Fluid Container:Sterile vial Temperature:Refrigerated (preferred) 28 days Place Of Service:MAYO	86255 x13 86341 83519
UHEALTH Labs	EPINE	EPINEPHRINE, PLASMA Container:Lavender-top (EDTA) tube or green-top (heparin) tube Temperature:Frozen Place Of Service:LABCORP	82384
UHEALTH Labs	EPIPL	EPINEPHRINE, PLASMA Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	EPORMD	EPOR Mutation Detection by Sequencing Container:Lavender (EDTA) Temperature:Refrigerated Place Of Service:ARUP	81479
UHEALTH Labs	EBVQT	EPSTEIN BARR VIRUS DNA, QN RT PCR Container:EDTA (lavender-top) tube Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	EPBAPC	EPSTEIN-BARR DNA QUANT, PCR Place Of Service:LABCORP	
UHEALTH Labs	EPSPCR	EPSTEIN-BARR VIRUS PCR Place Of Service:LABCORP	
UHEALTH Labs	PROTO	ERYTHROCYTE PROTOPORPHYRIN (EP) Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ESR	ERYTHROCYTE SEDIMENTATION RATE Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85652
UHEALTH Labs	EGSP	ERYTHROCYTOSIS GENE SEQUENCING PANEL Container:whole blood in lavender top (EDTA) OR Saliva (use Saliva collection kit) Temperature:Store at room temperature Place Of Service:CINCIN	81404 81479 x3
UHEALTH Labs	ERYP	ERYTHROPOIETIN LEVEL Container:Serum Separator Tube (SST) Temperature:Room Temperature Place Of Service:UHEALTH LABS	83668
UHEALTH Labs	ERYPO	ERYTHROPOIETIN LEVEL Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	ESTD	ESTRADIOL LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82670
UHEALTH Labs	ESFRSR	ESTRADIOL, FREE SERUM Place Of Service:LABCORP	
UHEALTH Labs	ESTRFR	ESTRADIOL, FREE, LC/MS/MS Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	ESTRIO	ESTRIOL, LC/MS/MS, SERUM Container:RED-TOP TUBE (NO-GEL); TRANSPORT CONTAINER: PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	ESTROL	ESTRIOL, SERUM Place Of Service:LABCORP	
UHEALTH Labs	ESGNF	ESTROGENS FRACTIONATED, S Place Of Service:LABCORP	
UHEALTH Labs	ESTRGT	ESTROGENS, TOTAL Place Of Service:LABCORP	
UHEALTH Labs	ELCMS	ESTRONE, LCMS, ENDO SCI Place Of Service:LABCORP	
UHEALTH Labs	ESTRNE	ESTRONE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	ETHOS	ETHOSUXIMIDE Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	ETHGYL	ETHYLENE GLYCOL, SERUM Place Of Service:LABCORP	
UHEALTH Labs	EVEROL	Everolimus Level Container:EDTA (LAVENDER-TOP) TUBE Temperature:Send Room Temperature	
UHEALTH Labs	AENA	Extractable Nuclear Ag Place Of Service:LabCorp	
UHEALTH Labs	EYECUL	EYE CULTURE Container:STERILE CUP, SCREW-CAPPED, STERILE TUBE OR CONTAINER Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	EGGWHE	F001-IGE EGG WHITE Place Of Service:LABCORP	
UHEALTH Labs	CODFE	F003-IGE CODFISH Place Of Service:LABCORP	
UHEALTH Labs	WHEATE	F004-IGE WHEAT Place Of Service:LABCORP	
UHEALTH Labs	RYEE	F005-IGE RYE Place Of Service:LABCORP	
UHEALTH Labs	BARLE	F006-IGE BARLEY Place Of Service:LABCORP	
UHEALTH Labs	CORNE	F008-IGE CORN Place Of Service:LABCORP	
UHEALTH Labs	RICEE	F009-IGE RICE Place Of Service:LABCORP	
UHEALTH Labs	SASAME	F010-IGE SESAME SEED Place Of Service:LABCORP	
UHEALTH Labs	SESAME	F010-IGE SESAME SEED Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	BUCKE	F011-IGE BUCKWHEAT Place Of Service:LABCORP	
UHEALTH Labs	GRPEE	F012-IGE GREEN PEA Place Of Service:LABCORP	
UHEALTH Labs	PNUTE	F013-IGE PEANUT Place Of Service:LABCORP	
UHEALTH Labs	SOYBE	F014-IGE SOYBEAN Place Of Service:LABCORP	
UHEALTH Labs	WHBNE	F015-IGE WHITE BEAN Place Of Service:LABCORP	
UHEALTH Labs	HAZELE	F017-IGE HAZELNUT (FILBERT) Place Of Service:LABCORP	
UHEALTH Labs	BRAZLE	F018-IGE BRAZIL NUT Place Of Service:LABCORP	
UHEALTH Labs	ALMNDE	F020-IGE ALMOND Place Of Service:LABCORP	
UHEALTH Labs	CRABE	F023-IGE CRAB Place Of Service:LABCORP	
UHEALTH Labs	SHRMPE	F024-IGE SHRIMP Place Of Service:LABCORP	
UHEALTH Labs	TOMATE	F025-IGE TOMATO Place Of Service:LABCORP	
UHEALTH Labs	PORKE	F026-IGE PORK Place Of Service:LABCORP	
UHEALTH Labs	CAROTE	F031-IGE CARROT Place Of Service:LABCORP	
UHEALTH Labs	ORNGE	F033-IGE ORANGE Place Of Service:LABCORP	
UHEALTH Labs	WHPOEB	F035-IGE POTATO, WHITE Place Of Service:LABCORP	
UHEALTH Labs	COCNTE	F036-IGE COCONUT Place Of Service:LABCORP	
UHEALTH Labs	MUSSLE	F037-IGE MUSSEL Place Of Service:LABCORP	
UHEALTH Labs	TUNAE	F040-IGE TUNA Place Of Service:LABCORP	
UHEALTH Labs	SALMOE	F041-IGE SALMON Place Of Service:LABCORP	
UHEALTH Labs	EGGYLK	F075-IGE EGG (YOLK) Place Of Service:LABCORP	
UHEALTH Labs	ALLAE	F076-IGE ALPHA LACTALBUMIN Place Of Service:LABCORP	
UHEALTH Labs	LACTOE	F077-IGE BETA LACTOGLOBULIN Place Of Service:LABCORP	
UHEALTH Labs	CASINE	F078-IGE CASEIN Place Of Service:LABCORP	
UHEALTH Labs	GLUIGE	F079-IGE GLUTEN Place Of Service:LABCORP	
UHEALTH Labs	LOBSE	F080-IGE LOBSTER Place Of Service:LABCORP	
UHEALTH Labs	CHICKE	F083-IGE CHICKEN Place Of Service:LABCORP	
UHEALTH Labs	CELRYE	F085-IGE CELERY Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	PARSYE	F086-IGE PARSLEY Place Of Service:LABCORP	
UHEALTH Labs	MANGOE	F091-IGE MANGO Place Of Service:LABCORP	
UHEALTH Labs	PECANE	F201-IGE PECAN NUT Place Of Service:LABCORP	
UHEALTH Labs	CASHE	F202-IGE CASHEW NUT Place Of Service:LABCORP	
UHEALTH Labs	PISTE	F203-IGE PISTACHIO NUT Place Of Service:LABCORP	
UHEALTH Labs	CLAME	F207-IGE CLAM Place Of Service:LABCORP	
UHEALTH Labs	WHEGGE	F245-IGE EGG, WHOLE Place Of Service:LABCORP	
UHEALTH Labs	WLNTB	F256-IGE WALNUT Place Of Service:LABCORP	
UHEALTH Labs	SQUIDE	F258-IGE SQUID Place Of Service:LABCORP	
UHEALTH Labs	OYSTRE	F290-IGE OYSTER Place Of Service:LABCORP	
UHEALTH Labs	SWCHE	F299-IGE SWEET CHESTNUT Place Of Service:LABCORP	
UHEALTH Labs	SCALPE	F338-IGE SCALLOP Place Of Service:LABCORP	
UHEALTH Labs	MACNTE	F345-IGE MACADAMIA NUT Place Of Service:LABCORP	
UHEALTH Labs	FACTH	Factor H (B1H) Level Place Of Service:LabCorp	
UHEALTH Labs	FII	FACTOR II ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85210
UHEALTH Labs	FIIDNA	FACTOR II, DNA ANALYSIS Place Of Service:LABCORP	
UHEALTH Labs	FIX	FACTOR IX ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85250
UHEALTH Labs	FIXAC	FACTOR IX ACTIVITY AND HUMAN INHIBITOR Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85250
UHEALTH Labs	F9INHB	FACTOR IX INHIBITOR Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	FV	FACTOR V ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	83891 83903 x2 83912 83892 83896 x2
UHEALTH Labs	F5ACIN	Factor V Activity and Human Inhibitor Container:Two plastic screw-cap vials Temperature:Frozen Place Of Service:UHEALTH LABS	85220 85335
UHEALTH Labs	F5INH	FACTOR V INHIBITOR Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	F5INHB	FACTOR V INHIBITOR Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	83891 83903 x2 83912 83892 83896 x2
UHEALTH Labs	F5LEID	FACTOR V LEIDEN Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85220
UHEALTH Labs	FVMAW	FACTOR V(LEIDEN) MUTATION ANALYSIS W/RFL HR2 MUT Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	FVII	FACTOR VII ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85230
UHEALTH Labs	F8SA	Factor VIII (F8) Sequence Analysis Container:Fetal: Amniotic Fluid, CVS, Cultured Amniocytes or Cultured CVS; Parental/Patient: EDTA Whole Blood (lavender top) Temperature:Room Temperature Place Of Service:UHEALTH LABS	81407
UHEALTH Labs	F8CHRG	FACTOR VIII ACT, CLOTTING W/RFL CHROMOGENIC Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	FVIII	FACTOR VIII ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85240
UHEALTH Labs	F8CH	Factor VIII Activity, Chromogenic Container:Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	F8INH	FACTOR VIII INHIBITOR, HUMAN Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	FX	FACTOR X ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85260
UHEALTH Labs	F10CHR	FACTOR X CHROMOGENIC Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	FXI	FACTOR XI ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85270
UHEALTH Labs	FXIAC	FACTOR XI ACTIVITY AND HUMAN INHIBITOR Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85270
UHEALTH Labs	F11IN	FACTOR XI INHIBITOR Container:BLUE Temperature:FROZEN Place Of Service:UHealth Labs	
UHEALTH Labs	FXII	FACTOR XII ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85280
UHEALTH Labs	FXIIIF	FACTOR XIII FUNCTIONAL Temperature:FROZEN Place Of Service:QUEST	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	F13	Factor XIII Semi-quantitative Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85291
UHEALTH Labs	FMFMA	Familial Mediterranean Fever Mutation Analysis Container:EDTA (lavender or royal blue-top),ACD (yellow-top), sodium or lithium heparin (green-top) tube Temperature:Room temperature Place Of Service:QUEST	81402
UHEALTH Labs	FANAN	Fanconi Anemia DNA Mutation Analysis Container:EDTA (lavender-top) tube Temperature:Room temperature Place Of Service:QUEST	81242
UHEALTH Labs	FANAGS	Fanconi Anemia Gene Sequencing Panel Container:EDTA (lavender top) for whole blood. Saliva kit for saliva. Temperature:Room temperature Place Of Service:UHEALTH LAB	81443
UHEALTH Labs	FATPRO	FATTY ACID PROFILE, ESSENTIAL Place Of Service:LABCORP	
UHEALTH Labs	FEBABP	FEBRILE ANTIBODY PROFILE Place Of Service:LABCORP	
UHEALTH Labs	FECQU	FECAL FAT, QUALITATIVE Place Of Service:LABCORP	
UHEALTH Labs	FFQUA	FECAL FAT, QUANTITATIVE Place Of Service:LABCORP	
UHEALTH Labs	FECGIM	FECAL GLOBIN BY IMMUNOCHEM Container:FOBT (STOOL COLLECTION KIT) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82710
UHEALTH Labs	FGIMM	FECAL GLOBIN BY IMMUNOCHEM. (MEDICARE) Container:FOBT (STOOL COLLECTION KIT) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82710
UHEALTH Labs	FECWBC	FECAL WBC Container:STERILE CUP Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	89055
UHEALTH Labs	FELB	FELBAMATE (FELBATOL) Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	URFENT	FENTANYL AND NORFENTANYL, URINE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	FNTNOR	FENTANYL/NORFENTANYL, CONFIRM Place Of Service:LABCORP	
UHEALTH Labs	FER	FERRITIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82728
UHEALTH Labs	FGF23B	FGF-23 Place Of Service:LABCORP	
UHEALTH Labs	FIB	FIBRINOGEN ACTIVITY Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85384
UHEALTH Labs	FIBP	Fibrinogen Evaluation Profile Container:Blue-top (sodium citrate) tubes Temperature:Frozen Place Of Service:UHealth Labs	85384 85385 85635

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	FDP	FIBRIN DEGRADATION PRODUCT Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85362
UHEALTH Labs	FILAB	FILARIA ANTIBODY (IGG4) Place Of Service:LABCORP	
UHEALTH Labs	FIREAN	Fire Ant (Invicta), IgE, I070 Place Of Service:LabCorp	
UHEALTH Labs	FCEP12	FISH (CLL PANEL) Container:GreenNa, Lav Temperature:Room Temperature Place Of Service:CYTO	
UHEALTH Labs	FDUP	FISH, Duplication 15q11q13 Container:Sodium heparin (green-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	FAML	FISH,AML M3,PML/RARA,TRANS 15,17 Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	FCA	FLD CALCIUM Container:RED top tube (No gel) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FLEC	FLECAINIDE (TAMBOCOR) LEVEL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	FLECSP	Flecainide, Serum or Plasma Place Of Service:LABCORP	
UHEALTH Labs	FLOUN	Flounder IgE (F147) Place Of Service:LABCORP	
UHEALTH Labs	FWCYTO	FLOW CYTOMETRY PNH Place Of Service:LABCORP	
UHEALTH Labs	FLT3	FLT3 MUTATION ANALYSIS ITD/TKD Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	FALB	FLUID ALBUMIN Container:Sterile Transport tube Temperature:Frozen Place Of Service:UHealth Labs	82042
UHEALTH Labs	FCL	FLUID CHLORIDE Container:RED top tube (No gel) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FCR	FLUID CRYSTALS Container:RED TOP TUBE- NO ADDITIVE Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FGLUC	FLUID GLUCOSE Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FHCT	FLUID HEMATOCRIT Container:RED TOP TUBE- NO ADDITIVE Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FLDH	FLUID LDH Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	FMAG	FLUID MAGNESIUM Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FPHOS	FLUID PHOSPHORUS Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FK	FLUID POTASSIUM Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FNA	FLUID SODIUM Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FTP	FLUID TOTAL PROTEIN Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FUN	FLUID UREA NITROGEN Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	FTRAB	Fluorescent treponema Ab,CSF Place Of Service:LabCorp	
UHEALTH Labs	FLUOX	FLUOXETINE LEVEL Container:PLASTIC SCREW-CAP VIAL Temperature:SHIP REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	FLUP	FLUPHENAZINE,SERUM/PLASMA Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	FOL	FOLATE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82746
UHEALTH Labs	FLRBC	FOLATE, RBC Place Of Service:LABCORP	
UHEALTH Labs	FSH	FOLLICLE STIMULATING HORMONE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83001
UHEALTH Labs	FOOD	FOOD ALLERGY PROFILE Place Of Service:LABCORP	
UHEALTH Labs	FOODRX	FOOD PROF W/COMPONENT RFLX Place Of Service:LABCORP	
UHEALTH Labs	FORM	FORMIC ACID Container:Specimen Container: PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	FXSD	Fragile X Syndrome, Diagnostic Place Of Service:LABCORP	
UHEALTH Labs	FRNTUL	FRANCISELLA TULARENSIS ABS Place Of Service:LABCORP	
UHEALTH Labs	INSTF	FREE AND TOTAL INSULIN Place Of Service:LabCorp	
UHEALTH Labs	FKLCHN	FREE K+L LT CHAINS, QN, S Place Of Service:LabCorp	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	FREKLT	FREE K+L LT CHAINS, QN, UR Place Of Service:LABCORP	
UHEALTH Labs	FREEKL	FREE KAPPA/LAMBDA LIGHT CHAINS QNT Container:Goldtop, Red Temperature:Refrigerated Place Of Service:JMH	
UHEALTH Labs	KLFREE	FREE KAPPA/LAMBDA LIGHT CHAINS QNT Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86334
UHEALTH Labs	FREET3	FREE T3 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84481
UHEALTH Labs	FREET4	FREE T4 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84439
UHEALTH Labs	FT4DMS	FREE T4 BY DIALYSIS/MASS SPEC Place Of Service:LABCORP	
UHEALTH Labs	FRUCTO	FRUCTOSAMINE Place Of Service:LABCORP	
UHEALTH Labs	FSHLH	FSH AND LH Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83001 83002
UHEALTH Labs	FSHKID	FSH, PEDIATRIC Place Of Service:LABCORP	
UHEALTH Labs	FTAABS	FTA-ABS Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	FNGBAL	FUNGITELL (BAL) Container:sterile, screw top tube. Temperature:Frozen Place Of Service:VIRACOR	
UHEALTH Labs	YFNGIT	FUNGITELL(SERUM) Container:sterile, screw top tube. Temperature:Specimen can be stored refrigerated at 2 to 8°C or preferably stored frozen and shipped for overnight delivery. Place Of Service:VIRACOR	
UHEALTH Labs	F7INHB	FVII INHIBITOR-HUMAN Container:Plastic screw-cap vials (x3) Temperature:Plastic screw-cap vials (x3)-Frozen Place Of Service:UHEALTH LABS	85240 85335
UHEALTH Labs	F8INHB	FVIII INHIB HUMAN Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85335 85240 85732 x3 85730
UHEALTH Labs	FXIIP	FXII ACTIVITY AND INHIBITOR PANEL Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85280
UHEALTH Labs	FYAB	Fya/Fyb RBC Antigen Typing Place Of Service:LabCorp	
UHEALTH Labs	KNTKE	G008-IGE BLUEGRASS, KENTUCKY Place Of Service:LABCORP	
UHEALTH Labs	JHNGRS	G010-IGE JOHNSON GRASS Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	LG6PD	G-6-PD, QUANT, BLOOD AND HGB Place Of Service:LABCORP	
UHEALTH Labs	G6PDQ	G-6-PD, QUANT, BLOOD AND RBC Place Of Service:LABCORP	
UHEALTH Labs	GABPNT	GABAPENTIN (NEURONTIN), SERUM Place Of Service:LABCORP	
UHEALTH Labs	GAD65L	GAD-65 AUTOANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	GADIIA	GAD65, IA-2, and Insulin Autoantibody Container:Transport tube Temperature:Room temperature Place Of Service:Quest	86337 86341 x2
UHEALTH Labs	GALAC	Galactose-alpha-1,3-galactose (Alpha-Gal) IgE Container:Gold or Red Temperature:Room Temperature Place Of Service:Quest	86003
UHEALTH Labs	GALEC3	Galectin-3 Place Of Service:LabCorp	
UHEALTH Labs	GM1	GANGLIOSIDE ANTIBODY PANEL 6 Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	GGM1	GANGLIOSIDE GM 1 ABS IGG, IGM Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	GARLIC	Garlic, IgE, F047 Place Of Service:LabCorp	
UHEALTH Labs	GASTNB	GASTRIN, SERUM Place Of Service:LABCORP	
UHEALTH Labs	GCNAAP	GC NAA, PHARYNGEAL Place Of Service:LABCORP	
UHEALTH Labs	GELPOR	Gelatin Porcine IgE Place Of Service:LabCorp	
UHEALTH Labs	GENP	GeneSeq PLUS, HEXA Place Of Service:LABCORP	
UHEALTH Labs	GENCUL	GENITAL CULTURE WITH GRAM STAIN Container:CULTURETTE SWAB Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	GNOSAR	GENOSURE ARCHIVE (SM) Place Of Service:LABCORP	
UHEALTH Labs	GSRMGB	GENOSURE(R) MG Place Of Service:LABCORP	
UHEALTH Labs	GENOQN	GENOTRACE BCR-ABL1 BY QUANTITATIVE RT-PCR Place Of Service:GENOPT	
UHEALTH Labs	GENTP	GENTAMICIN PEAK LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERRED), GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80170
UHEALTH Labs	GENT	GENTAMICIN RANDOM LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERRED), GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80170
UHEALTH Labs	GENTT	GENTAMICIN TROUGH LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERRED), GoldTop / Tiger top (Serum Separator Tube) (SST®)	80170

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	
UHEALTH Labs	GGTP	GGTP Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82977
UHEALTH Labs	GIPR	GI Profile, Stool, PCR Place Of Service:LabCorp	
UHEALTH Labs	GIARD	GIARDIA LAMBLIA AB, IFA Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	GIACRY	Giardia lamblia, Direct Detection EIA and Cryptosporidium, Direct Detection EIA Place Of Service:LABCORP	
UHEALTH Labs	GLASEQ	GLA Sequencing Place Of Service:LabCorp	
UHEALTH Labs	GLIGA	GLIADIN IGG/IGA AB PROF, EIA Place Of Service:LABCORP	
UHEALTH Labs	GTT2	GLU TOLERANCE LIMIT 2HR Container:3-GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	
UHEALTH Labs	GTT4	GLU TOLERANCE LIMIT 4HR Container:5-GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	
UHEALTH Labs	GLUCA	GLUCAGON, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	GLU1PP	GLUCOSE 1HR POST PRANDIAL Container:Grey Top Tube fluoride/oxalate Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	GLU2PP	GLUCOSE 2HR POST PRANDIAL Container:GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82947 82950
UHEALTH Labs	GLUPP	GLUCOSE 2HR POST PRANDIAL Container:Fluoride/oxalate (gray-top) tubes Temperature:Room temperature Place Of Service:UHEALTH LABS	82947 x2
UHEALTH Labs	GLUF	GLUCOSE FASTING Container:GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82947
UHEALTH Labs	GLU	GLUCOSE LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82951
UHEALTH Labs	GLUTF	GLUCOSE TOL. FASTING Container:Grey Top Tube fluoride/oxalate Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	GTT1	GLUCOSE TOLERANCE 1HR Container:Grey Top Tube fluoride/oxalate Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	GTT3	GLUCOSE TOLERANCE 3HR Container:4-GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	GTT5	GLUCOSE TOLERANCE 5HR Container:Grey Top Tube fluoride/oxalate Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	GTT6	GLUCOSE TOLERANCE 6HR Container:Grey Top Tube fluoride/oxalate Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	GTT2S	Glucose Tolerance Test, 2 Specimens (75g) Container:Fluoride/oxalate (gray-top) tubes Temperature:Room Temperature Place Of Service:UHEALTH LABS	82947 82950
UHEALTH Labs	GCCUL	GONOCOCCAL CULTURE Container:CULTURETTE SWAB Temperature:STORE AT RT, DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	GOFTAL	GOOSE FEATHERS (E70) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	GSO	GRAM STAIN ONLY Container:STERILE CUP OR SWAB Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87205
UHEALTH Labs	GRAPES	Grape IgE, F529 Place Of Service:LabCorp	
UHEALTH Labs	STRPA	GROUP A STREPTOCOCCUS DETECTION, NAA Place Of Service:LABCORP	
UHEALTH Labs	GBSTRP	GROUP B STREPTOCOCCUS CULTURE Container:CULTURETTE SWAB Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87081
UHEALTH Labs	GRHRAB	GROWTH HORMONE ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	GHSER	GROWTH HORMONE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	GRHRSR	GROWTH HORMONE, SERUM (5 SPEC) Place Of Service:LABCORP	
UHEALTH Labs	G360	Guardant 360 Container:KIT - BLOOD Place Of Service:UHEALTH LABS	
UHEALTH Labs	HPYBT	H. PYLORI BREATH TEST Place Of Service:LABCORP	
UHEALTH Labs	DSTGRE	H001-IGE HOUSE DUST (GREER) Place Of Service:LABCORP	
UHEALTH Labs	DSTHOL	H002-IGE HOUSE DUST HOLLISTER Place Of Service:LABCORP	
UHEALTH Labs	FLUBIG	HAEMOPHILUS INFLUENZAE B IGG Place Of Service:LABCORP	
UHEALTH Labs	HALIBF	Halibut IgE Ab, Serum Place Of Service:LabCorp	
UHEALTH Labs	HNTVAB	HANTAVIRUS ANTIBODIES, ELISA Place Of Service:LABCORP	
UHEALTH Labs	HPT	HAPTOGLOBIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83010
UHEALTH Labs	HZNIGE	HAZEL NUT TREE (T4) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HBVDP	HBV DNA QNT, RT-PCR W/RFLX Drug Resistance and PreCore Mut Container:PEARL Temperature:Frozen Place Of Service:UHEALTH LABS	87517
UHEALTH Labs	HBVRG	HBV DNA QNT, RT-PCR W/RFLX GENOTYPE Container:PEARL Temperature:Frozen Place Of Service:UHEALTH LABS	87517
UHEALTH Labs	HBVGDP	HBV DNA QNT,RT-PCR W/RFLX Genotype,Drug Resistance and PreCore Mut Container:PEARL Temperature:Frozen Place Of Service:UHEALTH LABS	87517
UHEALTH Labs	HBPCR	HBV DNA QUANT (REAL TIME PCR) Container:LAVENDER (EDTA) Temperature:FROZEN Place Of Service:UHealth Labs	87517
UHEALTH Labs	HBVDRU	HBV DRUG RESISTANCE Container:GoldTop;Lavendr;PEARL;RED Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	HBVGEN	HBV GENOTYPE Container:GoldTop;Lavendr;PEARL;RED Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	HBVPRE	HBV PRECORE Container:GoldTop;Lavendr;PEARL;RED Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	HCGCSF	hCG, Beta, CSF Container:sterile, plastic leak-proof container Temperature:Room temperature Place Of Service:Quest	84702
UHEALTH Labs	HCGT	HCG, TOTAL, CSF Container:STERILE TRANSPORT TUBE Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	FIBROS	HCV FIBROSURE Place Of Service:LABCORP	
UHEALTH Labs	HCVGS	HCV GENOSURE(R) NS3/4A Place Of Service:LABCORP	
UHEALTH Labs	HCVGNS	HCV GENOTYPE REFLEX NS3/RA Place Of Service:LABCORP	
UHEALTH Labs	HCVGEN	HCV GENOTYPING (PCR/LIPA) Container:LAVENDER TOP TUBE (EDTA) PREFERED, GoldTop TUBE, PEARL TOP TUBE, RED TOP TUBE Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	87902
UHEALTH Labs	HCVRFX	HCV RNA BY PCR, QN RFX GENO Place Of Service:LABCORP	
UHEALTH Labs	HCVRD	HCV RNA DIAGNOSIS, NAA Place Of Service:LABCORP	
UHEALTH Labs	HCVPCR	HCV RNA QUANT (REAL TIME PCR) Container:PEARL TOP OR SST Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86803
UHEALTH Labs	HCVRGN	HCV RNA, QN PCR W/REFL TO GENOTYPE, LIPA(R) Container:LAVENDER TOP TUBE (EDTA) PREFERED GoldTop TUBE, PEARL TOP TUBE Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86803

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HE4OVA	HE4, OVARIAN CANCER MONITORING Container:Transport tube Temperature:Frozen Place Of Service:QUEST	
UHEALTH Labs	HMC24U	HEAVY METALS 24H UR WITH CADMIUM Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HMCPU	Heavy Metals Comprehensive Panel, Random Urine Container:Acid-washed plastic urine container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HMPU	HEAVY METALS PANEL, 24 HR URINE Container:ACID-WASHED PLASTIC URINE CONTAINER Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HMPIBL	HEAVY METALS PROFILE I, BLOOD Place Of Service:LABCORP	
UHEALTH Labs	HVYMT2	HEAVY METALS PROFILE II, URINE Place Of Service:LABCORP	
UHEALTH Labs	HVYMET	HEAVY METALS PROFILE, URINE Place Of Service:LABCORP	
UHEALTH Labs	HPYGST	HELICOBACTER PYLORI ANTIGEN, STOOL Container:PARA PAK STERILE CONTAINER Temperature:STORE REFRIGERATED AT 2-8 C IF DELAY TRANSPORT ANTICIPATED DELIVER IMMEDIATELY TO MICROBIOLOGY Place Of Service:UHealth Labs	87338
UHEALTH Labs	HPYL	Helicobacter pylori, Culture with Reflex to Susceptibility Container:Brucella broth or equivalent with 10-20% glycerol Temperature:Frozen -70° C Place Of Service:QUEST	
UHEALTH Labs	HTLCD4	HELPER T-LYMPH-CD4 Place Of Service:LABCORP	
UHEALTH Labs	HCT	HEMATOCRIT Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85018
UHEALTH Labs	HGB	HEMOGLOBIN Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85014
UHEALTH Labs	GLYCOH	HEMOGLOBIN A1C Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83036
UHEALTH Labs	HH	HEMOGLOBIN AND HEMATOCRIT Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85014 85018
UHEALTH Labs	HGELPI	HEMOGLOBIN ELECTROPHORESIS QUANT,W/INTERPRETATION Container:LAVENDER TOP Temperature:ROOM TEMP Place Of Service:UHEALTH LAB	83020
UHEALTH Labs	HGELP	HEMOGLOBIN ELECTROPHORESIS QUANTITATION Container:LAVENDER TOP Temperature:ROOM TEMP Place Of Service:UHEALTH LAB	83020
UHEALTH Labs	HGBFRE	HEMOGLOBIN, FREE, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	HGBP	HEMOGLOBINOPATHY EVALUATION Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:QUEST	
UHEALTH Labs	HMOFER	HEMOGLOBATHY+FER W/A THAL RFX Place Of Service:LABCORP	
UHEALTH Labs	HLGSP	Hemophagocytic Lymphohistiocytosis Gene Sequencing Pnl Container:EDTA (lavender top) for whole blood. Saliva kit for saliva. Temperature:Room temperature Place Of Service:UHEALTH LAB	81443
UHEALTH Labs	HMOFEC	HEMOQUANT, FECES FECAL HGB Place Of Service:LABCORP	
UHEALTH Labs	HEMOUR	HEMOSIDERIN, URINE Place Of Service:LABCORP	
UHEALTH Labs	HEPBCT	HEP B CORE AB, TOT Place Of Service:LABCORP	
UHEALTH Labs	HBCRBM	HEP B CORE ANTIBODY TOTAL, WITH RFLX TO IGM Container:Gold top Temperature:Room Temp Place Of Service:UHealth Labs	
UHEALTH Labs	HBCRM	HEP B CORE ANTIBODY TOTAL, WITH RFLX TO IGM Container:Serum Separator Tube (SST) Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	HBAGCF	HEP B SURFACE Ag,CONFIRMATION Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room Temperature Place Of Service:UHealth Labs	
UHEALTH Labs	HEPPLT	HEPARIN INDUCED PLATELET Place Of Service:LABCORP	
UHEALTH Labs	HITPP	HEPARIN INDUCED THROMBOCYTOPENIA Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	HFP	HEPATIC FUNCTION PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80076
UHEALTH Labs	HAABM	HEPATITIS A AB IGM Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86709
UHEALTH Labs	HAABT	HEPATITIS A AB TOTAL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86708
UHEALTH Labs	HAABTR	HEPATITIS A ANTIBODY, TOTAL W RFLX Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86708
UHEALTH Labs	HAVAB	HEPATITIS A VIRUS (HAV) AB, TOTAL WITH RFLX TO IGM Place Of Service:LABCORP	
UHEALTH Labs	HBCT	HEPATITIS B CORE Ab TOTAL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86704
UHEALTH Labs	HBCMAB	HEPATITIS B CORE IGM ANTIBODY Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86705
UHEALTH Labs	HPBVG	HEPATITIS B GENOTYPE BCP/PRECORE Container:PPT potassium EDTA (white-top) tube Temperature:FROZEN Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HBVPAN	HEPATITIS B IMMUNITY PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	86317
UHEALTH Labs	HBSBB	HEPATITIS B SURFACE Ab Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86317
UHEALTH Labs	HBAG	HEPATITIS B SURFACE Ag Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	87340
UHEALTH Labs	HEPBSA	HEPATITIS B SURFACE ANTIBODY, QUALITATIVE Place Of Service:LABCORP	
UHEALTH Labs	HBVDR	Hepatitis B Virus (HBV) Genotyping Plus Drug Resistance Place Of Service:LabCorp	
UHEALTH Labs	HBVSD	HEPATITIS B VIRUS (HBV) SCREENING AND DIAGNOSIS Place Of Service:LABCORP	
UHEALTH Labs	HBPRVI	Hepatitis B Virus (Profile VI) Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	
UHEALTH Labs	HBVW6	HEPATITIS B VIRUS DNA, QN PCR W/REFL HBV GENOTYPE Container:Plastic screw-cap vial(s) Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	HEPBB	Hepatitis Be Antibody Place Of Service:LabCorp	
UHEALTH Labs	HEPBG	Hepatitis Be Antigen Place Of Service:LabCorp	
UHEALTH Labs	HBEPAN	Hepatitis Be Panel Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	86707
UHEALTH Labs	HEPC	HEPATITIS C ANTIBODY SCREEN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86803
UHEALTH Labs	HCVRQN	HEPATITIS C ANTIBODY W/RFLX HCV RNA QNT, RT PCR Container:Serum Separator Tube (SST) Temperature:Room Temperature Place Of Service:Uhealth Labs	86803
UHEALTH Labs	HRNAQL	HEPATITIS C VIRAL RNA, QUALITATIVE TMA Container:PPT potassium EDTA (white-top) or plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	HCVABR	HEPATITIS C VIRUS (HCV) AB WITH REFLEX TO QN RT PCR Place Of Service:LABCORP	
UHEALTH Labs	HEPAD	HEPATITIS D QRT-PCR (SERUM) Place Of Service:LABCORP	
UHEALTH Labs	HEPD	HEPATITIS D VIRUS (HDV) IGM AB Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HEPEAB	Hepatitis E IgG Antibodies Place Of Service:LABCORP	
UHEALTH Labs	HEPEQN	Hepatitis E Virus by Quantitative PCR Container:Lavender (EDTA), pink (K2EDTA) or Serum separator tube (SST) Temperature:Frozen Place Of Service:ARUP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HEPPAN	HEPATITIS PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	86704
UHEALTH Labs	HERELA	HER2, ELISA Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HER2Q	HER-2/NEU QUANTITATIVE ELISA Place Of Service:LABCORP	
UHEALTH Labs	HHDNA	HERED.HEMOCHROMATOSIS, DNA Place Of Service:LABCORP	
UHEALTH Labs	HSVCT	HERPES SIMPLEX VIR CX RFL TO TYP Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HSV2AB	Herpes Simplex Virus (HSV) Type 2-Specific Antibodies, IgG Place Of Service:LABCORP	
UHEALTH Labs	HSVAB	Herpes Simplex Virus (HSV) Types 1 and 2-Specific Antibodies, IgG Place Of Service:LABCORP	
UHEALTH Labs	HSVIGG	Herpes Simplex Virus 1/2 (IgG) Type Antibody, CSF Container:Transport tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	HSVIGM	Herpes Simplex Virus 1/2 Antibody (IgM), IFA with Reflex to Titer, CSF Place Of Service:QUEST	
UHEALTH Labs	HSV7	HERPES VIRUS 7 DNA, QN RT PCR Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HEXPP	HEXAGONAL PHASE PHOSPHOLIPID Place Of Service:LABCORP	
UHEALTH Labs	HGBA2Q	HGB A2, QUANT Place Of Service:LABCORP	
UHEALTH Labs	HGBFLQ	HGB,FETAL,QUANT Place Of Service:LABCORP	
UHEALTH Labs	HHV6PS	HHV-6 DNA, QUANT RT PCR (PLASMA) Place Of Service:VIRACOR	
UHEALTH Labs	HHV6AB	HHV-6 IGG/M ANTIBODIES, QUANT Container:Red-top tube, gel-barrier tube, or serum transfer tube Temperature:Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	HHV6CF	HHV-6 QN Real-time PCR, CSF Container:sterile, screw top tube Temperature:Frozen Place Of Service:Viracor	87533
UHEALTH Labs	HHV6DA	HHV-6 QUANT DNA PCR Place Of Service:LABCORP	
UHEALTH Labs	HHV8WB	HHV8 QPCR (WHOLE BLOOD) Place Of Service:LABCORP	
UHEALTH Labs	YHHV8Q	HHV-8DNA,QUANT RT PCR(WB) Place Of Service:VIRACOR	
UHEALTH Labs	HICKAL	HICKORY/PECAN TREE (T22) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HDL	HIGH DENSITY LIPOPROTEIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83695
UHEALTH Labs	HISDEU	HISTAMINE DETERMINATION, URINE Place Of Service:LABCORP	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HRCU	Histamine Release (Chronic Urticaria) Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HIST	HISTAMINE, PLASMA Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	HISTPL	HISTAMINE, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	HISAG	HISTOPLASMA AG, QN, EIA, SERUM Container:Collect serum specimens in serum separator or red top tube. Temperature:Refrigerated/Frozen/Room Temperature Place Of Service:UHEALTH LABS	87385
UHEALTH Labs	HISAGU	HISTOPLASMA AG, QN, EIA, URINE Container:sterile screw cap container Temperature:Room Temperature Place Of Service:UHEALTH LABS	87385
UHEALTH Labs	HISIG	Histoplasma Antibody IgG, IgM EIA Container:Serum separated from the clot in a leak-proof container. Temperature:Frozen Place Of Service:MiraVista	86635 x2
UHEALTH Labs	HISAB	Histoplasma Antibody,Immunodiffusion Container:Serum separated from the clot in a leak-proof container. Temperature:Frozen Place Of Service:MiraVista	86698
UHEALTH Labs	HISTUR	HISTOPLASMA ANTIGEN, URINE Place Of Service:LABCORP	
UHEALTH Labs	HIVDNA	HIV 1 DNA, QUALITATIVE PCR Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HIV1QN	HIV 1 RNA, QUANTITATIVE RT PCR Container:PEARL TOP TUBE Temperature:ROOM TEMPERATURE Place Of Service:UHealth Labs	87536
UHEALTH Labs	HIVCON	HIV 1/2 AB DIFFERENTIATION W/REFLEX TO HIV1 RNA QL Place Of Service:LabCorp	
UHEALTH Labs	HIVSCR	HIV SCREEN 4TH GENERATION W/REFLEX Container:Gold top SST Temperature:Room Temp Place Of Service:UHEALTH LABS	87389
UHEALTH Labs	HIVRGE	HIV-1 RNA, Quantitative, PCR with Reflex to Genotype Container:Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	RNAPCR	HIV1,RNA PCR (NonGraph) rfx/Geno Place Of Service:LABCORP	
UHEALTH Labs	HLAAD	HLA A 02:01 Determination Container:Preferred: EDTA (lavender-top) tube *** EDTA (lavender-top) tube cannot be used for other testing prior to HLA testing *** Alternative: ACD (yellow-top) tube 4 buccal swabs Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	LHLABC	HLA A,B,C (IR) Place Of Service:LABCORP	
UHEALTH Labs	HLADRB	HLA A,B,C,DRB1,3,4,5,DQB1 (IR) Place Of Service:LABCORP	
UHEALTH Labs	HLA27B	HLA B 27 DISEASE ASSOCIATION Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HLABLR	HLA B LOW RESOLUTION Container:N EDTA (LAVENDER-TOP) TUB Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HLA27D	HLA B27 DETERMINATION Container:N EDTA (LAVENDER-TOP) TUB Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	B5701	HLA B5701 TEST Place Of Service:LABCORP	
UHEALTH Labs	HLAABS	HLA Class 1,2 Antibody Screen Place Of Service:LABCORP	
UHEALTH Labs	HLADQB	HLA DQB High Resolution Place Of Service:LabCorp	
UHEALTH Labs	HLABIR	HLA-B (IR) Place Of Service:LABCORP	
UHEALTH Labs	HLAB51	HLA-B51 Determination Container:EDTA (lavender-top) tube Temperature:Room temperature Place Of Service:QUEST	81374
UHEALTH Labs	HLALR	HLA-DR/DQ Low Resolution Typing Container:N EDTA (LAVENDER-TOP) TUB Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HOMMET	HOMOCY+METHYL Place Of Service:LABCORP	
UHEALTH Labs	HOMOB	HOMOCYST(E)INE, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	UHVA	HOMOVANILLIC ACID, 24HR URINE Container:PLASTIC URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HOMOUR	HOMOVANILLIC ACID, RANDOM UR Place Of Service:LABCORP	
UHEALTH Labs	HRSDAL	HORSE DANDER (E3) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HPVHRR	HPV DNA, High Risk with Reflex to Genotypes 16,18 Container:PreservCyt or SurePath liquid-based pap smear collection vial Cell pellet with 2 mL SurePath fluid in a 15 mL conical centrifuge tube. Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HPVGN	HPV GENOTYPES 16 AND 18 Container:PreservCyt, liquid-based pap smear collection vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HPVG45	HPV Genotypes 16, 18/45 Container:APTIMA transport tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HSTDTU	HPV, 16/18,45 Place Of Service:LABCORP	
UHEALTH Labs	HSVHS	HSV 1/2 IGG, W/REFL HSV2 INHIB Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HSV12R	HSV 1/2 PCR Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	HSVG2R	HSV 2 (IgG),HerpeSelect w/ Rflx to HSV-2 Inhibition Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HSVAR	HSV Acyclovir Resistance Place Of Service:LabCorp	
UHEALTH Labs	HSVCUL	HSV CULTURE AND TYPING Place Of Service:LABCORP	
UHEALTH Labs	HSVNT	HSV CULTURE WITHOUT TYPING Place Of Service:LABCORP	
UHEALTH Labs	HSV12B	HSV TYPE 1&2 RT QUANT (BAL) Container:STERILE SCREW-CAP TUBE Temperature:ROOM TEMPERATURE Place Of Service:VIRACOR	
UHEALTH Labs	HSV1IG	HSV TYPE 1-SPECIFIC AB, IGG Place Of Service:LABCORP	
UHEALTH Labs	HSV12	HSV, TYPE 1 AND 2 DNA, QN RT PCR Container:STERILE SCREW-CAP CONTAINER Temperature:REFRIGERATED: WHOLE BLOOD; FROZEN: ALL OTHER SPECIMEN TYPES Place Of Service:QUEST	
UHEALTH Labs	HTLVPC	HTLV 1/2 QUAL RT PCR-WHOLE BLD Place Of Service:LABCORP	
UHEALTH Labs	HTLV2B	HTLV-I/II ANTIBODIES, QUAL Place Of Service:LABCORP	
UHEALTH Labs	HTLVIM	HTLV-I/II IMMUNOBLOT Place Of Service:LABCORP	
UHEALTH Labs	HTLVCF	HTLV-I-II Abs, CSF Place Of Service:LabCorp	
UHEALTH Labs	HUMGLP	HUMAN ANTI GIARDIA LAMBLIA PNL Place Of Service:LABCORP	
UHEALTH Labs	HUMSE	HUMAN ANTI-MOUSE ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	HUMEPI	HUMAN EPIDIDYMIS PROT 4, SERIAL Place Of Service:LABCORP	
UHEALTH Labs	HGE	Human Granulocytic Ehrlichiosis (HGE), IgG, IgM, IFA Place Of Service:LabCorp	
UHEALTH Labs	HIVIGA	Human Immunodeficiency Virus 1 (HIV-1) GenoSure Archive® Plus Trofile® DNA Place Of Service:LABCORP	
UHEALTH Labs	HPVHY	HUMAN PAPILLOMA VIRUS Container:DIGENE CERVICAL SAMPLER (BRUSH) OR CYTYC THINPREP SOLUTION OR DIGENE SPECIMEN COLLECTION KIT Temperature:CERVICAL BRUSH OR BROOM IN DIGENE SPECIMEN: ROOM TEMPERATURE CERVICAL BIOPSY-FROZEN CYTYC PRESERVCYT -THINPREP ROOM TEMPERATURE CYTORICH FRACTION SUREPATH -ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HDMA	HUNTINGTON DISEASE MUTATION ANALYSIS Container:EDTA (lavender-top) or ACD solution B (yellow-top) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HYALA	HYALURONIC ACID Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	HYDPLN	HYDROXYPROLINE, LC/MS, PLASMA Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	HYSPP	Hypersensitivity Pneumonitis Profile Place Of Service:LabCorp	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CKROCH	I006-IGE COCKROACH, GERMAN Place Of Service:LABCORP	
UHEALTH Labs	IA2AB	IA-2 AUTOANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	IBDXP	IBD EXPANDED PANEL Place Of Service:LABCORP	
UHEALTH Labs	IDH12M	IDH1 & IDH2 MUTATION ANL,EXON4 Temperature:BONE MARROW: REFRIGERATED FFPE TUMOR TISSUE: ROOM TEMPERATURE OR REFRIGERATED. SHIP IN COOLED CONTAINERS DURING SUMMER MONTHS. Place Of Service:ARUP	
UHEALTH Labs	UKAPPB	IFE U+KAPPA/LAMBDA QN U Place Of Service:LABCORP	
UHEALTH Labs	IFE24U	IFE+PROTEIN ELECTRO, 24-HR UR Place Of Service:LABCORP	
UHEALTH Labs	SIGLIT	Ig LIGHT CHAINS, TOTAL IN SERUM Place Of Service:JMH	
UHEALTH Labs	IGA	IGA Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82787
UHEALTH Labs	IGEAB	IGE ANTIBODY (ANTI IGE IGG) Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	EGGWHP	IGE EGG WHITE COMPONENT PROF Place Of Service:LABCORP	
UHEALTH Labs	MILKCP	IGE MILK COMPONENT PROFILE Place Of Service:LABCORP	
UHEALTH Labs	PEANCP	IGE PEANUT COMPONENT PROFILE Place Of Service:LABCORP	
UHEALTH Labs	IGE	IgE, TOTAL, SERUM Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IGF1	IGF-1 Place Of Service:LABCORP	
UHEALTH Labs	IGFBP1	IGF-BP1 Place Of Service:LABCORP	
UHEALTH Labs	IGF2BP	IGFBP-2 Place Of Service:LABCORP	
UHEALTH Labs	IGFB3P	IGF-BP3 Place Of Service:LABCORP	
UHEALTH Labs	IGG	IGG Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82784
UHEALTH Labs	IGGSR	IgG SYNTHESIS RATE Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	INMIN	IgG, IgA, Indirect Immunofluorescence Container:Serum Separator (red/black tiger-top) Tube (SST) Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	IGG3SB	IGG, SUBCLASS 3 Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	IGG4SB	IGG, SUBCLASS 4 Place Of Service:LABCORP	
UHEALTH Labs	IGGSU4	IGG, SUBCLASSES(1-4) Place Of Service:LABCORP	
UHEALTH Labs	IMMFLU	IGG,IGA INDIRECT IMMUNOFLUORESCEN Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	IGG4F1	IGG4 FOOD PANEL I Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IGG4F8	IGG4 FOOD PANEL VIII Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IGM	IGM Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82784
UHEALTH Labs	IGVHHA	IgVH Hypermutation Analysis Container:EDTA (purple-top) tube Temperature:Refrigerated (cold packs) Place Of Service:Genoptix	81263
UHEALTH Labs	IL1BS	IL-1b (serum)* Place Of Service:LABCORP	
UHEALTH Labs	IMIP	IMIPRAMINE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IMMKNW	IMMUKNOW (R) Place Of Service:LABCORP	
UHEALTH Labs	YIMMUK	ImmuKnow Immune Cell Function (WB) Container:sodium heparin (green top) tube Temperature:ROOM TEMPERATURE Place Of Service:VIRACOR	
UHEALTH Labs	IMMCF	IMMUNE CELL FUNCTION Container:Sodium heparin (green-top) tube Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IFFE	IMMUNOFIXATION ELECTROPHORESIS Container:SST Temperature:ROOM TEMPERATURE Place Of Service:Uhealth Labs	
UHEALTH Labs	IFE	IMMUNOFIXATION ELECTROPHORESIS Container:Transport tube Temperature:Refrigerate Place Of Service:UHealth Labs	
UHEALTH Labs	IMMUPR	IMMUNOFIXATION IGA,IGG,IGM QT Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IMMUN	IMMUNOFIXATION, CSF Container:STERILE PLASTIC LEAK-PROOF CONTAINER Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	IMMUR	IMMUNOFIXATION, URINE Container:PLASTIC URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	IMMOF	IMMUNOFIXATION, URINE Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	IMMDQS	IMMUNOGLOBULIN D, QUANT Place Of Service:LABCORP	
UHEALTH Labs	IMMGET	IMMUNOGLOBULIN E, TOTAL Place Of Service:LABCORP	
UHEALTH Labs	IGGIND	Immunoglobulin G Index Place Of Service:LabCorp	
UHEALTH Labs	IGGS1S	IMMUNOGLOBULIN G, SUBCLASS 1 Place Of Service:LABCORP	
UHEALTH Labs	IGGS2S	IMMUNOGLOBULIN G, SUBCLASS 2 Place Of Service:LABCORP	
UHEALTH Labs	IMGCSF	IMMUNOGLOBULIN G, SYN RATE, CSF Place Of Service:LABCORP	
UHEALTH Labs	IMQ	IMMUNOGLOBULINS QUANTITATION Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82784
UHEALTH Labs	IMT	IMMUNOTYPING Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86334
UHEALTH Labs	MSPT	INFECTIOUS MONONUCLEOSIS Container:RED TOP TUBE (NO SERUM SEPARATOR) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86308
UHEALTH Labs	INFXAB	Infliximab+Ab (Serial Monitor) Place Of Service:LabCorp	
UHEALTH Labs	FLUAGS	INFLUENZA A & B ANTIGEN DETECTION Container:BD universal viral transport Temperature:Refrigerate Place Of Service:JMH Virology	
UHEALTH Labs	H1N1	INFLUENZA A H1N1(2009)RTPCR(BL) Place Of Service:VIRACOR	
UHEALTH Labs	CFSP	Inheritest (R) CF / SMA Panel Place Of Service:LABCORP	
UHEALTH Labs	INAULT	INHIBIN A, ULTRASENSITIVE Place Of Service:LABCORP	
UHEALTH Labs	IHBNAB	INHIBIN A+B Place Of Service:LABCORP	
UHEALTH Labs	INHBB	INHIBIN B Place Of Service:LABCORP	
UHEALTH Labs	INSULN	INSULIN Place Of Service:LABCORP	
UHEALTH Labs	INSABS	INSULIN ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	INTPR	INTELLIGENT FLOW PROFILE Container:sodium heparin (green-top) tube,EDTA is acceptable, RPMI Temperature:Use refrigerated cold pack for transport. Place Of Service:UHEALTH LABS	88184 88185 88189
UHEALTH Labs	IL2R	Interleukin 2 Receptor (CD25), Soluble Container:Serum separator tube, plain red, or green (lithium heparin) Temperature:CRITICAL FROZEN Place Of Service:ARUP	83520
UHEALTH Labs	IL10P	Interleukin-10, Plasma Place Of Service:LABCORP	
UHEALTH Labs	IL2RA	Interleukin-2 Soluble Receptor Alpha Container:Red-top tube or gel-barrier tube Temperature:frozen Place Of Service:LABCORP	83520



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	IL2S	Interleukin-2, Serum Place Of Service:LABCORP	
UHEALTH Labs	INTR6B	INTERLEUKIN-6, SERUM Place Of Service:LABCORP	
UHEALTH Labs	IFASER	INTRINSIC FACTOR ABS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	IDN24U	IODINE, 24 HR URINE Place Of Service:LABCORP	
UHEALTH Labs	IODINU	IODINE, RANDOM URINE Place Of Service:LABCORP	
UHEALTH Labs	IODINE	IODINE, SERUM OR PLASMA Place Of Service:LABCORP	
UHEALTH Labs	CAI	IONIZED CALCIUM Container:DARK GREEN TOP TUBE. (Lithium Heparin) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82330
UHEALTH Labs	IRON	IRON Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83540
UHEALTH Labs	IRONBC	IRON AND TOTAL IRON BINDING CAPACITY Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83550
UHEALTH Labs	IRONPR	IRON PROFILE Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	ISAV	Isavuconazole (CRESEMBA) LC-MS/MS Container:Red top tube Temperature:Frozen Place Of Service:Viracor	80299
UHEALTH Labs	ITRMTB	Itraconazole and Mtb s/p Place Of Service:LABCORP	
UHEALTH Labs	JAK2QN	JAK2 MUTATION (V617F) ANALYSIS,QN Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	JAK2MT	JAK2 MUTATION ANALYSIS, QUAL Place Of Service:LABCORP	
UHEALTH Labs	JAK2QB	JAK2 V617F, QUANTITATIVE Place Of Service:LABCORP	
UHEALTH Labs	JAK2WR	JAK2 V617F,QL,LEUMETA(R) W/REFL EXONS 12 AND 13 Container:EDTA (lavender-top) Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	JCPV	JC POLYOMA VIRUS DNA,Qual RT PCR Container:EDTA (lavender-top) Temperature:SHIP FROZEN      SHIPPING ROOM TEMP UNACCEPTABLE; REFRIGERATED ACCEPTABLE Place Of Service:QUEST	SHIPPING
UHEALTH Labs	JCVCF	JC Virus (JCV) Qn Real-time PCR, CSF Container:Sterile, screw top tube Temperature:Store frozen and ship on dry ice for overnight delivery. Place Of Service:Viracor	87799
UHEALTH Labs	JCPV	JC Virus (JCV) Qn Real-time PCR, Plasma Container:EDTA or ACD tube Temperature:ambient or frozen Place Of Service:Viracor	87799

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	JCVWB	JC VIRUS DNA, PCR (WHOLE BLOOD) Place Of Service:LABCORP	
UHEALTH Labs	JO1AB	JO-1 ANTIBODY Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	JKBG	JUNE (KENTUCKY BLUE) GRASS IgE Container:Serum Separator Tube (SST®) Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	LATEXE	K082-IGE LATEX Place Of Service:LABCORP	
UHEALTH Labs	LATXE	K082-IGG LATEX Place Of Service:LABCORP	
UHEALTH Labs	SNFLWE	K084-IGE SUNFLOWER SEED Place Of Service:LABCORP	
UHEALTH Labs	KAPLLC	KAPPA/LAMBDA LIGHT CHAINS TOTAL, SERUM Place Of Service:QUEST	
UHEALTH Labs	KPLRU	Kappa/Lambda Light Chains, Total, Random Urine Container:Plastics leak-proof container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	D816V	KIT (D816V) Digital PCR Place Of Service:LABCORP	
UHEALTH Labs	KITPCR	KIT (D816V) MUTATION BY PCR Place Of Service:ARUP	
UHEALTH Labs	KOH	KOH PREP Container:STERILE CUP Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87220
UHEALTH Labs	KRASB	KRAS EXTENDED ANALYSIS Place Of Service:LABCORP	
UHEALTH Labs	KRASMA	KRAS Mutation Analysis Container:Formalin-fixed, paraffin embedded tissue block Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LEGAB	L. PNEUMOPHILA IGG & IGM ABS 1-6 Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	LACOS	LACOSAMIDE Place Of Service:LABCORP	
UHEALTH Labs	LACISO	LACTATE DEHYDROGENASE ISOENZYMES Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LAC	LACTIC ACID Container:GREY TOP TUBE (Fluoride/Oxalate) Temperature:TRANSPORT ON ICE Place Of Service:UHealth Labs	83605
UHEALTH Labs	LACCSF	Lactic Acid, CSF Place Of Service:LabCorp	
UHEALTH Labs	LCTOFQ	LACTOFERRIN, FECAL, QUANT Place Of Service:LABCORP	
UHEALTH Labs	LACTS	LACTOFERRIN, QN, STOOL Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LACTOL	LACTOSE TOLERANCE TEST Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	LACT1	LACTOSE TOLERANCE TEST 5 SPEC 50G Container:Fluoride/oxalate (gray-top) tubes Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	LAMOT	LAMOTRIGINE (LAMICTAL), SERUM Place Of Service:LABCORP	
UHEALTH Labs	LDISOB	LD ISOENZYMES Place Of Service:LABCORP	
UHEALTH Labs	LDSF	LD, SYNOVIAL FLUID Container:Sterile transport tube Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LDH	LDH Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83615
UHEALTH Labs	LDLPN	LDL-P Place Of Service:LabCorp	
UHEALTH Labs	ULEAD	LEAD, 24HR URINE Container:ACID-WASHED PLASTIC URINE CONTAINER Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	ALEAD	LEAD, BLOOD (ADULT) Place Of Service:LABCORP	
UHEALTH Labs	LEADPV	Lead, Blood (Pediatric), Venous Place Of Service:LabCorp	
UHEALTH Labs	URNLEG	LEGIONELLA ANTIGEN, URINE Container:Sterile urine container Temperature:Refrigerate Place Of Service:UHEALTH LABS	
UHEALTH Labs	LGPNAB	Legionella Pneumophila Antibody Place Of Service:LABCORP	
UHEALTH Labs	LEGPNE	LEGIONELLA PNEUMOPHILA DFA Place Of Service:LABCORP	
UHEALTH Labs	LEGPM	LEGIONELLA PNEUMOPHILA, IGM Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ULEGAG	LEGIONELLA URINARY ANTIGEN Container:PLASTIC URINE CONTAINER Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	LEISH	LEISHMANIA ANTIBODY, IGG Place Of Service:LABCORP	
UHEALTH Labs	LEMON	Lemon F208-IgE Place Of Service:LABCORP	
UHEALTH Labs	LENTIL	Lentil (F235) IgE Place Of Service:LABCORP	
UHEALTH Labs	LEPTIN	Leptin, Serum or Plasma Place Of Service:LABCORP	
UHEALTH Labs	LEPDNA	LEPTOSPIRA DNA,QUALITATIVE REAL-Time PCR Container:Plastic leak-proof container Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LPTOMB	LEPTOSPIRA IGM Place Of Service:LABCORP	
UHEALTH Labs	LETTU	LETTUCE IgE ALLERGEN Container:Serum Separator Tube (SST)	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LAPQ	LEUCINE AMINOPEPTIDASE Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	LEVKEP	LEVETIRACETAM (KEPPRA), S Place Of Service:LABCORP	
UHEALTH Labs	LIDOS	LIDOCAINE (XYLOCAINE), SERUM Place Of Service:LABCORP	
UHEALTH Labs	LIPASE	LIPASE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83690
UHEALTH Labs	FLIP	LIPASE FLUID Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	LIPID	LIPID PROFILE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80061
UHEALTH Labs	LIPIDD	LIPID PROFILE W/ REFLEX DLDL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80061
UHEALTH Labs	LASA	LIPID-ASSOCIATED SIALIC ACID Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LPOPOA	LIPOPROTEIN (A) Place Of Service:LABCORP	
UHEALTH Labs	LDLLIP	LIPOPROTEIN FRACTIONATIONION MOB Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LPOPHE	LIPOPROTEIN PHENOTYPING Place Of Service:LABCORP	
UHEALTH Labs	LI	LITHIUM LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80178
UHEALTH Labs	LVCTAB	LIVER CYTOSOL AUTOANTIBODIES Container:Transport tube Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	LFPNL	LIVER FIBROSIS PANEL (HEPASCORE) Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LVKM	LIVER-KIDNEY MICROSOMAL AB Place Of Service:LABCORP	
UHEALTH Labs	LOBST	LOBSTER IgE ALLERGEN Place Of Service:QUEST	
UHEALTH Labs	LORQ	LORAZEPAM Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	LRP4	LRP4 Autoantibody Test Container:PLASTIC TRANSPORT TUBE Temperature:REFRIGERATE	86255

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:QUEST	
UHEALTH Labs	LSDSR	LSD (LYSERGIDE), SERUM Place Of Service:LABCORP	
UHEALTH Labs	LSNKP	LS-NK Container:EDTA (LAVENDER-TOP) TUBE Place Of Service:UHEALTH LABS	
UHEALTH Labs	LPSLE	Lupus (SLE) Panel Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	LUPEVL	LUPUS ANTICOAGULANT EVALUATION WITH REFLEX Container:Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LUPACP	Lupus Anticoagulant Panel Place Of Service:LABCORP	
UHEALTH Labs	SLECDP	LUPUS ERYTHEMATOSUS (SLE) PANEL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LUPUSA	Lupus Inhibitor Profile Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85613 85730
UHEALTH Labs	LH	LUTEINIZING HORMONE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83002
UHEALTH Labs	PLH	Lutenizing Hormone, Pediatric Place Of Service:LabCorp	
UHEALTH Labs	LYMIG	LYME DISEASE AB IGG, IMMUNOBLOT Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LYMEQ	LYME DISEASE AB W/REFL TO BLOT Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LYMETA	Lyme Disease Total Antibody W/Reflex to Immunoassay Place Of Service:LABCORP	
UHEALTH Labs	LYMED	Lyme Disease, Borrelia burgdorferi, RT PCR Place Of Service:LABCORP	
UHEALTH Labs	LYMCSF	LYME, WESTERN BLOT, CSF Place Of Service:LABCORP	
UHEALTH Labs	LWBSER	LYME, WESTERN BLOT, SERUM Place Of Service:LABCORP	
UHEALTH Labs	LPA3P	Lymphocyte Proliferation to Anti-CD3/Anti-CD28 & Anti-CD 3/IL-2, Flow Cytometry Container:Green top (sodium heparin) Temperature:Room Temperature Place Of Service:MAYO	86353
UHEALTH Labs	LSUBS	LYMPHOCYTE SUBSET ANALYSIS Container:LAVENDER (EDTA) Top Tube for peripheral blood / Yellow Top tube for body fluid Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86355
UHEALTH Labs	LSUB2	LYMPHOCYTE SUBSET PANEL 2 Container:TS-Lavender Temperature:Room Temperature Place Of Service:UHEALTH LABS	86355 86359 86360
UHEALTH Labs	LSUB4	LYMPHOCYTE SUBSET PANEL 4 Container:LAVENDER (EDTA) Top Tube for peripheral blood / Yellow Top tube for body fluid Temperature:TRANSPORT AT ROOM TEMPERATURE	86360

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:UHealth Labs	
UHEALTH Labs	CSFLCM	LYMPHOCYTIC CHORIOMENIN VIRUS AB Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	LYMMP	LYMPHOID MOLECULAR PROFILE Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube. Bone marrow: 2-3 mL in EDTA (purple-top) tube. Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	LYSO	LYSOZYME, SERUM Place Of Service:LABCORP	
UHEALTH Labs	MTUBC	M.TUBERCULOSIS COMPLEX, PCR, NON RESPIRATORY Container:Sterile, leak-proof container Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	CLADOE	M002-IGE CLADOSPORIUM HERBARUM Place Of Service:LABCORP	
UHEALTH Labs	ASFUME	M003-IGE ASPERGILLUS FUMIGATUS Place Of Service:LABCORP	
UHEALTH Labs	CNALBE	M005-IGE CANDIDA ALBICANS Place Of Service:LABCORP	
UHEALTH Labs	ALTERE	M006-IGE ALTERNARIA ALTERNATA Place Of Service:LABCORP	
UHEALTH Labs	STCHBS	M024-IGE STACHYBOTRYS ATRA Place Of Service:LABCORP	
UHEALTH Labs	ASPERE	M207-IGG ASPERGILLUS NIGER Place Of Service:LABCORP	
UHEALTH Labs	MACIGE	MACADAMIA NUT (RF345) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MACKF	Mackerel IgE Ab, Serum Place Of Service:LabCorp	
UHEALTH Labs	MAGIGG	MAG Antibody, IgG Place Of Service:LABCORP	
UHEALTH Labs	MGSG	MAG SGPG AB (IGM),EIA Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	MAG	MAGNESIUM LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83735
UHEALTH Labs	MAGRBC	MAGNESIUM, RBC Place Of Service:LABCORP	
UHEALTH Labs	MANGBD	MANGANESE, BLOOD Place Of Service:LABCORP	
UHEALTH Labs	MANGAP	MANGANESE, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	MANGO	MANGO FRUIT IgE ALLERGEN Container:Serum Separator Tube (SST®) Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	MBCHCP	Mature B-Cell/Hairy Cell Panel (Add-On) Container:sodium heparin (green-top) tube Temperature:Refrigerated (cold packs) Place Of Service:Genoptix	88184 88185 x7
UHEALTH Labs	MDSFPE	MDS FISH PROFILE (EXTENDED) Container:Peripheral Blood or Bone Marrow: sodium heparin (green-top) tube; Bodily fluids: Bodily fluids	88368 88369 x3 88377 x5



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		for diagnostic testing of hematological cancers (subject to specimen viability). Preferred transport in equal parts RPMI. Lymphoid tissue: Transport in RPMI. FFPE tissue: no decalcification Temperature: Use refrigerated cold pack for transport Place Of Service: GENOPTIX	
UHEALTH Labs	MEASAB	Measles (Rubeola) (IgG, IgM) Antibody Panel, IFA, CSF Container: Plastic screw-cap vial Temperature: Room temperature Place Of Service: QUEST	
UHEALTH Labs	MEASG	MEASLES ANTIBODY (IGG) Container: PLASTIC SCREW-CAP VIAL Temperature: ROOM TEMPERATURE Place Of Service: QUEST	
UHEALTH Labs	MEASM	MEASLES ANTIBODY (IgM) Container: GOLD TOP/RED Temperature: ROOM TEMPERATURE Place Of Service: QUEST	
UHEALTH Labs	MMRUIM	MEASLES/MUMPS/RUBELLA IMMUNITY Place Of Service: LABCORP	
UHEALTH Labs	MELAL	Melaleuca IgE (T021) Place Of Service: LabCorp	
UHEALTH Labs	MEP	MEPERIDINE Container: PLASTIC SCREW-CAP VIAL Temperature: REFRIGERATED (COLD PACKS) Place Of Service: QUEST	
UHEALTH Labs	MEPHP	MEPHENYTOIN AND METABOLITE Container: PLASTIC SCREW-CAP VIAL Temperature: ROOM TEMPERATURE Place Of Service: QUEST	
UHEALTH Labs	MERCB	MERCURY, BLOOD Place Of Service: LABCORP	
UHEALTH Labs	MERCUR	Mercury, Urine Place Of Service: LABCORP	
UHEALTH Labs	MERSNP	MERS (NASOPHARYNGEAL) Place Of Service: DOHM	
UHEALTH Labs	MERSOP	MERS (OROPHARYNGEAL) Place Of Service: DOHM	
UHEALTH Labs	MERSER	MERS (SERUM) Place Of Service: DOHM	
UHEALTH Labs	MERSPU	MERS (SPUTUM) Place Of Service: DOHM	
UHEALTH Labs	MERSTO	MERS (STOOL) Place Of Service: DOHM	
UHEALTH Labs	MESQ	MESORIDAZINE Container: PLASTIC SCREW-CAP VIAL Temperature: ROOM TEMPERATURE Place Of Service: QUEST	
UHEALTH Labs	METSCU	METABOLIC SCREEN, URINE QL Place Of Service: GENE	
UHEALTH Labs	META24	METANEPHRINES, FRAC, QN, 24-HR Place Of Service: LABCORP	
UHEALTH Labs	METANE	METANEPHRINES, FRAC., PL. FREE Place Of Service: LABCORP	
UHEALTH Labs	METPH	METANEPHRINES, PHEOCHROMOCYT Place Of Service: LABCORP	
UHEALTH Labs	MMETAU	METHADONE METABOLITE QUAL, URINE Container: URINE random collection CUP/ Yellow top tube no additive Temperature: Room temperature Place Of Service: UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	MRSAP	METHICILLIN RESIST STAPH AUR PCR Container:COLLECT IN: BBL(TM) CULTURE SWAB(TM) LIQUID AMIES, DOUBLE SWAB OR COPAN VENTURI TRANSYSTEM(TM) LIQUID AMIES, DOUBLE SWAB. Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	METHX	METHOTREXATE Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	MTXSER	METHOTREXATE (MTX), SERUM Place Of Service:LABCORP	
UHEALTH Labs	METHO	METHOTREXATE LEVEL Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	METHYM	METHYLMALONIC ACID AND HOMOCYSTEINE Place Of Service:QUEST	
UHEALTH Labs	METHY	METHYLMALONIC ACID, SERUM Place Of Service:LABCORP	
UHEALTH Labs	RITAL	METHYLPHENIDATE (U) Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	MTHMTB	METHYLPHENIDATE + MTB, URINE Place Of Service:LABCORP	
UHEALTH Labs	MEXI	MEXILETINE, SERUM/PLASMA Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MI2AB	MI-2 AUTOANTIBODIES Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	UMALB	MICROALBUMIN, 24 HR URINE Container:Transport Container: Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MCROST	MICROSPORIDIA STAIN Place Of Service:LABCORP	
UHEALTH Labs	MILLET	Millet IgE (F55) Place Of Service:LABCORP	
UHEALTH Labs	M2AB	MITOCHONDRIAL (M2) ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	MITOAB	MITOCHONDRIAL AB W/REFL TITER Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MITOD	Mitochondrial Disorders (mtDNA Del/Dup Analysis) Container:Lavender Temperature:Room Temperature Place Of Service:CINCINNATI	81465
UHEALTH Labs	MILPP	MITOGEN INDUCED LYMPHOCYTE PROLIF Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MITOT	MITOTANE (LYSODREN) Place Of Service:LABCORP	
UHEALTH Labs	MIXSTD	MIXING STUDY Container:Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	MLLPTD	MLL-PTD MUTATION ANALYSIS Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	MOGS	MOG ANTIBODY RFLX TO TITER,SERUM Container:2 mL serum collected in red-top tube (no gel) Place Of Service:QUEST	86255
UHEALTH Labs	LMOL	MOLECULAR DIAGNOSTICS & FLOW CYTOMETRY REPORT Place Of Service:UHEALTH LABS	
UHEALTH Labs	DRGPRU	MONITOR DRUG PROFILE 10(MW) Place Of Service:LABCORP	
UHEALTH Labs	NEUROP	Motor Neuropathy Complete Antibody Panel Container:Serum: Plastic screw-cap vial; Urine: Plastic urine container Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	MPLMTA	MPL MUTATION ANALYSIS Place Of Service:LABCORP	
UHEALTH Labs	MPNTP	MPN TARGETED PROFILE Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	MRSA	MRSA CULTURE Container:CULTURETTE SWAB Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87081
UHEALTH Labs	MTBRNA	MTB-RIF NAA W Cult, Nonsputum Place Of Service:LABCORP	
UHEALTH Labs	MTHFR	MTHFR Mutation Container:Citrated Light Blue Top Tube Place Of Service:Uhealth Labs	
UHEALTH Labs	MUMPGM	MUMPS (IGG/M) Place Of Service:LABCORP	
UHEALTH Labs	CSFMUM	MUMPS AB PANEL, IFA (CSF) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MUIGG	MUMPS ANTIBODIES, IGG Place Of Service:LABCORP	
UHEALTH Labs	MMPSIM	MUMPS ANTIBODIES, IGM Place Of Service:LABCORP	
UHEALTH Labs	MUMPAB	MUMPS VIRUS ANTIBODIES (IGG,IGM) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MUMGQ	MUMPS VIRUS ANTIBODY (IGG) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MUMPCR	Mumps Virus RNA, Qualitative Real-Time PCR Container:Swab in sterile, leak proof container in 3 mL M4 media,or V-C-M medium (green-cap) tube or equivalent (UTM). Temperature:Refrigerated (cold packs) Place Of Service:QUEST	87798
UHEALTH Labs	MUSK	Muscle Specific Kinase (MuSK) Antibodies Place Of Service:LabCorp	
UHEALTH Labs	MUSKLR	MuSK and LRP4 Antibodies Panel Container:PLASTIC TRANSPORT TUBE Temperature:REFRIGERATE Place Of Service:QUEST	83519 86255
UHEALTH Labs	MUSTAR	Mustard IgE, F089 Place Of Service:LabCorp	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	MYA2	MYASTHENIA GRAVIS PANEL 2 Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	MYCOSM	MYCOBACTERIA, CULTURE, WITH FLUOROCHROME SMEAR Container:Refer to Preferred Specimen(s) and Instructions Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	MICPNL	MYCOBACTERIUM AVIUM COMPLEX MIC PANEL Container:Double-walled safety container Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	MTBSUS	MYCOBACTERIUM TUBERCULOSIS SUSCEPTIBILITY Place Of Service:DHLTH	
UHEALTH Labs	AFBC	Mycobacterium, Culture, Blood Container:BACTEC Myco/F Lytic aerobic blood bottle Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MYBC	MYCOLOGY BLOOD CULTURE Container:MYCO/F LYTIC BOTTLE (RED TOP) Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87103
UHEALTH Labs	MYCC	MYCOLOGY CULTURE Container:STERILE CUP OR CULTURETTE Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87102
UHEALTH Labs	MYCD	MYCOLOGY CULTURE SKIN, HAIR OR NAILS Container:STERILE CUP Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87101
UHEALTH Labs	MYCOAM	MYCOPHENOLIC ACID AND METABO. Place Of Service:LABCORP	
UHEALTH Labs	MCOGM	MYCOPLASMA PNEU. IGG/IGM ABS Place Of Service:LABCORP	
UHEALTH Labs	MYPNCL	MYCOPLASMA PNEUMONIAE CULTURE Place Of Service:LABCORP	
UHEALTH Labs	MYCIGG	MYCOPLASMA PNEUMONIAE, IGG AB Place Of Service:LABCORP	
UHEALTH Labs	MYCIGM	MYCOPLASMA PNEUMONIAE, IGM AB Place Of Service:LABCORP	
UHEALTH Labs	MYECSE	MYELIN BASIC PROTEIN, CSF Place Of Service:LABCORP	
UHEALTH Labs	MAGP	MYELIN-ASSOC. GLYCOPROTEIN Ab(IgM) Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	MYMP	MYELOID MOLECULAR PROFILE Container:Peripheral blood: 2-3 mL in EDTA (purple-top) tube Place Of Service:UHEALTH LABS	81450
UHEALTH Labs	MYEMRD	Myeloma MRD Profile Container:sodium heparin (green-top) tube Temperature:Refrigerated (cold packs) Place Of Service:Genoptix	88184 88185 x9 88188
UHEALTH Labs	MPNMT	MYELOPROLIF NEOPLASM M WO B/A/J/M Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	MYOG	MYOGLOBIN, SERUM Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	MYOS	MYOGLOBIN, SERUM Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83874
UHEALTH Labs	MYOU	MYOGLOBIN, URINE Place Of Service:LABCORP	
UHEALTH Labs	MYOM	MyoMarker 3 Profile (RDL) Place Of Service:LabCorp	
UHEALTH Labs	MYOSAB	Myositis Specific 11 Antibodies Panel Container:Transport tube Temperature:Room temperature Place Of Service:QUEST	84182 x6 86235 x5
UHEALTH Labs	NMEN	N.Meningitidis A/C/W-135/Y IgG Place Of Service:LABCORP	
UHEALTH Labs	NAPA	N-ACETYLPROCAINAMIDE LEVEL Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	NASALC	NASAL CULTURE Container:CULTURETTE SWAB Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	NATAB	Natalizumab and Anti-Natalizumab Ab,DoseASSURE NTZ Place Of Service:LABCORP	
UHEALTH Labs	YNKF	NATURAL KILLER CELL (NK) FUNCTION Place Of Service:ARUP	
UHEALTH Labs	NSIGON	NEISSERIA GONORRHOEA AB, TOTAL Place Of Service:LABCORP	
UHEALTH Labs	NGRNAT	Neisseria gonorrhoeae RNA, TMA Place Of Service:QUEST	
UHEALTH Labs	NEOPR	NeoCerebellar Degeneration Paraneoplastic Profile with Recombx Container:Transport tube Temperature:Refrigerated (cold packs) Place Of Service:Quest	83520 x7 86341
UHEALTH Labs	NPPE	NEOENCEPHALITIS PARANEOPLS RECOMB Container:A RED-TOP TUBE (NO GEL) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	NCLLPP	NeoTYPE CLL Prognostic Profile Container:Peripheral blood: 5 mL in EDTA tube. Bone marrow: 2 mL in EDTA tube. Fresh tissue: 0.5 - 1 cm3 in RPMI. Note: Please exclude biopsy needles, blades, and other foreign objects from transport tubes. These can compromise specimen viability and yield, and create hazards for employees. Temperature:Refrigerate specimen. Use cold pack for transport, making sure cold pack is not in direct contact with specimen. Ship same day as drawn whenever possible.	88374 x4 81445
UHEALTH Labs	NEURLC	Neurofilament Light Chain, Serum Place Of Service:LABCORP	
UHEALTH Labs	NEUCSF	Neuron Specific Enolase,CSF Place Of Service:LabCorp	
UHEALTH Labs	ANTIHU	NEURONAL NUCLEAR AB(Hu) SCREEN Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	NEUENO	NEURON-SPECIFIC ENOLASE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	NEUAB	NEUTROPHIL AB FLOW CYTOMETRY Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	NEUHL	Neutrophil Ab ID and HLA Ab Screen Container:Serum Temperature:Refrigerated Place Of Service:UHEALTH LAB	86021 86828
UHEALTH Labs	NEUABS	Neutrophil Ab Screen RFX to 5113 (Neut ID & HLA Scrn) Container:Serum- not specified Temperature:Refrigerate Place Of Service:UHEALTH LABS	86021
UHEALTH Labs	NEUSWR	Neutrophil Ab Screen w/ RFX HLA Ab Screen Container:Serum- not specified Temperature:Refrigerate Place Of Service:UHEALTH LABS	86021
UHEALTH Labs	OXBUR	NEUTROPHIL FUNCTION Container:SODIUM HEPARIN (GREEN-TOP) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	NICMET	NICOINE AND METABOLITE, QUANT Place Of Service:LABCORP	
UHEALTH Labs	NICMU	NICOTINE METABOLITE,URINE Place Of Service:LABCORP	
UHEALTH Labs	NMDARA	NMDA Receptor (NR1-subunit) Autoantibody Test Container:Transport tube or CSF collected in: Sterile leak-proof container Temperature:Refrigerated (cold packs) Place Of Service:Quest	86255
UHEALTH Labs	NMOIGA	NMO IGG AUTOANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	NMRLIP	NMR LIPOPROFILE Place Of Service:LABCORP	
UHEALTH Labs	NOCBAL	Nocardia Real-time PCR (BAL) Container:sterile, screw top tube Temperature:ambient or frozen temperature Monday through Friday Place Of Service:Viracor	87798
UHEALTH Labs	FAFS	NONESTERIFIED FATTY ACIDS Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	LCF	NON-GYN, FINE NEEDLE ASPIRATION SPECIMEN(S) Container:Slide(s), Coplin jar(s) Temperature:Refrigerate Place Of Service:UHEALTH LABS	
UHEALTH Labs	NOREPI	NOREPINEPHRINE, PLASMA Container:sodium heparin (green-top) tube Temperature:Frozen Place Of Service:QUEST	
UHEALTH Labs	NORWAL	NOROVIRUS, EIA (STOOL) Container:STERILE CONTAINER Temperature:FROZEN Place Of Service:VIRACOR	
UHEALTH Labs	NRVPCR	NOROVIRUS, RT-PCR Place Of Service:LABCORP	
UHEALTH Labs	NORTRI	NORTRIPTYLINE (AVENTYL), SERUM Place Of Service:LABCORP	
UHEALTH Labs	NMP1F	NPM1 MUTATION BY PCR& FRAGMENT ANAL Temperature:REFRIDGERATED Place Of Service:ARUP	
UHEALTH Labs	NTELOS	N-Telopeptide, Serum Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	NTELU	N-TELOPEPTIDE, URINE Place Of Service:LABCORP	
UHEALTH Labs	BNP	NT-proBNP Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83880
UHEALTH Labs	NUCMP	Nuclear Matrix Proteins (NMP), Random Urine Container:Matritech NMP-22 collection kit Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	OAKW	OAK, WHITE Place Of Service:LABCORP	
UHEALTH Labs	OAKLV	Oak,Live/Virginia IgE (T218) Place Of Service:LabCorp	
UHEALTH Labs	OAT	OAT (F7) IgE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	OBGAST	Occult Blood, Gastric Container:ARUP Standard Transport Tube (ARUP Supply #43115) Temperature:Refrigerated Place Of Service:ARUP	
UHEALTH Labs	OLIBND	OLIGOCLONAL BANDING Place Of Service:LABCORP	
UHEALTH Labs	OLIVE	Olive Tree (T009) IgE Place Of Service:LabCorp	
UHEALTH Labs	COCCAB	OMEGA Coccidioides Antibody EIA Container:Serum Separator Tube Temperature:Frozen Place Of Service:Viracor	86335
UHEALTH Labs	CODE	OPIATES BY GC/MS (U) Container:PLASTIC URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ORGACU	ORGANIC ACID ANALYSIS, URINE Place Of Service:LABCORP	
UHEALTH Labs	ORTHOC	ORTHOPEDIC CULTURE Temperature:STORE AT RT DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	OSM	OSMOLALITY Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83930
UHEALTH Labs	FEOSMO	OSMOLALITY, FECAL Place Of Service:LABCORP	
UHEALTH Labs	OSMOFE	OSMOTIC FRAGILITY ERYTHROCYTE Place Of Service:LABCORP	
UHEALTH Labs	OSTEO	OSTEOCALCIN, HUMAN (BGP) Container:No additive (red-top) Temperature:REFRIGERATED PREFERRED ROOM TEMPERATURE UNACCEPTABLE FROZEN ACCEPTABLE Place Of Service:QUEST	
UHEALTH Labs	OVAPAR	OVA AND PARASITES Container:PARA PAK ULTRA ECOFIX SYSTEM (GREEN) STOOL OR STERILE CUP Temperature:STORE PRESERVE SPECIMENS AT ROOM TEMPERATURE, IF UNPRESERVED REFRIGERATE WITHIN 1 HOUR OF COLLECTION. Place Of Service:UHealth Labs	87177
UHEALTH Labs	OVAPS2	Ova and Parasites, Stool, Concentrate and Permanent Smear, 2 Specimens Container:10% formalin and Polyvinyl Alcohol Transport or Total-Fix Transport vial; Urine submitted in a sterile screw-cap container; Sputum submitted in a sterile screw-cap container or with 10% formalin	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Temperature:STOOL AND SPUTUM (PRESERVED): ROOM TEMPERATURE URINE AND SPUTUM (UNPRESERVED): REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	OVAB	OVARIAN ANTIBODY SCREEN W/REFL Container:PLASTIC SCREW CAP VIAL Temperature:SHIP REFRIGERATED; SHIPPING ROOM TEMP ACCEPTABLE; SHIPPING FROZEN ACCEPTABLE Place Of Service:QUEST	
UHEALTH Labs	ROMAB	OVARIAN MALIGNANCY RISK-ROMA Place Of Service:LABCORP	
UHEALTH Labs	OXARU	OXALIC ACID, RANDOM URINE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	OXCARB	OXCARBAZEPINE (TRILEPTAL), S Place Of Service:LABCORP	
UHEALTH Labs	OXYCOD	OXYCODONE & METABOLITE, QN, SER Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	OXCOMO	OXYCODONE/OXYMORPHONE CONFIRM Place Of Service:LABCORP	
UHEALTH Labs	P50B	Oxygen Dissociation, P50, Erythrocytes Container:Patient : Green top (sodium heparin) and Normal Shipping Control:Green top (sodium heparin) Temperature:Immediately refrigerate specimens after collection. Place Of Service:Mayo	82820
UHEALTH Labs	GREAR	P190 BCR ABL1 Container:EDTA (lavender-top) tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	GREAP	P210 BCR ABL1 Container:EDTA (lavender-top) tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	PMBARB	PAIN BARBITURATES CONF medMATCH Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PMFNT	PAIN FENTANYL, QN,W/medMATCH,U Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PMGABP	PAIN GABAPENTIN QN medMATCH U Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PMPMED	PAIN MANAGEMENT PROFILE 1 W/ CONFIRMATION, URINE Container:Clinical drug test transport vial Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	PMPTX	PAIN MANAGEMENT PROFILE 1 W/CONFIRMATION,URINE Container:Urine Collection Container Temperature:Room Temperature Place Of Service:UHealth Labs	80307
UHEALTH Labs	PMFQU	Pain Management, Fentanyl, Quantitative, Urine Container:Clinical drug test transport vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PMMETU	Pain Management, Methadone Metabolite, QN, Urine Container:Clinical drug test transport vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PMTAPU	Pain Management, Tapentadol, Quantitative, Urine Container:Clinical drug test transport vial Temperature:ROOM TEMPERATURE	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:QUEST	
UHEALTH Labs	PMOXC	PAIN OXYCODONE QN,W/medMATCH,U Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PALMQ	Palm, Queen IgE (T072) Place Of Service:LabCorp	
UHEALTH Labs	PNCTFE	PANCREATIC ELASTASE, FECAL Place Of Service:LABCORP	
UHEALTH Labs	PANPLY	PANCREATIC POLYPEPTIDE Place Of Service:LABCORP	
UHEALTH Labs	PAPRIK	Paprika (F218) IgE Place Of Service:LabCorp	
UHEALTH Labs	PABRAS	Paracoccid. brasiliensis Ab Place Of Service:LabCorp	
UHEALTH Labs	PARARC	Paraneoplastic Ab Eval. w/Rflx Titer and LB,Basic,CSF Container:Sterile leak-proof container Temperature:Refrigerated (cold packs) Place Of Service:Quest	86255 x10 86341
UHEALTH Labs	PARAEV	Paraneoplastic Antibody Evaluation w/Reflex to Titer and Western Blot,Basic Container:Transport tube Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	PARACF	Paraneoplastic Autoantibody Eval. w/Recombx, CSF Container:Sterile leak-proof container Temperature:Refrigerated (cold packs) Place Of Service:Quest	86255 x3 83520 x8
UHEALTH Labs	QPTHI	PARATHYROID HORMONE,INTACT Container:Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	PARCAB	PARIETAL CELL AB W/REFL TITER Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PARO	PAROXETINE, SERUM/PLASMA Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PRSIGE	PARSLEY (F86) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PARIGG	PARVOVIRUS B19 ANTIBODY (IGG) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PARIGM	PARVOVIRUS B19 ANTIBODY (IGM) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PB19Q	PARVOVIRUS B19 DNA, QN RT PCR Container:PLASTIC LEAK-PROOF CONTAINER Temperature:FROZEN. ROOM TEMPERATURE: WHOLE BLOOD OR BONE MARROW Place Of Service:QUEST	
UHEALTH Labs	PRVPCR	PARVOVIRUS B19 PCR Place Of Service:LABCORP	
UHEALTH Labs	PARVGM	PARVOVIRUS B19, HUMAN, IGG/IGM Place Of Service:LABCORP	
UHEALTH Labs	PB19IG	PARVOVIRUS B19, IGG Place Of Service:LABCORP	
UHEALTH Labs	PB19IM	PARVOVIRUS B19, IGM Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	THCAPX	PAX8/PPAR TRANSLOCATION THYRD CNC Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PDGIST	PDGFRA ANALYSIS IN GISTS Place Of Service:LABCORP	
UHEALTH Labs	PEIRFX	PE+INTERP(RFX IFE),24-HR UR Place Of Service:LABCORP	
UHEALTH Labs	PEAAL	PEA (F12) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PEMAB	Pemphigus Antibody Panel Container:Plain red or serum separator tube. Temperature:Refrigerate Place Of Service:ARUP	88346 88350 83516
UHEALTH Labs	PENTA	PENTA SCREEN Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PENTOQ	PENTOBARBITAL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SCANRV	PERIPHERAL SMEAR SCAN Container:Lavender Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	PERPH	PERPHENAZINE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PHFLD	PH, BODY FLUID Place Of Service:LABCORP	
UHEALTH Labs	PHSTO	PH, STOOL Place Of Service:LABCORP	
UHEALTH Labs	PHENCU	PHENCYCLIDINE, CONFIRM, URINE Place Of Service:LABCORP	
UHEALTH Labs	PHNO	PHENOBARBITAL LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80184
UHEALTH Labs	STPHEN	PHENOLPHTHALEIN, FECES Container:STERILE, PLASTIC, LEAK-PROOF, FECES CONTAINER Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	PHENGT	PHENOSENSE GT (R) Place Of Service:LABCORP	
UHEALTH Labs	PHOSRU	PHOSPHATE, RAND UR (W/CREATININE) Container:PLASTIC URINE CONTAINER Temperature:PRESERVED: ROOM TEMPERATURE UNPRESERVED: REFRIGERATED (COLD PACK) Place Of Service:QUEST	
UHEALTH Labs	PHSPGM	PHOSPHATIDYL SERINE (IGG/M) Place Of Service:LABCORP	
UHEALTH Labs	PHLIP	PHOSPHOLIPIDS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	PHOS	PHOSPHOROUS Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84100

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	PIGNAL	PIGEON DROPPINGS (E7) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PINWRM	Pinworm Preparation Place Of Service:LabCorp	
UHEALTH Labs	PJIS	PJI Detection (Synovasure) ARUP 2013008 Container:Red tube Temperature:Refrigerated. Place Of Service:ARUP	83516 84311 86140
UHEALTH Labs	PAI1P	Plasminogen Activator Inhibitor 1 4G/5G Polymorphism Place Of Service:LabCorp	
UHEALTH Labs	PLAS	PLASMINOGEN ACTIVITY Container:Citrated Light Blue Top Tube Place Of Service:Uhealth Labs	85420
UHEALTH Labs	PLSMAI	PLASMINOGEN ACTVINHBT PAI 1 4G/5G Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PLT	PLATELET Container:LAVENDER TOP TUBE Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85049
UHEALTH Labs	WBAGG	PLATELET AGGREGATION Container:BLUE Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	PLTPRO	PLATELET ANTIBODY PROFILE Place Of Service:LABCORP	
UHEALTH Labs	PLTIGG	PLATELET ANTIBODY, DIRECT (IGG) Container:ACD SOLUTION A (YELLOW-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PLTIND	Platelet Antibody, Indirect (IgG) Container:Transport tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	PFA	PLATELET FUNCTION ANALYSIS Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE IS PREFERRED Place Of Service:UHealth Labs	86022
UHEALTH Labs	PSGP	Platelet Surface Glycoprotein by Flow Cytometry, Blood Container:Yellow (ACD) Temperature:Room Temperature Place Of Service:MAYO	88184 88185 88187
UHEALTH Labs	PTEMS	Platelet Transmission Electron Microscopic Study MAYO PTEM Container:Whole Blood ACD Temperature:Refrigerate Place Of Service:MAYO	85390 88348
UHEALTH Labs	PLTBLU	PLATELET, CITRATED BLOOD Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85025
UHEALTH Labs	PLATAB	PLT ANTIBODY COMPREHENSIVE Place Of Service:LABCORP	
UHEALTH Labs	PMSCL	PM SCL ANTIBODY Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PNE23S	PNEUMOCOCCAL AB (23 SEROTYPE) Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	PNE14B	PNEUMOCOCCAL IM (14 SEROTYPE) Place Of Service:LABCORP	
UHEALTH Labs	DFPCP	PNEUMOCYSTIS CARINII, DFA Container: Sterile, leak-proof container Place Of Service: JMH	
UHEALTH Labs	PNEDNA	Pneumocystis DNA detection by PCR Container: Sterile screw cap container Place Of Service: University of Washington	
UHEALTH Labs	PJPCR	PNEUMOCYSTIS JIROVECI CARINI(SER) Container: sterile, screw top tube. Temperature: Can be shipped at ambient or frozen temperature Place Of Service: VIRACOR	
UHEALTH Labs	PNEJBW	Pneumocystis jiroveci Qn RT PCR, BW Container: Sterile, screw top tube Temperature: Specimen can be shipped at ambient or frozen temperature Monday through Friday. Place Of Service: Viracor	87799
UHEALTH Labs	POLIOV	POLIOVIRUS ANTIBODIES Place Of Service: LABCORP	
UHEALTH Labs	POLMAB	POLYMYOSITIS DERMATOMYOSITIS AB Container: RED-TOP TUBE (NO GEL) Temperature: ROOM TEMPERATURE Place Of Service: QUEST	
UHEALTH Labs	POPPY	Poppy Seed IgE, F224 Place Of Service: LabCorp	
UHEALTH Labs	POR24U	PORPHOBILINOGEN, QN, 24-HR UR Place Of Service: LABCORP	
UHEALTH Labs	PORURB	PORPHOBILINOGEN, QN, RANDOM, UR Place Of Service: LABCORP	
UHEALTH Labs	PORPHF	PORPHYRINS FRACTIONATED QN 24H U Temperature: REFRIGERATED (COLD PACKS) Place Of Service: QUEST	
UHEALTH Labs	PRPHQN	PORPHYRINS FRACTIONATED QN RND UR Temperature: REFRIGERATED (COLD PACKS) Place Of Service: QUEST	
UHEALTH Labs	PRP24U	PORPHYRINS, QN, 24 HR UR. Place Of Service: LABCORP	
UHEALTH Labs	PORQNU	PORPHYRINS, QN, RANDOM U Place Of Service: LABCORP	
UHEALTH Labs	PORTOT	PORPHYRINS, TOTAL PLASMA Place Of Service: LABCORP	
UHEALTH Labs	POSLEV	Posaconazole LC-MS/MS (serum)* Container: Red Temperature: Frozen Place Of Service: Virocore	80187
UHEALTH Labs	POSAC	POSACONAZOLE LEVEL Container: Red top Temperature: Frozen Place Of Service: FTL	
UHEALTH Labs	PVASA	POST VAS SEMEN ANALYSIS, AUA Place Of Service: LABCORP	
UHEALTH Labs	K	POTASSIUM Container: GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature: TRANSPORT AT ROOM TEMPERATURE Place Of Service: UHealth Labs	84132
UHEALTH Labs	UKCR	POTASSIUM, 24 HOUR URINE (W/ CREAT) Container: JUG- PLASTIC URINE CONTAINER Temperature: TRANSPORT AT ROOM TEMPERATURE Place Of Service: UHealth Labs	84133



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	POHEP	POTASSIUM, HEPARIN PLASMA Place Of Service:LABCORP	
UHEALTH Labs	KPLAS	POTASSIUM, PLASMA Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PREALB	PREALBUMIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84134
UHEALTH Labs	PREGNE	PREGNENOLONE, LC/MS/MS Container:R RED-TOP TUBE (NO GEL) Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	PREGMS	PREGNENOLONE, MS Place Of Service:LABCORP	
UHEALTH Labs	PREKLA	PREKALLIKREIN ASSAY Place Of Service:LABCORP	
UHEALTH Labs	PRENAB	Prenatal Infectious Disease Antibodies,Qualitative, IgG Place Of Service:LABCORP	
UHEALTH Labs	PRISER	PRIMIDONE (MYSOLINE), SERUM Place Of Service:LABCORP	
UHEALTH Labs	PROCAN	PROCAINAMIDE LEVEL Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	PANAPA	PROCAINAMIDE LEVEL (INCLUDES NAPA) Container:RED top tube (No gel) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	PRCALT	PROCALCITONIN Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	PROG	PROGESTERONE LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84144
UHEALTH Labs	PRINS	PROINSULIN Place Of Service:LABCORP	
UHEALTH Labs	PROL	PROLACTIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84146
UHEALTH Labs	PROLAC	PROLACTIN Container:SERUM SEPARATOR TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PROLPM	Prolactin, Pituitary Macroadenoma, Serum Place Of Service:LABCORP	
UHEALTH Labs	PROMDG	PROMETHEUS IBD sgi Diagnostic Container:red-top (no gel) tube --AND-- EDTA (lavender-top) tube Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	PROPFN	PROPAFENONE Container:RED-TOP TUBE (NO GEL) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RYTHM	PROPAFENONE (RYTHMOL) Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	PRPTDI	PROPEPTIDE TYPE I COLLEGEN Place Of Service:LABCORP	
UHEALTH Labs	PRAPH	PROSTATIC ACID PHOS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	PAPH	PROSTATIC ACID PHOSPHATASE (PAP) Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	PSA	PROSTATIC SPECIFIC ANTIGEN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84153
UHEALTH Labs	PSAFT	PROSTATIC SPECIFIC ANTIGEN,FREE/TOTAL Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84153 84154
UHEALTH Labs	PCACT	PROTEIN C ACTIVITY Container:Plastic screw cap transfer vial Temperature:Frozen Place Of Service:UHealth Labs	85303
UHEALTH Labs	PCACAG	PROTEIN C ACTIVITY AND ANTIGEN Container:Plastic screw cap transfer vial Temperature:Frozen Place Of Service:UHealth Labs	85303
UHEALTH Labs	PCRAG	Protein C Activity w/ Reflex to Protein C Antigen Container:Blue-top (sodium citrate) tube Temperature:Freeze; three freeze/thaw cycles are acceptable. Stable at room temperature for eight hours. Do not refrigerate. Place Of Service:UHealth	85303 85302
UHEALTH Labs	PCSACT	PROTEIN C AND PROTEIN S, FUNCTIONAL Container:Plastic screw cap transfer vial Temperature:Frozen Place Of Service:UHealth Labs	85302
UHEALTH Labs	PCSWR	PROTEIN C AND S ACTIVITY W/REFLEX TO ANTIGEN Container:Plastic screw-cap vial (x2) Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	PCPSFN	Protein C and S Activity w/Reflex to Antigen Container:Blue-top (sodium citrate) tube Temperature:Freeze; three freeze/thaw cycles are acceptable. Stable at room temperature for eight hours. Do not refrigerate. Place Of Service:LabCorp	85303 85306
UHEALTH Labs	PCANT	PROTEIN C ANTIGEN Container:Plastic screw cap transfer vial Temperature:Frozen Place Of Service:UHealth Labs	85306
UHEALTH Labs	PRCFNC	PROTEIN C FUNCTIONAL W/RFLX Place Of Service:LABCORP	
UHEALTH Labs	PREL24	PROTEIN ELECTRO, 24-HOUR URINE Place Of Service:LABCORP	
UHEALTH Labs	URPEP	Protein Electrophoresis, 24-Hour Urine (UPEP) Container:PLASTIC URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PSFUN	Protein S Activity Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85306
UHEALTH Labs	PSARAG	PROTEIN S ACTIVITY W/REFL PROTEIN S ANTIGEN TOTAL Container:Plastic screw-cap vial Temperature:FROZEN	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
		Place Of Service:QUEST	
UHEALTH Labs	PSFREE	Protein S Free Antigen Container:Citrated Light Blue Top Tube Place Of Service:Uhealth Labs	85306
UHEALTH Labs	PSPANL	PROTEIN S PANEL Container:Blue top (Citrated plasma) Temperature:Room Temperature Place Of Service:Uhealth Labs	85305 85306
UHEALTH Labs	PSTOT	Protein S Total Antigen Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85305
UHEALTH Labs	PEP	PROTEIN, TOTAL AND PROT ELECTRO Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	PTPEWR	PROTEIN, TOTAL AND PROTEIN ELECTROPHORESIS, URINE W/REFL Container:Plastic urine container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	UTPCR	Protein, Total, with Creat.,24 hr URN Container:JUG- Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84155
UHEALTH Labs	PTPE24	PROTEIN,TOTAL AND PROTEIN ELECTROPHORESIS,24 HOUR URINE AND IMMUNOFIXATION Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	PTFRAG	PROTHROMBIN 1+2 Container:BLUE Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	P20210	PROTHROMBIN 20210 Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	81240
UHEALTH Labs	PTABGM	PROTHROMBIN ANTIBODIES IGG IGM Container:3.2% sodium citrate (light blue-top) tube Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	PTPTTM	PT AND PTT MIXING STUDIES Container:collect in blue top tube (sodium citrate) and transfer to transfer vial Temperature:TRANSPORT FROZEN Place Of Service:UHealth Labs	85610 85730
UHEALTH Labs	PTINH	PT Mixing Study Container:Citrated Light Blue Top Tube Place Of Service:UHealth Labs	85611 x2
UHEALTH Labs	PTH	PTH, INTACT Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83970
UHEALTH Labs	PTHNCA	PTH, Intact (ICMA) and Ionized Calcium Container:PTH Intact: Transport tube. Ionized Calcium: Unopened gel barrier tube.Transport Temperature-PTH Intact: Frozen.Ionized Calcium: Room temperature Temperature:Room Temperature AND FROZEN Place Of Service:QUEST	82310 82330 83970
UHEALTH Labs	PTHICA	PTH, INTACT + CALCIUM Container:2 TUBES ( One LAVENDER top tube and one GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82310 83970
UHEALTH Labs	PTHRPB	PTHRP (PTH-RELATED PEPTIDE) Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	PTINR	PT-INR Container:Light blue Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	PTINRL	PT-INR Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85610
UHEALTH Labs	PTTINH	PTT Mixing Study Container:Citrated Light Blue Top Tube Temperature:FROZEN Place Of Service:UHealth Labs	85732 x5
UHEALTH Labs	PTTLAB	PTT-LA Container:Blue-top (sodium citrate) tube Temperature:Frozen Place Of Service:LABCORP	
UHEALTH Labs	PURK	PURKINJE CELL (Yo) Ab SCREEN Container:Transport tube Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	FPURK	PURKINJE CELL (Yo) Ab SCREEN,CSF Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	PYKINB	PYRUVATE KINASE Place Of Service:LABCORP	
UHEALTH Labs	PYRACI	PYRUVIC ACID, BLOOD Place Of Service:LABCORP	
UHEALTH Labs	QFEVER	Q FEVER ANTIBODIES (IGG,IGM) Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	QFVIGG	Q FEVER ANTIBODIES, IGG Place Of Service:LABCORP	
UHEALTH Labs	QTFTBG	QuantiFERON-TB GOLD PLUS Container:QTF5 QTFTBG (Quantiferon kit 4 tubes) Temperature:Room Temperature. Place Of Service:Uhealth Labs	86480
UHEALTH Labs	QUINI	QUINIDINE LEVEL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	QUIIGE	Quinoa (f347) IgE Container:Serum Separator Ttube (SST®) Temperature:Room temperature Place Of Service:QUEST	86003
UHEALTH Labs	RABBIT	Rabbit Epithelium IgE (E082) Place Of Service:LabCorp	
UHEALTH Labs	CRW	RAGWEED(W1),COMMON IgE ALLERGEN Container:SST/Red Top Temperature:Room Temp Place Of Service:QUEST	
UHEALTH Labs	SRAGW	Ragweed, Short IgE (W001) Place Of Service:LabCorp	
UHEALTH Labs	RGBSUS	RAPID GROWER BROTH SUSCEP. Place Of Service:LABCORP	
UHEALTH Labs	HIVRAP	RAPID HIV Container:Lav Top-EDTA Temperature:Room Temp Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	RPDQN	RAPID PLASMA REAGIN, QUANT Place Of Service:LABCORP	
UHEALTH Labs	NRAPDA	RAPID STREP GROUP A Container:regular swab Temperature:Room Temperature Place Of Service:JMH Virology	
UHEALTH Labs	RASTSH	RAS Mutation Analysis, Thyroid Cancer Container:Fixative, slides or paraffin embedded tissue block Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RAT	Rat Epithelium IgE (E073) Place Of Service:LABCORP	
UHEALTH Labs	RBCMDG	RBC Membrane Disorders Gene Sequencing Panel Container:EDTA (lavender top) for whole blood. Saliva kit for saliva. Temperature:Room temperature Place Of Service:UHEALTH LAB	81443
UHEALTH Labs	STREDS	REDUCING SUBSTANCES, FECES Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	RFP	RENAL FUNCTION PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80069
UHEALTH Labs	RENINP	RENIN ACTIVITY, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	REPT	REPTILASE TIME Container:Citrated Light Blue Top Tube Place Of Service:Uhealth Labs	85635
UHEALTH Labs	RASTFL	RESP ALLERGY PANEL REGION IV Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RSPBAC	Respiratory Bacterial Panel TEM-PCR (BAL) Container:Collect BAL in a sterile, screw top tube. Temperature:Can be shipped at ambient or frozen temperature Place Of Service:Viracor	87798 x9 87541 87486 87581 87641 87651
UHEALTH Labs	RESPC	RESPIRATORY CULTURE WITH GRAM STAIN Container:STERILE CUP Temperature:REFRIGERATE IF DELAYED DELIVERY IS ANTICIPATED Place Of Service:UHealth Labs	87070
UHEALTH Labs	RPP	RESPIRATORY PATHOGEN PANEL,NUC ACID PROBE Container:1 SWAB IN UNIVERSAL TRANSPORT MEDIA (UTM) Temperature:STORE AT RT DO NOT REFRIGERATE Place Of Service:UHealth Labs	
UHEALTH Labs	RSVIA	Respiratory Syncytial Virus (RSV), Immunoassay Place Of Service:LabCorp	
UHEALTH Labs	RVH1N1	Respiratory Virus PCR Panel with 2009 H1N1 Container:VCM medium (green-cap) tube Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	RETSDD	RET Sequencing and Deletion/Duplication Container:EDTA (lavender-top) tube or ACD (yellow-top) tube Temperature:Room Temperature Place Of Service:QUEST	
UHEALTH Labs	THYCA	RET/PTC REARRANGEMENT, THYRD CANC Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RETIGA	RETICULIN IGA ANTIBODIES Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	RETC	RETICULOCYTE PANEL Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85045
UHEALTH Labs	RT3SER	REVERSE T3, SERUM Place Of Service:LABCORP	
UHEALTH Labs	RADPI	RHEUMATOID ARTH IDENTRA W1433 ETA Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RHEPRO	RHEUMATOID ARTHRITIS PROFILE Place Of Service:LABCORP	
UHEALTH Labs	RA	RHEUMATOID FACTOR Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86431
UHEALTH Labs	RISPPR	RISPERIDONE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RISTOC	Ristocetin Cofactor Activity of VWF Container:Citrated Light Blue Top Tube Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85245
UHEALTH Labs	RKYGM	RKY (IGG/M) Place Of Service:LABCORP	
UHEALTH Labs	RNAP3	RNA POLYMERASE III AB Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RNPAB	RNP Antibodies Container:Red-top tube or gel-barrier tube Temperature:Room temperature Place Of Service:LABCORP	
UHEALTH Labs	ROTOVRB	ROTAVIRUS AG, EIA Place Of Service:LABCORP	
UHEALTH Labs	RPRDX	RPR (DX) W/REFL TITER AND FTA Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RPRQNR	RPR, RFX QN RPR/CONFIRM TP Place Of Service:LABCORP	
UHEALTH Labs	RUBAG	RUBELLA AB IgG Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	86762
UHEALTH Labs	RBLAIM	RUBELLA ANTIBODIES, IGM Place Of Service:LABCORP	
UHEALTH Labs	RUBIGG	RUBELLA IMMUNE STATUS Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RUBIGB	RUBEOLA ANTIBODIES, IGG Place Of Service:LABCORP	
UHEALTH Labs	BANZEL	RUFINAMIDE (BANZEL) Place Of Service:LABCORP	
UHEALTH Labs	RYEF	RYE (F5) IgE Container:Serum Separator Tube (SST®) Temperature:Room temperature Place Of Service:QUEST	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	SPNAG	S. PNEUMONIAE ANTIGEN DETECTION,SERUM OR CSF Container:Transport tube Temperature:Frozen Place Of Service:QUEST	86403
UHEALTH Labs	SACAB	SACCHAROMYCES CEREVISIAE AB Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ASCAIG	Saccharomyces cerevisiae Antibodies (ASCA) (IgG) Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SACCCP	SACCHAROMYCES CEREVISIAE PANEL Place Of Service:LABCORP	
UHEALTH Labs	SALI	SALICYLATE LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80329
UHEALTH Labs	SLCORT	SALIVARY CORTISOL, MS Place Of Service:LABCORP	
UHEALTH Labs	SALMON	SALMON (F41) IgE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	COV2AB	Sars-CoV-2 Semi-Quantitative Total Antibody, Spike Place Of Service:LABCORP	
UHEALTH Labs	SCHAB	SCHISTOSOMA IGG AB FMI (SERUM) Container:RED TOP OR GOLD TOP Temperature:REFRIGERATED Place Of Service:QUEST	86682
UHEALTH Labs	SCLERA	SCL-70 ANTIBODY Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SDP	Scleroderma Diagnostic Profile Place Of Service:LABCORP	
UHEALTH Labs	SELE	SELENIUM, SERUM/PLASMA Place Of Service:LABCORP	
UHEALTH Labs	PATHRV	SEND FOR PATHOLOGIST REVIEW Container:Lavender Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	SERREL	SEROTONIN RELEASE ASSAY Place Of Service:LABCORP	
UHEALTH Labs	SEROTS	SEROTONIN, SERUM Place Of Service:LABCORP	
UHEALTH Labs	SROTNB	SEROTONIN, WHOLE BLOOD Place Of Service:LABCORP	
UHEALTH Labs	ZOLOFT	SERTRALINE/DESMETHYLSERTRALINE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SPREG	SERUM PREGNANCY (QUAL) Container:GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	SPEPSC	SERUM PROTEIN ELECTROPHORESIS Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84155

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	SHBGR	SEX HORM BINDING GLOB, SERUM Place Of Service:LABCORP	
UHEALTH Labs	SEZAC	SEZARY CELL COUNT Container:EDTA (lavender-top) tube Temperature:Room temperature Place Of Service:UHEALTH LABS	
UHEALTH Labs	SHRIMP	SHRIMP (F24) IgE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SICKLE	SICKLE CELL SCREEN Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85660
UHEALTH Labs	SILICA	SILICON Container:ACID WASHED PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SINEME	SINEMET(R) Container:Specimen Container: PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	SIROB	SIROLIMUS (RAPAMUNE), BLOOD Place Of Service:LABCORP	
UHEALTH Labs	YSIROL	SIROLIMUS (RAPAMYCIN) LEVEL Container:Lavender-top tube Temperature:SHIP AND STORE REFRIDGERATED.SHIPPING AT AMBIENT TEMPERATURE (<25 C) IS ACCEPTABLE. Place Of Service:IML	
UHEALTH Labs	SSASSB	SJOGREN'S AB, ANTI-SS-A/-SS-B Place Of Service:LABCORP	
UHEALTH Labs	SJOG	SJOGREN'S ANTIBODIES (SS-A,SS-B) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SASAA	SJOGREN'S ANTI-SS-A Place Of Service:LABCORP	
UHEALTH Labs	SASAB	SJOGREN'S ANTI-SS-B Place Of Service:LABCORP	
UHEALTH Labs	SMADTL	SMA Diagnostic Test Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SMITH	SMITH ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	SMIRNP	SMITH/RNP ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	NA	SODIUM Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84295
UHEALTH Labs	SLLVRG	SOLUBLE LIVER AG (IGG AB) Place Of Service:LABCORP	
UHEALTH Labs	STRB	SOLUBLE TRANSFERRIN RECEPTOR Place Of Service:LABCORP	
UHEALTH Labs	USPG	SPECIFIC GRAVITY Container:URINE random collection / tiger top UA container with preservative Temperature:Room temperature Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	LSMA	Spinal Muscular Atrophy (SMA) Place Of Service:LABCORP	
UHEALTH Labs	SSIGG	SSDNA IGG ANTIBODY Container:Plastic screw-cap vial Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	STCIGA	Stachybotrys chartarum/atra (RGm24) IgA Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	STCIGE	Stachybotrys chartarum/atra (RGm24) IgE Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	STCIGG	STACHYBOTRYS CHARTARUM/ATRA (RGM24) IGG Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	STONEB	STONE ANALYSIS Place Of Service:LABCORP	
UHEALTH Labs	NKSTN	STONE ANALYSIS Container:STERILE SCREW-CAP CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	SRDIAG	STONERISK(R) DIAGNOSTIC PROFILE Container:URINE COLLECTION CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	STOCSC	STOOL CULTURE WITH SHIGATOXIN/CAMPY SCREEN Container:STERILE CUP, CULTURETTE SWAB, CARY BLAIR TRANSPORT MEDIA Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87045
UHEALTH Labs	SJCVA	STRATIFY JCV(TM) AB (INDEX,RFL) Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	STRAWB	STRAWBERRY (F44) IgE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	STCPNE	STREPTOCOCCUS PNEUMONIAE AG, URN/CSF Place Of Service:LABCORP	
UHEALTH Labs	STRIGG	STRONGYLOIDES IGG ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	MYCOTS	SUCEPTIBILITY, MYCOBACTERIUM TUBERCULOSIS Container:Double-walled container with biohazard label Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	SULFU	Sulfonylurea Analysis, Urine Place Of Service:LabCorp	
UHEALTH Labs	SULFQT	SULFONYLUREA SCREEN, QT Place Of Service:LABCORP	
UHEALTH Labs	SWNDC	SURFACE WOUND CULTURE WITH GRAM STAIN Container:STERILE CONTAINER Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	SURGC	SURGICAL SWAB CULTURE WITH GS Container:STERILE CUP OR CULTURETTE Temperature:STORE AT RT DO NOT REFRIGERATE Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	SWORDF	Swordfish IgE Ab, Serum Place Of Service:LabCorp	
UHEALTH Labs	SLEPA	Systemic Lupus Profile A Place Of Service:LabCorp	
UHEALTH Labs	TBLYMP	T + B-LYMPHOCYTE DIFFERENTIAL Place Of Service:LABCORP	
UHEALTH Labs	TBCTOT	T AND B CELLS, TOTAL Container:EDTA (LAVENDER-TOP) TUB Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TCELLF	T Cell Fx-Complete Place Of Service:LABCORP	
UHEALTH Labs	TCRGGR	T CELL RECEPT TCR GAMMA RERNG,PCR Container:EDTA (lavender-top) tube or an EDTA (lavender-top) tube Temperature:Room temperature Place Of Service:QUEST	81342 84999
UHEALTH Labs	TPSCD	T PALLIDUM SCREENING CASCADE Place Of Service:LABCORP	
UHEALTH Labs	TUPT	T UPTAKE Container:SST Temperature:Room Temp Place Of Service:UHealth Labs	
UHEALTH Labs	FTAGM	T. PALLIDUM IGG, IGM AB PANEL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	HZLNTR	T004-IGE HAZELNUT TREE Place Of Service:LABCORP	
UHEALTH Labs	WLNTE	T010-IGE WALNUT Place Of Service:LABCORP	
UHEALTH Labs	WPINE	T-016 IgE Pine, White Place Of Service:LABCORP	
UHEALTH Labs	EUCLTS	T018-IGE EUCALYPTUS Place Of Service:LABCORP	
UHEALTH Labs	PCNHIE	T022-IGE PECAN, HICKORY Place Of Service:LABCORP	
UHEALTH Labs	T3	T3 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84480
UHEALTH Labs	T3UP	T3 UPTAKE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	T3UPB	T3 UPTAKE Place Of Service:LABCORP	
UHEALTH Labs	FRT3	T3, FREE, NON-DIALYSIS Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	REVT3	T3, REVERSE Temperature:REFRIGERATED Place Of Service:QUEST	
UHEALTH Labs	T4	T4 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84436
UHEALTH Labs	LCT4AB	T4 (THYROXINE) ANTIBODY Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	T4FRD	T4 Free, Direct Dialysis Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	T4FDD	T4, FREE, DIR DIAL AND T4, TOTAL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	YTACRO	TACROLIMUS (FK506) Container:EDTA (Lavender-Top) Temperature:ROOM TEMPERATURE PREFERRED Place Of Service:IML	
UHEALTH Labs	FK506	TACROLIMUS (FK506), BLOOD Place Of Service:LABCORP	
UHEALTH Labs	TACRHS	TACROLIMUS, HIGHLY SENSITIVE Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TFVA	Targeted Familial Variant Analysis Container:2-5 mL EDTA Whole Blood or Bone Marrow (lavendar top), 3-4 Buccal Swabs, 1 ug DNA at 50 ng/ul, 7-15 mL Amniotic Fluid, 5-10 mg CVS,2-T25 flasks Cultured Amniocytes or CVS (2x10 <sup>6</sup> minimum) Temperature:Room Temperature Place Of Service:Versiti	81479
UHEALTH Labs	TAYSAC	TAY-SACHS, BIOCHEMICAL, SERUM Place Of Service:LABCORP	
UHEALTH Labs	TCCAB	T-CELL CLONALITY ASSESSMENT (BETA) Container:EDTA (purple-top) tube Temperature:Use refrigerated cold pack for transport Place Of Service:GENOPTIX	81340
UHEALTH Labs	TCCAG	T-CELL CLONALITY ASSESSMENT (GAMMA) Container:EDTA (purple-top) tube Temperature:Use refrigerated cold pack for transport Place Of Service:GENOPTIX	81342
UHEALTH Labs	TCRBGC	TCRB GENE CLONALITY Place Of Service:LABCORP	
UHEALTH Labs	TEICHA	TEICHOIC ACID AB SCR W/RFL TITER Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	TSTFTB	TESTOSTERON, F/WKLYBD+T LC/MS Place Of Service:LABCORP	
UHEALTH Labs	TESTO	TESTOSTERONE LEVEL, TOTAL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84403
UHEALTH Labs	TESTFQ	TESTOSTERONE, FREE & TOTAL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TSTFR	TESTOSTERONE, FREE, LC/MS/MS Container:A RED-TOP TUBE (NO GEL) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TSTFT	TESTOSTERONE, FREE+TOTAL LC/MS Place Of Service:LABCORP	
UHEALTH Labs	URTEST	TESTOSTERONE, URINE Container:PLASTIC SCREW-CAP VIAL Temperature:SHIP FROZEN; SHIPPING ROOM TEMP UNACCEPTABLE; SHIPPING REFRIGERATED UNACCEPTABLE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	TETOID	TETANUS ANTITOXOID IGG AB Place Of Service:LABCORP	
UHEALTH Labs	THYABH	TgAb+Thyroglobulin,IMA or LCMS Place Of Service:LABCORP	
UHEALTH Labs	THYABR	TGAB+THYROGLOBULIN,IMA OR RIA Place Of Service:LABCORP	
UHEALTH Labs	THALM	THALLIUM Place Of Service:LABCORP	
UHEALTH Labs	THEO	THEOPHYLINE LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80198
UHEALTH Labs	LCG	THINPREP PAP REFLEX HPV MRNA E6/E7 Container:ThinPrep® vial Temperature:Room Temperature Place Of Service:UHEALTH LABS	88175
UHEALTH Labs	THIOSR	THIOCYANATE, SERUM Place Of Service:LABCORP	
UHEALTH Labs	THMETA	THIOPURINE METABOLITES Place Of Service:LABCORP	
UHEALTH Labs	THMTB	THIOPURINE METHYLTRANSFERASE Place Of Service:LABCORP	
UHEALTH Labs	THIOR	THIORIDAZINE AND METABOLITE Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	THRC	THROAT CULTURE SCREEN Container:CULTURETTE SWAB Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87070
UHEALTH Labs	TT	THROMBIN TIME Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT FROZEN Place Of Service:Uhealth Labs	85670
UHEALTH Labs	THROMB	Thrombopoietin Container:Gold top Temperature:frozen Place Of Service:CINCINNATI	83520
UHEALTH Labs	THYAB	THYROGLOBULIN ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	TLCMSB	THYROGLOBULIN BY LCMS Place Of Service:LABCORP	
UHEALTH Labs	THYROG	THYROGLOBULIN PANEL Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	THYLNA	THYROGLOBULIN, LN ASPIRATE Place Of Service:LABCORP	
UHEALTH Labs	THYRAB	THYROID ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	THYCAN	Thyroid Cancer (Thyroglobulin) Monitor Container:Plastic screw-cap vial Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	THYCMP	THYROID CANCER MUTATION PANEL (BRAF, RAS, RET/PTC, PAX8/PPAR) Container:Fixative, slides or paraffin embedded tissue block Temperature:Room temperature Place Of Service:QUEST	



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	FTHY	THYROID PANEL Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	THYPNL	THYROID PANEL (NO TSH) Container:Serum Separator Tube (SST®) Temperature:Room temperature Place Of Service:UHEALTH LABS	84436
UHEALTH Labs	THYT	THYROID PANEL WITH TSH Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TPOAB	THYROID PEROXIDASE (TPO) AB Place Of Service:LABCORP	
UHEALTH Labs	THYROI	THYROID STIM IMMUNOGLOBULIN Place Of Service:LABCORP	
UHEALTH Labs	TDFNA	Thyroid, FNA Cytomorphology with Molecular Reflex Container:Plastic leak-proof container and slide container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TYROAB	THYROTROPIN RECEPTOR AB, SERUM Place Of Service:LABCORP	
UHEALTH Labs	THYBGL	THYROXINE BINDING GLOBULIN Place Of Service:LABCORP	
UHEALTH Labs	THYESE	Thyroxine-Binding Protein Electrophoresis, Serum Container:Plastic screw-cap vial Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	TIAGA	Tiagabine (Gabitril), Serum or Plasma Place Of Service:LabCorp	
UHEALTH Labs	TILIGE	Tilapia, IgE (F414) Place Of Service:LABCORP	
UHEALTH Labs	TIMGR	Timothy Grass IgE (G006) Place Of Service:LabCorp	
UHEALTH Labs	TTGIGG	TISSUE TRANSGLUTAMINASE AB, IGG Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TLYHS	T-LYMPHOCYTE HELPER/SUPPRESSOR Place Of Service:LABCORP	
UHEALTH Labs	TOBRAP	TOBRAMYCIN PEAK LEVEL Container:Red top Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	TOBRA	TOBRAMYCIN RANDOM LEVEL Container:Red top Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	TOBRAT	TOBRAMYCIN TROUGH LEVEL Container:Red top Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	TOPIR	TOPIRAMATE (TOPAMAX), SERUM Place Of Service:LABCORP	
UHEALTH Labs	TBIL	TOTAL BILIRUBIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:PROTECTED FROM LIGHT, ROOM TEMPERATURE. Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	CO2	TOTAL CO2 CONTENT Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82374
UHEALTH Labs	TP	TOTAL PROTEIN Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84155
UHEALTH Labs	TOXUR	TOXASSURE SELECT 13 (MW) Place Of Service:LABCORP	
UHEALTH Labs	TXIGG	Toxocara Ab, IgG, S Place Of Service:LABCORP	
UHEALTH Labs	TXOGM	TOXOPLASMA ABS IGG/IGM Place Of Service:LABCORP	
UHEALTH Labs	TOXGM	Toxoplasma Abs, IgG/IgM, CSF Place Of Service:LabCorp	
UHEALTH Labs	TOXIGG	TOXOPLASMA GONDII AB, IGG, QN Place Of Service:LABCORP	
UHEALTH Labs	TXGON	TOXOPLASMA GONDII AB, IGM, QN Place Of Service:LABCORP	
UHEALTH Labs	TOXGON	Toxoplasma gondii DNA, Quantitative Real-Time PCR Container:Plastic leak-proof container Temperature:Refrigerated (cold packs) Place Of Service:Quest	87799
UHEALTH Labs	TXOPCR	TOXOPLASMA GONDII PCR Place Of Service:LABCORP	
UHEALTH Labs	TOXDNA	TOXOPLASMA GONDII,DNA,QUAL,PCR Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	TXPIGA	TOXOPLASMOSIS ADULT PANEL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TPNUV	TPMT and NUDT15 Genotyping, Varies Container:Lavender top (EDTA) Place Of Service:MAYO	0034U
UHEALTH Labs	TRANSF	TRANSFERRIN Container:SERUM SEPARATOR TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TRNSFE	TRANSFERRIN Place Of Service:LABCORP	
UHEALTH Labs	TRAZ	TRAZODONE Container:RED-TOP TUBE (NO GEL) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TRPPAL	TREPONEMA PALLIDUM ANTIBODIES Place Of Service:LABCORP	
UHEALTH Labs	TRPDNA	TREPONEMA PALLIDUM DNA, QUALITATIVE REAL-TIME PCR Container:EDTA or ACD Temperature:Transport Temperature Refrigerated (cold packs) Place Of Service:UHealth Labs	87798
UHEALTH Labs	TCHVAG	TRICH VAG BY NAA Place Of Service:LABCORP	
UHEALTH Labs	TRLAB	TRICHINELLA IgG ANTIBODY, ELISA Temperature:ROOM TEMPERATURE Place Of Service:QUEST	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	THVRNA	Trichomonas vaginalis RNA, Qualitative, TMA, Males Container:APTIMA® Unisex Swab Specimen Collection kit Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	PMTCA	TRICYCLIC ANTIDEPRESS QN medMATCH Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TRIG	TRIGLYCERIDE Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84478
UHEALTH Labs	TROWHI	Tropheryma whipplei (Whipple's) DNA, QL, RT PCR Container:EDTA (lavender-top) or ACD (yellow-top) tube Temperature:Refrigerated (cold packs) Place Of Service:Quest	87798
UHEALTH Labs	TROPT	TROPONIN-T Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84484
UHEALTH Labs	TROUTF	Trout IgE Ab, Serum Place Of Service:LabCorp	
UHEALTH Labs	TRYIGG	Trypanosoma Cruzi Ab IgG Place Of Service:LABCORP	
UHEALTH Labs	TRYPCT	TRYPANOSOMA CRUZI AB, TOTAL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TYPSIN	TRYPSIN Place Of Service:LABCORP	
UHEALTH Labs	TRYPTB	TRYPTASE Place Of Service:LABCORP	
UHEALTH Labs	TSH	TSH Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84443
UHEALTH Labs	TSHHAM	TSH with HAMA Treatment Container:Red-top tube (no gel), Transport container Temperature:Room temperature Place Of Service:QUEST	84443 x2
UHEALTH Labs	THYFS	TSH WITH REFLEX TO FREET4 Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84439
UHEALTH Labs	TTGAG	TTG IGA/G Place Of Service:LABCORP	
UHEALTH Labs	TTGIGA	T-TRANSGLUTAMINASE (TTG) IGA Place Of Service:LABCORP	
UHEALTH Labs	TTIGGG	T-TRANSGLUTAMINASE (TTG) IGG Place Of Service:LABCORP	
UHEALTH Labs	TNECFA	TUMOR NECROSIS FACTOR-ALPHA Place Of Service:LABCORP	
UHEALTH Labs	TUNA	TUNA (F40) IgE Container:Transport Container: Serum Separator Tube (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	TURKF	Turkey Feathers IgE (E089) Container:Red-top tube or gel-barrier tube Temperature:Room temperature Place Of Service:LABCORP	86003

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	UGT1A1	UGT1A1 Irinotecan Toxicity Container:Lavender-top (EDTA) tube, yellow-top (ACD) tube, or LC buccal swab kit Temperature:Room Temp or Refrig Place Of Service:LABCORP	81350
UHEALTH Labs	URMYHO	UREAPLASMA/MYCOPLASMA HOMINIS Place Of Service:LABCORP	
UHEALTH Labs	URIC	URIC ACID Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84550
UHEALTH Labs	URICR	Uric Acid, 24-Hour Urine with Creatinine Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	82570
UHEALTH Labs	UABFLD	URIC ACID, BODY FLUID Place Of Service:LABCORP	
UHEALTH Labs	URICRR	URIC ACID, RANDOM URINE (W/CREATININE) Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84560
UHEALTH Labs	UASF	URIC ACID, SYNOVIAL FLUID Container:A STERILE PLASTIC TRANSPORT TUBE. Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	UMACUC	URINALYSIS WITH REFLEX TO URINE CULTURE Container:URINE random collection / tiger top UA container with preservative/ GRAY Top Tube urine with preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	81001
UHEALTH Labs	UAMYR	URINE AMYLASE CONC.,RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82150
UHEALTH Labs	UAMY	URINE AMYLASE, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	82150
UHEALTH Labs	UCAR	URINE CALCIUM CONC., RANDOM Container:URINE CUP / Urine yellow top with no preservative Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	RUCA	URINE CALCIUM WITH CREATININE, RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82340 82570
UHEALTH Labs	UCA	URINE CALCIUM, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	82340
UHEALTH Labs	UCLR	URINE CHLORIDE CONC., RANDOM Container:Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82436
UHEALTH Labs	RUCL	URINE CHLORIDE CONC., RANDOM Container:Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82436
UHEALTH Labs	UCL	URINE CHLORIDE, 24HR Container:JUG- Plastic urine container Temperature:Refrigerated Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	UCRCL	URINE CREATININE CLEARANCE, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82575
UHEALTH Labs	UCRR	URINE CREATININE CONC., RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82570
UHEALTH Labs	UCREAT	URINE CREATININE, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	82570
UHEALTH Labs	UC	URINE CULTURE Container:STERILE CONTAINER OR BORIC ACID GRAY TOP TUBE Temperature:TRANSPORT AT ROOM TEMP. DO NOT REFRIGERATE. Place Of Service:UHealth Labs	87086
UHEALTH Labs	URNSCR	URINE DRUG SCREEN 1 Container:Plastic leak-proof container Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ULYTES	URINE ELECTROLYTES, 24HR Container:JUG- Plastic urine container Temperature:Refrigerated Place Of Service:UHealth Labs	
UHEALTH Labs	ULYTER	URINE ELECTROLYTES, RANDOM Container:URINE CUP / Urine yellow top with no preservative Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	UEOCT	URINE EOSINOPHIL SMEAR Container:Random urine cup Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	UGLUR	URINE GLUCOSE CONC., RANDOM Container:URINE CUP / Urine yellow top with no preservative Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	UGLU	URINE GLUCOSE, 24HR Container:JUG- Plastic urine container Temperature:Refrigerated Place Of Service:UHealth Labs	
UHEALTH Labs	UPREG	URINE HCG, QUALITATIVE Container:URINE CUP / Urine yellow top with no preservative Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	UIMQ	URINE IMMUNOGLOBULINS Container:sterile urinalysis container then transfer to sterile, screw top tube for shipment Temperature:Refrigerate Place Of Service:UHealth Labs	
UHEALTH Labs	UMAC	URINE MACROSCOPIC Container:Urine Cup Urine Transport Tube Temperature:Urine cup: Room temperature up to 2 hours and refrigerated up to 48 hours. Urine transport tube: Refrigerated or room temperature for up to 72 hours. Place Of Service:UHealth Labs	81003
UHEALTH Labs	UMAGR	URINE MAGNESIUM CONC., RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	83735
UHEALTH Labs	UMAG	URINE MAGNESIUM, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	83735

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	UMAR24	URINE MICROALBUMIN LEVEL 24HRS Container:JUG - Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82232
UHEALTH Labs	UMAR	URINE MICROALBUMIN, RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82232
UHEALTH Labs	UMIC	URINE MICROSCOPIC Container:URINE CUP / Urine tiger top with preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	81001
UHEALTH Labs	UOSMR	URINE OSMOLALITY, RANDOM Container:Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83935
UHEALTH Labs	UPH	URINE pH Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83986
UHEALTH Labs	UPHOSR	URINE PHOSPHORUS CONC., RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84105
UHEALTH Labs	UPHOS	URINE PHOSPHORUS, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	84105
UHEALTH Labs	UKR	URINE POTASSIUM CONC., RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84133
UHEALTH Labs	RUK	URINE POTASSIUM CONC., RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84133
UHEALTH Labs	UK	URINE POTASSIUM, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	84133
UHEALTH Labs	UPCONC	URINE PROTEIN CONC. Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84156
UHEALTH Labs	RUPRO	URINE PROTEIN CONC.,RANDOM Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84156
UHEALTH Labs	CREPEU	Urine Protein Electrophoresis and Creatinine, (24 Hour) Container:12-mL plastic transport tube (LabCorp N° 23597) without preservative Temperature:Refrigerate at 2°C to 8°C. Place Of Service:LABCORP	
UHEALTH Labs	CRPEUU	Urine Protein Electrophoresis and Creatinine, Random Container:12-mL plastic transport tube (LabCorp N° 23597) without preservative Temperature:Refrigerate at 2°C to 8°C. Place Of Service:LABCORP	
UHEALTH Labs	UELP	URINE PROTEIN ELP Place Of Service:UHealth Labs	
UHEALTH Labs	UNAR	URINE SODIUM CONC., RANDOM Container:JUG- Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84300



# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	RUNA	URINE SODIUM CONC., RANDOM Container:Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84300
UHEALTH Labs	UNA	URINE SODIUM, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	84300
UHEALTH Labs	UTP	URINE TOTAL PROTEIN, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84156
UHEALTH Labs	UUNR	URINE UREA NITROGEN CONC.,RANDOM Container:URINE CUP / Urine yellow top with no preservative Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	UUREA	URINE UREA NITROGEN, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT REFRIGERATED Place Of Service:UHealth Labs	84540
UHEALTH Labs	UURICR	URINE URIC ACID CONC.,RANDOM Container:URINE CUP / Urine yellow top with no preservative Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	UURIC	URINE URIC ACID, 24HR Container:JUG- Plastic urine container Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	84560
UHEALTH Labs	UMACR	URN MICROALBUMIN/CREATININE, RND Container:URINE CUP / Urine Yellow Top Tube with no preservative Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82043
UHEALTH Labs	UPGD	Uroporphyrinogen Decarboxylase (UPG D), Whole Blood Container:Preferred: Green top (sodium heparin). Preferred: Green top (sodium heparin) Temperature:Immediately place specimen on wet ice. Place Of Service:Mayo	82657
UHEALTH Labs	URORSK	URO-RISK® DIAGNOSTIC PROFILE Container:See Instructions Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	VAGWM	VAGINAL WET MOUNT Container:STERILE TUBE OR AIMS TRANSPORT MEDIA Temperature:STORE AT RT. Place Of Service:UHealth Labs	
UHEALTH Labs	VAGHSV	Vaginitis Plus (VG+), Herpes Simplex Virus (HSV), NuSwab Place Of Service:LABCORP	
UHEALTH Labs	VALAC	VALPROIC ACID (TOTAL+FREE) Place Of Service:LABCORP	
UHEALTH Labs	DEPAK	VALPROIC ACID LEVEL Container:RED TOP TUBE (NO SERUM SEPARATOR PREFERED), GoldTop / Tiger top (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80164
UHEALTH Labs	VANCOP	VANCOMYCIN PEAK LEVEL Container:RED TOP TUBE Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80202
UHEALTH Labs	VANCO	VANCOMYCIN RANDOM LEVEL Container:RED TOP TUBE Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80202

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	VANCOT	VANCOMYCIN TROUGH LEVEL Container:RED TOP TUBE Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	80202
UHEALTH Labs	VAN24	VANILLYLMADELIC ACID, 24-HR U Place Of Service:LABCORP	
UHEALTH Labs	VMA24	VANILLYLMADELIC ACID,24 HR UR Container:24-HOUR URINE CONTAINER Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	VZABS	VARICELLA ZOSTER ABS, IGG/IGM Container:Red-top tube or gel-barrier tube Temperature:Room Temperature or Refrigerated Place Of Service:LABCORP	
UHEALTH Labs	VZVPCR	VARICELLA ZOSTER DNA, PCR QUAL Container:EDTA (lavender-top) tube Temperature:REFRIGERATED (COLD PACKS) Place Of Service:QUEST	
UHEALTH Labs	VZGQ	VARICELLA ZOSTER VIRUS AB (IGG) Container:Transport tube (Gold tube) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	CSFVZ	VARICELLA ZOSTER VIRUS AB, CSF Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	VARZOS	VARICELLA-ZOSTER AB, IGM Place Of Service:LABCORP	
UHEALTH Labs	VZIGG	VARICELLA-ZOSTER V AB, IGG Place Of Service:LABCORP	
UHEALTH Labs	VZVABT	Varicella-Zoster Virus (VZV) Antibodies (Total, IgM), ACIF/ Container:sterile, leak-proof container Temperature:Room temperature Place Of Service:QUEST	
UHEALTH Labs	VEGFP	Vascular Endothelial Growth Factor, Plasma Place Of Service:LabCorp	
UHEALTH Labs	VDRLCF	VDRL, CSF Place Of Service:LABCORP	
UHEALTH Labs	VDRLS	VDRL, SERUM Container:Transport tube Temperature:Refrigerated Place Of Service:QUEST	
UHEALTH Labs	VEGFSB	VEGF, SERUM Place Of Service:LABCORP	
UHEALTH Labs	EFFXR	VENLAFAXINE (EFFEXOR) Place Of Service:LABCORP	
UHEALTH Labs	VRAPAM	VERAPAMIL, SERUM/PLASMA Place Of Service:LABCORP	
UHEALTH Labs	VGCCAB	VGCC ANTIBODY Place Of Service:LABCORP	
UHEALTH Labs	VIBC	VIBRIO CULTURE,STOOL Container:STERILE CUP, CULTURETTE SWAB, CARY BLAIR TRANSPORT MEDIA Temperature:STORE AT RT, DO NOT REFRIGERATE. Place Of Service:UHealth Labs	
UHEALTH Labs	VIPPL	VIP, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	VIRCUL	VIRAL CULTURE, GENERAL Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	VRCLRL	VIRAL CULTURE, RAPID, LESION Place Of Service:LABCORP	
UHEALTH Labs	VHEPS3	Viral Hepatitis Screening and Diagnosis (HAV and HBV HCV) Place Of Service:LABCORP	
UHEALTH Labs	VHEPSD	Viral Hepatitis Screening and Diagnosis (HAV and HBV) Place Of Service:LABCORP	
UHEALTH Labs	VISCO	VISCOSITY, SERUM Place Of Service:LABCORP	
UHEALTH Labs	VB12B	VIT B12 BINDING CAPACITY, UNSAT Container:PLASTIC SCREW-CAP VIAL Temperature:REFRIGERATED (COLD-PACKS) Place Of Service:QUEST	
UHEALTH Labs	VITAS	VITAMIN A, SERUM Place Of Service:LABCORP	
UHEALTH Labs	VITB1B	VITAMIN B1 (THIAMINE), BLOOD Place Of Service:LABCORP	
UHEALTH Labs	VITB1	VITAMIN B1 LEVEL (THIAMINE), BLOOD Container:PLASTIC SCREW-CAP VIAL PROTECTED FROM LIGHT Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	VB12DC	Vitamin B12 Deficiency Panel Place Of Service:LABCORP	
UHEALTH Labs	B12	VITAMIN B12 LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82607
UHEALTH Labs	B12FOL	VITAMIN B12/FOLATE,SERUM PANEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82607
UHEALTH Labs	VITB2	VITAMIN B2 LEVEL Container:PLASTIC SCREW-CAP VIAL PROTECTED FROM LIGHT Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	VB2WB	VITAMIN B2, WHOLE BLOOD Place Of Service:LABCORP	
UHEALTH Labs	VTB3	VITAMIN B3 (NIACIN+METABOLITE) Place Of Service:LABCORP	
UHEALTH Labs	VITB6P	VITAMIN B6, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	VITEC	VITAMIN C Place Of Service:LABCORP	
UHEALTH Labs	VTD25	VITAMIN D 25-HYDROXY LEVEL Container:GOLD TOP/ TIGER TOP (Serum Separator Tube) (SST®) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	82306
UHEALTH Labs	VITD	VITAMIN D, 1,25-DIHYDROXY LEVEL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	VITD25	VITAMIN D, 25-HYDROXY LEVEL Container:PLASTIC SCREW-CAP VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	VITES	VITAMIN E, SERUM Place Of Service:LABCORP	
UHEALTH Labs	VITK1	VITAMIN K1 Place Of Service:LABCORP	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	VMAUR	VMA, RANDOM URINE Place Of Service:LABCORP	
UHEALTH Labs	VOLA	Volatiles, Whole Blood Place Of Service:LabCorp	
UHEALTH Labs	VONWCP	von Willebrand Comprehensive Panel Container:Light Blue Tube Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	VWD2N	von Willebrand Disease Type 2N(Factor VIII Binding Activity) Container:3.2% sodium citrate (light blue-top) tube Temperature:Frozen Place Of Service:QUEST	83520
UHEALTH Labs	VWFAG	VON WILLEBRAND FACTOR ANTIGEN Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85670
UHEALTH Labs	VORIP	VORICONAZOLE LC-MS/MS, PLASMA Place Of Service:LABCORP	
UHEALTH Labs	VORIS	VORICONAZOLE LC-MS/MS, SERUM Place Of Service:LABCORP	
UHEALTH Labs	VORCZ	VORICONAZOLE, HPLC Container:PLASTIC SCREW-CAP VIAL Temperature:FROZEN Place Of Service:QUEST	
UHEALTH Labs	VRE	VRE CULTURE Container:CULTURETTE SWAB Temperature:STORE AT RT Place Of Service:UHealth Labs	
UHEALTH Labs	VWFMUL	VWF MULTIMER PLASMA Container:BLUE TOP TUBE (SODIUM CITRATE) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	83520
UHEALTH Labs	VZVRT	VZV REAL TIME PCR Place Of Service:LABCORP	
UHEALTH Labs	PLNTAN	W009-IGE PLANTAIN, ENGLISH Place Of Service:LABCORP	
UHEALTH Labs	LAMQRT	W010-IgE Lamb's Quarters Place Of Service:LABCORP	
UHEALTH Labs	FENDOG	W046-IgE Fennel, Dog Place Of Service:LABCORP	
UHEALTH Labs	WALIGE	WALNUT (F256) IGE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	WALNUT	WALNUT TREE (T10) IGE Container:SERUM SEPARATOR TUBE (SST) Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	WARSEN	Warfarin Sensitivity (CYP2C8, CYP2C9, CYP4F2, VKORC1) Genotyping Container:Whole Blood: Lavender (EDTA), Pink (K2EDTA), or Yellow (ACD Solution A or B). Saliva: Collection Device by DNA Genotek (OCD-100, ARUP Supply #49295) available online through eSupply using ARUP Connect™ or by contacting ARUP Client Services at (800) 522-2787. Temperature:Whole Blood: Refrigerated. Saliva: Room temperature. Place Of Service:ARUP	81227 81355
UHEALTH Labs	WARF	Warfarin, Serum/Plasma Place Of Service:LABCORP	
UHEALTH Labs	WATERC	WATER CULTURE Container:STERILE CUP Temperature:REFRIGERATE Place Of Service:UHealth Labs	

# Online Test Menu

## LABORATORY CATALOG

LAB	CODE	TEST NAME	CPT4 CODE
UHEALTH Labs	WBC	WBC COUNT Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85048
UHEALTH Labs	DIFF	WBC DIFFERENTIAL Container:Lavender Temperature:Room temperature Place Of Service:UHealth Labs	
UHEALTH Labs	WNV	WEST NILE VIRUS AB Place Of Service:DOHM	
UHEALTH Labs	WNVCSF	WEST NILE VIRUS ANTIBODY, CSF Place Of Service:LABCORP	
UHEALTH Labs	WNVABS	WEST NILE VIRUS ANTIBODY, SERUM Place Of Service:LABCORP	
UHEALTH Labs	WHTBAL	WHITE BEAN (F15) IGE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	WBCDI	WHITE CELL COUNT + DIFFERENTIAL Container:LAVENDER TOP TUBE (EDTA) Temperature:TRANSPORT AT ROOM TEMPERATURE Place Of Service:UHealth Labs	85004
UHEALTH Labs	WASP0	Wiskott-Aldrich Syndrome Protein (WASP) Container:Sodium Heparin (DARK Green) Temperature:Room Temperature Place Of Service:CINCINNATI	88184 88185 x2 88187
UHEALTH Labs	DNAFRG	XSENSE(R), FRAGILE X W/REFLEX Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	YERSC	YERSINIA CULTURE, STOOL Container:STERILE CUP, CULTURETTE SWAB, CARY BLAIR TRANSPORT MEDIA Temperature:STORE AT RT, DO NOT REFRIGERATE. Place Of Service:UHealth Labs	
UHEALTH Labs	ZN	ZINC Container:ACID-WASHED OR METAL-FREE VIAL Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	RBCPRO	ZINC PROTOPORPHYRIN Container:EDTA (LAVENDER-TOP) TUBE Temperature:ROOM TEMPERATURE Place Of Service:QUEST	
UHEALTH Labs	ZINC	ZINC, PLASMA OR SERUM Place Of Service:LABCORP	
UHEALTH Labs	ZNRBC	Zinc, RBC Place Of Service:LABCORP	
UHEALTH Labs	ZONIB	ZONISAMIDE(ZONEGRAN), SERUM Place Of Service:LABCORP	
LABCORP	789120	789120 10+Alc-Bund Container:Urine container(1) Temperature:Refrigerate	80307