



LAB	CODE	TEST NAME	CPT4 CODE
ACL	MYCSER	M PNEUMONIAE AB IgA Container :One plain red 4.0 MI (ALSO ACCEPTABLE: GOLD GEL 3.5 ML) Stability :Frozen: 120 Days	86738
ACL	MACPRO	Macroprolactin Container :One gold gel 5.0 mL (also acceptable: one green lithium heparin 6.0 mL or one lavender (EDTA) 3.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 7 Days / Frozen: 1 Month	84146 x2
ACL	MAGRBC	Magnesium RBC Container :One royal blue (EDTA) 6.0 mL Stability :Ambient: 2 Weeks / Refrigerated: 2 Weeks / Frozen: Unacceptable	83735
PHC	MG	Magnesium, Plasma or Serum Stability :Room temperature: 7 days Refrigerated: 7 days Frozen: 1 year	83735
ACL	UMGT	Magnesium, Urine 24 Hour Container :One 24 hour urine Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	83735
ACL	UMGR	Magnesium, Urine Random Container :One random urine Stability :Ambient: Unacceptable / Refrigerated: 7 Days / Frozen: 14 Days	83735
ACL	RMALP	MALARIA RAPID SCREEN Container :One lavender (EDTA) 3.0 mL blood Stability :Ambient: 24 Hours / Refrigerated: 24 Hours / Frozen: Unacceptable	
PHC	VR01713	Measles Virus PCR Container :UTM Stability :Refrigerated: 48 hours	87798
ACL	DRGSCM	Meconium Drug Screen Container :All meconium (blackish material) excreted until milk/formula-based stool (yellow-green) appears. Stability :Ambient: 1 Week / Refrigerated: 3 Weeks / Frozen: 1 Year	80326 80347 80364 80355 80323
ACL	MSHA	Melanocyte Stimulating Hormone Container :One lavender (EDTA) 3.0 mL Stability :Ambient: Unacceptable / Refrigerated: 24 Hours / Frozen: 1 Month	83519
ACL	MELAT	Melatonin Container :One plain red 6.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: 2 Weeks / Refrigerated: 1 Month / Frozen: 1 Year	80375
ACL	MERCB	MERCURY, BLOOD Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	83825
ACL	URHMS	METALS PNL, UR/OCCU Stability :Refrigerated: 5 Days / Frozen: 14 Days Stability based upon creatinine requirements.	82175 82300 82570 83655 83825
ACL	HMSB	METALS PNL1, BLOOD Stability :Ambient: 14 Days / Refrigerated: 14 Days (Samples older than 14 days may be accepted at the discretion of the department supervisor.)	82175 82300 83655 83825
ACL	PMET	Metanephrines, Plasma Container :One lavender (EDTA) 6.0 mL (also acceptable: One pink K2EDTA and Green Sodium or Lithium heparin, NO GEL) Stability :Ambient: Unacceptable / Refrigerated: Unacceptable / Frozen: 1 Month	83835