



LAB	CODE	TEST NAME	CPT4 CODE
ACL	T3REV	T3 (Triiodothyronine), Reverse Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL or two lavender (EDTA) 3.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 3 Months	84482
ACL	TT3	T3, Blood Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 5 Days / Frozen: 14 Days	84480
ACL	FT3B	T3, Free, Blood Container :One gold gel 3.5 mL (also acceptable: one gold gel microtainer 600 mcL (clear or amber) or one plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	84481
ACL	T3U	T3, Uptake, Blood Container :One gold gel 5.0 mL (also acceptable: one plain red 6.0 mL) Stability :Ambient: 8 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	84479
ACL	T4	T4 (Thyroxine), Blood Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	84436
PHC	FT4	T4 (Thyroxine), Free, Plasma or Serum Stability :Refrigerated: 7 days	84439
ACL	FK506S	Tacrolimus Level, Blood Container :One lavender (EDTA) 3.0 mL Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 14 Days	80197
ACL	TCELLB	TCell T4/T8 Ratio Container :One green (sodium heparin no gel) 4.0 mL (also acceptable: one lavender (EDTA) 3.0 mL) Stability :Ambient: 72 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable	86355 86359 86360 86357
ACL	TSFFR	Testosterone Free, Females or Children Container :One gold gel 3.5 mL (also acceptable: one green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 6 Months	84402
ACL	TSTFPB	Testosterone Free/Total/Bioavailable/SHBG, Adult Male Container :One gold gel 3.5 mL blood (also acceptable: one plain red 4.0 mL) Stability :Ambient: 4 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	84270 84403
ACL	TESTG	Testosterone Total & Free, Males Only Container :One gold gel 3.5 mL (also acceptable: one green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 1 Week / Frozen: 2 Months	84402 84403
ACL	TESTO	Testosterone Total, Males Only Container :One gold gel 3.5 mL (also acceptable: one plain red 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 48 Hours / Frozen: 14 Days	84403
ACL	TSSHFR	Testosterone, Bioavailable & Sex Hormone Binding Globulin (w/ Total), Female &/or Child Container :One gold gel 5.0 mL (also acceptable: one green (sodium or lithium heparin, no gel) 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 6 Months	84402 84403 84270
ACL	TESTFR	Testosterone, Females or Children Container :One gold gel 3.5 mL (also acceptable: one green (sodium or lithium heparin) 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 6 Months	84403
ACL	TSTFFR	Testosterone, Free and Total, Includes Sex Hormone-Binding Globulin, Female or Child Container :One gold gel 3.5 mL (also acceptable: one green (sodium or lithium heparin, no gel) 4.0 mL) Stability :Ambient: 24 Hours / Refrigerated: 7 Days / Frozen: 6 Months	84402 84403 84270