



LAB	CODE	TEST NAME	CPT4 CODE
ACL	YSTC	<p>Yeast Culture</p> <p>Container :Collect One of the Following:</p> <p>One eSwab or aerobic culture swab (acceptable sources: mouth, throat, vagina, or penis)</p> <p>One random urine in a sterile container</p> <p>Stability :eSwab or Aerobic Culture Swab: Ambient: 72 Hours / Refrigerated: Unacceptable / Frozen: Unacceptable</p> <p>Urine in Gray Boric Acid Tube: Ambient: Unacceptable / Refrigerated: 48 Hours / Frozen: Unacceptable</p> <p>Urine in Sterile Container: Ambient: Unacceptable / Refrigerated: 24 Hours / Frozen: Unacceptable</p>	87102
ACL	YERC	<p>Yersinia Culture, Stool</p> <p>Container :One stool specimen in a clean leak-proof container</p> <p>IMPORTANT: Stool must be transferred to Cary-Blair media within one hour of collection.</p> <p>Stability :Unpreserved Stool: Ambient: 1 Hour / Refrigerated: Unacceptable / Frozen: Unacceptable</p> <p>Stool in Cary-Blair: Ambient: 4 Days / Refrigerated: Unacceptable / Frozen: Unacceptable</p>	87046